


Azenta Life Sciences Sample Consumables and Instruments Catalog



AZENTA
LIFE SCIENCES

azenta.com

Introduction

Unrivaled Sample Solutions to Accelerate Discovery, Development, and Delivery

Azenta Life Sciences provides unrivaled sample exploration and management solutions to help our customers accelerate discovery, development, and delivery to bring impactful breakthroughs and therapies to market faster. We are the global leader in automated compound management for drug discovery, biological storage, and sample processing solutions. Azenta understands the importance of sample integrity and provides a comprehensive range of solutions across our leading capabilities of multiomic services, sample repository services (SRS), consumables and instruments, data management and informatics, and automated ultra-cold storage.

Offer Categories

GENEWIZ Multiomics & Synthesis Solutions

Sequencing


- Next Generation Sequencing
- Sanger Sequencing
- PCR + Sanger Services

DNA & RNA Synthesis

- Gene Synthesis
- Oligo Synthesis
- Gene-to-Virus Services
- Plasmid DNA Preparation

Other Services

- Gene Editing
- Biofluid Processing



Sample Management LIMS

Powerful software tools for sample and product lifecycle management



Consultative Services

Expert guidance and planning for genomics, cold-chain logistics, automated storage, and clinical trial management



Automated Storage & Cryogenics

Automated Storage Systems

- SampleArc™ Ambient to -20°C range
- BioArc™ -80°C range
- CryoArc™ LN2-based range

Clinical Trial Management



- Study Initiation
- Sample Preparation & Lab Services
- Short- & Long-Term Sample Storage
- Sample Support Services

Drug Products & Therapies

- Ultra-Cold Storage
- Cold-Chain Logistics
- Inventory Management

Transportation

- Biological Transport
- Laboratory Moving



Consumables & Instruments

Sample Tube Solutions


- Tri- and Dual-Coded Sample Tubes
- Tube Code Reading Systems
- Tube Capping & Sealing Systems
- Tube Management Systems

PCR & Microplate Solutions

- PCR Plates, Strips & Individual Tubes
- Assay Microplates
- Adhesive & Heat Seals
- Automated Sealing Instrumentation

Sample Cooling And Heating

- Ice-Free Cooling Workstations
- Thermoconductive Tube Racks
- Alcohol Free Cell Freezing Containers



Accelerate Discovery, Development & Delivery

Research & Development

Expedite target identification, hit discovery, and lead optimization with comprehensive solutions for genomics research and high-throughput screening.



Pre-Clinical & Clinical

Streamline sample management and exploration in regulated environments with state-of-the-art automation and unparalleled expertise.



Manufacturing & Distribution

Store and deploy therapies with the utmost security, sophisticated informatics, and full compliance—from the industry leader in cold-chain logistics.



Our Customers

20 of 20

Top pharma/biotech companies are served by Azenta

13 of 15

Top pharmas trust Azenta with their samples

Top 5

Pharma products (by sales) have their clinical samples managed by Azenta

1 in 3

US molecular biologists use Azenta genomic services

20,000

Citations in scientific journals mention Azenta and our brands


33

Nobel laureate labs use Azenta



AZENTA
LIFE SCIENCES

Our History




Over 10 Years of Building Out a Leading-Edge Life Science Company

Since expanding into life sciences in 2011, we have dedicated ourselves to supporting our customers' ambitions. But with advancements in biotech, genomics, and AI, and a rapidly changing landscape, we saw an opportunity to do more.

To help our customers achieve their goals with greater agility and insight, to bring more breakthroughs and therapies to market—sooner. To do so, we brought together leading capabilities across genomics, cryogenic storage, automation, and informatics.

Our Global Reach

Azenta Life Sciences is the global leader in automated compound management for drug discovery, biological storage, and sample processing solutions with 30 locations in seven countries.



AZENTA
LIFE SCIENCES

Contents

Sample Tube Consumables and Instruments

Introduction

Guide to FluidX Sample Tubes	8
Glossary of Terms	8
Introduction	9
Sample Storage Consumables and Devices	10-13
Anatomy of a Tube	14-18
Customization Options	19
1. Sample Storage Tube Range	
Tube Range Table	20-21

Consumables

2. Tri-Coded Sample Storage Tubes	
Introduction to External Thread Tri-Coded Tubes	23
Introduction to Internal Thread Tri-Coded Tubes	29
0.5ml external thread	25
0.8ml external thread	26
1.0ml external thread	27
1.5ml external thread	28
1.6ml external thread	29
1.9ml external thread	30
3.8ml external thread	31
7.6ml external thread	32
0.48ml internal thread	33
0.65ml internal thread	34
0.9ml internal thread	35
3. Dual-Coded Sample Storage Tubes	
Introduction to 96-Format Dual-Coded Tubes	36
0.26ml external thread	37
0.3ml internal thread	38
0.7ml internal thread	39
0.9ml internal thread	40
0.5ml external thread	41
0.9ml external thread	42
4. 2D-Coded Sample Storage Tubes	
External Thread 2D-Coded Tubes	43
Internal Thread 2D-Coded Tubes	43
Acoustic Sample Tube Echo® Qualified Consumable	44
0.2ml external thread	45
5. Non-Coded Sample Storage Tubes	
Non-Coded Tubes	46-47

6. Capping and Sealing Options	
Tube Screw caps	49-51
TPE septum caps	52
SBS cap carriers	53
7. Treatment Services	
Treatment Services	56-57
8. Recommended Temperature Range for Sample Storage Tubes	
Recommended Temperature Range for Sample Storage Tubes	58
How Safe Are Your Samples?: Leachables, Working Volume and Pressure Testing	59

Instruments


9. Ziath Tube and Rack Code Reading Systems	
Uno and Duo Single tube readers	62-63
Mirage Rack Reader	64-65
Express Rack Reader	66-67
Cube Rack Reader	68-69
Mirage RITrack RFID Reader	70-71
Handheld 2 Tube Reader	72
Cyclops 1D Rack Code Scanner	73
DP5 Decoding Software	74
10. Sample Tube Capping & Sealing Systems	
Semi-Automated Screw Cap Single Tube Decapper/Recapper	76-78
Semi-Automated Handheld Screw Cap Decapper, 8-channel	79
Semi-Automated Screw Cap Decapper/Recapper, Single Channel	80-81
IntelliXcap™ Screw Cap Decapper/Recappers	82-86
IntelliXcap Automated Screw Cap Decapper/Recapper 96-format, Extended Height	84
IntelliXcap Automated Septum Cap	85
Semi-Automated Septum Cap Capper	86
11. Sample Tube Management Systems	
Manual Tube Pickers and Manual Decappers	88
Mohawk Semi-Automated Tube Picker	89-90
Direct Tube Marker	91-92
Tube Auditor™	93
FreezerPro® Sample Management System	94
Automated Sample Storage Ambient to -190°C	95-96
Cryogenic Storage and Carriers	97-98

4titude PCR/Microplate Consumables and Instruments

4titude Manufacturing & Quality Standards	100	Heat Sealing Consumables Compatibility Table	206-207
PCR Plates: Clear, Frosted or White		Thermosensitive Colour Forming Film	208
Wells & Low DNA Binding Properties	102-103	Automated Roll Heat Sealers	
FrameStar® 2-Component PCR Plates	104-122	Automated Roll Heat Sealer	209-210
FrameStar Breakable Plates	123-126	Automated Individual Access Heat Sealer	211
Individual Access 96 Well PCR Plates	127-135	Semi-Automated Sheet Heat Sealer	213-214
8 Well PCR Tube Strip with PC Frame	137-140	Automated Plate Seal Remover	215
96 Well PCR Plates with 8 Well Removeable Tube Strips	141-142	Adhesive Sealing Consumables & Accessories	217-231
Standard PCR Plates, Strips & Tubes	143-157	Adhesive Sealing Consumables Comparison Table	232-233
96 Well Non-Skirted PCR Plate		Plate Lids, Caps & Mats	234-241
Horizontally or Vertically Breakable	158-159	Custom Capabilities	
Instrument Compatibility Table	161-165	Custom Capabilities for OEM Diagnostic Kit and	
Storage Microplates	166-175	Medical Device Manufacturing	243
Assay Microplates	176-185	Customizable PCR Plates, Microplates and Seals	244
Heat Sealing Consumables	186-205		

Sample Cooling and Heating

Sample Cooling and Heating Standardization	246	Accessories	
CoolCell Alcohol-Free Cell Freezing Containers	247-243	Ice Pans	270
CoolBox Ice-Free Cooling Workstations	254-260	Hinged CryoBoxes	271
Thermoconductive Tube Racks	261-266	Thermoconductive Tube Module Temperature Strip	272
Thermoconductive Sinks	267-268	Thermoconductive Tube Module Sleeves	272
		Cryo Tube Grippers	272
		Cryo Tube Locking Racks	272
		1D-coded Cryo Tubes	273
		Cryo Tubes Cap Inserts	273



AZENTA
LIFE SCIENCES

FluidX Sample Tube Consumables and Instruments



Guide to FluidX™ from Azenta Sample Tubes

Tri-coded Tubes

Each tube features a permanent 2D-code laser etched in high-contrast on the tube base, a permanent 1D (linear barcode) and Human-Readable Number laser etched in high-contrast on the tube side.

The tube is manufactured using an advanced manufacturing process which results in a one-piece tube with a clear window and black side for coding. All three codes are identical and auditing processes guarantee all three codes match. Our standard 2D option is black on white.

Dual-coded Tubes

Each tube features a permanent 2D-code and Human Readable Number laser etched in high-contrast on the tube base or tube side.

The Dual-coded tube is manufactured using an advanced manufacturing technique used to integrate 2 resin colors into the same tube for high resolution coding.

Both codes are identical and auditing processes guarantee both codes match.

Our standard 2D option is white on black.

2D-coded Tubes

2D-code laser etched on the base only. Allows users to add additional information or a second identifier in high-contrast to the tube side. Our standard 2D option is black on white.

Non-coded Tubes

Simply come as they are.

Glossary of Terms

2D Datamatrix code

A 2D code made of black and white cells, for Azenta tubes these are arranged in a square. The L-shape found on the border is its finder pattern, which is used by tube readers and scanners to recognize and read the code.

1D-code / Linear barcode

Unique barcode represented by parallel lines of different widths and spacings printed on the side of a tube, SBS rack or CryoBox. Our standard code is code 128. Code 128 is able to encode alpha-numeric data.

2D-datamatrix

2D Datamatrix code on tube base.

2D4-coded

Quad Code on Acoustic Sample Tube base.

Dual-coded

2D code and Human Readable Number on tube base or tube side.

Tri-coded

2D-code on tube base, 1D (linear barcode) and Human-Readable Number on tube side.

External Thread

Thread is on the outside of the tube so there is no loss in working volume when a cap is added, saving valuable freezer space. Azenta external thread caps feature a double start thread. The thread is fully engaged after a maximum rotation of 180°, making the cap easier to use especially in automated environments.

Internal Thread

Thread is on the inside of the tube.

Working Volume

The maximum sample volume that will still allow space between the sample and underside of the cap for ice expansion during freezing.

Fill Volume

The total capacity of the tube at 21°C.

All About the Sample!

As a global leader in innovative sample management solutions Azenta Life Sciences is all about the sample.

Working across a wide range of industries Azenta Life Sciences offers unparalleled knowledge and experience of 2D-coded sample storage tubes, readers and sample management systems. As part of the team that developed the original 2D-coded sample tubes in 1999, we have been at the global forefront of developing sample storage consumables and instruments for over 20 years and we continue our philosophy of innovation to this day.

In addition to the innovation behind our products, we are committed to providing the highest levels of customer service, support and quality. Our Technical Support Team provides expert assistance, making sure our products give the reliable and optimal performance you expect.

We believe that a quality sample is the cornerstone in the generation of reliable, reproducible and quantifiable data – which is why Azenta Life Sciences is all about the sample.



FluidX Sample Storage Consumables and Devices



Azenta tubes are available in a range of formats, including racked, bulk, capped and uncapped. Our robust code management system ensures each and every code is permanently affixed, unique and never duplicated. Tubes are suitable for applications ranging from +121°C to -196°C.

All Azenta tubes are developed with broad compatibility in mind, performing without compromise in conjunction with automated code reading, capping and sample management systems from Azenta and all other industry-recognized manufacturers.

Key Features

Our coding systems are designed to deliver the highest levels of sample security and labeling flexibility. We offer tubes with a choice of code, either 2D-code, 1D (linear barcode) or Human Readable Number (HRN). Our Tri-coded tubes have all three options whilst our Dual-coded tubes, can have a 2D-code and Human Readable Number on the tube base or side. A robust code management system ensures each and every code is unique and never duplicated, preventing any possibility of misidentification involving your samples. All Azenta polypropylene sample storage tubes sealed with a screw cap are suitable for use in cryogenic storage conditions.

Superior Datamatrix code quality:

Azenta sample storage tubes are easy to read even in harsh conditions or when damaged. Our 2D-codes adhere fully to the ECC200 standard. All codes are generated using the latest and most sophisticated error correction methods and high quality, permanent laser etching provides sharp detail.


Automation Friendly:

Screw capped tubes are compatible with our range of cappers and de-cappers, including our IntelliXcap range. Available in 24, 48, and 96 format to automatically remove and re-cap a complete rack of tubes. The TwistLock feature prevents the tubes rotating in the rack during capping and de-capping.

Every tube is quality checked to ensure they meet our exacting standards for readability.

Excellent Sealing: Tubes work equally well with either TPE septum or screw caps. A double start thread engages in a maximum rotation of 180°, facilitating automation.

Ideal for Cold Storage: Temperature range from -196°C (with screw cap) to 121°C. All our polypropylene tubes are suitable for cryogenic storage, but not for submersion in liquid phase nitrogen.



Integration Friendly



AZENTA
LIFE SCIENCES

Choice of Coding Options

Tri-coded:

Each tube features a permanent 2D-code laser etched in high-contrast on the tube base, a permanent 1D (linear barcode) and Human-Readable Number laser etched in high-contrast on the tube side.

Each tube is manufactured using an advanced process which results in a one-piece tube with a clear window and black side for coding. All three codes are identical and auditing processes guarantee all three codes match. Our standard 2D option is black on white.

Dual-coded:

Each tube features a permanent 2D-code and Human-Readable Number laser etched in high-contrast on the tube base or side.

The Dual-coded tube is created using an advanced manufacturing technique used to integrate 2 resin colors into the same tube for high resolution coding.

Both codes are identical and auditing processes guarantee both codes match.

Our standard 2D option is white on black.

Non-coded – simply come as they are.

2D datamatrix code – a 2D code made of black and white cells, for Azenta tubes these are arranged in a square. The L-shape found on the border is its finder pattern, which is used by tube readers and scanners to recognize and read the code.

1D-coded / Linear barcode – unique barcode represented by parallel lines of different widths and spacings printed on the side of a tube, SBS rack or CryoBox. Our standard code is code 128. Code 128 is able to encode alpha-numeric data.

2D-coded – 2D-code on tube base.

Dual-coded – 2D-code and Human-Readable Number on tube base or side.

Tri-coded – 2D-code on tube base, 1D (linear barcode) and Human-Readable Number on tube side.




*Dual-Coded
(white on black)*

*2D-Orientation Coded
(black on white)*



*Unique 2D4 Code allows code
reading without interfering with
Acoustic dispensing window*



*2D orientation coded rack base allows automatic
detection of orientation by instrumentation*




SBS Storage Racks

SBS Storage Rack Options

Azenta tubes are compatible with industry standard SBS sized racks in either: 24, 48, 96, 240 or 384 format. Dependent upon tube type, the following rack options are available.


TwistLock:

Prevents tubes rotating within the rack to enable automated capping and de-capping of screw caps. TwistLock is provided as standard with the option available to remove.




TubeLock:

Tubes can be locked into the rack to prevent falling out, even when there is no lid present. Tubes can be placed in either locked or unlocked positions. TubeLock is used in manual workflows and is activated by applying pressure to the tube top, clicking the tube into place. Racked tubes can be ordered either pre-locked or non-locked.



LidLock:

Racks fitted with a LidLock latch are designed to withstand a 1m drop for added sample security.




Automatic Rack Orientation:

Racks can be supplied with a unique 2D-code identifier which can be read at the same time as the tube 2D-code. This provides automatic rack orientation and more secure sample tracking.

Direct Laser Etching:

A cutout window on the rack sides allow the linear barcode to be read more easily; linear barcodes can be laser etched directly onto racks.



Cryo Boxes

In addition to industry standard SBS boxes, a range of tube and application specific cryo boxes are available. Each cryo box incorporates:

- Open bottom for 2D-code decoding on camera-based readers
- Direct laser etched 1D linear barcode and Human-Readable Number on rack side and 2D box orientation ID on base of rack.

9 x 9 Cryo Storage Boxes

- 136.2mm x 136.2mm polycarbonate cryo box option for cryogenic sample storage
- Holds 81 tubes in 9 x 9 array



9 x 9 Box

10 x 10 Cryo Storage Boxes


- 136.2mm x 136.2mm polycarbonate cryo box option for cryogenic sample storage
- Holds 100 tubes in 10 x 10 array



10 x 10 Box

14 x 14 Cryo Storage Boxes

- 136.2mm x 136.2mm polycarbonate cryobox box option for cryogenic sample storage
- Holds 196 tubes in 14 x 14 array



14 x 14 Box

Ordering Information

66-1801	Cryo Box 9x9, black, polycarbonate, 10 boxes per case, suitable for 1.5ml, 1.6ml and 1.9ml External Thread Tubes
66-1802	Cryo Box 9x9, black, polycarbonate, 10 boxes per case, suitable for Taller Cryo Tubes
66-1800	Cryo Box 10x10, black, polycarbonate, 10 boxes per case, suitable for 1.9ml External Thread Tubes
66-1803	Cryo Box 10x10, black, polycarbonate, 10 boxes per case, suitable for Taller Cryo Tubes
66-0196-01	Cryo Box 14x14, black, polycarbonate, 10 boxes per case, suitable for 0.3ml and 0.48ml Internal Thread Tubes and 0.5ml External Thread Tubes
66-0196-02	Cryo Box 14x14, black, polycarbonate, 10 boxes per case, suitable for 0.8ml External Thread Tubes
66-0196-03	Cryo Box 14x14, black, polycarbonate, 10 boxes per case, suitable for 1.0ml External Thread Tubes

WARNING

Do not store tubes in liquid phase nitrogen

Ingress of nitrogen into the tube can occur causing the tube to rupture when taken out of storage

May cause injury and loss of tube contents.



AZENTA
LIFE SCIENCES

Anatomy of a FluidX Tube – Internal Thread, Dual-Coded

96-format, Internal Thread, Dual-coded tubes have a range of features that are only possible with advanced manufacturing techniques. Co-molded caps offer a far superior seal over traditional O-ring caps, which can 'pop' when a sample is frozen or the cap is overtightened.


Internal Threaded Cap:

- Automation friendly
- Co-molding prevents o-ring 'popping'
- Our most secure internal threaded cap




High Quality Virgin Polypropylene:

- No detectable leachables or extractables
- Manufactured in ISO Class 8 clean-room environment
- Endotoxin, DNase/RNase, heavy metals & animal free



Clear Window:

- Allows visual inspection
- Optional space to print direct due to high quality manufacturing



Dual-coded:

- 2D-code and Human-Readable Number (HRN) on tube base
- Enables whole rack or single tube reading
- High contrast enabling reliable reading
- Permanent laser etching



AZENTA
LIFE SCIENCES

Anatomy of a FluidX Tube – External Thread, Tri-Coded

The External Thread Tri-coded tube offers our most secure seal. This highly secure design also offers significant benefits over internal thread caps including, a higher working volume.

External Threaded Cap:

- Automation friendly
- Our most secure cap
- Enables greater working volume
- Designed to prevent over tightening

High Quality Virgin Polypropylene:

- No detectable leachables or extractables
- Manufactured in ISO Class 8 clean-room environment
- Endotoxin, DNase/RNase, heavy metals & animal free

Side Coding:

- Linear Barcode
- Human-Readable Number
- Permanent laser etching

Clear Window:

- Allows visual inspection
- Optional space to print direct due to high quality manufacturing

2D-coded:

- Enables whole rack or single tube reading
- High contrast enabling reliable reading
- Permanent laser etching




Anatomy of a FluidX Tube – External Thread, Dual-Coded

The Dual-Coded tube with external thread cap offers our most secure seal. This highly secure design also offers significant benefits over internal thread caps including, a higher working volume.


External Threaded Cap:

- Automation friendly
- Our most secure cap
- Enables greater working volume
- Designed to prevent over tightening



High Quality Virgin Polypropylene:

- No detectable leachables or extractables
- Manufactured in ISO Class 8 clean-room environment
- Endotoxin, DNase/RNase, heavy metals & animal free



Clear Window:

- Allows visual inspection
- Optional space to print direct due to high quality manufacturing



Dual-Coded:

- 2D-code and Human-Readable Number (HRN) on tube base
- Enables whole rack or single tube reading
- High contrast enabling reliable reading
- Permanent laser etching



AZENTA
LIFE SCIENCES

Anatomy of a FluidX Acoustic Sample Tube – Echo[®] Qualified Consumable

Internal Threaded Cap:

- New internal co-molded design
- Increases capping and de-capping performance to over 500 cycles

Optimized Geometry:

- Optimized tube geometry for Acoustic dispensing technology
- 70µl working volume

High Quality Virgin Polypropylene:

- No detectable leachables or extractables
- Manufactured in ISO Class 8 clean-room environment
- Endotoxin, DNase/RNase, heavy metals & animal free

2D4 Coded:

- Unique 2D4 Code allows reading of the code without interfering with Acoustic dispensing window
- Linear barcoded rack with Tube Retention for sample handling with Acoustic Dispenser and traditional liquid handlers
- Permanent laser etching



AZENTA
LIFE SCIENCES

Anatomy of a FluidX Maximum Recovery Tube

Featuring a novel design, the 1.6ml Maximum Recovery tube is ideally suited for automated workflows using liquid handling, where more sample can be retrieved due to the unique shaped bottom.

External Threaded Cap:

- Automation friendly
- Our most secure cap
- Enables greater working volume
- Designed to prevent over tightening



High Quality Virgin Polypropylene:

- No detectable leachables or extractables
- Manufactured in ISO Class 8 clean-room environment
- Endotoxin, DNase/RNase, heavy metals & animal free



Maximum Recovery:

- Unique tapered design with combined 'v' and 'u' shaped bottom
- Enables higher amount of sample retrieval
- Design creates a reduction in sample wastage resulting in cost savings



Tri-Coded:

- 1D code and Human-Readable Number (HRN) on tube base
- 2D datamatrix code on tube base
- Enables whole rack or single tube reading
- High contrast enabling reliable reading
- Permanent laser etching



Customization Options

In addition to standard catalogue products the Azenta range can also be customized to provide tubes and racks ideally matched to your particular workflow. Available customization options are detailed below.

Tube Coding

Custom prefixing with user defined two letter prefix followed by an 8 digit number sequence determined by Azenta. e.g. AB12345678.

Full custom coding with a user generated 10 digit alphanumerical sequence, e.g. ABC1234567.

Choice of Black on White or White on Black 2D datamatrix coding.

Rack Coding

Code 128 barcodes applied by Laser etching or labelling.

Barcode position selectable on any or multiple sides.

Cap Colors

In addition to Orange, caps are available in a variety of user selectable colors.

96 format caps available in: Natural, Clear/Light Blue, Dark Blue, Green, Purple, Red, White, Yellow, Amber and Black.

48 format caps available in: Dark Blue, Green, Red & Yellow for externally threaded tubes and in Natural for Internally threaded tubes.

24 format caps (automation friendly) for externally threaded tubes available in Natural.

TPE Septum caps available in Dark Blue, Green, Red, and Yellow.

Treatment Options

Gamma Irradiation, Ethylene Oxide, Dual Ethylene Oxide, or Electron Beam treatment treatments are available.









Special Options

Further options may be available depending upon product selection and details of requirements including custom product packaging and user definable tube/rack combinations.












Please note that for all custom products a minimum order quantity will be applied.

For further information on how to order please contact your local Azenta Life Sciences representative

FluidX Sample Storage Tubes from Azenta Life Sciences

								
	24-format, 7.6ml External Thread, Tri-coded	48-format, 3.8ml External Thread, Tri-coded	48-format, 1.9ml External Thread, Tri-coded	48-format, 1.6ml External Thread, Maximum Recovery	48-format, 1.5ml External Thread, Tri-coded	96-format 1ml External Thread, Tri-coded	96-format 0.9ml External Thread, Dual-coded	96-format, 0.9ml Internal Thread, Tri-coded
Further details see Page:	32	31	30	29	28	27	40	35
Max Fill Volume 21°C (ml) Screw Cap	9.2	4.6	2.3	2.0	1.8	1.2	1.0	1.1
Max Working Volume (ml) Screw Cap Frozen	7.6	3.8	1.9	1.6 Approx.	1.5	1.0	0.9	0.9
Max Working Volume (µl) Screw Cap	7660	3830	1910	1600	1500	1000	911	916
Max Working Volume (µl) Septum Cap	-	-	-	-	-	916	887	999
Tube Height (mm)	77.4	75	38.2	38.2	30.6	46.2	42.3	44.2
Tube Height with Cap (mm)	83.6	80.7	43.9	43.9	36.3	49.6	45.7	52.5
Tube Height with Septum Cap (mm)	-	-	-	-	-	47.4	43.5	45.4
Inner Diameter (mm)	13	9.6	9.6	9.6	9.6	6.5	6.5	6.8
Outer Diameter with Cap (mm)	17	12.8	12.8	12.8	12.6	8.7	8.7	8.6
Center to Center (mm)	18	13.5	13.5	13.5	13.5	9	9	9
Min Temperature °C Screw Cap	-196	-196	-196	-196	-196	-196	-196	-196
Min Temperature °C Septum Cap	-	-	-	-	-	-80	-80	-80
2D-coded	Base	Base	Base	Base	Base	Base	Base	Base
Human Readable Number	Side	Side	Side	Side	Side	Side	Base	Side
Linear Barcode	Side	Side	Side	Side	Side	Side	-	Side
Product Codes								
Bulk, Uncapped	65-9303	65-7516	65-7640	65-7650	65-7660	68-1003-00	68-1001-00	67-0757-00
Bulk, Capped	66-9302	65-7517	65-7641	65-7651	65-7661	68-1003-10	68-1001-10	67-0757-10
Racked, Uncapped	-	65-7514	65-7642	65-7652	65-7662	68-1003-01	68-1001-01	67-0757-01
Racked, Capped	-	65-7515	65-7643	65-7653	65-7663	68-1003-11	68-1001-11	67-0757-11



										
96-format, 0.9ml Internal Thread, Dual-coded	96-format, 0.8ml External Thread, Tri-coded	96-Format 0.7ml Internal Thread, Dual-coded	96-format 0.65ml Internal Thread, Tri-coded	96-format, 0.5ml External Thread, Tri-coded	96-format, 0.5ml External Thread, Dual-coded	96-format, 0.48ml Internal Thread, Tri-coded	96-format, 0.3ml Internal Thread, Dual-coded	96-format, 0.2ml External Thread, 2D-coded	96-format, 0.26ml External Thread, Dual-coded	Acoustic Sample Tube, Echo® Qualified Consumable
40	26	39	34	25	41	33	38	45	37	44
1.1	0.96	0.88	0.80	0.66	0.66	0.58	0.40	0.24	0.31	0.15
0.9	0.8	0.7	0.65	0.5	0.5	0.48	0.3	0.2	0.26	0.07
929	800	731	666	552	550	482	336	204	261	70
1018	-	821	749	525	525	572	425	-	238	-
43.5	36.9	36.2	36.8	26.4	26.4	26.4	21	24.1	15.2	13.4
51.8	40.3	44.5	45.1	29.8	29.8	34.7	29.3	27.5	18.6	14.7
44.7	-	37.4	38	27.6	27.6	27.6	22.1	-	16.4	
6.8	6.5	6.8	6.8	6.5	6.5	6.8	6.8	3.9	6.5	6.5
8.6	8.7	8.6	8.6	8.7	8.5	8.3	8.7	5.8	8.7	7.6
9	9	9	9	9	9	9	9	9 (96)	6 (240)	9
-196	-196	-196	-196	-196	-196	-196	-196	-196	-196	-80
-80	-	-80	-80	-80	-80	-80	-80	-	-80	-
Base	Base	Base	Base	Base	Base	Base	Base	Base	Base & Side	2D4 Quad Code
Base	Side	Base	Side	Side	Base	Side	Base	-	Side	-
-	Side	-	Side	Side	-	Side	-	-	-	-
66-62345	68-0801-00	66-62318	67-0755-00	68-0703-00	68-0701-00	67-0753-00	66-62326	67-0203-01	68-0303-00	-
66-62345-Y6	68-0801-10	66-62318-Y6	67-0755-10	68-0703-10	68-0701-10	67-0753-10	66-62326-Y6	67-0203-10	68-0303-10	-
66-62330	68-0801-01	66-62319	67-0755-01	68-0703-02	68-0701-02	67-0753-02	66-62325	67-0203-02	68-0303-01	-
66-62330-Y6	68-0801-11	66-62319-Y6	67-0755-11	68-0703-12	68-0701-12	67-0753-12	66-62325-Y6	67-0203-11	68-0303-11	69-0200-11



Azenta Sample Storage Tubes



Introduction to FluidX External Thread Tri-Coded Tubes

Overview

External Thread Tri-coded tubes have been developed to exceed the demands of sample security, management and tracking in modern high-density storage applications. Each tube features a permanent 2D-code laser etched in high-contrast on the tube base, a permanent 1D (linear barcode) and Human-Readable Number laser etched in high-contrast on the tube side. The tube is manufactured using an advanced process which results in a one-piece jacket tube, therefore, you never lose the code. All three codes are identical and auditing processes guarantee all three codes match. Our standard 2D coding option for these tubes is black on white.

The tubes provide a lifelong and secure chain of custody for samples in biobanks, compound libraries and a broad range of biological and chemical stores, including cryogenic storage.

Designed and developed with broad compatibility in mind, the tubes perform without compromise in conjunction with automated barcode reading, capping and sample management systems from Azenta and all other industry-recognized manufacturers.

Key Features

- Permanently laser etched, 2D-code on base, 1D (linear barcode) and Human-Readable Number on the side
- Tri-coded tubes offer unequalled sample audit traceability, enabling sample tracking and data sharing between multiple users, labs, locations and automation capabilities
- 100% Quality Control, each tube is tested to ensure both readability and uniqueness
- Manufactured using an advanced process which results in a one-piece jacket tube from high-quality virgin polypropylene
- External thread tubes provide greater working volume than internal thread tubes
- Suitable for cryogenic storage as low as -196°C*
- Secure sample storage and tracking
- Available bulk uncapped or capped
- Available pre-racked and capped in 96, 48 and 24 well format SBS racks
- 2D-codes readable without removing tubes from racks
- Suitable for sealing with either screw caps or TPE septum caps

**not for use in liquid phase Nitrogen*



*External Thread
Tri-coded tube*

Capping options

Screw Caps

- A deforming compression seal more effective than a silicone alternative
- A non-silicone seal means the cap can never be over-tightened
- Caps and tubes are manufactured from the same material, preventing differential expansion during freeze-thaw cycles
- A double-start thread engages in a maximum rotation of 180°, thereby facilitating automation

Introduction to FluidX Internal Thread Tri-Coded Tubes

Overview

Internal Thread Tri-coded tubes have been developed to exceed the demands of sample security, management and tracking in modern high-density storage applications and comply with ISBER standards. Each tube features a permanent 2D-code laser etched in high-contrast on the tube base, a permanent 1D (linear barcode) and Human-Readable Number laser etched in high-contrast on the tube side. The tube is manufactured using an advanced process which results in a one-piece jacket tube, therefore, you never lose the code. All three codes are identical and auditing processes guarantee all three codes match. Our standard 2D coding option for these tubes is black on white.

The tubes provide a lifelong and secure chain of custody for samples in biobanks, compound libraries and a broad range of biological and chemical stores, including cryogenic storage.

Designed and developed with broad compatibility in mind, the tubes perform without compromise in conjunction with automated barcode reading, capping and sample management systems from Azenta and all other industry-recognized manufacturers.

Key Features

Secure Sample Storage and Tracking

- Permanently laser etched, 2D-code on base, 1D (linear barcode) and Human-Readable Number on the side
- Tri-coded offers unequalled sample audit traceability, enabling sample tracking and data sharing between multiple users, labs, locations and automation capabilities
- 100% Quality Control, each tube is tested to ensure both readability and uniqueness
- Manufactured using an advanced process which results in a one-piece jacket tube from high-quality virgin polypropylene
- Suitable for cryogenic storage as low as -196°C^*
- Secure sample storage and tracking
- Available bulk uncapped or capped
- Available pre-racked and capped in 96 well format SBS racks
- 2D-codes readable without removing tubes from racks
- Suitable for sealing with either screw caps or TPE septum caps

**not for use in liquid phase Nitrogen*



*Internal Thread
Tri-coded tube*

Capping options


Screw Caps

- Automation friendly co-molded screw caps for internal thread tubes eliminate the possibility of over tightening
- Co-molded caps eliminate the failures in caps using a silicone O-ring
- A double-start thread engages in a maximum rotation of 180° , thereby facilitating automation

TPE Septum Caps

- Cost-effective sealing option for samples that are only accessed occasionally
- Septum caps are supplied in 96-format back mats to facilitate automation

0.5ml Tri-coded Tube, 96-format, External Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed using screw caps or TPE septum caps

	0.5ml screw cap	0.5ml TPE septum cap
		
	0.5ml Tri-coded Tube, 96-format, External Thread with screw cap	0.5ml Tri-coded Tube, 96-format, External Thread with TPE septum cap
Max Working Volume (ml)	0.52	0.50
Tube Height (mm)	26.4	26.4
Tube Height with Cap (mm)	29.8	27.6
Inner Diameter (mm)	6.5	6.5
Outer Diameter with Cap (mm)	8.7	8.7
Center to Center	9.0	9.0
Minimum Temperature (°C)	-196	-80
Tube Height in Rack (mm)	30.3	28.1
Overall Rack Height including lid (mm)	32.9	32.9

Ordering Information

68-0703-00	0.5ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk, 960 tubes per case
68-0703-10	0.5ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, bulk, 960 tubes per case
68-0703-02	0.5ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-51026
68-0703-11	0.5ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case, LowBase rack, empty rack part number: 66-51004
68-0703-12	0.5ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case, HighBase rack, empty rack part number: 66-51026
68-0704-00	0.5ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, uncapped, bulk, 960 tubes per case
68-0704-10	0.5ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, capped, bulk, 960 tubes per case
68-0704-02	0.5ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-51026
68-0704-12	0.5ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, capped, 10 racks per case, HighBase rack, empty rack part number: 66-51026



0.8ml Tri-coded Tube, 96-format, External Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed using screw caps or TPE septum caps

	0.8ml screw cap	0.8ml TPE septum cap
		
	0.8ml Tri-coded Tube, 96-format, External Thread with screw cap	0.8ml Tri-coded Tube, 96-format, External Thread with TPE septum cap
Max Working Volume (ml)	0.8	0.7
Tube Height (mm)	36.9	36.9
Tube Height with Cap (mm)	40.3	38.1
Inner Diameter (mm)	6.5	6.5
Outer Diameter with Cap (mm)	8.7	8.7
Center to Center	9.0	9.0
Minimum Temperature (°C)	-196	-80
Tube Height in Rack (mm)	41.2	39
Overall Rack Height including lid (mm)	43.9	43.9

Ordering Information

68-0801-00	0.8ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk, 960 tubes per case
68-0801-10	0.8ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, bulk, 960 tubes per case
68-0801-01	0.8ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-51021
68-0801-11	0.8ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case, HighBase rack, empty rack part number: 66-51021
68-0802-00	0.8ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, uncapped, bulk, 960 tubes per case
68-0802-10	0.8ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, capped, bulk, 960 tubes per case
68-0802-01	0.8ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-51021
68-0802-11	0.8ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, capped, 10 racks per case, HighBase rack, empty rack part number: 66-51021




AZENTA
LIFE SCIENCES

1.0ml Tri-coded Tube, 96-format, External Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed using screw caps or TPE septum caps

	1.0ml screw cap	1.0ml TPE septum cap
		
	1.0ml Tri-coded Tube, 96-format, External Thread with screw cap	1.0ml Tri-coded Tube, 96-format, External Thread with TPE septum cap
Max Working Volume (ml)	1.0	0.9
Tube Height (mm)	46.2	46.2
Tube Height with Cap (mm)	49.6	47.4
Inner Diameter (mm)	6.5	6.5
Outer Diameter with Cap (mm)	8.7	8.7
Center to Center	9.0	9.0
Minimum Temperature (°C)	-196	-80
Tube Height in Rack (mm)	50.5	47.1
Overall Rack Height including lid (mm)	53.2	53.2

Ordering Information

68-1003-00	1.0ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk , 960 tubes per case
68-1003-10	1.0ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, bulk , 960 tubes per case
68-1003-01	1.0ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, 10 racks per case , HighBase rack, empty rack part number: 66-51020
68-1003-11	1.0ml Tri-coded Tube, 96-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case , HighBase rack, empty rack part number: 66-51020
68-1004-00	1.0ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, uncapped, bulk , 960 tubes per case
68-1004-10	1.0ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, capped, bulk , 960 tubes per case
68-1004-01	1.0ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, uncapped, 10 racks per case , HighBase rack, empty rack part number: 66-51020
68-1004-11	1.0ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, capped, 10 racks per case , HighBase rack, empty rack part number: 66-51020



AZENTA
LIFE SCIENCES

1.5ml Tri-coded Tube, 48-format, External Thread



- Supplied in 48-well format SBS racks or bulk, empty 9x9 cryo rack also available (see page 13)
- Securely sealed using screw caps

1.5ml screw cap



1.5ml Tri-coded Tube, 48-format, External Thread with screw cap

Max Working Volume (ml)	1.5
Tube Height (mm)	30.6
Tube Height with Cap (mm)	36.3
Inner Diameter (mm)	9.6
Outer Diameter with Cap (mm)	12.8
Center to Center	13.5
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	45.4
Overall Rack Height including lid (mm)	49.4

Ordering Information

65-7660	1.5ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk, 480 tubes per case
65-7661	1.5ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, bulk, 480 tubes per case
65-7662	1.5ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D linear barcode and Human Readable Number on side, uncapped, 10 racks per case, 48-format rack (1 piece rack base), empty rack part number: 65-9451
65-7663	1.5ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case, 48-format rack (1 piece rack base), empty rack part number: 65-9451
65-7664	1.5ml 2D-coded Tube, 48-format, External Thread, 2D Code on base, uncapped, bulk, 480 tubes per case
65-7665	1.5ml 2D-coded Tube, 48-format, External Thread, 2D Code on base, capped, bulk, 480 tubes per case
65-7666	1.5ml 2D-coded Tube, 48-format, External Thread, 2D Code on base, uncapped, 10 racks per case, 48-format rack (1 piece rack base), empty rack part number: 65-9451
65-7667	1.5ml 2D-coded Tube, 48-format, External Thread, 2D Code on base, capped, 10 racks per case, 48-format rack (1 piece rack base), empty rack part number: 65-9451
66-1801	Cryo Rack 9x9, black, polycarbonate, 10 racks per case, suitable for 1.5ml and 1.9ml External Thread Tubes



1.6ml Tri-coded Tube, Maximum Recovery



- Supplied in 48-well format SBS racks or bulk, empty 9x9 cryo rack also available (see page 13)
- Securely sealed using screw caps

1.6ml screw cap



1.6ml Tri-coded Tube, 48-format, External Thread with screw cap

Approx. Working Volume (ml)	1.6
Tube Height (mm)	38.2
Tube Height with Cap (mm)	43.9
Inner Diameter (mm)	9.6
Outer Diameter with Cap (mm)	12.8
Center to Center	13.5
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	45.0
Overall Rack Height including lid (mm)	49.4

Ordering Information

65-7650	1.6ml Tri-coded Tube, 48-format, External Thread, Maximum Recovery, 48-format, external thread, 2D code on base, 1D linear barcode and human readable code on side, uncapped, bulk; 480 tubes
65-7651	1.6ml Tri-coded Tube, 48-format, External Thread, Maximum Recovery, 48-format, external thread, 2D code on base, 1D linear barcode and human readable code on side, capped, bulk, 480 tubes per case
65-7652	1.6ml Tri-coded Tube, 48-format, External Thread, Maximum Recovery, 48-format, external thread, 2D code on base, 1D linear barcode and human readable code on side, uncapped, racked (65-9451); 10 racks/480 tubes
65-7653	1.6ml Tri-coded Tube, 48-format, External Thread, Maximum Recovery, 48-format, external thread, 2D code on base, 1D linear barcode and human readable code on side, capped, racked (65-9451); 10 racks/480 tubes



AZENTA
LIFE SCIENCES

1.9ml Tri-coded Tube, 48-format, External Thread



- Supplied in 48-well format SBS racks or bulk, empty 9x9 cryo rack also available (see page 13)
- Securely sealed using screw caps

1.9ml screw cap



1.9ml Tri-coded Tube, 48-format, External Thread with screw cap

Max Working Volume (ml)	1.9
Tube Height (mm)	38.2
Tube Height with Cap (mm)	43.9
Inner Diameter (mm)	9.6
Outer Diameter with Cap (mm)	12.8
Center to Center	13.5
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	45
Overall Rack Height including lid (mm)	49.4

Ordering Information

65-7640	1.9ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk, 480 tubes per case
65-7641	1.9ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, bulk, 480 tubes per case
65-7642	1.9ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, 10 racks per case, 48-format rack (1 piece rack base), empty rack part number: 65-9451
65-7643	1.9ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case, 48-format rack (1 piece rack base), empty rack part number: 65-9451
65-7644	1.9ml 2D-coded Tube, 48-format, External Thread, 2D Code on base, uncapped, bulk, 480 tubes per case
65-7645	1.9ml 2D-coded Tube, 48-format, External Thread, 2D Code on base, capped, bulk, 480 tubes per case
65-7646	1.9ml 2D-coded Tube, 48-format, External Thread, 2D Code on base, uncapped, 10 racks per case, 48-format rack (1 piece rack base), empty rack part number: 65-9451
65-7647	1.9ml 2D-coded Tube, 48-format, External Thread, 2D Code on base, capped, 10 racks per case, 48-format rack (1 piece rack base), empty rack part number: 65-9451
66-1800	Cryo Rack 10x10, black, polycarbonate, 10 racks per case, suitable for 1.9ml External Thread Tubes



3.8ml Tri-coded Tube, 48-format, External Thread



- Supplied in 48-well format SBS racks or bulk
- Securely sealed with screw caps

3.8ml screw cap



3.8ml Tri-coded Tube, 48-format, External Thread with screw cap

Max Working Volume (ml)	3.8
Tube Height (mm)	75.0
Tube Height with Cap (mm)	80.7
Inner Diameter (mm)	9.6
Outer Diameter with Cap (mm)	12.8
Center to Center	13.5
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	81.7
Overall Rack Height including lid (mm)	86.2

Ordering Information

65-7516	3.8ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk, 480 tubes per case
65-7517	3.8ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, bulk, 480 tubes per case
65-7514	3.8ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, 10 racks per case, 48-format rack (2 piece rack base), empty rack part number: 65-9460
65-7515	3.8ml Tri-coded Tube, 48-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case, 48-format rack (2 piece rack base), empty rack part number: 65-9460



7.6ml Tri-coded Tube, 24-format, External Thread



- Supplied bulk, empty
24-format SBS racks
available separately
- Securely sealed
with screw caps

7.6ml screw cap



7.6ml Tri-coded Tube, 24-format, External Thread with screw cap

Max Working Volume (ml)	7.6
Tube Height (mm)	77.4
Tube Height with Cap (mm)	83.6
Inner Diameter (mm)	13.0
Outer Diameter with Cap (mm)	17.0
Center to Center	18.0
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	84.7
Overall Rack Height including lid (mm)	88.5

Ordering Information

65-9303	7.6ml Tri-coded Tube, 24-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk, 240 tubes per case
66-9302	7.6ml Tri-coded Tube, 24-format, External Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped with Automation Friendly Screw Cap, bulk, 240 tubes per case

Large format Caps

66-9401	Screw Cap, 24-format, External Thread, orange, automation friendly, bulk, 240 caps per case, suitable for 7.6ml Tubes
---------	---

Large format Racks

66-9455	Rack, 24-format, 1 piece rack base, with open bottom for reading on rack readers, 10 racks per case, suitable for 7.6ml External Thread, Tri-coded Tubes (part number 66-9302)
---------	--



0.48ml Tri-coded Tube, 96-format, Internal Thread




- Supplied in 96-well format SBS racks or bulk
- Securely sealed with screw caps or TPE septum caps

	0.48ml screw cap	0.48ml TPE septum cap
		
	0.48ml Tri-coded Tube, 96-format Internal Thread with screw cap	0.48ml Tri-coded Tube, 96-format Internal Thread with TPE septum cap
Max Working Volume (ml)	0.48	0.57
Tube Height (mm)	26.4	26.4
Tube Height with Cap (mm)	34.7	27.6
Inner Diameter (mm)	6.8	6.8
Outer Diameter with Cap (mm)	8.6	8.6
Center to Center	9.0	9.0
Minimum Temperature (°C)	-196	-80
Tube Height in Rack (mm)	36.2	35.2
Overall Rack Height including lid (mm)	44.9	43.9

Ordering Information

67-0753-00	0.48ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk, 960 tubes per case
67-0753-10	0.48ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, bulk, 960 tubes per case
67-0753-02	0.48ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-51025
67-0753-12	0.48ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case, HighBase rack, empty rack part number: 66-51025



0.65ml Tri-coded Tube, 96-format, Internal Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed with screw caps or TPE septum caps

	0.65ml screw cap	0.65ml TPE septum cap
		
	0.65ml Tri-coded Tube, 96-format, Internal Thread with screw cap	0.65ml Tri-coded Tube, 96-format, Internal Thread with TPE septum cap
Max Working Volume (ml)	0.65	0.75
Tube Height (mm)	36.8	36.8
Tube Height with Cap (mm)	45.1	38.0
Inner Diameter (mm)	6.8	6.8
Outer Diameter with Cap (mm)	8.4	8.4
Center to Center	9.0	9.0
Minimum Temperature (°C)	-196	-80
Tube Height in Rack (mm)	46	38.9
Overall Rack Height including lid (mm)	50.8	43.9

Ordering Information


67-0755-00	0.65ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk, 960 tubes per case
67-0755-01	0.65ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-51022
67-0755-10	0.65ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, bulk, 960 tubes per case
67-0755-11	0.65ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case, HighBase rack, empty rack part number: 66-51022



0.9ml Tri-coded Tube, 96-format, Internal Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed with screw caps

	0.9ml screw cap	0.9ml TPE septum cap
		
	0.9ml Tri-coded Tube, 96-format Internal Thread with screw cap	0.9ml Tri-coded Tube, 96-format, Internal Thread with TPE septum cap
Max Working Volume (ml)	0.9	0.99
Tube Height (mm)	44.2	44.2
Tube Height with Cap (mm)	52.5	45.4
Inner Diameter (mm)	6.8	6.8
Outer Diameter with Cap (mm)	8.4	8.6
Center to Center	9.0	9.0
Minimum Temperature (°C)	-196	-80
Tube Height in Rack (mm)	53.4	46.3
Overall Rack Height including lid (mm)	61.8	50.8

Ordering Information

67-0757-00	0.9ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, bulk, 960 tubes per case
67-0757-10	0.9ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, bulk, 960 tubes per case
67-0757-01	0.9ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-51023
67-0757-11	0.9ml Tri-coded Tube, 96-format, Internal Thread, 2D Code on base, 1D Linear Barcode and Human Readable Number on side, capped, 10 racks per case, HighBase rack, empty rack part number: 66-51023



Introduction to 96-Format Dual-Coded Tubes

Overview

The Dual-Coded Tube features a 2D-code and Human Readable Number (HRN) on the tube base, allowing compatibility with low throughput manual workflows, semi-automated workflows or fully automated workflows on integrated platforms.

The tubes provide a lifelong and secure chain of custody for samples in biobanks, compound libraries and a broad range of biological and chemical stores, including cryogenic storage.

Designed and developed with broad compatibility in mind, the tubes perform without compromise in conjunction with automated barcode reading, capping and sample management systems from Azenta and all other industry-recognized manufacturers.



*Dual-Coded tube
internal thread*

Key Features

- Permanently laser etched, 2D-code and a Human Readable Number (HRN) on the tube base
- Developed to exceed the demands of sample security, management and tracking in modern high-density storage
- 100% Quality Control, each tube is tested to ensure both readability and uniqueness
- Manufactured using an advanced process which results in a one-piece tube from high-quality virgin polypropylene
- 2D-code and HRN ensure a permanent link between sample and data
- High-contrast 2D-codes are more reliably readable in frost or condensation conditions
- 2D-codes can be scanned and decoded without removing tubes from storage racks, enabling data to be associated with individual tubes
- Equally suitable for sealing with either screw caps or TPE septum caps
- Azenta tubes have been leak tested to ensure sample security
- Suitable for cryogenic storage
- Manufactured from high-quality virgin polypropylene

Screw Caps

- Automation friendly co-molded screw caps for internal thread tubes eliminate the possibility of over tightening
- Co-molded caps eliminate the failures in caps using a silicone O-ring
- A double-start thread engages in a maximum rotation of 180°, thereby facilitating automation

96-Well Format SBS Racks


- **TwistLock:** prevents tubes rotating within the rack to enable automated capping and de-capping of screw caps, provided as standard
- **LidLock:** racks fitted with a LidLock latch are designed to withstand a 1m drop test for added sample security
- **TubeLock:** tubes can be locked in position in a rack, even without a lid, preventing them from falling out even if the rack is inverted. Lock or unlock simply by pushing the tube downwards or upwards
- **Automatic Rack Orientation:** racks are supplied with a unique 2D code identifier which can be read at the same time as the tube 2D code, to provide automatic rack orientation and more secure sample tracking
- **Direct Laser Etching:** linear barcodes are permanently etched directly onto the rack

14 x 14 Cryo Storage Boxes

- 136.2mm x 136.2mm polycarbonate cryobox option available for cryogenic sample storage
- Holds 196 tubes in 14 x 14 array
- Open bottom for 2D code decoding on Azenta Camera-Based Reader for SBS Racks and Cryo Boxes
- Cryo boxes can be supplied with a unique 2D code identifier which can be read at the same time as the tube 2D barcode, to provide more secure sample tracking



0.26ml Dual-coded Tube, 96-format, External Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed using screw caps or TPE septum caps

	0.26ml screw cap	0.26ml septum cap
		
	0.26ml Dual-coded Tube, 96-format, External Thread with screw cap	0.26ml Dual-coded Tube, 96-format, External Thread with septum cap
Max Working Volume (ml)	0.26	0.23
Tube Height (mm)	15.2	15.2
Tube Height with Cap (mm)	18.6	16.4
Inner Diameter (mm)	6.5	6.5
Outer Diameter with Cap (mm)	8.7	8.7
Center to Center	9.0	9.0
Minimum Temperature (°C)	-196	-80
Tube Height in Rack (mm)	19.13	15.73
Overall Rack Height including lid (mm)	22	22

Ordering Information

68-0303-00	0.26ml Dual-coded Tube, 96-format, External Thread, 2D Code on base, 2D and Human Readable Number on side, uncapped, bulk, 960 tubes per case
68-0303-10	0.26ml Dual-coded Tube, 96-format, External Thread, 2D Code on base, 2D and Human Readable Number on side, capped, bulk, 960 tubes per case
68-0303-01	0.26ml Dual-coded Tube, 96-format, External Thread, 2D Code on base, 2D and Human Readable Number on side, uncapped, 10 racks per case, LowBase rack, 2D rack ID position H12, empty rack part number: 68-0300-20
68-0303-11	0.26ml Dual-coded Tube, 96-format, External Thread, 2D Code on base, 2D and Human Readable Number on side, capped, 10 racks per case, LowBase rack, 2D rack ID position H12, empty rack part number: 68-0300-20

Note: available with 2D-code only on request.



AZENTA
LIFE SCIENCES

0.3ml Dual-coded Tube, 96-format, Internal Thread



- Supplied in 96-well format SBS racks, or bulk compatible with 14x14 cryo storage racks
- Securely sealed with screw caps or TPE septum caps

0.3ml screw cap




0.3ml Dual-coded Tube, 96-format, Internal Thread with screw cap

Max Working Volume (ml)	0.33
Tube Height (mm)	21
Tube Height with Cap (mm)	29.3
Inner Diameter (mm)	6.8
Outer Diameter with Cap (mm)	8.7
Center to Center	9
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	29.5
Overall Rack Height including lid (mm)	33.1

Ordering Information

66-62326	0.3ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, uncapped, bulk, 960 tubes per case
66-62326-Y6	0.3ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, capped, bulk, 960 tubes per case
66-62325	0.3ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, uncapped, 10 racks per case, LowBase rack, empty rack part number: 66-51004
66-62325-Y6	0.3ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, capped, 10 racks per case, LowBase rack, empty rack part number: 66-51004



AZENTA
LIFE SCIENCES

0.7ml Dual-coded Tube, 96-format, Internal Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed with screw caps or TPE septum caps

0.7ml screw cap



0.7ml Dual-coded Tube, 96-format, Internal Thread with screw cap

Max Working Volume (ml)	0.73
Tube Height (mm)	36.2
Tube Height with Cap (mm)	44.5
Inner Diameter (mm)	6.8
Outer Diameter with Cap (mm)	8.6
Center to Center	9
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	46.4
Overall Rack Height including lid (mm)	50.8

Ordering Information

66-62318	0.7ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, uncapped , bulk, 960 tubes per case
66-62318-Y6	0.7ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, capped , bulk, 960 tubes per case
66-62319	0.7ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, uncapped , 10 racks per case, HighBase rack, empty rack part number: 66-61002
66-62319-Y6	0.7ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, capped , 10 racks per case, HighBase rack, empty rack part number: 66-61002



AZENTA
LIFE SCIENCES

0.9ml Dual-coded Tube, 96-format, Internal Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed with screw caps or TPE septum caps

0.9ml screw cap




0.9ml Dual-coded Tube, 96-format, Internal Thread with screw cap

Max Working Volume (ml)	0.92
Tube Height (mm)	43.5
Tube Height with Cap (mm)	51.8
Inner Diameter (mm)	6.8
Outer Diameter with Cap (mm)	8.6
Center to Center	9
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	53.8
Overall Rack Height including lid (mm)	61.8

Ordering Information

66-62345	0.9ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, uncapped, bulk, 960 tubes per case
66-62345-Y6	0.9ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, capped, bulk, 960 tubes per case
66-62330	0.9ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, uncapped, 10 racks per case, HighBase rack, lid suitable for Screw Caps and TPE Caps, empty rack part number: 66-61003
66-62330-Y6	0.9ml Dual-coded Tube, 96-format, Internal Thread, 2D Code and Human Readable Number on base, capped, 10 racks per case, HighBase rack, lid suitable for Screw Caps and TPE Caps, empty rack part number: 66-61003



AZENTA
LIFE SCIENCES

0.5ml Dual-coded Tube, 96-format, External Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed with screw caps or TPE septum caps

0.5ml screw cap



0.5ml Dual-coded Tube, 96-format, External Thread with screw cap

Max Working Volume (ml)	0.55
Tube Height (mm)	26.4
Tube Height with Cap (mm)	29.8
Inner Diameter (mm)	6.5
Outer Diameter with Cap (mm)	8.5
Center to Center	9
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	30.0
Overall Rack Height including lid (mm)	33.1

Ordering Information

68-0701-00	0.5ml Dual-coded Tube, 96-format, External Thread, 2D Code and Human Readable Number on base, uncapped, bulk, 960 tubes per case
68-0701-10	0.5ml Dual-coded Tube, 96-format, External Thread, 2D Code and Human Readable Number on base, capped, bulk, 960 tubes per case
68-0701-02	0.5ml Dual-coded Tube, 96-format, External Thread, 2D Code and Human Readable Number on base, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-51026
68-0701-12	0.5ml Dual-coded Tube, 96-format, External Thread, 2D Code and Human Readable Number on base, capped, 10 racks per case, HighBase rack, empty rack part number: 66-51026
68-0701-11	0.5ml Dual-coded Tube, 96-format, External Thread, 2D Code and Human Readable Number on base, capped, 10 racks per case, LowBase rack, empty rack part number: 66-51004



AZENTA
LIFE SCIENCES

0.9ml Dual-coded Tube, 96-format, External Thread



- Supplied in 96-well format SBS racks or bulk
- Securely sealed with screw caps or TPE septum caps

0.9ml screw cap




0.9ml Dual-coded Tube, 96-format, External Thread with screw cap

Max Working Volume (ml)	0.9
Tube Height (mm)	42.3
Tube Height with Cap (mm)	45.7
Inner Diameter (mm)	6.5
Outer Diameter with Cap (mm)	8.7
Center to Center	9
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	47.8
Overall Rack Height including lid (mm)	50.8

Ordering Information

68-1001-00	0.9ml Dual-coded Tube, 96-format, External Thread, 2D Code and Human Readable Number on base, uncapped, bulk, 960 tubes per case
68-1001-10	0.9ml Dual-coded Tube, 96-format, External Thread, 2D Code and Human Readable Number on base, capped, bulk, 960 tubes per case
68-1001-01	0.9ml Dual-coded Tube, 96-format, External Thread, 2D Code and Human Readable Number on base, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-61002
68-1001-11	0.9ml Dual-coded Tube, 96-format, External Thread, 2D Code and Human Readable Number on base, capped, 10 racks per case, HighBase rack, empty rack part number: 66-61002
66-61002	Rack, 96-format, HighBase, 10 racks per case, suitable for 0.7ml Internal Thread Tubes with Screw Caps, 0.9ml Internal Thread, Dual-Coded Tubes with Septum TPE Caps or Individual Thermal Seals and 0.9ml External Thread Tubes with Screw Caps
66-51016	Rack, 96-format, HighBase, with TubeLock, 10 racks per case, suitable for 0.9ml External Thread Tubes with Screw Caps or Septum TPE Caps, 0.7ml Internal Thread Tubes with Screw Caps or Septum TPE Caps and 0.9ml Internal Thread Tubes with Septum TPE Caps
66-0196-01	Cryo Rack 14x14, black, polycarbonate, 10 racks per case, suitable for 0.3ml and 0.48ml Internal Thread Tubes and 0.5ml External Thread Tubes



AZENTA
LIFE SCIENCES

External Thread 2D-Coded Tubes

External Thread 2D-coded tubes carry a unique and permanent high-contrast 2D-code tube identifier on the base of the tube readable in frost or condensation conditions, or when damaged. The tubes provide a lifelong and secure chain of custody for samples in biobanks, compound libraries and a broad range of biological and chemical stores, including cryogenic storage.


Key Features

- Permanently laser etched, 2D-code on base
- 100% Quality Control, each tube is tested to ensure both readability and uniqueness
- Manufactured from high-quality virgin polypropylene
- External thread tubes provide greater working volume than internal thread tubes
- Suitable for cryogenic storage as low as -196°C
- Secure Sample Storage and Tracking
- Available bulk uncapped or capped

- Available pre-racked and capped in 96, 48 and 24 well format SBS racks
- 2D-codes readable without removing tubes from racks

Screw Caps

- A deforming compression seal is more effective than a silicone alternative
- Non-silicone seal means the cap can never be over-tightened
- Caps and tubes are manufactured from the same material, preventing differential expansion during freeze-thaw cycles
- A double-start thread engages in a maximum rotation of 180° thereby facilitating automation



Internal Thread 2D-Coded Tubes

2D-coded internal thread tubes carry a unique and permanent high-contrast 2D-code tube identifier on the base of the tube providing a lifelong and secure chain of custody for samples in biobanks, compound libraries and a broad range of biological and chemical stores, including cryogenic storage.

Key Features

- Permanently laser etched, 2D-code on base
- 100% Quality Control, each tube is tested to ensure both readability and uniqueness
- Manufactured from high-quality virgin polypropylene
- Suitable for cryogenic storage as low as -196°C
- Secure sample storage and tracking
- Available bulk uncapped or capped
- Available pre-racked and capped in 96, 48 and 24 well format SBS racks
- 2D-codes readable without removing tubes from racks

Screw Caps


- Co-molded screw caps for internal thread tubes eliminate the possibility of over-tightening
- Co-molded caps eliminate the failures common in caps using silicon O-rings
- Improved internal thread provides a more secure seal when using screw caps

TPE Septum Caps

- Cost-effective sealing option for samples that are only accessed occasionally
- Septum caps are supplied in 96-format back mats to facilitate automation



Acoustic Sample Tube – Echo® Qualified Consumable



- Acoustic dispense direct from the tube
- Assay robustness based on a higher sample integrity
- Compound saving vs. Conventional storage methods
- Faster turnaround time from request for compound to compound ready for assay
- Established benefits of 2D-coded Tubes vs. Plates which have to be fully thawed and unsealed to access one sample
- Combined technologies of 2D-coded Tubes and Acoustic Transfer
- Unique and permanent 2D4 Quad Code laser-etched on tube base
- Secure Sample Storage and Tracking - barcode ensures a permanent link between sample and data
- Scan and decode without removing tubes from storage racks, enables data to be associated with individual tubes
- Leak tested to ensure sample security
- Screw Caps - non-silicone seal means the cap can never be overtightened
- Double-start thread engages in a maximum rotation of 180°, facilitating automation

Acoustic Sample Tube – Echo® Qualified Consumable



Acoustic Sample Tube – Echo® Qualified Consumable

Max Working Volume	70 µL @5mm
Dead Volume	Approx. 15µL (with DMSO)
Total Volume	85µL
Rack Option	SBS format 96-way rack
Tube Height (mm)	13.4
Tube Height with Cap (mm)	14.7
Outer Diameter with Cap (mm)	7.9
Rack Height with lid (mm)	22.1
Coding	2D4 Quad Code on base
Capping Option	Internal Thread


Ordering Information

Please contact your local Azenta Life Sciences representative.



AZENTA
LIFE SCIENCES

0.2ml 2D-coded Tube, 96-format, External Thread



- Supplied in 96-well format SBS racks
- Securely sealed using screw caps
- Available in 240-well format on request

0.2ml screw cap



0.2ml 2D-coded Tube, 96-format, External Thread with screw cap

Max Working Volume (ml)	0.2
Tube Height (mm)	24.1
Tube Height with Cap (mm)	27.5
Inner Diameter (mm)	3.9
Outer Diameter with Cap (mm)	5.8
Center to Center	9.0
Minimum Temperature (°C)	-196
Tube Height in Rack (mm)	28
Overall Rack Height including lid (mm)	31

Ordering Information

67-0203-01	0.2ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, uncapped, bulk, 960 tubes per case
67-0203-10	0.2ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, capped, bulk, 960 tubes per case
67-0203-02	0.2ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, uncapped, 10 racks per case, empty rack part number: 67-0203-00
67-0203-11	0.2ml 2D-coded Tube, 96-format, External Thread, 2D Code on base, capped, 10 racks per case, empty rack part number: 67-0203-00
67-0203-51	Screw Cap, 96-format, External Thread, bulk, 960 caps per case, suitable for 0.2ml 2D-coded Tube, 96-format, External Thread
67-0203-00	Rack, 96-format, 10 racks per case, suitable for 0.2ml Tubes
67-0200-00	Rack, 240-format, 10 racks per case, suitable for 0.2ml Tubes



Non-Coded Tubes

In addition to our fully traceable coded tubes offering sample security, management and tracking in modern high-density storage applications, Azenta tubes are also available non-coded or with alpha numeric coding.

Ordering Information

0.5ml Non-Coded External Thread Screw Cap Tubes

66-0700-01	0.5ml Non-coded Tube, 96-format, External Thread, uncapped, 10 racks per case, LowBase rack, empty rack part number: 66-51004
66-0700-11	0.5ml Non-coded Tube, 96-format, External Thread, capped, 10 racks per case, LowBase rack, empty rack part number: 66-51004
66-0700-02	0.5ml Non-coded Tube, 96-format, External Thread, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-51026
66-0700-12	0.5ml Non-coded Tube, 96-format, External Thread, capped, 10 racks per case, HighBase rack, empty rack part number: 66-51026
66-0700-00	0.5ml Non-coded Tube, 96-format, External Thread, uncapped, bulk, 960 tubes per case
66-0700-10	0.5ml Non-coded Tube, 96-format, External Thread, capped, bulk, 960 tubes per case

0.9ml Non-Coded External Thread Screw Cap Tubes

66-1000-01	0.9ml Non-coded Tube, 96-format, External Thread, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-61002
66-1000-11	0.9ml Non-coded Tube, 96-format, External Thread, capped, 10 racks per case, HighBase rack, empty rack part number: 66-61002
66-1000-02	0.9ml Non-coded Tube, 96-format, External Thread, uncapped, 10 racks per case, HighBase rack, with TubeLock, empty rack part number: 66-51016
66-1000-12	0.9ml Non-coded Tube, 96-format, External Thread, capped, 10 racks per case, HighBase rack, with TubeLock, empty rack part number: 66-51016
66-1000-00	0.9ml Non-coded Tube, 96-format, External Thread, uncapped, bulk, 960 tubes per case
66-1000-10	0.9ml Non-coded Tube, 96-format, External Thread, capped, bulk, 960 tubes per case



Non-Coded Tubes

0.3ml Non-Coded Internal Thread Tubes

66-32041	0.3ml Non-coded Tube, 96-format, Internal Thread, V-bottom, uncapped, 10 racks per case, LowBase rack, suitable for use with Screw Cap Tubes, empty rack part number: 66-51004
66-32041-Y6	0.3ml Non-coded Tube, 96-format, Internal Thread, V-bottom, capped, 10 racks per case, LowBase rack, empty rack part number: 66-51004
66-32041-Y6-L	0.3ml Non-coded Tube, 96-format, Internal Thread, V-bottom, capped, 10 racks per case, LowBase rack, with TubeLock, includes standard profile non-locking lid, empty rack part number: 66-51014
66-32141	0.3ml Non-coded Tube, 96-format, Internal Thread, uncapped, 10 racks per case, LowBase rack, lid suitable for use with TPE Caps/Thermal Individual Tube Seal only, empty rack part number: 66-51003
66-32040	0.3ml Non-coded Tube, 96-format, Internal Thread, V-bottom, uncapped, bulk, 960 tubes per case
66-32040-Y6	0.3ml Non-coded Tube, 96-format, Internal Thread, V-bottom, capped, bulk, 960 tubes per case

0.7ml Non-Coded Internal Thread Tubes

66-32033	0.7ml Non-coded Tube, 96-format, Internal Thread, V-bottom, uncapped, bulk, 960 tubes per case
66-32033-Y6	0.7ml Non-coded Tube, 96-format, Internal Thread, V-bottom, capped, bulk, 960 tubes per case
66-32034	0.7ml Non-coded Tube, 96-format, Internal Thread, V-bottom, uncapped, 10 racks per case, HighBase rack, empty rack part number: 66-61002
66-32034-Y6	0.7ml Non-coded Tube, 96-format, Internal Thread, V-bottom, capped, 10 racks per case, HighBase rack, empty rack part number: 66-61002
66-32034-L	0.7ml Non-coded Tube, 96-format, Internal Thread, V-bottom, uncapped, 10 racks per case, HighBase rack, with TubeLock, empty rack part number: 66-51016
66-32034-Y6-L	0.7ml Non-coded Tube, 96-format, Internal Thread, V-bottom, capped, 10 racks per case, HighBase rack, with TubeLock, empty rack part number: 66-51016

0.9ml Non-Coded Internal Thread Screw Cap Tubes

66-32062	0.9ml Non-coded Tube, 96-format, Internal Thread, V-bottom, uncapped, bulk, 960 tubes per case
66-32062-Y6	0.9ml Non-coded Tube, 96-format, Internal Thread, V-bottom, capped, bulk, 960 tubes per case
66-32042	0.9ml Non-coded Tube, 96-format, Internal Thread, uncapped, 10 racks per case, HighBase rack, lid suitable for TPE caps only, empty rack part number: 66-61002
66-32043	0.9ml Non-coded Tube, 96-format, Internal Thread, uncapped, 10 racks per case, HighBase rack, lid suitable for Screw Caps and TPE Caps, empty rack part number: 66-61003
66-32043-Y6	0.9ml Non-coded Tube, 96-format, Internal Thread, capped, 10 racks per case, HighBase rack, lid suitable for Screw Caps and TPE Caps, empty rack part number: 66-61003
66-32042-L	0.9ml Non-coded Tube, 96-format, Internal Thread, uncapped, 10 racks per case, HighBase rack, with TubeLock, lid suitable for TPE caps only, empty rack part number: 66-51016
66-32043-L	0.9ml Non-coded Tube, 96-format, Internal Thread, uncapped, 10 racks per case, HighBase rack, with TubeLock, lid suitable for Screw Caps and TPE Caps, empty rack part number: 66-51017
66-32043-Y6-L	0.9ml Non-coded Tube, 96-format, Internal Thread, capped, 10 racks per case, HighBase rack with TubeLock, empty rack part number: 66-51017



Capping and Sealing Options



Tube Screw Caps



Developed to exceed the demands of sample security, management and tracking in modern high-density storage applications, screw caps are manufactured from high-quality virgin polypropylene and are designed for optimal seal quality and sample security.

Compatible with all Azenta 96-format sample storage tubes with a screw top, caps are available for use with either external or internal thread Azenta screw top tubes and are supplied in bags of 960 caps.

Practical Design Based on Experience of Applications

- High chemical resistance
- Broad operating temperature range -196°C to +121°C, do not use in liquid phase nitrogen
- Automation friendly, available in Azenta Cap Carrier for use with automated capping and de-capping systems
- Autoclavable
- Available in up to 10 different colors to aid sample identification
- Manufactured from high-quality virgin polypropylene

Screw Caps for External Thread Tubes

- Deforming compression seal is more effective than a silicone alternative
- Non-silicone seal means the cap can never be over-tightened
- Cap and tube manufactured from the same material, preventing differential expansion during freeze-thaw cycles
- Double-start thread engages in a maximum rotation of 180°, facilitating automation

Screw Caps for Internal Thread Tubes

- Co-molded screw caps for internal thread tubes eliminate the possibility of over-tightening
- Co-molded caps eliminate the failures common in caps using silicon O-rings
- Improved internal thread provides a more secure seal when using screw caps

Designed for Reduced Sample Loss

- Significant amounts of liquid can become trapped within a standard design screw cap, resulting in potential loss of valuable sample
- This hollow area on the cap can catch sample, and surface tension then makes cap removal difficult
- Studies demonstrate that Azenta cap design, with a reduced hollow, reduces this effect by lowering liquid retention whilst maintaining seal quality



AZENTA
LIFE SCIENCES

Ordering Information

96 Format External Thread Screw Caps

68-53111-10N	Cap Carrier (External Thread) , with orange caps, 96-format, SBS stackable, re-usable, 10 carriers/960 caps per case, suitable for all 96-format External Thread Screw Cap Tubes
68-53111-50N	Cap Carrier (External Thread) , with orange caps, 96-format, SBS stackable, re-usable, 50 carriers/4,800 caps per case, suitable for all 96-format External Thread Screw Cap Tubes
68-53100-Z1N	Screw Cap, 96-format, External Thread , white, bulk, 960 caps per case
68-53100-Z2N	Screw Cap, 96-format, External Thread , red, bulk, 960 caps per case
68-53100-Z3N	Screw Cap, 96-format, External Thread , yellow, bulk, 960 caps per case
68-53100-Z4N	Screw Cap, 96-format, External Thread , blue, bulk, 960 caps per case
68-53100-Z5N	Screw Cap, 96-format, External Thread , transparent blue, bulk, 960 caps per case
68-53100-Z6N	Screw Cap, 96-format, External Thread , orange, bulk, 960 caps per case
68-53100-Z8N	Screw Cap, 96-format, External Thread , green, bulk, 960 caps per case
68-53100-Z10N	Screw Cap, 96-format, External Thread , amber, bulk, 960 caps per case
68-53100-Z11N	Screw Cap, 96-format, External Thread , purple, bulk, 960 caps per case
68-53100-Z12N	Screw Cap, 96-format, External Thread , natural, bulk, 960 caps per case
68-53100-Z13N	Screw Cap, 96-format, External Thread , black, bulk, 960 caps per case

48 Format External Thread Screw Caps

65-7572	Screw Cap, 48-format, External Thread , orange, bulk, 480 caps per case, suitable for 48-format Cryo Tubes
65-7573	Screw Cap, 48-format, External Thread , red, bulk, 480 caps per case, suitable for 48-format Cryo Tubes
65-7574	Screw Cap, 48-format, External Thread , blue, bulk, 480 caps per case, suitable for 48-format Cryo Tube
65-7575	Screw Cap, 48-format, External Thread , green, bulk, 480 caps per case, suitable for 48-format Cryo Tubes
65-7576	Screw Cap, 48-format, External Thread , yellow, bulk, 480 caps per case, suitable for 48-format Cryo Tubes
65-7577	Screw Cap, 48-format, External Thread , purple, bulk, 480 caps per case, suitable for 48-format Cryo Tubes



Tube Screw Caps

96 Format Internal Thread Screw Caps

66-63100-Y1	Screw Cap, 96-format, Internal Thread, white, bulk, 960 caps per case
66-63100-Y2	Screw Cap, 96-format, Internal Thread, red, bulk, 960 caps per case
66-63100-Y3	Screw Cap, 96-format, Internal Thread, yellow, bulk, 960 caps per case
66-63100-Y4	Screw Cap, 96-format, Internal Thread, blue, bulk, 960 caps per case
66-63100-Y5	Screw Cap, 96-format, Internal Thread, transparent blue, bulk, 960 caps per case
66-63100-Y6	Screw Cap, 96-format, Internal Thread, orange, bulk, 960 caps per case
66-63100-Y8	Screw Cap, 96-format, Internal Thread, green, bulk, 960 caps per case
66-63100-Y10	Screw Cap, 96-format, Internal Thread, amber, bulk, 960 caps per case
66-63100-Y11	Screw Cap, 96-format, Internal Thread, purple, bulk, 960 caps per case
66-63100-Y12	Screw Cap, 96-format, Internal Thread, natural, bulk, 960 caps per case
66-63100-Y13	Screw Cap, 96-format, Internal Thread, black, bulk, 960 caps per case

24 Format External Thread Screw Caps

66-9401	Screw Cap, 24-format, External Thread, orange, automation friendly, bulk, 240 caps per case, suitable for 7.6ml PP Tubes
66-9402	Screw Cap, 24-format, External Thread, clear, automation friendly, bulk, 240 caps per case, suitable for 7.6ml PP Tubes



TPE Septum Caps



Developed to meet the needs of sample security, management and tracking in modern high-density storage applications, TPE septum caps are a disposable, thermo plastic elastomer (TPE) cap designed for optimal seal quality. Compatible with 96-format tubes, TPE caps are available for use with internal thread tubes and are manufactured from high-quality TPE, supplied as 96-cap mats or in bulk.



Practical Design Based on Experience of Applications

- Piercable cap, for use with any 96-format internal thread tube
- Broad operating temperature range, suitable for use -80°C to +110°C
- Ideal solution for -20°C storage with occasional access
- Available in a choice of five colors to aid sample identification
- Suitable for Automatic and Semi-Automatic Cappers and De-cappers
- Natural color only recommended for automatic systems
- Supplied in 96-format back mats to facilitate automation

Septum Cap Compatibility


- Suitable for 96-format Tri-coded Sample Storage Tubes, internal thread 0.65ml, 0.9ml
- Suitable for 96-format Dual-coded Tubes internal thread 0.3ml, 0.7ml, 0.9ml



Ordering Information


Internal Thread TPE Septum Seals

65-73000	TPE Septum Cap, natural, 96-format, on backing mat, 50 mats/4,800 caps per case, suitable for all Internal Thread Tubes
65-73001	TPE Septum Cap, blue, 96-format, on backing mat, 50 mats/4,800 caps per case, suitable for all Internal Thread Tubes
65-73002	TPE Septum Cap, green, 96-format, on backing mat, 50 mats/4,800 caps per case, suitable for all Internal Thread Tubes
65-73003	TPE Septum Cap, red, 96-format, on backing mat, 50 mats/4,800 caps per case, suitable for all Internal Thread Tubes
65-73004	TPE Septum Cap, yellow, 96-format, on backing mat, 50 mats/4,800 caps per case, suitable for all Internal Thread Tubes
65-74000	TPE Septum Cap, natural, 96-format, bulk, 960 caps per case, suitable for all Internal Thread Tubes



AZENTA
LIFE SCIENCES

SBS Cap Carriers



SBS-format screw cap carriers enable automated capping of 96-format sample storage tubes. Cap carriers are suitable for manual use or for use with Semi-Automated Screw Cap Decapper/Recapper, Single Channel systems and Automated Single Column Screw Cap Capper/Decapper systems.

Cap carriers are supplied in cases of 10 or 50, pre-filled with low-retention screw caps that can reduce sample loss (960 caps per case for 10 cap carriers, 4800 caps per case for 50 cap carriers).

Improved Process Efficiency

- Cap Carriers enable uncapped tubes to be used immediately and then capped using an automated capping and de-capping system
- Leads to improved process efficiency and allows batch filling of tubes prior to capping
- When used manually, by transferring caps from the carrier to the tube using a single tube transfer device, the operation becomes faster and the risk of contamination is reduced as cap handling is eliminated
- In alignment with our ESG commitments, we now manufacture Cap Carriers using at least 50% recycled polypropylene, reducing waste and reliance on virgin materials. While maintaining the same high-quality standards, this sustainable improvement is reflected in a color change from blue to black.

Ordering Information

68-53111-50X	Cap Carrier (External and Internal Thread) , empty, 96-format, SBS stackable, re-usable, 50 carriers per case, suitable for all 96-format External and Internal Thread Screw Cap Tubes
67-63111-10	Cap Carrier (Internal Thread) , with orange caps, 96-format, SBS stackable, re-usable, 10 carriers/960 caps per case, suitable for all Internal Thread Screw Cap Tubes
67-63111-50	Cap Carrier (Internal Thread) , with orange caps, 96-format, SBS stackable, re-usable, 50 carriers/4,800 caps per case, suitable for all Internal Thread Screw Cap Tubes
68-53111-10N	Cap Carrier (External Thread) , with orange caps, 96-format, SBS stackable, re-usable, 10 carriers/960 caps per case, suitable for all 96-format External Thread Screw Cap Tubes
68-53111-50N	Cap Carrier (External Thread) , with orange caps, 96-format, SBS stackable, re-usable, 50 carriers/4,800 caps per case, suitable for all 96-format External Thread Screw Cap Tubes
68-53111-10X	Cap Carrier (External and Internal Thread) , empty, 96-format, SBS stackable, re-usable, 10 carriers per case, suitable for all 96-format External and Internal Thread Screw Cap Tubes
66-9951	Cap Carrier, 48-format , empty, SBS, for use with Semi-Automated and Automated Screw Cap Decapper series, 10 carriers per case
65-9801	Cap Carrier, 24-format , empty, SBS stackable, 10 carriers per case, suitable for 24-format Automation Friendly External Thread Screw Caps (part number 66-9401)





AZENTA
LIFE SCIENCES

Treatment Services



Treatment Services

We understand how important it is to provide labware in a contamination-free condition and ideally suited for your research. This is why we provide alternative methods for treatment of products.

Standard Product

Prevention is better than a cure, so we manufacture all sample storage consumables in an ISO class 8 clean-room environment. To gain access to this area, everyone must wear a gown, gloves, face mask, hair net and overshoes. In addition, they must pass through a double door with an “air-shower” designed to remove any particles. No material that could cause contamination is allowed in the clean-room.

The products are treated at the point of molding, which takes place at temperatures that melt the resin – typically between 140°C to 150°C for copolymer and 165°C to 175°C for homopolymer.

From the injection molding machine, our consumables are processed and bagged in the clean-room. Only once they have left the clean-room are they placed in the transport box for shipment.

The efficacy of this process is tested periodically, which allows us to have a high degree of confidence that our products meet the standards outlined in the table below.

Contaminant	Testing Standard
Endotoxin (Pyrogen)	Product(s) are tested on a periodic basis and found to be below the acceptance level ≤ 0.05 EU/ml
DNA/RNA DNase/RNase	Product(s) are free of contaminations based on 3rd party electrophoresis evaluation of degradation
Heavy Metals	No heavy metal is contained in any of our product(s) that meet CONEG requirements of 100 ppm/weight
Animal (TSE/BSE)	No product(s) are manufactured from or come into contact with, animal materials
PCR Inhibitors	Product(s) are PCR-inhibitor free

Gamma Irradiation

Gamma irradiation is not guaranteed to destroy DNA contamination to the point where it doesn't amplify or interfere with subsequent analysis. Ensuring that no amplifiable DNA is present relies on (i) the manufacturing process and testing described above and (ii) fulfilling requirements for sensitive applications such as using Ethylene Oxide (EtO).

Gamma irradiation can negatively affect polypropylene, making it more brittle over time and increasing the binding characteristics of the plastic. For this reason, Azenta Life Sciences recommends that wherever possible, the standard product is relied on to deliver the contamination-free product required for research and long-term storage.

Ethylene Oxide (EtO)

EtO is just as effective as gamma irradiation and in addition, it ensures there is no DNA that can be amplified to a level where there is risk of subsequent interference. Therefore, EtO is the treatment method of choice in critical areas such as forensics.

Dual Ethylene Oxide (EtO2)

The newly published ISO 18385:2016 standard centers on minimizing the risk of human DNA contamination in consumable products used in collection, storage, or analysis of biological material for forensic DNA purposes. Dual EtO treatment is being recommended to further minimize the presence of PCR-amplifiable DNA

Ordering Information

Contaminant	Standard Product	Gamma Irradiation	Ethylene Oxide	Dual Ethylene Oxide
General Description	Product is made in an ISO 8 (class 100K) clean room. Product is endotoxin (pyrogen), DNAase/RNAase, heavy metals & animal-free.	Ensures no viable micro-organisms but cannot guarantee destroyed DNA contamination. Affects PP, making it more brittle over time & increasing binding of plastic.	Effective at ensuring no DNA can be amplified to a level where there is any risk of interference.	Recommended for consumables used in the collection, storage, or analysis of biological material for forensic DNA.
Recommended for standard use	Yes	No	No	No
Recommended for critical forensic applications (DNA amplification)	No	No	Yes	Yes
Modifies Polypropylene	NO	Yes	No	No
Product types	All	Capped or uncapped tubes	Uncapped tubes only, bulk caps, cap carriers	Uncapped tubes only, bulk caps, cap carriers
Packaging	Bulk, case of 10 racks	Sleeve of five racks as standard – individually wrapped as option	Individually wrapped	Individually wrapped
Min Quantity	None	10 cases	25 cases	25 cases
Part number	N/A	-GS	-S	-DS
Example part number	68-0703-02	68-0703-02-GS	68-0703-02-S	68-0703-02-DS



Recommended Temperature Range for Sample Storage Tubes

Cryogenic storage at very cold temperatures is designed to provide an indefinite, if not nearly infinite, longevity to biological materials.

By reducing sample temperatures to below the glass transition phase of water, all metabolic activity comes to a halt. Storage below this temperature therefore offers the most secure form of long term cryopreservation.

Storage in gas phase liquid nitrogen provides a convenient way of storing samples at temperatures below this transition temperature.

Screw cap storage tubes are designed for long term sample storage at cryogenic temperatures, including storage in the gas phase of liquid nitrogen.

Gas phase liquid nitrogen temperatures have been typically quoted as between -150°C and -178°C , depending on the location and the distance away from the liquid nitrogen reservoir surface. However, the latest generation of cryogenic storage systems are powered by liquid nitrogen in such a way to ensure a consistent temperature within the storage chamber, and can be as low as -195°C .

All screw cap storage tubes are fully compatible with this storage temperature and can be safely used for long term storage in this environment. However, immersion in liquid nitrogen is not recommended for two reasons:

1. Liquid nitrogen is extremely pervasive.

If a tube cap is not completely and correctly tightened nitrogen may seep inside the tube. When the tube is then removed from storage this liquid nitrogen will instantly boil, expanding to over 700 times its liquid volume. This creates a significant risk of explosion and/or biohazard risk from aerosols.

2. Many researchers have noted that immersion in liquid nitrogen carries with it a significant risk of cross contamination from bacteria, viruses and DNA that can be present within the liquid nitrogen.

Therefore, if required, snap-freezing is best achieved by immersing the tube in liquid nitrogen to a depth that avoids the nitrogen coming into direct contact with the tube cap.



CryoArc™ Deca with CryoPod® Carrier



CryoArc™ Pico Automated Sample Storage System



How Safe Are Your Samples? - Leachables and Extractables, Working Volume and Pressure Testing

Part I: Evaluation of External Thread Tubes for Potential Leachable Compounds

Most plastics are supplied pyrogen- and DNase/RNase-free, it is generally accepted that this guarantees product integrity. Few researchers ever consider that despite these apparent guarantees, plastic ware can still provide a potential source of error. Evidence shows that bioactive compounds can diffuse into solutions that come into contact with the surface of the plastic. These compounds, typically referred to as “leachables” or “extractables”, are used during the manufacturing process to improve product stability and durability.

The aim of this evaluation was to determine if “extractables” could be detected in solvent solutions stored in Azenta external thread tubes and to compare the performance against competitor tubes.

Part II: Working Volume

Knowledge of the working volume is a key criterion to consider when selecting the most appropriate tube and, is dependent on a range of factors including: Fill volume of the tube; Accuracy of the volume dispensed; Freezing conditions; Thawing conditions; Cap selected; Burst pressure of the cap.

Manufacturers of sample storage tubes specify the size of tubes in a variety of ways, but rarely state the sample type or storage conditions required to achieve the volumes stated. This can lead to catastrophic results if the working volume of the tube is exceeded.

Part III: Sample Integrity Through Secure Capping

When using sample tubes, secure capping is paramount to protect sample integrity during the storage and handling process for a range of applications. These applications could include; Biobanking, Compound Management, Cell Therapy, Benchtop Research or applied industry. Two major hazards that can be mitigated through safe and secure sample tube capping are Evaporation and Cross contamination.



White Papers Available



AZENTA
LIFE SCIENCES



AZENTA
LIFE SCIENCES

Ziath 2D-Coded Tube Reading and Handling Systems for Effective Sample Management



Uno and Duo Single Tube Readers

**1D and 2D-code reader
for individual tubes**



Key Features

- Easy set-up, straight out of the box
- Reads 2D datamatrix coded tubes, including Azenta Acoustic tubes (with Acoustic add-on)
- Uno features higher resolution camera that can often decode older, poorly printed codes
- Duo combines the datamatrix reader with a linear barcode reader for rack codes or 1D linear codes along the edge of some tubes
- Super-fast electronics capable of reading 2D coded tubes in multiple formats including 384-well tubes
- Large imaging windows make it easy to present the tubes

Uno and Duo single tube readers

Single tube readers offer easy set-up, straight out of the box, with only a few minutes needed to install the application software. The instruments are provided already calibrated, able to read 2D datamatrix coded tubes, including Azenta Acoustic tubes. The new Uno single tube reader features a higher resolution camera than the previous model that can often decode older, poorly printed codes easier.

The Duo reader combines the datamatrix reader with a linear code reader to decode rack or 1D-codes.

Specifications

- Dimensions (LxWxH) 66 mm x 53.5 mm x 51.5 mm
- Single USB connection is all that is required for power and data communication



AZENTA
LIFE SCIENCES

Easy to set up and easy to use

With the Uno measuring just 66 mm x 53.5 mm x 51.5 mm (LxWxH), these rugged new tube readers feature a high resolution camera with super-fast electronics capable of reading 2D-coded tubes in multiple formats including 384-well tubes.



Simply hold your tube to the imaging window and the code is instantly displayed. The readers can be supplied with or without optional cryoprotection to enable trouble-free reading of tubes taken straight from cold storage. These USB 3 devices act as a Microsoft® Keyboard Wedge and insert the code into any open programme on the desktop. The enhanced drivers allow the Uno or Duo to be ready to read immediately with no lag time and even allow adjustment of the camera for changes in contrast between different manufacturers’ tubes. The software for Uno and Duo offers a keyboard wedge mode which allows you to enter the decoded data directly into your application. The large imaging windows make it easy to present the tubes, even when wearing cryoprotective or anti-contamination clothing. A single USB connection is all that these single tube readers need for power and data communication.

The Cryoprotection™ option reduces condensation caused by cold racks. It uses no heat or fans and will not damage sensitive cold samples being read. The formation of condensation on the instrument window whilst reading tubes straight from cryogenic storage is a well-known problem.

Cryoprotection is a special coating that eliminates the condensation and ensures uninterrupted reading. This method for eliminating condensation uses no extra electrical components or heating, resulting in no risk of accidentally thawing your samples, increased reliability and ensuring that the cost of this option is at acceptable levels.

Ordering Information

ZTS-UNO-80	Uno 2D Single tube reader with Cryoprotection
ZTS-DUO-80	Duo 2D Single tube reader with 1D linear code reader with Cryoprotection

Mirage Rack Reader

**Setting new standards for
2D-code SBS rack readers**



Key Features

- The most compact camera-based 2D rack reader available
- Low profile allows use on robotic liquid handling systems
- Innovative and unique patented reading technology
- Reads commonly used SBS format racks
- Economically priced

Included with the Mirage

- Mirage instrument
- Power adaptor and country relevant power cable
- DP5 software (download link supplied)
- USB lead to connect instrument to host computer
- User guide (download link supplied)

Easy to set-up and easy to use


The Mirage 2D reader offers easy, out-of-the-box set up and is delivered pre-calibrated and ready to read.

This instrument can read commonly used racks in the SBS format, and its low profile means that it can be easily added into liquid handling platforms and other automation systems.

The low profile and economic price makes the Mirage the perfect replacement for a traditional flat bed scanner.



Options



Cryoprotection

Prevents the window from misting up when a cold rack is read. Our passive technology does not use heat or blown air.



Cyclops 1D Linear Reader

Decodes the 1D-code on the side/edge of your racks.



AZENTA
LIFE SCIENCES

Mirage Rack Reader



Technical Specification

Dimensions (D x W x H)	208 x 135 x 80 mm
Weight	1.35 kg
Power adapter	Input 100-240VAC; output 5Vdc 4A
Power consumption	≤10 watts
Labware compatibility	24, 48, 96, 240 position racks
Imaging technology	CMOS camera
Read time	Approx. 2 seconds
Code compatibility	Datamatrix 2D codes
Working temperature	5 to 35 °C
Data export	Export data as .csv, .text, Excel®, .XML, .JSON, Python, .DB, .PNG or email
Operating system on host PC	Windows® 10, 11
Warranty	2 years
Software	Includes remote control capability for integration applications

Ordering Information

DP5-M-STD	Mirage Standalone and DP5 Standard; includes Mirage instrument, power adaptor and country relevant power cable, DP5 software (download link supplied), USB lead to connect instrument to host computer, and user guide (download link supplied); 1 reader
DP5-M-STD-80	Mirage Standalone and DP5 Standard; includes Mirage instrument, power adaptor and country relevant power cable, DP5 software (download link supplied), Cryoprotection, USB lead to connect instrument to host computer, and user guide (download link supplied); 1 reader
DP5-M-NET	Mirage LAN and DP5 Network; includes Mirage LAN instrument, power adaptor and country relevant power cable, DP5 Network software (download link supplied), USB lead to connect instrument to host computer, and user guide (download link supplied); 1 reader
DP5-M-NET-80	Mirage LAN with Cryoprotection and DP5 Network; includes Mirage LAN instrument, power adaptor and country relevant power cable, DP5 Network software (download link supplied), Cryoprotection, USB lead to connect instrument to host computer, and user guide (download link supplied); 1 reader
DP5-CYS-M	Cyclops 1D linear reader for Mirage



Express Rack Reader

Small, automation-friendly tube and rack code reader

Key Features

- Smallest footprint 2D-code rack reader
- Fast camera-based reader
- Ideal for integration into large automated systems
- Available with either USB or RJ45 connector
- Capable of reading SBS Rack formats
- Patented technology allows imaging under difficult lighting conditions
- Cryoprotection™ and 1D rack code reading options

Easy to set-up and easy to use

The Express 2D reader offers easy, out-of-the-box set up and is delivered pre-calibrated and ready to read. This instrument can read racks in the SBS format, regardless of the configuration or number of tubes, and is supplied with all the features needed for robotic integration.


Our new automatic rack type detection feature means that you no longer have to configure the Express for all rack types. Simply place your rack on the imaging window and the software will work out what type of rack it is and read it - easy!



Included with the Express


- Express instrument
- Power adaptor, power cable
- DP5 software (download link supplied)
- USB3 or RJ45 lead
- User guide (download link supplied)

Options



Cryoprotection

Prevents the window from misting up when a cold rack is read. Our passive technology does not use heat or blown air.



Cyclops 1D Linear Reader

Decodes the 1D-code on the side/edge of your racks.



AZENTA
LIFE SCIENCES

Express Rack Reader



Technical Specification

Dimensions (D x W x H)	Imager: 85 x 185 x 95 mm Controller: 95 x 195 x 55 mm
Weight	1.4 kg
Power adapter	Input = 110-220V +/- 10V AC; output = 5V DC, 4A
Labware compatibility	24, 48, 96, 240, 384 position SBS racks
Code compatibility	Datamatrix 2D-coded tubes
Data export	Export data as .csv, .text, Excel®, .XML, .JSON, Python, .DB, .PNG or email
Operating system on host PC	Windows® 10, 11
Warranty	2 years
Software	Includes remote control capability for integration applications

Ordering Information

DP5-E-STD	Express Standalone and DP5 Standard; includes Express instrument, power adaptor, power cable, DP5 software (download link supplied), USB lead to connect instrument to host computer, user guide (download link supplied); 1 reader
DP5-E-STD-80	Express Standalone with Cryoprotection and DP5 Standard; includes Express instrument, power adaptor, power cable, DP5 software (download link supplied), Cryoprotection, USB lead to connect instrument to host computer, user guide (download link supplied); 1 reader
DP5-E-NET	Express LAN and DP5 Network; includes Express LAN instrument, power adaptor, power cable, DP5 Network software (download link supplied), USB lead to connect instrument to host computer, user guide (download link supplied); 1 reader
DP5-E-NET-80	Express LAN with Cryoprotection and DP5 Express; includes Express LAN instrument, power adaptor, power cable, DP5 Network software (download link supplied), Cryoprotection, USB lead to connect instrument to host computer, user guide (download link supplied); 1 reader
DP5-CYS-E	Cyclops 1D linear code reader for Express; 1 side reader



Cube Rack Reader

**2D tube and rack
readers for cryoboxes
and SBS racks**



Key Features

- Fast camera-based 2D reader
- Capable of reading SBS racks and cryobox formats
- Patented technology allows imaging under difficult lighting conditions
- Cryoprotection and 1D rack/box decoding option


Included with the Cube

- Cube instrument complete with masks for racks and boxes
- Power adaptor, power cable
- DP5 software (download link supplied)
- USB lead to connect instrument to host computer
- User guide (download link supplied)

Easy to set-up and easy to use


The Cube offers fast identification of racked, coded sample storage tubes, without the need to remove tubes from racks. With easy, out-of-the-box set up, the Cube is delivered pre-calibrated and ready to read. Able to read both cryoboxes and SBS format racks (using a simple adaptor, included as standard), regardless of the configuration or number of tubes, the new automatic rack type detection feature means that you no longer have to configure the Cube for all rack types. Simply place your rack on the reading window and the software will work out what type of rack it is and read it - easy!

Options



Cryoprotection

Prevents the window from misting up when a cold rack is read. Our passive technology does not use heat or blown air.



Cyclops 1D Linear Reader

Decodes the 1D-code on the side/edge of your racks.



AZENTA
LIFE SCIENCES

Cube Rack Reader



Technical Specification

Dimensions (D x W x H)	152 x 174 x 146 mm
Weight	1.2 kg
Power adapter	Input = 110-220V +/- 10V AC; output = 5V DC, 4A
Labware compatibility	24, 48, 96, 240, 384 position SBS racks 81, 100, 196 position cryoboxes
Code compatibility	Datamatrix 2D-coded tubes
Data export	Export data as .csv, .text, Excel®, .XML, .JSON, Python, .DB, .PNG or email
Operating system on host PC	Windows® 10, 11
Warranty	2 years
Software	Includes remote control capability for integration applications


Ordering Information

DP5-C-STD	Cube Standalone and DP5 Standard; includes Cube instrument complete with masks for racks and boxes, power adaptor, power cable, DP5 software (download link supplied), USB lead to connect instrument to host computer, user guide (download link supplied); 1 reader
DP5-C-STD-80	Cube Standalone with Cryoprotection and DP5 Standard; includes Cube instrument complete with masks for racks and boxes, power adaptor, power cable, DP5 software (download link supplied), Cryoprotection, USB lead to connect instrument to host computer, user guide (download link supplied); 1 reader
DP5-CYS-C	Cyclops 1D linear code scanner for Cube; 1 side reader

Mirage RITrack RFID Reader

A revolution in cryogenic sample management – locate tubes with RFID even if the rack is covered in ice

Azenta Life Sciences is proud to present an era-defining change in the way that frozen samples stored in 2D-coded tubes can be managed. Overcoming the problems of ice on the base of tubes and racks using novel RITrack radio-frequency tagging eliminates the frustration of racks and tubes that cannot be read due to excessive frosting.




Key Features

- Stores individual tube identities, up to 12 characters, for up to 96 tubes in a rack
- Works with racks covered in ice and with racks stored in vapor-phase LN2 – no warming or thawing necessary – protects thermally sensitive samples
- No need to buy special racks; Tags fit most commercially available SBS racks – lowers your cost of operation
- Tags last and function for the whole life of your racks – increase your ROI and lower your cost of ownership
- Quick and easy read out with any smartphone through a simple App – improves productivity in the lab
- Write and re-write unlimited times to the same tag with a RITrack Mirage – quick and convenient whilst keeping costs under control
- Unique RFID technology exclusive to Azenta Life Sciences readers – take advantage of the smarter solution to cryogenic sample management

Extensive testing has demonstrated that RITrack RFID technology can endure over 10,000 cycles of immersion in vapor-phase liquid nitrogen without data loss or degradation. This makes RITrack ideal for all forms of long-term cryogenic storage.

The RITrack system is compatible with standard SBS format laboratory consumables, eliminating the need for costly dedicated racks. Simply insert a RITrack tag into your existing racks using our convenient tag insertion tool. The tags fit securely within the internal space of most SBS tube racks and are firmly held in place, ensuring they remain secure, and providing added peace of mind.

Standard laboratory consumables – in the form of SBS format tube racks – can be used with the RITrack system. There is usually no need to buy expensive dedicated racks – just insert a RITrack tag into your existing racks using our handy tag insertion tool. RITrack tags will fit securely into the internal space of most SBS tube racks. Once fitted, they are securely held in place and almost impossible to remove giving added peace of mind that the tag will stay in position for the whole lifetime of the rack.



AZENTA
LIFE SCIENCES

Mirage RITrack RFID Reader

Revolutionizing the handling of frozen samples stored in 2D-coded tubes, RITrack utilizes innovative radiofrequency tagging to eliminate the challenges posed when no line of sight is clear for 2D coded tubes, such as ice formation on tube and rack bases. As RITrack tags use standard RFID technology, they can also be read by any smartphone device with Near Field Communication (NFC). An easy-to-use app can be downloaded free from the App Store or Google Play to any compatible smartphone. The app allows the tag information to be read and the data displayed in either graphical or tabular format. This makes it easy to identify any individual tube in the frozen rack together with its location. Simply open the freezer, withdraw a rack and offer it up to the phone to see all the available tubes in that rack, eliminating the inefficiency of taking one or more cold racks back to the lab.

Using the Pick List function of the app, tubes from several different racks can be checked-off right in front of the freezer. Simply open the freezer, withdraw a rack and offer it up to the phone to see all the available tubes in that rack - again saving valuable technician time. With RITrack, data for every individual tube as well as the rack barcode can be read or stored to the tag in a matter of seconds. This makes it very quick and easy to manage frozen samples without thawing or removing ice, thus saving valuable time and protecting delicate thermally sensitive materials in the stored tubes.



Ordering Information

DP5-M-RIT-80	Mirage RITrack reader with special RFID sensor, DP5 Software, Cryoprotection
DP5-NET2	Upgrade DP5 Standard to a full DP5 Network license, for Mirage RITrack
DP5-RIT-TAG	RITrack RFID Rack Tags, with magnetic RAM for unlimited re-writes, includes 1 tag applicator and 12 tags per case
DP5-RIT-APL	RITrack tag applicator

For full details of DP5 software and Mirage code readers, please see page 74.

Handheld 2 Tube Reader

The Ziath Handheld 2 reader is a truly portable reader which can be held and operated in one hand; leaving the other free for tube handling.

Key Features

- Compatible with 2D-coded tubes
- Portable, independent, and compact
- No external computer or tablet needed

Easy to Set-Up and Easy to Use

The Handheld 2 Tube Reader requires no external computer and is easily operated with one hand; leaving your other free for precise tube handling. The software is pre-loaded and boots up immediately with a battery life of up to a week with casual use. This device allows you to read tubes standing in the basement at the freezer or liquid nitrogen tank. Even out in the field acquiring samples.

Operation Modes

The reader has four operation modes:

- **Code reading** – Read a tube code to view it on the screen
- **Pick list** – Enter a pick list using the associated Client and read the tubes with the handheld device to ensure you have the correct tubes while at the freezer
- **Data View** – Read a code to see the details of the sample in the tube
- **Data Entry** – Read a code and enter the information on the same to store



Client App

To assist with the transfer of data on and off the handheld so that it can be used in environments with no computer network available, the handheld comes with an associated client application to help define the data to store and load/unload it onto the handheld. This allows you to go into the deepest darkest corners of your sample storage and still be able to access your data!

Technical Specification

Dimensions (D x W x H)	155 x 65 x 38 mm
Weight	150 g
Battery charger input	100-240VAC 60/50Hz, 260mA
Battery charger output	5VDC, 2.4A
Battery life	Typically 8 hours between charges
Data storage on instrument	Typically up to 60,000 sample

Ordering Information

ZTS-HND2	Handheld 2 Single Tube 2D Barcode Reader; 1 reader
----------	--



AZENTA
LIFE SCIENCES



Cyclops 1D Linear Code Reader

Key Features

- Use with Ziath 2D-code rack readers – Mirage, Express, and Cube
- Reads in milliseconds
- Wide range of placement options on each instrument
- Quick & easy installation, USB connection to operating laptop or PC
- Attaches to the rack reader using dedicated robust bracket
- Compatible with Windows® 10 & 11



Essential Tool for Sample Tracking

The Cyclops reader is designed to illuminate and read linear 1D-codes printed on the end of an SBS rack of 2D-coded tubes. Once the rack is read, the Cyclops reports the rack identity alongside the individual tube identities when the full rack is read in conjunction with a Ziath Tube and Rack Reading Instrument. For use in conjunction with Ziath Tube and Rack Readers - Mirage, Express, or Cube - each instrument has a dedicated, easy-to-install bracket that attaches securely to the compatible reader using the instrument's feet, while the Cyclops device slides into position and is held firmly in place by a thumb-wheel.



Easy to set-up and easy to use

The Cyclops Linear Rack Barcode Scanner is very simple to attach to existing 2D readers, and our software makes it simple to export the linear (rack) barcode data and the tube data to Excel®, text, csv, XML, JSON or Python. The linear rack code takes less than one second to read/decode, and the Cyclops is powered and controlled by USB cable, so needs no additional power cables.

Technical Specification

Dimensions (W x D x H)	30 x 25 x 80 mm
Weight	100 g (without bracket)
Bracket	Requires an extra 8 cm ahead of each instrument
Micro-USB to USB cable	included
Warranty	2 years

Ordering Information

DP5-CYS-M	Cyclops Advanced Rack Code 1D Scanner for Mirage); 1 side reader
DP5-CYS-E	Cyclops Advanced Rack Code 1D Scanner for Express; 1 side reader
DP5-CYS-C	Cyclops Advanced Rack Code 1D Scanner for Cube; 1 side reader




AZENTA
LIFE SCIENCES

DP5 Decoding Software

DP5 represents the very latest high-speed communications and data processing technology applied to laboratory sample tracking and identification in an easy to use and scalable format that can adjust to the needs of your individual laboratory. Combined with the unrivalled accuracy and speed of Ziath camera-based readers, tube identification and sample tracking has never been easier or more convenient.


- Runs from a laptop, desktop, phone, tablet
- Offers data encryption and multi-level user authority management



Key Features

The Very Latest Code Reading Technology

- DP5 is the latest version of Ziath's market leading 2D rack decoding software
- Using the latest technologies DP5 can be run from a laptop, desktop, phone, tablet or any other network connected device
- The same robustness, reliability and speed are there but now you can work unchained from your lab bench
- Using high-speed, web-based, communications protocols allows DP5 to run just one scanner in your lab, or multiple scanners connected over a local area network
- In the advanced DP5 Network application, scanners can be controlled directly from Wi-Fi connected smart phones, tablets or laptops
- Ziath's DP5 Secure application offers data encryption and multi-level user authority management as part of the package
- Also included along with a full audit trail is the revolutionary ScanStore feature in which every scan is stored with both image and tabular data for later review or audit



Versions

DP5 Standard

DP5 Standard simplifies day-to-day code reading operation through an intuitive and uncluttered user interface that is ready to go in just minutes, straight out of the box.

- Runs like any other desktop application from a connected PC or laptop
- Autoscaling allows DP5 Standard to work with the latest 4K displays
- Runs all Ziath camera based readers: Mirage, Cube & Express
- Export data directly to a database, text, PDF, XML, JSON or Excel
- Remote control function allows direct integration with popular liquid handling robots, compound stores and other external compatible devices
- Reader and export settings can be assigned individually for each type of rack in use
- Automatic rack type detection and rack template creation as standard

DP5 Network

DP5 Network unchains you from the bench, allowing a scanner to operate anywhere on your LAN from your own tablet or phone, even right alongside your freezers.

- Advanced operation via a LAN allows scanners to be controlled directly from a Wi-Fi connected phone, tablet, desktop PC or laptop
- Control multiple scanners on one server from one hand-held device, laptop, tablet or desktop computer
- Scripting feature allows high level of post-run data processing automation – send files to LIMS, e-mail, Excel or any other compatible program
- Take a reader anywhere on site, run a scan on a hand-held device and send that file back to any computer in a variety of formats



AZENTA
LIFE SCIENCES

Sample Tube Capping & Sealing Systems



Tube Capping and Decapping Systems

Azenta provides dependable manual, semi- and fully automated capping and de-capping solutions. The range incorporates products from single-tube manual de-cappers, through column cappers / de-cappers to fully automated and whole-rack solutions. All products are the result of over 15 years of continuous development and innovation, bringing the best solutions, with the most useful features, to the sample storage market.

Semi-Automated Screw Cap Single Tube Decapper/Recapper

The Semi-Automated Screw Cap Single Tube Decapper/Recapper is the first-of-its-kind, single tube capper/de-capper in the industry. It provides a low-cost, reliable solution for managing the automated capping and de-capping of single and larger sized master sample tubes such as: centrifuge, cell culture and even glass tubes. Unlike many products on the market, where an engineering change of the gripper is required in order to tighten different tube types, the Semi-Automated Screw Cap Single Tube Decapper/Recapper allows you to manage a variety of tube types with a simple, interchangeable cap-driver system which can be changed in 5 seconds.

A flexible, faster and lower cost option for managing the capping/de-capping of larger sample tubes in the laboratory, this Semi-Automated Screw Cap Single Tube Decapper/Recapper is compatible with Azenta sample storage tubes, as well as a wide variety of tubes, with internal and external threads, from almost all major manufacturers.

Our Semi-Automated Screw Cap Single Tube Decapper/Recapper is designed to improve productivity while reducing the risk of repetitive strain injury.



Key Features

- Fast and consistent single tube capping/decapping
- Perfect for automating the handling of larger sized master sample tubes such as: centrifuge, cell culture or even glass tubes
- Interchangeable cap driver kit for any tube format (9mm-40mm)
- Change the cap driver in 5 seconds
- Easy, single button, multi-torque control
- Suitable for use inside a laminar flow cabinet
- Low cost, automated solution
- Please contact us to discuss your specific tube requirements

Automation for Fast and Consistent Capping Performance

- The same torque is applied to each cap, according to tube manufacturer specification, in order to ensure a consistent seal quality
- Automation reduces sample handling time, increasing sample throughput and workflow
- Screw caps can be held within the Semi-Automated Screw Cap Single Tube Decapper/Recapper unit whilst the user fills, or accesses, the sample tube
- Caps can subsequently be re-applied to the same tubes, eliminating the risk of cross-contamination

Broad Compatibility for Every Laboratory

- The Semi-Automated Screw Cap Single Tube Decapper/Recapper is compatible with sample storage tubes from almost every manufacturer, as long as the screw cap can be gripped with cap drivers and the cap has a diameter from 9mm up to 40mm
- By adjusting the tube pitch required for each tube manufacturer, the Semi-Automated Screw Cap Single Tube Decapper/Recapper ensures that caps are not cross threaded during the capping cycle



Multiple Torque Control and Interchangeable Cap Driver

- Quick, user-friendly, interchangeable cap driver
- Cap driver can be changed in less than 5 seconds
- Torque is adjusted to fit the majority of tubes, allowing users to determine the optimal setting for their tubes
- 3 different torque settings: Low, Medium and High
- Easy to Use

Technical Specification

Dimensions	333(L) x 110(W) x 300(H) mm
Weight	6.0 kg
Power Requirements	AC 100 -240V, 50/60 Hz, 26W
Cycle Time De-cap	< 3 sec per tube
Cycle Time Cap	< 3 sec per tube
Compatibility	Standard sized tubes with screw cap (verify gripper list)

Ordering Information

46-6001	Semi-Automated Screw Cap Decapper, 1-channel, with interchangeable cap drivers and multiple torque settings
---------	---

Ordering Information - Grippers

46-6002-1	Gripper kit, TRP tube 50ml Cap (Fujifilm)
46-6002-2	Gripper kit, BD Falcon tube 50ml Cap
46-6002-3	Gripper kit, BD Greiner tube Cap
46-6002-4	Gripper kit, Sarstedt 0.5 & 1.5 & 2.0ml Cap
46-6002-5	Gripper kit, Azenta 6 & 10ml Automation-Friendly Cap
46-6002-6	Gripper kit, Corning 21mm Cap
46-6002-7	Gripper kit, Azenta 6.0 & 10ml External Cap
46-6002-8	Gripper kit, Azenta 4.0ml Glass Jacket Tube Cap
46-6002-9	Gripper kit, Ø 50mm Cap for glass container
46-6002-10	Gripper kit, Azenta / Greiner Ext & Sarstedt Internal Cap
46-6002-11	Gripper kit, Greiner/Corning 1.2ml – 5.0ml Cap
46-6002-12	Gripper kit, Fisher Custom Tube Cap dia. 16mm
46-6002-13	Gripper kit, Falcon 15mm dia. centrifuge tube
46-6002-14	Gripper kit, Azenta 96-format External & Internal Co-Molded Caps
46-6002-15	Gripper kit, Sarstedt 12ml screwcap tube (60.9922.937); dia. 16mm
46-6002-16	Gripper kit, 3.5ml 14mm dia. False Bottom MarketLab Tube (Inpeco tube, EXT)
46-6002-17	Gripper kit, Glass Vials, 4ml (dia. 15mm,length 48mm, GNF)
46-6002-18	Gripper kit, 50ml, green cap, dia. 35.8mm, customized tube for optimum processing



Semi-Automated Handheld Screw Cap Decapper, 8-channel

The Semi-Automated Handheld Screw Cap Decapper, 8-channel provides users with a flexible solution for capping and de-capping 96 format screw caps for low-to-mid throughput environments, with easily interchangeable cartridges to switch between tube types. The 8-channel hand-held capper and de-capper offers consistent sealing at an affordable price.



Key Features

Designed with the User in Mind

- Lightweight, hand-held, semi-automatic single column-based capper and de-capper
- Compatible with 96-SBS rack format screw capped tubes
- Caps a single column in under 4 seconds; caps/de-caps a complete rack of 96 tubes in less than 90 seconds
- All caps are tightened to the same torque to create a secure seal
- Battery operation allows capping/de-capping for more than 40 complete racks
- Docking station for instrument placement/storage and charging also provides contamination-free handling of caps
- Single 'action' button operation to de-cap, re-cap and eject caps
- Option to eject caps onto cap carrier using second 'eject' button



Flexibility for Ease of Use

- Designed to be used by both left and right-handed users
- Unit can be operated by battery charge or powered through AC plug in
- Easily transportable
- Interchangeable cartridges, easy to switch between cap types

Get in touch with your local representative to discuss your labware requirements with us.

Ordering Information

46-9012	Semi-Automated Handheld Screw Cap Decapper, 8-channel , Includes decapper (46-9008), docking station (46-9001), set of power cables, operation manual, cartridges for Azenta Internal Thread (48-9013-01) and External Thread (48-9013-02) for 96-SBS format racks
Accessories	
46-9001	Docking Station with set of power cables
Compatible cartridges	
48-9013-01	Cartridge for use with Azenta internal thread , screw capped tubes
48-9013-02	Cartridge for use with Azenta external thread , screw capped tubes
48-9013-03	Cartridge for use with Azenta Acoustic tubes
48-9013-04	Cartridge for use with Micronic external thread , screw capped tubes
48-9013-05	Cartridge for use with Matrix internal thread , screw capped tubes
48-9013-06	Cartridge for use with Rhinostics Swabs Internal thread
48-9013-07	Cartridge for use with Micronic external thread IG , screw capped tubes



AZENTA
LIFE SCIENCES

Semi-Automated Screw Cap Decapper/Recapper

The Semi-Automated Screw Cap Decapper/Recapper systems are compact, bench top units designed for efficient tube capping in labs with medium throughput. Offering the consistency of an automated de-capping system, but at much lower cost, the Semi-Automated Screw Cap Decapper/Recapper can cap a single column of tubes, from a cap carrier, in under 10 seconds and will cap, or de-cap, a complete rack of 96 tubes in under 2 minutes.

Key Features

Flexible Product Options to Suit a Range of Tube Types



*Semi-Automated Screw Cap
Decapper, 4-channel*



*Semi-Automated Screw Cap
Decapper, 6-channel*



*Semi-Automated Screw Cap
Decapper, 8-channel*

- **4-channel capping and de-capping system** compatible with screw-capped tubes in 24-format SBS racks
- **6-channel capping and de-capping system** compatible with screw-capped tubes in 48-format SBS racks
- **8-channel capping and de-capping system** compatible with screw-capped tubes in 96-format SBS racks

Your choice of model can be configured for use with one of the following tube types:

- Azenta internal and external thread screw-capped tubes
- Thermo-Matrix internal thread screw-capped tubes
- Micronic internal and external thread screw-capped tubes
- Thermo-Nunc Bank-IT tubes
- LVL external thread screw capped tubes

Please contact your local representative for the latest list of supported tubes

Fast and Consistent Sealing Performance

- Has a set-down position that allows the simple insertion of screw caps using an Azenta Cap Carrier
- Caps a single tube in under 10 seconds
- Caps or de-caps a complete rack of 96 tubes in less than 2 minutes
- Individually spring-loaded cap drivers mean easy and consistent capping
- All caps tightened to the same torque, delivering a secure seal and peace of mind
- Designed for ease of use by both left and right handed people



AZENTA
LIFE SCIENCES

Semi-Automated Screw Cap Decapper/Recapper

Space Saving Design

With a small footprint, the semi-automated screw capper decapper/recapper fits easily onto the bench top and is easy to install in a laminar flow cabinet.

Weight: 6 kg, Dimensions: 310mm (w) x 345mm (l) x 345mm (h)

Ordering Information

Semi-Automated Screw Cap Decapper, 8-channel for use with 96-format racked tubes	
46-6501	Semi-Automated Screw Cap Decapper, 8-channel, compatible with Azenta Life Sciences internal thread tubes, in 96-SBS format rack
46-6502	Semi-Automated Screw Cap Decapper, 8-channel, compatible with Azenta Life Sciences external thread tubes, in 96-SBS format rack
46-6601	Semi-Automated Screw Cap Decapper, 8-channel, compatible with Matrix and Thermo internal thread tubes, in 96-SBS format rack
46-6602	Semi-Automated Screw Cap Decapper, 8-channel, compatible with Micronic internal thread tubes, in 96-SBS format rack
46-6604	Semi-Automated Screw Cap Decapper, 8-channel, compatible with LVL external thread tubes, in 96-SBS format rack
46-6606	Semi-Automated Screw Cap Decapper, 8-channel, compatible with Micronic external thread tubes, in 96-SBS format rack
Semi-Automated Screw Cap Decapper, 6-channel for use with 48-format racked tube	
46-6511	Semi-Automated Screw Cap Decapper, 6-channel, compatible with Azenta Life Sciences external thread tubes, in 48-SBS format rack
46-6512	Semi-Automated Screw Cap Decapper, 6-channel, compatible with Azenta Life Sciences internal thread tubes, in 48-SBS format rack
46-6513	Semi-Automated Screw Cap Decapper, 6-channel, compatible with Nunc external thread tubes in Greiner racks, in 48-SBS format rack
46-6605	Semi-Automated Screw Cap Decapper, 6-channel, compatible with LVL external thread tubes, in 48-SBS format rack
Semi-Automated Screw Cap Decapper, 4-channel for use with 24-format racked tube	
46-6521	Semi-Automated Screw Cap Decapper, 4-channel, compatible with Azenta Life Sciences external thread tubes, in 24-SBS format rack



IntelliXcap™

IntelliXcap™ is the next-generation capper and de-capper, engineered for increased throughput, ease of use and high system reliability, making the system ideal for any laboratory managing compound libraries or biological sample stores.

Key Features

Advanced Automation Means Increased Throughput

- IntelliXcap is extremely fast and works with multiple sample storage tube types from a range of tube manufacturers in 24, 48 and 96 format including taller 96 format tubes
- Capable of de-capping a complete rack of 96 tubes in as little as 20 seconds
- Automation and high speed reduces sample handling time, increasing sample throughput workflow
- Ideal for medium to high-throughput laboratories managing compound libraries, biobanks or other biological samples
- Specific model available with verification camera, designed for use with Azenta Life Sciences Acoustic Sample Tube – Echo Qualified Consumable

Compact, Modular Design Means Ultimate Flexibility


- IntelliXcap features the unique, fully automated, interchangeable cartridge cap driver system
- Cartridges allow change between different pre-configured tube types in less than 2 minutes
- Simply by fitting the relevant cartridge, IntelliXcap is compatible with all Azenta sample storage tubes, as well as sample storage tubes from Thermo Matrix, Thermo Nunc, LVL and Micronic
- Cartridges are available for both internal and external thread screw-cap tubes
- Control is via a quick installation, easy-to-use touch panel
- Light Curtain System detects the height of tube rack on the stage reducing the risk of damage to tubes, samples or the instrument itself that could be caused by failed de-capping / capping or the use of incorrect consumables



Integration Friendly



Capper & De-capper for 96, 48 and 24 format Screw Capped Tubes



AZENTA
LIFE SCIENCES

Easy to Use, Easy to Integrate

- Quick installation, easy-to-use touch panels bring instant control
- Wide operating temperature range of 5°C to 40°C
- Electric torque control means less persistent wear on cap drivers
- Easily integrated, sample rack stage extension facilitates fully automated workflow integration

Driver and Cap Compatibility

- We provide a wide range of custom made and off the shelf driver cartridges to accommodate your specific labware needs
- The modular design of IntelliXcap™ allows you to quickly and easily change the cap driver cartridge and seamlessly swap between tube types, making the instrument the ideal capper/de-capper for your automated high throughput workflows
- Due to its flexible cartridge design, with one IntelliXcap instrument you can cap, de-cap and re-cap both internal and external thread tubes
- Please contact us to discuss your specific tube requirements

Weights and Dimensions

	Weight (kg)	Height (mm)	Width (mm)	Depth (mm)
IntelliXcap 96	26.80	320	256	468
IntelliXcap 48	28.74	386.2	256	464.4
IntelliXcap 24	27.96	386.2	256	464.4
IntelliXcap 96 Acoustic	27	316	256	634
IntelliXcap Extended Height	24	370	256	464

Ordering Information

46-8012	IntelliXcap Automated Screw Cap Decapper, 96-format, compatible with validated sample storage tubes and cartridges
46-8011	IntelliXcap Automated Screw Cap Decapper, 48-format, compatible with validated sample storage tubes and cartridges
46-8010	IntelliXcap Automated Screw Cap Decapper, 24-format, compatible with validated sample storage tubes and cartridges
46-8014	IntelliXcap Automated Screw Cap Decapper, Acoustic 96-format, compatible with validated sample storage tubes and cartridges including Echo Qualified Acoustic Tube
46-8112	IntelliXcap Extended Stage Kit



IntelliXcap Automated Screw Cap Decapper/Recapper 96-format, Extended Height

The IntelliXcap™ Automated Screw Cap Decapper/Recapper, Extended Height is specifically designed to work with taller 96-format tubes, including nasal swabs with integral caps. The IntelliXcap range is engineered for increased throughput, ease of use, and high system reliability — making it ideal for any laboratory managing a large number of samples.

The system is extremely fast and works with multiple 96-format sample storage tube types from a range of tube manufacturers. The modular design makes changing cap driver cartridges to accommodate various tube types quick and easy.

Key Features

Advanced Automation for Increased Throughput


- Extremely fast, works with multiple sample storage tube types from a range of tube manufacturers
- Suitable for use with taller 96-format tubes including those with swab caps
- Capable of decapping a complete rack of 96 tubes in as little as 20 seconds
- Ideal for medium-to-high throughput laboratories managing compound libraries, biobanks, or other biological samples

Compact, Modular Design for Ultimate Flexibility

- Features a unique, fully automated, interchangeable cartridge cap driver system
- Cartridges allow format change between different tube types in less than 2 minutes
- Simply by fitting the relevant cartridge, IntelliXcap Extended Height 96-format is compatible with all Azenta Life Sciences sample storage tubes, and supports other sample tube ranges available in the market
- Cartridges are available for both internal and external thread screw cap tubes
- No need for specialist intervention to change cartridge drivers
- Control is via a quick installation, easy-to-use touch panel

Easy to Use, Easy to Integrate

- Quick installation, easy-to-use touch panels bring instant control
- Wide operating temperature range of 5°C to 40°C
- Electric torque control means less persistent wear on cap drivers



- Sample rack stage extension facilitates fully automated workflow integration

Driver and Cap Compatibility: We provide a wide range of custom and off-the-shelf driver cartridges to accommodate your specific labware needs. The system's modular design allows you to quickly and easily change the cap driver cartridge and seamlessly swap between tube types — making it ideal for automated high throughput workflows. It can cap, decap, and recap both internal and external thread tubes.

Technical Specification

Dimensions	370(D) x 256(W) x 464(H) mm
Weight	24 kg
Ambient operating temperature	5°C to 40°C
Capping/decapping speed	20–40 seconds to decap or to recap a rack of tubes

Ordering Information

IntelliXcap Automated Screw Cap Decapper/Recapper 96-format, Extended Height	
46-8018	Suitable for all Azenta Life Sciences tubes and almost all other brands with relevant cartridges; 1 decapper/capper
IntelliXcap Accessories	
46-8112	IntelliXcap Extended Stage for Automation, suitable for IntelliXcap 24, IntelliXcap 48 and IntelliXcap 96
Compatible cartridges	
48-8013-18	IntelliXcap Automated Screw Cap Decapper Cartridge, 96-format, for Rhinostic Swab internal thread

Compatible with all other IntelliXcap Automated Screw Cap Decapper/Recapper 96-format (46-8012) cartridges.

Please specify when ordering.



AZENTA
LIFE SCIENCES

IntelliXcap Automated Septum Cap Decapper/Recapper 96-format

The IntelliXcap™ Automated Septum Cap Decapper/Recapper is ideal for use in high-throughput workflows that require capping and decapping of septum capped tubes. The instrument provides secure sealing – eliminating the need for manual intervention, while preserving vital sample integrity.

This high-speed whole rack tube decapper and capper automatically removes, recaps, or disposes of septum caps from a full 96-format SBS rack of sample storage tubes as well as allowing for the capping of Azena tubes from a cap mat. The system is ideal for use within high-throughput environments such as in biobanking, compound libraries, and other storage applications.

Compatible with multiple sample storage tube types from a range of tube manufacturers, the system provides flexibility in either fully automated robotic systems or as a standalone unit



Key Features

Advanced Automation for Increased Throughput

- Suitable for use with 96-format tubes from Azena as well as those from other manufacturers
- Automation and speed reduces sample handling time, increasing sample throughput workflow
- Ideal for medium to high-throughput laboratories managing compound libraries, biobanks, or other biological samples

Compact, Modular Design for Ultimate Flexibility

- No requirement for compressed air for flexible use within the lab
- Capping module available for Azena labware, enabling capping of tubes from cap mat
- Single column harpoon decapping and recapping available as add-on

Easy to Use, Easy to Integrate

- Quick installation and easy-to-use touch panels bring instant control
- Wide operating temperature range of 15°C to 35°C
- Control is via a quick installation, easy-to-use touch panel

Technical Specification

Dimensions	770(D) x 255(W) x 330(H) mm
Weight	18.5 kg
Weight (with Capping Module)	20.5 kg
Ambient operating temperature	15°C to 35°C
Power Supply Voltage	100-240 ±10% Vac

Ordering Information

IntelliXcap Automated Septum Cap Decapper/Recapper 96-format	
46-5030	Suitable for validated 96-format septum capped tubes, please specify tube type when ordering; 1 decapper/capper
IntelliXcap Accessories	
46-5061	Single Column Harpoon, 1 column, compatible with Azena consumables and other manufacturerd labware
46-5051	Adaptor Kit for Thermo Matrix tubes, compatible with 1.4 ml Matrix 2D Tubes (3791), Matrix 2D Tubes caps (4465); 1 kit



AZENTA
LIFE SCIENCES

Semi-Automated Septum Cap Capper

Semi-Automated Septum Cap Capper for TPE septum cap sealing of sample tubes in 96-well SBS format.

Offering secure storage tube sealing, this Semi-Automated Septum Cap Capper helps preserve sample integrity and audit trails in biobanks, compound libraries and other high-throughput storage applications.

Key Features

Flexible Sealing Performance, Broad Compatibility

- Semi-Automated Septum Cap Capper suitable for use with all 96-format Azenta sample storage tubes with TPE Septum Caps, as well as Micronic tubes with TPE Capmats, Matrix Technologies tubes with Sepraseals and Abgene Sealing Mats
- Depth of capping action is adjustable with the use of spacer plates, which ensures effective and efficient capping of all tube types
- Semi-Automated Septum Cap Cappers can be used to seal either full, or partially-full, racks

Fast and Easy to Use

- Cycle time for capping is 20-30 seconds, depending on tube type and on the spacer plate used
- Simply place a rack in the Semi-Automated Septum Cap Capper drawer, with a TPE septum cap mat fitted loosely on top and, if needed, a spacer plate (adapter). Push the drawer shut, and the capping action will start automatically. Once all caps have been inserted, the drawer will open automatically and the sealed rack of tubes can be removed



Semi-Automated Septum Cap Capper

Additional TubeLock Functionality

- Using the Semi-Automated Septum Cap Capper in conjunction with Azenta racks featuring TubeLock enables automatic locking and unlocking of tubes within the rack
- Azenta racks with TubeLock are dual position, allowing tubes to be locked in place for added sample security during handling
- Individual sample tubes can be pushed downwards to lock and pulled upwards to unlock in the rack
- For a faster approach, the Semi-Automated Septum Cap Capper, with an adapter set, can be used to lock, or unlock, all tubes in a 96-format rack simultaneously

Ordering Information

System	
46-2004-115V	Semi-Automated Septum Cap Capper, 115V
46-2004-230V	Semi-Automated Septum Cap Capper, 230V
Adapters	
42-1003	Adapter for Azenta tube rack [A] C:2mm H:4mm - Included
6.09.661	Adapter for Azenta tube rack [B] C:1.6mm H:14mm - Included
42-1001	Adapter for Azenta tube rack [C] C:2mm H:21mm - Included
6.09.663	Adapter for Azenta tube rack [D] C:0.4mm H:30mm - Optional
6.09.664	Adapter for Azenta tube rack [E] C:0.4mm H:37mm - Optional



AZENTA
LIFE SCIENCES

Sample Tube Management Systems



AZENTA
LIFE SCIENCES

Manual Tube Pickers and Manual Decappers

A manual device, with one-handed operation, the Manual Single Tube Picker incorporates an eject button enabling tube pick and place operations whilst avoiding contact between the hand and tube.

The Manual Single Tube Picker minimizes the risk of heat transfer to a frozen sample, as well as minimizing the risk of possible cross contamination.

Simple to use, even when wearing laboratory or cryogenic-handling gloves, the tube picker is available in either 96-rack or 48-rack format.



Key Features

- Single-handed operation, soft grip handle
- Easy to use whilst wearing gloves
- All components that potentially come into contact with samples are made from temperature-resistant and chemically-resistant plastics
- Tubes can be selected from any position in the rack
- Removes the need for manual contact with tubes reducing the risk of heat transfer to frozen samples and cross contamination
- Ideal for use with Azenta sample storage tubes and compatible with sample storage tubes from most other manufacturers

Specifications

Manual Single Tube Picker, 96-format	length 150mm	diameter 17mm	weight 20g
Manual Single Tube Picker, 48-format	length 165mm	diameter 17mm	weight 23g

Manual Decappers: The ideal range of accessories for the manual removal of TPE caps and Screw caps.

Our TPE septum cap options comprise: A 1-way decapper for the individual removal of TPE caps and an 8-way decapper for the simultaneous removal of 8 TPE caps. For screw caps, we offer our Azenta screw cap Capper/Decapper, designed specifically for external thread tubes.



Ordering Information

10-5010	Manual Single Tube Picker, 96-format, 1 picker per case Suitable for Azenta and most other manufacturers 96-format tubes
10-5020	Manual Single Tube Picker, 48-format, 1 picker per case Suitable for Azenta and most other manufacturers 48-format Cryo Tubes
65-54000	Manual Decapper-1, for removal of 1 TPE cap at a time, 1 decapper per case
65-54001	Manual Decapper-8, for removal of 8 TPE caps at a time, 1 decapper per case
65-54004	Manual Screw Cap Capper/Decapper, 1 capper/decapper per case. Suitable for Azenta 96-format External Thread tubes

Mohawk Semi-Automated Tube Picker

Key Features

- Fast and reliable tube picking
- Pick from 48 or 96 SBS racks
- Works with frozen and thawed tubes
- Create picklists in Excel®, csv, text, xml or json
- 1D scanner for rack barcode included as standard



The Mohawk Semi-Automatic Tube Selector enables you to pick individual or groups of tubes directly from 48- or 96-well racks with ease. Trying to extract a single tube from a rack can be difficult and doing so manually can lead to positional errors. By using software to control 96 individual solenoids, the Mohawk can lift one, two or up to sixteen individual tubes at a time, making it easy to remove them from the rack and place in a subsidiary rack or other tube device.

The pins to be selected can be set in a 'pick list' which can be generated manually from the on-screen graphical representation of the Mohawk deck or from an Excel spreadsheet. The Excel file can also be created manually or imported from the FreezerPro® inventory management system or from your LIMS or another sample tracking programme that can export Excel format files. The picklist can even include tubes that are contained in multiple racks. The integrated 1D barcode reader identifies individual racks and the system software then uses this ID to look up the pre-defined picklist. This reduces the risks of the wrong tubes being selected during picking.

The Mohawk works out of the box and needs no set up or calibration. In Automatic mode, The Mohawk allows for operation with minimal interaction of the control PC, while in Worklist mode, a file that describes the tubes that are to be selected is used to control the pins. It can be in Excel, csv, text, xml or json format. The Mohawk works with both frozen and thawed tubes. It can work in combination with a Ziath Mirage, Express or Cube reader to automatically create pick lists and then verify the contents of the subsidiary (or Daughter) racks and reupload their new positions to FreezerPro® or to your LIMS via an Excel file.



Mohawk Semi-Automated Tube Picker



Technical Specification

Dimensions (W x D x H)	300 x 250 x 180 mm
Total weight	5.7 kg
Power adapter input	110-220V AC +/- 10V
Power adapter output	12V DC, 5A
Labware compatibility	48 position SBS racks
Picklist import	Excel, text, csv, json and xml
Operating system on host PC	Windows® 10, 11
Compatibility	48 and 96 position SBS racks

Included with the Cube

- Mohawk tube picker instrument
- Power adaptor, power cable
- Mohawk software and license
- 1D rack code scanner and bracket
- User guide

Ordering Information

ZTS-MHK-48	Mohawk Tube Picker 48 position; includes Mohawk tube picker instrument, power adaptor, power cable, Mohawk software and license, 1D rack code scanner and bracket, user guide; 1 picker
ZTS-MHK-96	Mohawk Tube Picker 96 position; includes Mohawk tube picker instrument, power adaptor, power cable, Mohawk software and license, 1D rack code scanner and bracket, user guide; 1 picker
ZTS-MHK-LOW	Low Profile Tube Rack Adapter for Mohawk (if using tubes 0.5 ml or less); 1 adapter

Direct Tube Marker

Tube labels or stickers are no longer required for sample identification - Direct Tube Marker is designed to print text, 1D linear barcodes, 2D data matrix codes and / or graphics directly onto most polypropylene sample storage and general-purpose laboratory tubes. Using Thermal Pixel Printing technology, Direct Tube Marker prints directly onto the surface of tubes sized from 0.5ml to 50ml.

Clear and durable marking with consistently legible text and barcodes printed at high resolutions is easier and faster than marking by hand. Direct Tube Markers permanent printing is resistant to water, alcohols (methanol, ethanol, isopropanol), DMSO, haematoxylin, liquid nitrogen and mechanical abrasion.*

Printed information is stable over a wide temperature range from -196°C to 100°C and markings do not transfer when tubes are manually handled.

The Direct Tube Marker delivery package contains tube adapters with hole diameters 8.5mm, 11.5mm and 12.7mm. These are compatible with most commonly-used tubes. Other adapters are available to order separately, please specify at time of ordering.



Key Features

- Direct Tube Marker prints directly on the surface of plastic laboratory and sample storage tubes in seconds
- Most 0.5ml to 50ml tubes with a smooth surface can be printed
- Adjustable tube support platform allows printing on tubes without a lip or collar
- Tube profiles created in the Direct Tube Marker software can be saved for future reference and sample tracking
- Adjustable sled pivot point improves overall print quality for larger tubes and labels
- Wide range of adapters available for differing diameter tubes

Flexible Label Printing to Suit Your Needs

- Tubes can be marked with any orientation
- Direct Tube Marker can print text, 1D linear barcodes, 2D datamatrix codes and graphic files (monochrome .bmp, .gif, .tiff) such as logos
- Several lines of text can be printed, such as: sample name, your name, date, time etc.
- All TrueType fonts available on the connected PC can be printed
- Font size and type will determine the maximum available printable content


High-Quality and Durable Printing

- Thermal Pixel Printing technology produces clearer and more durable marking than marker pens
- Metallic blue, black and white ink ribbons are available for a range of colored tubes and colored contents
- Printing straight onto the tube surface is much easier than marking by hand or adhesive label
- Quick and easy to change ink ribbon
- Marking resistant to 70% ethanol, 70% isopropanol, >98% methanol, DMSO, pure haematoxylin, liquid nitrogen & mechanical abrasion*
- Printing is temperature resistant from -196°C (liquid nitrogen) to 100°C
- 300m ink ribbon will mark up to 100,000 tubes**
- Enables full integration within robotic automation systems
- Easy integration with SiLA compliant devices using SiLA driver
- Instrument reliability

*ribbon chemical resistance dependent on ribbon type

**dependent on size of label print

Direct Tube Marker



Specifications

Parameter	Value
Communication	USB port
Print Method	Thermal Pixel Printing
Label Detail	Print height: 6mm, Print width: 40mm, Any print orientation
Tube Compatibility	Most plastic laboratory tubes 0.5 to 50ml
Dimensions L x W x H	280mm x 270mm x122mm
Electrical	V in: AC 100-240V V out: DC 18V
Weight	4.8kg
Throughput	Up to 200 Plate Seal Removals per Hour

Ordering Information

4ti-0680-1	Direct Tube Marker, includes: 3 x tube adapters (4ti-0681, 4ti-0683 and 4ti-0684)
4ti-0681	Direct Tube Marker Adapter, tube adapter for 1.5 ml/2 ml tubes
4ti-0683	Direct Tube Marker Adapter, tube adapter for 2D cluster tubes
4ti-0684	Direct Tube Marker Adapter, tube adapter for cryo tubes
4ti-0685-1	Direct Tube Marker Adapter, tube adapter for 2 ml screw cap tubes (10 mm diameter)
4ti-0685-2	Direct Tube Marker Adapter, tube adapter for 15 ml tubes
4ti-0685-3	Direct Tube Marker Adapter, tube adapter for 50 ml tubes
4ti-0685-4	Direct Tube Marker Adapter, tube adapter for Matrix tubes
4ti-0685-5	Direct Tube Marker Adapter, tube adapter for 0.2 ml tubes
4ti-0686	Direct Tube Marker Ribbon, metallic blue
4ti-0688	Direct Tube Marker Ribbon, white
4ti-0689	Direct Tube Marker Ribbon, black
4ti-0689-1	Direct Tube Marker Ribbon, black, universal





Tube Auditor™



Tube Auditor™ is a fast, accurate and non-invasive device used to measure sample volume in microtubes. The instrument performs quick and easy Goods-In and Goods-Out QA checks to help ensure that your suppliers are providing what you ordered and that your customers are receiving what you promised. Tube Auditor gives confidence in sample quality (at both collection and output) and helps to reduce downstream costs and waste, which can arise from empty wells and precipitated samples

- High-speed measurement - accurate to better than +/- 10µl - minimizes downstream costs from the processing of empty plate wells
- Compatible with 96-SBS or 48-SBS format racks
- Precipitate detection - increases confidence in the concentration of delivered output samples
- Cap detection - helps avoid damage to liquid handling tips caused by failed de-cap operations
- No need to de-cap tubes - eliminates possibility of cross contamination
- Manual or remote operation - if integrated into an automated system, an OPC license is required to enable remote control
- Image storage and recall - allows audit trail and provides ability to re-assess or re-analyze the image
- Tube Auditor uses high-speed vision technology to measure sample volume and detect precipitate



Specifications

Parameter	Tube Auditor
Volume measurement	✓
Cap detection	✓
Manual and remote operation	✓
User interface software	✓
Image capture and recall	Available with "Pro" version enabled only
Precipitate detection	Available with "Pro" version enabled only
Dimensions (L x W x H)	872 x 433 x 433mm (34.33 x 17.05 x 17.05in)
Weight	35kg (77lbs)
Electrical	110-240 VAC 50/60Hz
PC	Microsoft Windows 10
Output data format	CSV or XML(user configurable)


Ordering Information

System	
96-0001	Tube Auditor, 96-format, for high-speed non-contact volume measurement
96-0002	Tube Auditor, 96-format, for high-speed non-contact volume measurement, with 2D tube reader
96-0003	Tube Auditor, 96-format, for high-speed non-contact volume measurement, with opc license
96-0004	Tube Auditor, 96-format, for high-speed non-contact volume measurement, with 2D tube reader and opc license
Tube Auditor compatible with 48-SBS racks	
98-0001	Tube Auditor, 48-format, for high-speed non-contact volume measurement
98-0002	Tube Auditor, 48-format, for high-speed non-contact volume measurement, with 2D tube reader
98-0003	Tube Auditor, 48-format, for high-speed non-contact volume measurement, with opc license
98-0004	Tube Auditor, 48-format, for high-speed non-contact volume measurement, with 2D tube reader and opc license
Software License for Pro version	
97-0001	Tube Auditor Precipitate Detection Software, used to activate the Precipitate Detection feature of Tube Auditor, provides a facility to save and recall audited tube images, suitable for Tube Auditor 48-format (part number 98-0001-98-0004) and Tube Auditor 96-format (part number 96-0001-96-0004), compatible with all existing Tube Auditor configurations

FreezerPro® Sample Management System


The FreezerPro system is a scalable web-based sample inventory management system ideal for users managing hundreds of collections in all types of organizations ranging from lower throughput labs to central biorepositories with millions of records; delivers secure management of samples and sample information.

- Complete sample management
- Track samples in and out of freezers
- Reporting and data export/import (multiple formats)
- Live and smart search
- Full audit trails



Using cross-region servers and industry-defining storage technologies, the system is designed to keep samples in the right conditions and make retrieving them for study easier. We manage security through sample data encryption, both during sample transit and while they are at rest. FreezerPro is designed and developed as a web-based sample management system that provides access to sample information from anywhere in the world.

- **No IT Requirements Whatsoever** - By not requiring installation, individual labs within larger organizations gain more autonomy. Start-ups and other small labs can completely avoid all IT staffing and investment
- **Regulatory Liability Protection** - Remove nearly all liability concerns related to data security. Using FreezerPro Cloud negates expensive information security upgrades for individual labs
- **Lower Capital Expenditure** - With FreezerPro Cloud, the total cost of deployment is significantly lower, compared to other sample management software options. FreezerPro is a class-leading frozen sample management solution that is indispensable to any modern scientific or pharmaceutical laboratory
- **Deployed in Hours, Not Weeks** - FreezerPro systems automate more daily sample management operations than expensive software solutions
- **Makes Daily Lab Operations Easier** - FreezerPro Cloud records information for sample check in and out, aliquotting and inheritance, plus storage location data with only a few clicks or a drag of the mouse
- **Enhances Accuracy and Availability of Lab Information** - Advanced data mining, automatic data validity checks, hundreds of alert settings and embedded integration with Microsoft Excel are just a few of the tools users and admins have to ensure that data is input properly and is immediately accessible for analysis and reporting.



AZENTA
LIFE SCIENCES

Automated Sample Storage Solutions

In addition to our extensive range of consumables and instruments, Azenta Life Sciences is also the global leader in automated sample storage systems.

Our deep application experience and proven, highly reliable technology protects sample integrity and improves sample visibility to ensure the continuity of your sample management. From compound management and high-throughput screening, to biobank sample management, to storage of cellular products at cryogenic temperatures; Azenta provides flexible, modular solutions that offer the security and reliability to optimize our customers workflow. As the preferred storage partner to the world's top biotechnology companies, Azenta reduces risk and delivers industry-leading uptime that maximizes the return of your sample management investment.

Flexible automated storage technology that fits your applications allowing you to focus on delivering your research goals.

Transforming workflows and productivity to improve research timelines, our stores accommodate collections from fewer than 7,000 to up to 100 million samples, and are adaptable to meet future needs.

- **SampleArc™ Space Efficient** – Ambient to -20°C
- **BioArc™ Space Efficient** – -80°C
- **SampleArc™** – Ambient to -20°C
- **BioArc™** – -80°C
- **CryoArc™ Pico** – -190°C


*Temperature variation is +/- 6°C due to thermocouple variation and thermal loading



*SampleArc™
Ambient to -20° C*



*BioArc™ Deca
Automated Cryogenic
Storage Solutions -190° C*



*CryoArc™ Pico -190°C
LN2-Based Automated
Storage System*



*BioArc™
-80° C*





AZENTA
LIFE SCIENCES

ARCON Control Software

ARCON™ control software (previously known as Strata) enables precise sample inventory management within the Azenta automated storage systems. Strata's intuitive, information rich interface makes automated storage flexible, easy and convenient.

- No special engineering expertise is required
- Efficiently processes and securely protects samples

- Sample data is accessible from any location using a laptop, tablet, or other web enabled devices
- Simplified integration with corporate IT and LIMS networks
- Intelligent diagnostics and error recovery ensure reliable performance



High Efficiency Cryogenic Freezers

Simple, Secure -190°C LN2 Vapor Storage

Azenta Life Sciences' range of high efficiency freezers is designed to offer innovative cryogenic solutions for life science research and therapy. Successful use of samples in the lab or clinic is the ultimate goal and cryopreservation is a critical link. Decades of experience across the cold chain, vertical integration and a commitment to quality make our unique freezers possible.

Increase capacity, improve ergonomics and stay connected with simple, secure cryo storage designed for samples by users. Azenta is committed to preserving your sample potential!



*Azenta A220 and
Extended Height
E264 Cryo Freezers*



*Azenta A440 and
Extended Height
E528 Cryo Freezers*



*Azenta A700 and
Extended Height
E840 Cryo Freezers*



*Azenta A1000 and
Extended Height E1200
Cryo Freezers*

Key Features

Performance


- Temperature, storage density, LN2 usage

Redundancy and security

- Redundant valves and sensors
- Secondary temperature & LN2 level monitoring
- Lifecycle testing

Ergonomics and ease of use

- Better sample access with reduced reach and lift
- Full sample visibility with cryo LED and fog clear
- Workspace to maintain cold chain



Remote monitoring and connectivity

- Simple controller with WiFi/LAN
- 24/7 real-time monitoring – monitor all freezers from one interface
- Text/Email alerts, Cloud backup, Web App
- View key parameter charts and detailed equipment temperature history
- Download reports for custom date ranges



AZENTA
LIFE SCIENCES

CryoPod® Carrier

Portable Liquid Nitrogen (LN2)-Based Cryogenic Transport.

CryoPod® Carrier provides a safe, reliable and portable < -150°C cryogenic environment for the handling and transport of biological specimens for over 3 hours. The instrument displays and logs temperature, date and time, and features audible and visual alarms, and integrates into an optional automated filling station for hands-free replenishing of the LN2 charge in less than 15 minutes.

Ensures operator safety

- Allows safe and quick transportation of cryogenic samples
- Hands-free auto-fill option

Maintains sample cold chain integrity

- Temperature display with audible and visual alarms
- Temperature logging and retrieval

Delivers reliable performance

- Over 3 hours <-150°C
- No direct sample contact with LN2

Portable

- Compact footprint; only ~9 lbs
- Built-in handle and bottom finger grips

Specifications

Hold Time	Over 3 hours < -150°C with lid closed
Capacity	One 2" cryobox, 2 SBS plates or some cassette sizes
Charge	Requires < 3L liquid nitrogen (LN2)
Alarm	Two settings with audible and visual signals
Lid	Magnetized foam lid for safer transport and insulation
Power	3 AA alkaline batteries (included)
Temperature Audit Trail	Downloadable temperature log data via USB port and CryoPod Data Log Software
Weight	4.1 kg (9 lbs) without LN2 5.9 kg (13 lbs) fully charged with 3L LN2 - no samples
External dimensions (L x W x H)	34.0 x 32.0 x 26.0 cm (13.4 x 12.6 x 10.2 in)
Cryogenic chamber basket dims (L x W x H)	17.4 x 18.8 x 7.8 cm (6.9 x 7.4 x 3.1 in)



Ordering Information

CryoPod® Carrier	
243354-001	Suitable for cryoboxes, includes orange lid and manual fill kit; 1 carrier
284574	Suitable for cassettes, includes purple lid and manual fill kit; 1 carrier
Lids & Accessories	
252888-001	Orange; 1 lid
252888-002	Green; 1 lid
252888-003	Grey; 1 lid
252888-004	Pink; 1 lid
252888-005	Purple; 1 lid
252885	CryoPod Carrier for Cryoboxes manual fill kit (compatible with 243354-001); 1 kit
321706	CryoPod Carrier for Cassettes manual fill kit (compatible with 284574); 1 kit
252886	CryoPod® LN2 Absorbent Pads, 4pc, 2/pk



AZENTA
LIFE SCIENCES

4titude PCR/Microplate Consumables and Instruments



Attitude Manufacturing & Quality Standards

Azenta is ISO 13485:2016 certified to manufacture and supply microplates and associated consumables for the life sciences sector. We also provide complete custom design solutions from prototyping to tool design and contract manufacturing.



Manufacturing Standard for Microplate and Associated Consumables

- ISO 13485:2016 certified
- Process validation & mapping
- Fully document controlled manufacturing processes
- Statistical analysis of production processes
- Continuous improvement programs
- Injection molding in ISO class 7 cleanrooms
- Virgin, medical grade polymers

Quality Standard

Azenta performs visual, physical and biological tests to ensure the integrity of our PCR/Microplate consumables and that they are contamination free at all times.

- Consumables are certified free from human genomic DNA, nucleases and pyrogens
- Skirted microplates and PCR plates meet the SBS standard footprint
- PCR inhibition tests are performed on polymers used
- Leak tests are performed on every well of every PCR plate
- White-well plates are checked for background fluorescence



AZENTA
LIFE SCIENCES

4titude PCR Plates: Clear, Frosted or White Wells & Low DNA Binding Properties





4titude Clear, Frosted or White Wells & Low DNA Binding Properties

Selection of the right plastic material of your PCR consumables has a measurable effect on your (q)PCR results. Azenta has carried out extensive research and development into our PCR consumables to offer customers a range of products suitable for the diverse applications and instruments they are required for.

Well color options - A question of sensitivity

Well color is often not considered when choosing a PCR plate, but can in fact have a significant impact on results. PCR plates are available in three color options: clear, frosted and white, each of which has specific advantages and disadvantages.



Greater sensitivity in a qPCR reaction enables earlier Ct values and higher fluorescence readings.

Identical qPCR assays in plates with clear, frosted and white wells (4ti-0770/C, 4ti-0772 and 4ti-0771, respectively). Clear and frosted wells perform similarly whereas white wells gave earlier Ct values and higher fluorescence intensity.

Azenta recommends the use of white wells, where possible, to achieve the greatest sensitivity and consistency of qPCR reactions. Another point to consider are well color recommendations from (q)PCR instrument manufacturers, for example, Roche recommend the use of white wells on their instruments and, ABI® recommend the use of frosted wells. In addition to using the recommended well color for your instrumentation, this must be combined with thin walled tubes for optimal heat transfer and optimal sealing to prevent evaporation.

Low DNA binding - Smarter plastics for advanced applications

Polypropylene (PP) is the best material for PCR tubes as PP is chemically inert, resistant to solvents and well suited for injection molding, allowing for production of thin-walled tubes for optimum PCR results.

DNA has been shown to bind PP, especially at high ionic strength, despite the very hydrophobic nature of this material. This has typically not been an issue but due to progressing miniaturisation of reaction volumes and the introduction of new technologies such as NGS, ultra-low DNA binding consumables have become essential for use in sensitive assays.

Please see our dedicated application note on low DNA binding products which can be found on our website for information outlining the characteristics of our selected low bind polymer and for the range of low-binding PCR consumables available.




Coding Options & Ethylene Oxide Treatment

Linear and 2D datamatrix coding – Sample tracking made easy

All skirted and semi-skirted plates are available with linear Code 128 format barcodes for identification and traceability of your samples.

The labels are highly scratch-resistant and can withstand cold storage (-80°C), temperatures of up to 100°C, and solvents such as DMSO. Single, double, triple or quadruple barcodes are available, and a variety of custom options.



Position of standard barcode labels



8 Well PCR Tube Strip with PC Frame with off-the-shelf 2D code

2D datamatrix coding uses a defined number of fields to encode alphanumeric information. The code uses data redundancy so even if codes become partly destroyed, the information will be retained. Azenta offers 2D coding on several products including 8 Well PCR Tube Strips with PC Frame and PCR tubes with flat caps.

Should additional customizations not covered by our standard barcoding and 2D coding service be required, further information and our custom linear barcode request form can be found on our website.

Ethylene Oxide Treatment – Reliable consumables for forensic applications



Azenta has a stringently controlled clean-room production facility, for production of PCR consumables free from DNA and RNA contamination. However, some applications require the absolute highest quality of consumables such as forensic workflows and tissue culture.

For these applications, Azenta offer treatment of selected products with Ethylene Oxide, a technique proven to reduce traces of amplifiable DNA, for peace of mind in your reactions. Additional plate types can be treated on request.



AZENTA
LIFE SCIENCES

4titude FrameStar 2-Component PCR Plates

FrameStar® PCR plates prevent sample loss by minimizing thermal expansion during PCR, enabling reductions in PCR volumes and cost savings on reagents.

The 2-component design combines the advantages of thin walled polypropylene (PP) tubes, for optimum PCR results, with a rigid polycarbonate (PC) frame for highest thermal stability and rigidity, making them the plates of choice for any robotic workflows.

- **Multiple frame color options with clear, frosted, white or black tubes are available**

Flexible solutions for every application

- **No warping due to stable polycarbonate frame**

Reliable use with stackers and liquid handlers

- **Minimizing thermal expansion**

Better sealing properties & reduced evaporation for improved PCR consistency

- **Downscaling of reaction volumes possible**

Cost saving

- **Standard and custom barcoding options available**

Error-free sample tracking

- **Plates with ultra-low DNA binding properties and processing options such as ethylene oxide treatment available**


Tailor-made solutions




Evaporation from standard PCR plates vs. FrameStar plates

Thermal expansion of polypropylene (PP) plates leads to greater risk of evaporation from outer wells.

PP is the optimum material for PCR tubes. It provides the most efficient heat transfer, as well as an inert surface with low binding affinity to nucleic acids, proteins and other molecules. However, the material is not thermally stable in plate format, causing it to expand and contract during each PCR cycle. Such thermal expansion will weaken the plate seal and lead to sample evaporation, mainly from corner and outer wells.



Side-on view of a PCR plate in a thermal cycler. The sealed plate is sandwiched between the cycler block and the heated lid but it is only partly fixed in position at the bottom of tubes, allowing the plate to expand horizontally at the top.



Standard polypropylene plates expand by up to 2 mm during thermal cycling which leads to movement of wells away from the plate centre. This movement is most significant in corner and outer wells.



AZENTA
LIFE SCIENCES

FrameStar 2-Component PCR Plates

FrameStar 2-component technology allows for reduction of assay volumes and cost.


Due to the much improved seal integrity, reaction volumes can often be reduced when using FrameStar plates. Such downscaling of experiments can be successfully implemented without any loss of assay sensitivity or consistency and, reagent savings can be considerable.

Evaporation from standard PP plates is highest in the outer wells

Since thermal expansion and movement of wells in standard PP plates is greatest around the edges of the plates, evaporation is highest from the two outer rows of wells. The adjacent figure illustrates the level of risk of sample evaporation from different areas of PP plates. The inner 32 wells of a standard 96 well plate have low risk of evaporation whilst the risk of sample loss is much higher in the outer two rows which contain 65 per cent of the wells.

“FrameStar plates led to significantly better results and reduced evaporation compared to standard PCR plates.”

Dr. Andreas Dahl, MPI f. Molekulare Genetik, Berlin, Germany




Risk of evaporation from the outer rows (red) of a standard PP PCR plate is highest, medium level evaporation occurs in the second row (yellow) and sample loss from the inner 32 wells (green) is lowest (the dotted line represents an expanded standard PP plate).

FrameStar 2-component plates improve consistency of PCR results

We have compared the degree of evaporation from different areas of standard PP and FrameStar PCR plates. First, the 64 outer wells (two outer rows - above, red and yellow area) of both plate types were filled with 10 µl H₂O. Plates were then sealed with a qPCR adhesive seal (code 4ti-0560) and their total weight determined before and after PCR. The experiment was repeated with a set of plates of which the inner 32 wells (green area) were filled. Table 1 shows that evaporation from outer wells of standard PP plates was 65 higher than from inner wells. As a result, evaporation causes varying changes in reaction volume across standard PP plates.

The results below show that reaction volumes remain consistent across the 96 wells (or 384 wells, data not shown) in FrameStar plates. In contrast, the reaction volumes in standard plates differ significantly between wells during PCR. Reagent concentrations in outer rows will increase dramatically, resulting in sub-optimal reaction efficiency. In extreme cases samples may fully evaporate.



The polycarbonate frame of FrameStar plates is more heat resistant than standard polypropylene plates which reduces thermal expansion to a minimum. For this reason seal integrity remains intact even at elevated temperatures during PCR.

FrameStar minimizes sample loss across the plate

Plate Type	Well position	Starting weight (g)	Weight post PCR (g)	Weight loss (g)	Volume loss	
					Total	Per well
FrameStar 4ti-0710	outer 64 wells	26.230	26.193	0.037	37 µl	0.57 µl
Standard PP	outer 64 wells	17.299	17.118	0.181	181 µl	2.8 µl
FrameStar 4ti-0710	inner 32 wells	25.841	25.824	0.017	17 µl	0.53 µl
Standard PP	inner 32 wells	17.132	17.078	0.054	54 µl	1.69 µl

Table 1: Weight and volume loss from different sections of 96 well PCR plates. Results shown are averages from 5 plates of each plate type. Volume loss from the outer wells of standard PP plates was 5-times higher than from FrameStar plates.

Evaporation has a significant effect on the reaction conditions resulting in noticeable effects, especially for qPCR. Identical samples can exhibit significant differences in their Ct values, depending on their position on the plate.

A solution to the problem of evaporation related qPCR inaccuracies is the use of 2-component plates.



AZENTA
LIFE SCIENCES

4titude FrameStar 384 Well Skirted PCR Plate

Polypropylene wells, polycarbonate frame, cut corner A24; working volume: <30 µl, total well capacity: 55 µl; designed for use on standard thermal cyclers

- Our FrameStar 384 Well Skirted PCR Plates are designed for high-throughput PCR
- Compatible with the majority of 384 well block PCR, qPCR and sequencing instruments
- Rigid 2-component design eliminates warping and distortion during PCR making it ideal for use with robotic systems
- The skirt allows for labeling or barcoding
- Ultra-low DNA binding option available

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Stackable

Wells


- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- <30 µl working volume, 55 µl total well capacity

Frame

- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of the frame
- Alphanumeric grid reference to aid well and sample identification

Use

- Recommended for low volume PCR
- Compatible with standard multichannel pipettes
- Ideal for use with robotic systems
- Compatible with the majority of 384 well block PCR, qPCR and sequencing instruments



Options

- Available with the following frame color options for the clear well variety: purple, blue, clear, green, red, and black
- Also available as a black frame with white wells (for optimum signal-to-noise ratio when using fluorescent based assays) and as a black frame with black wells (for minimal light diffusion and interference)
- Also available as a clear frame with frosted wells for use with ABI®/LifeTechnologies® qPCR instruments
- Available barcoded upon request

4titude FrameStar 384 Well Skirted PCR Plate

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	10.30 ± 0.05 mm
Well depth	9.20 ± 0.10 mm
Well diameter	3.00 ± 0.10 mm
Distance to center of A1 from top edge	8.99 ± 0.25 mm
Distance to center of A1 from left edge	12.13 ± 0.25 mm
Pitch (distance between A1 and A2)	4.50 mm

Ordering Information

4ti-0384	FrameStar 384 Well Skirted PCR Plate, clear PP wells, purple PC frame, cut corner A24, 50 plates per case
4ti-0384/B	FrameStar 384 Well Skirted PCR Plate, clear PP wells, blue PC frame, cut corner A24, 50 plates per case
4ti-0384/C	FrameStar 384 Well Skirted PCR Plate, clear PP wells, clear PC frame, cut corner A24, 50 plates per case
4ti-LB0384/C	FrameStar 384 Well Skirted PCR Plate, clear PP wells, clear PC frame, low binding, cut corner at A24, 50 plates per case
4ti-0384/G	FrameStar 384 Well Skirted PCR Plate, clear PP wells, green PC frame, cut corner A24, 50 plates per case
4ti-0384/R	FrameStar 384 Well Skirted PCR Plate, clear PP wells, red PC frame, cut corner A24, 50 plates per case
4ti-0384/X	FrameStar 384 Well Skirted PCR Plate, clear PP wells, black PC frame, cut corner A24, 50 plates per case
4ti-0384/RIG	FrameStar 384 Well Skirted PCR Plate, clear PP wells, black PC frame, extra rigid, low profile, cut corner A24, 50 plates per case
4ti-LB0384/RIG	FrameStar 384 Well Skirted PCR Plate, clear PP wells, black PC frame, extra rigid, low binding, low profile, cut corner A24, 50 plates per case
4ti-0385	FrameStar 384 Well Skirted PCR Plate, white PP wells, black PC frame, cut corner A24, 50 plates per case
4ti-0386	FrameStar 384 Well Skirted PCR Plate, black PP wells, black PC frame, cut corner A24, 50 plates per case
4ti-0387	FrameStar 384 Well Skirted PCR Plate, frosted PP wells, clear PC frame, cut corner A24, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude FrameStar 384 Well Skirted PCR Plate, Roche Style

Polypropylene wells, polycarbonate frame, cut corners A24 and P24; working volume: <30 µl, total well capacity: 55 µl; designed for use on the Roche LightCycler® 480 with 384 well block

- The dimensions of these plates are designed for optimum compatibility with the Roche LightCycler® 480, and are in a 384 well format for reaction volumes of up to 30 µl
- The rigid two-component design eliminates warping and distortion during the PCR process, making it ideal for use with robotic systems

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Stackable

Wells


- Ultra-smooth, uniform, thin-walled polypropylene wells for optimum PCR and real-time (RT-qPCR) results
- <30 µl working volume, 55 µl total well capacity

Frame

- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification

Use

- Compatible with Roche LightCycler® 480
- Recommended for low volume PCR
- Compatible with standard multichannel pipettes
- Ideal for use with robotic systems



Options

- Available as a clear polycarbonate frame with clear polypropylene wells
- Also available with white wells for optimum signal-to-noise ratio when using fluorescent based assays
- Available barcoded upon request
- Combi packs available with qPCR Seal (4ti-0560)



AZENTA
LIFE SCIENCES

4titude FrameStar 384 Well Skirted PCR Plate, Roche Style

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	10.60 ± 0.25 mm
Well depth	9.20 ± 0.10 mm
Well diameter	3.00 ± 0.10 mm
Distance to center of A1 from top edge	8.99 ± 0.25 mm
Distance to center of A1 from left edge	12.13 ± 0.25 mm
Pitch (distance between A1 and A2)	4.50 mm

Ordering Information

4ti-0380/C	FrameStar 384 Well Skirted PCR Plate, Roche style, clear PP wells, clear PC frame, cut corner A24/P24, 50 plates per case
4ti-0381	FrameStar 384 Well Skirted PCR Plate, Roche style, white PP wells, clear PC frame, cut corner A24/P24, 50 plates per case
Combi Pack	
4ti-0382	FrameStar 384 Well Skirted PCR Plate, Roche style, plus qPCR Seal, 4ti-0381 plus 4ti-0560, combi pack, 50 plates and seals per case
4ti-0383	FrameStar 384 Well Skirted PCR Plate, Roche style, plus qPCR Seal, 4ti-0380/C plus 4ti-0560, combi pack, 50 plates and seals per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude FrameStar 96 Well Skirted PCR Plate

Low profile, 0.1 ml polypropylene wells, polycarbonate frame, cut corner H1; working volume: <100 µl, total well capacity: 200 µl

- The low profile wells of this plate are shorter than “standard” profile wells, which decrease the “dead space” between the heated lid of the thermal cycler and the sample in the well
- Eliminates condensation forming on the side wall of the tubes, preventing reduction in PCR volume and increasing the efficiency of the reaction
- Especially recommended for reaction volumes below 20 µl
- The rigid polycarbonate skirt of this design eliminates warping and distortion during PCR, and makes it ideal for use with robotic systems
- Ultra-low DNA binding option available

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Stackable

Wells


- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- <0.1 ml (100 µl) working volume, 0.2 ml (200 µl) total well capacity

Frame

- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification

Use

- Recommended for low volume PCR
- Ideal for use with robotic systems
- Compatible with standard multichannel pipettes



Options

- Available with the following frame color options for the clear well variety: purple, blue, clear, green, red, black, and white
- Also available with a black frame and with a clear frame for the white well variety, for optimum signal-to-noise ratio when using fluorescent based assays
- Also available as a black frame with black wells for non-PCR fluorescent applications
- Extra rigid skirt option (4ti-0960/RIG) for use with Perkin-Elmer® Sciclone, Beckman, Hamilton, and other automation systems: eliminates the robotic grip picking up more than one plate at a time
- Ethylene oxide treated option available (4ti-0X960C/SBC) for forensic use
- Ultra-low DNA binding option available (4ti-LB0960/RIG) for sensitive applications with ultra-low DNA input and for maximum DNA recovery after low temperature storage and high temperature incubation
- Available barcoded on request
- Combi packs available (for 4ti-0960) with Optically Clear Windowed qPCR Seal (4ti-0565)
- Clear Polystyrene Lid (4ti-0287) compatible with PCR plates & robotics

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azena Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

4titude FrameStar 96 Well Skirted PCR Plate

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	16.10 ± 0.05 mm
Well depth	15.10 ± 0.10 mm
Well diameter	5.50 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0960	FrameStar 96 Well Skirted PCR Plate, clear PP wells, purple PC frame, low profile, cut corner H1, 50 plates per case
4ti-0960/B	FrameStar 96 Well Skirted PCR Plate, clear PP wells, blue PC frame, low profile, cut corner H1, 50 plates per case
4ti-0960/C	FrameStar 96 Well Skirted PCR Plate, clear PP wells, clear PC frame, low profile, cut corner H1, 50 plates per case
4ti-0960/G	FrameStar 96 Well Skirted PCR Plate, clear PP wells, green PC frame, low profile, cut corner H1, 50 plates per case
4ti-0960/R	FrameStar 96 Well Skirted PCR Plate, clear PP wells, red PC frame, low profile, cut corner H1, 50 plates per case
4ti-0960/X	FrameStar 96 Well Skirted PCR Plate, clear PP wells, black PC frame, low profile, cut corner H1, 50 plates per case
4ti-0960/W	FrameStar 96 Well Skirted PCR Plate, clear PP wells, white PC frame, low profile, cut corner H1, 50 plates per case
4ti-0961	FrameStar 96 Well Skirted PCR Plate, white PP wells, black PC frame, low profile, cut corner H1, 50 plates per case
4ti-0961/C	FrameStar 96 Well Skirted PCR Plate, white PP wells, clear PC frame, low profile, cut corner H1, 50 plates per case
4ti-0966	FrameStar 96 Well Skirted PCR Plate, black PP wells, black PC frame, low profile, cut corner H1, 50 plates per case
4ti-0X960C/SBC	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, clear PC frame, with upstand, ethylene oxide treated, single barcoded, low profile, cut corner H1, 20 plates per case
4ti-0960/RIG	FrameStar 96 Well Skirted PCR Plate, clear PP wells, black PC frame, extra rigid, low profile, cut corner H1, 50 plates per case
4ti-LB0960/RIG	FrameStar 96 Well Skirted PCR Plate, clear PP wells, black PC frame, extra rigid, low binding, low profile, cut corner H1, 50 plates per case
Combi Pack	
4ti-0960/0565	FrameStar 96 Well Skirted PCR Plate, plus Optically clear windowed qPCR seal, 4TI-0960 plus 4TI-0565, combi pack, 50 plates and seals per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

4titude FrameStar 96 Well Skirted Optical Bottom PCR Plate

Low profile, flat optical bottom, 0.1ml clear polypropylene wells, clear polycarbonate frame, cut corner H1, working volume: <100µl, total well capacity: 180µl

- The FrameStar 96 Well Skirted Optical Bottom PCR Plates (next to their counterpart within our Individual Access range, 4ti-0970/RA) are currently unique in the market, being suitable for use in both microscopy and PCR
- Optical bottom plates are ideal for applications requiring single cell sorting followed by molecular biology techniques such as (q)PCR and sequencing
- The flat bottoms enable excellent stackability, making these plates well suited for small sample volume storage such as compound libraries, with no risk of damaging the seal of the plate below
- Additionally, the small well volume enables excellent sample recovery
- The low profile wells of this plate are shorter than “standard” profile wells, which decreases the “dead space” between the heated lid of the thermal cycler and the reagents in the well
- This reduces condensation forming on the side wall of the tubes, moderating changes in PCR volume and increasing the efficiency of the reaction
- This is especially recommended for reaction volumes below 20µl
- The rigid polycarbonate skirt of this design eliminates warping and distortion during PCR, and makes it ideal for use with robotic systems

Key Features


- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Stackable

Well

- Ultra-smooth, uniform, thin-walled polypropylene wells for optimum PCR and real-time (RT-qPCR) results
- <0.1ml (100µl) working volume, 0.18ml (180µl) total well capacity

Frame

- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification



Use

- Suitable for microscopy and small volume sample storage
- Recommended for low volume PCR
- Compatible with standard multichannel pipettes
- Ideal for use with robotic systems

Options

- Also available as an Individual Access plate (4ti-0970/RA) with individually detachable and sealable wells
- Available barcoded upon request

Specifications

Parameter	Value
Plate length	127.70 ± 0.25mm
Plate width	85.48 ± 0.25mm
Plate height	16.10 ± 0.25mm
Well depth	12 ± 0.10mm
Well diameter	5.46 ± 0.10mm
Distance to center of A1 from top edge	11.24 ± 0.25mm
Distance to center of A1 from left edge	14.38 ± 0.25mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0970	FrameStar 96 Well Skirted Optical Bottom PCR Plate, clear PP wells, clear PC frame, low profile, cut corner H1, 50 plates per case
-----------------	---

4titude FrameStar 96 Well Semi-Skirted PCR Plate, Roche Style

Low profile, 0.1ml polypropylene wells, polycarbonate frame, cut corner H12; working volume: <100µl, total well capacity: 200µl; designed for use on Roche LightCycler® 96 and 480 (with 96 well block)

- Our FrameStar Roche Style plates are designed to achieve optimized assay conditions on the Roche LightCycler® 96 and 480 (with 96 well block)
- This particular style of plate is in a low profile 96-well format; perfect for reaction volumes of 10-100µl
- The wells are shorter than “standard” profile wells, which decrease the “dead space” between the heated lid of the thermal cycler and the sample in the well
- This eliminates condensation forming on the side wall of the tubes, preventing reduction in PCR volume and increasing the efficiency of the reaction
- This is especially recommended for reaction volumes below 20µl
- The semi-skirted frame allows for labeling or barcoding

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Well

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- <0.1ml (100µl) working volume, 0.2ml (200µl) total well capacity

Frame


- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification

Use

- Compatible with the Roche LightCycler® 96 and 480 (with 96 well block)
- Recommended for low volume PCR
- Ideal for use with robotic systems
- Compatible with standard multichannel pipettes

Options

- Standard plate comes with clear polypropylene wells on a clear polycarbonate frame
- Also available with white wells for optimum signal-to-noise ratio when using fluorescent based assays



- Available barcoded upon request
- Combi packs available with qPCR adhesive seal (4ti-0560) and (for 4ti-0951) with Optically Clear Windowed qPCR Seal (4ti-0564)

Specifications

Parameter	Value
Plate length	127.70 ± 0.25mm
Plate width	85.48 ± 0.25mm
Plate height	15.60 ± 0.25mm
Well depth	15.10 ± 0.10mm
Well diameter	5.50 ± 0.10mm
Distance to center of A1 from top edge	11.24 ± 0.25mm
Distance to center of A1 from left edge	14.38 ± 0.25mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0950/C	FrameStar 96 Well Semi-Skirted PCR Plate, Roche style, clear PP wells, clear PC frame, low profile, cut corner H12, 50 plates per case
4ti-0951	FrameStar 96 Well Semi-Skirted PCR Plate, Roche style, white PP wells, clear PC frame, low profile, cut corner H12, 50 plates per case
Combi Packs	
4ti-0952	FrameStar 96 Well Semi-Skirted PCR Plate, Roche style, plus qPCR Seal, 4TI-0951 plus 4TI-0560, combi pack, 50 plates and seals per case
4ti-0953	FrameStar 96 Well Semi-Skirted PCR Plate, Roche style, plus qPCR Seal, 4TI-0950/C plus 4TI-0560, combi pack, 50 plates and seals per case
4ti-0951/0565	FrameStar 96 Well Semi-Skirted PCR Plate, Roche style, plus Optically clear windowed qPCR seal, 4TI-0951 plus 4TI-0565, combi pack, 50 plates and seals per case



AZENTA
LIFE SCIENCES

4titude FrameStar 96 Well Semi-Skirted PCR Plate, Roche Style, High Sensitivity

Low profile, 0.1 ml polypropylene wells, polycarbonate frame, cut corner H12, working volume: <100 µl, total well capacity: 200 µl; designed for use on Roche LightCycler® 96 and 480 (with 96 well block); extra white wells for improved sensitivity in fluorescent assays

- Our FrameStar Roche Style plates are designed to achieve optimized assay conditions on the Roche LightCycler® 96 and 480 (with 96 well block)
- This particular style of plate is in a low profile 96-well format; perfect for reaction volumes of 10-100 µl
- The wells are shorter than “standard” profile wells, which decrease the “dead space” between the heated lid of the thermal cycler and the sample in the well. This eliminates condensation forming on the side wall of the tubes, preventing reduction in PCR volume and increasing the efficiency of the reaction
- Especially recommended for reaction volumes below 20 µl
- The extra white wells allow for increased sensitivity in fluorescent assays
- The semi-skirted frame allows for labeling or barcoding

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- <0.1 ml (100 µl) working volume, 0.2 ml (200 µl) total well capacity

Frame


- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification

Use

- Compatible with the Roche LightCycler® 96 and 480 (with 96 well block)
- Recommended for low volume PCR
- Compatible with standard multichannel pipettes
- Ideal for use with robotic systems

Options

- Plates are available with a clear polycarbonate frame and extra white wells for improved sensitivity when using fluorescent based assays
- Combi packs available with qPCR Adhesive Seal (4ti-0560)
- Available barcoded upon request



Specifications

Parameter	Value
Plate length	127.70 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	15.60 ± 0.05 mm
Well depth	15.10 ± 0.10 mm
Well diameter	5.50 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ct values Comparison

FrameStar Roche Style, High Sensitivity Plate	
Ct value	17.23
Δ Ct	-3.05
Competitor R equivalent plates	
Ct value	20.28
Δ Ct	0

Ordering Information

4ti-0954	FrameStar 96 Well Semi-Skirted PCR Plate, Roche style, High Sensitivity, extra white PP wells, clear PC frame, low profile, cut corner H12, 50 plates per case
Combi Pack	
4ti-0954/0560	FrameStar 96 Well Semi-Skirted PCR Plate, Roche style, High Sensitivity, extra white PP wells, plus qPCR Seal, 4TI-0954 plus 4TI-0560, combi pack, 50 plates and seals per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude FrameStar 96 Well Semi-Skirted PCR Plate, ABI® FastPlate Style

Low profile, 0.1 ml polypropylene wells, polycarbonate frame with upstand, cut corner A1, working volume: <100 µl, total well capacity: 200 µl; designed for use on ABI® Fast Block cyclers

- This semi-skirted low profile plate is recommended for use with ABI® Fast block thermal cyclers
- The low profile wells of this plate are shorter than “standard” profile wells, which decrease the “dead space” between the heated lid of the thermal cycler and the sample in the well
- This eliminates condensation forming on the side wall of the tubes, preventing reduction in PCR volume and increasing the efficiency of the reaction
- Especially recommended for reaction volumes below 20 µl
- The rigid polycarbonate skirt of this design eliminates warping and distortion during PCR, and makes it ideal for use with robotic systems. Its skirt also allows for labeling or barcoding

Key Features



- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Stackable

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- <0.1 ml (100 µl) working volume, 0.2 ml (200 µl) total well capacity

Frame

- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results



- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification

Use

- Compatible with ABI® Fast block thermal cyclers
- Recommended for low volume PCR
- Compatible with standard multichannel pipettes
- Ideal for use with robotic systems

Options

- Standard plate comes with clear polypropylene wells with a clear polypropylene frame
- Also available with white wells for optimum signal-to-noise ratio when using fluorescent based assays
- Also available as a clear frame with frosted wells for use with ABI®/LifeTechnologies® qPCR instruments
- Combi packs available (for 4ti-0910/C) with Optically Clear Windowed qPCR Seal (4ti-0565)
- Available barcoded upon request

Ordering Information

4ti-0910/C	FrameStar 96 Well Semi-Skirted PCR Plate, ABI Fastplate style, clear PP wells, clear PC frame, with upstand, low profile, cut corner A1, 50 plates per case
4ti-0910/C/10P	FrameStar 96 Well Semi-Skirted PCR Plate, ABI Fastplate style, clear PP wells, clear PC frame, with upstand, low profile, cut corner A1, 10 plates per case
4ti-0911	FrameStar 96 Well Semi-Skirted PCR Plate, ABI Fastplate style, white PP wells, clear PC frame, with upstand, low profile, cut corner A1, 50 plates per case
4ti-0912	FrameStar 96 Well Semi-Skirted PCR Plate, ABI Fastplate style, frosted PP wells, clear PC frame, with upstand, low profile, cut corner A1, 50 plates per case
4ti-0912/10P	FrameStar 96 Well Semi-Skirted PCR Plate, ABI Fastplate style, frosted PP wells, clear PC frame, with upstand, low profile, cut corner A1, 10 plates per case
Combi Pack	
4ti-0910/C/0565	FrameStar 96 Well Semi-Skirted PCR Plate, ABI Fastplate style, plus Optically clear windowed qPCR seal, 4TI-0910/C plus 4TI-0656, combi pack, 50 plates and seals per case

Specifications

Parameter	Value
Plate length	127.70 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	16.70 ± 0.10 mm
Well depth	15.10 ± 0.10 mm
Well diameter	5.50 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

4titude FrameStar 96 Well Semi-Skirted PCR Plate With Upright, ABI® Style

High profile, 0.2 ml polypropylene wells, polycarbonate frame with upstand, cut corner A12, working volume: <200 µl, total well capacity: 300 µl; designed for use on ABI® instruments

- This ABI® Style PCR plate offers the benefits of our 2-component design to ABI® standard block users
- This design combines the advantages of ultra-thin wall polypropylene tubes for optimum PCR results with a rigid polycarbonate skirt and deck for the highest thermal stability
- We recommend this semi-skirted plate for use with ABI® thermal cyclers and sequencers; it can be used directly with ABI® instruments with no adapters necessary
- The only case where this is not true is with the ABI Fast Block thermal cyclers; in this case, using our FrameStar Fast Plate is recommended instead

Key Features



- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time quantitative PCR (RT-qPCR) results
- <0.2 ml (200 µl) working volume, 0.3 ml (300 µl) total well capacity

Frame

- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification



Use

- Ideal for use with ABI® thermal cyclers & sequencers
- Recommended for low volume PCR
- Compatible with standard multichannel pipettes
- Ideal for use with robotic systems

Options

- Available with the following frame color options for the clear well variety: purple, blue, clear, green, red, and black
- Also available as a black frame with white wells for optimum signal-to-noise ratio when using fluorescent based assays
- Ethylene oxide treated option available for forensic use
- Available barcoded upon request
- Combi packs available (for 4ti-0730/C) with Optically Clear Windowed qPCR Seal (4ti-0565)



AZENTA
LIFE SCIENCES

4titude FrameStar 96 Well Semi-Skirted PCR Plate With Upstand, ABI® Style

Specifications

Parameter	Value
Plate length	124.26 ± 0.25 mm
Plate width	83.97 ± 0.25 mm
Plate height	23.20 ± 0.05 mm
Well depth	20.20 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	10.495 ± 0.25 mm
Distance to center of A1 from left edge	12.63 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0730	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, purple PC frame, with upstand, high profile, cut corner A12, 50 plates per case
4ti-0730/B	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, blue PC frame, with upstand, high profile, cut corner A12, 50 plates per case
4ti-0730/C	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, clear PC frame, with upstand, high profile, cut corner A12, 50 plates per case
4ti-0730/C/10P	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, clear PC frame, with upstand, high profile, cut corner A12, 10 plates per case
4ti-0730/G	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, green PC frame, with upstand, high profile, cut corner A12, 50 plates per case
4ti-0730/R	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, red PC frame, with upstand, high profile, cut corner A12, 50 plates per case
4ti-0730/X	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, black PC frame, with upstand, high profile, cut corner A12, 50 plates per case
4ti-0X730C/SBC	FrameStar 96 Well Semi-Skirted PCR Plate, ABI Style, clear PP wells, clear PC frame, with upstand, ethylene oxide treated, single barcoded, high profile, cut corner A12, 20 plates per case
Combi Pack	
4ti-0730/C/0565	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, clear PC frame, with upstand, high profile, cut corner A12, 50 plates per case plus 4ti-0565

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

4titude FrameStar 96 Well Semi-Skirted PCR Plate, ABI® Style

High profile, 0.2 ml polypropylene wells, polycarbonate frame, cut corner A12; working volume: <200 µl, total well capacity: 300 µl; designed for use on all major cyclers, including ABI® instruments with standard 96 well blocks

- Specifically designed to be directly compatible with all major thermal cyclers, this plate can be used directly in ABI® 96 well instruments without the need for any adapters
- The rigid FrameStar 2-component design eliminates warping and distortion during PCR, making it ideal for use with robotic systems
- The semi-skirt allows for labeling or barcoding for sample tracking

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells



- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- <0.2 ml (200 µl) working volume, 0.3 ml (300 µl) total well capacity

Frame

- Cut corner at A12
- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification


Use

- Perfect for ABI® thermal cyclers & sequencers
- Recommended for low volume PCR
- Ideal for use with robotic systems
- Compatible with standard multichannel pipettes



Options

- Available with the following frame color options for the clear well variety: purple, blue, clear, green, red, and black
- Also available as a black frame with white wells for optimum signal-to-noise ratio when using fluorescent based assays
- Also available as a clear frame with frosted wells for use with ABI®/LifeTechnologies® qPCR instruments
- Ultra-low DNA binding option available (4ti-LB0770/C) for sensitive applications with ultra-low DNA input and for maximum DNA recovery after low temperature storage and high temperature incubation; learn more about our low binding range
- Ethylene oxide treated option available (4ti-OX770C/SBC) for forensic use
- Available barcoded upon request
- FrameStar 96 Lid (4ti-0289) available



AZENTA
LIFE SCIENCES

Specifications

Parameter	Value
Plate length	124.26 ± 0.25 mm
Plate width	83.97 ± 0.25 mm
Plate height	20.70 ± 0.10 mm
Well depth	20.20 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	10.495 ± 0.25 mm
Distance to center of A1 from left edge	12.63 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0770	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, purple PC frame, high profile, cut corner A12, 50 plates per case
4ti-0770/B	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, blue PC frame, high profile, cut corner A12, 50 plates per case
4ti-0770/C	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, clear PC frame, high profile, cut corner A12, 50 plates per case
4ti-0770/C/10P	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, clear PC frame, high profile, cut corner A12, 10 plates per case
4ti-0770/G	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, green PC frame, high profile, cut corner A12, 50 plates per case
4ti-0770/R	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, red PC frame, high profile, cut corner A12, 50 plates per case
4ti-0770/X	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, black PC frame, high profile, cut corner A12, 50 plates per case
4ti-0771	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, white PP wells, black PC frame, high profile, cut corner A12, 50 plates per case
4ti-0772	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, frosted PP wells, clear PC frame, high profile, cut corner A12, 50 plates per case
4ti-0772/10P	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, frosted PP wells, clear PC frame, high profile, cut corner A12, 10 plates per case
4ti-0X770C/SBC	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, clear PC frame, with upstand, ethylene oxide treated, single barcoded, high profile, cut corner A12, 20 plates per case
4ti-LB0770/C	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, clear PC frame, low binding, high profile, cut corner A12, 50 plates per case



4titude FrameStar 96 Well Semi-Skirted PCR Plate

High profile, 0.2 ml polypropylene wells, polycarbonate frame, cut corner H1, working volume: <200 µl, total well capacity: 300 µl; universal semi-skirted plate designed for use on standard thermal cyclers

- The rigid FrameStar 2-component design eliminates warping and distortion during PCR, making it ideal for use with robotic systems
- The semi-skirted allows for labeling or barcoding for sample tracking

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- <0.2 ml (200 µl) working volume, 0.3 ml (300 µl) total well capacity

Frame


- Cut corner at H1
- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification

Use

- Recommended for low volume PCR
- Compatible with standard multichannel pipettes
- Ideal for use with robotic systems

Options

- Also available with white wells for optimum signal-to-noise ratio when using fluorescent based assays
- Available barcoded upon request
- Similar plate with a cut corner at A12 for use with ABI® thermal cyclers and sequencers available: FrameStar 96 Well Semi-Skirted PCR Plate, ABI® Style



Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	20.70 ± 0.25 mm
Well depth	20.20 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0900/C	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, clear PC frame, high profile, cut corner H1, 50 plates per case
4ti-0901	FrameStar 96 Well Semi-Skirted PCR Plate, white PP wells, clear PC frame, high profile, cut corner H1, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude FrameStar 96 Well Non-Skirted PCR Plate

High profile, 0.2 ml polypropylene wells, polycarbonate frame, cut corner H1, working volume: <200 µl, total well capacity: 300 µl; universal non-skirted plate designed for use on all major thermal cyclers

- The rigid FrameStar 2-component design eliminates warping and distortion during PCR, meaning seal integrity is not compromised and less of your sample is lost through evaporation

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Greatest compatibility with different thermal cyclers

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- <0.2 ml (200 µl) working volume, 0.3 ml (300 µl) total well capacity

Frame



- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification

Use

- Recommended for low volume PCR
- Compatible with standard multichannel pipettes

Options

- Available with the following frame color options for the clear well variety: purple, blue, clear, green, red, and black
- Also available as a black frame with white wells for optimum signal-to-noise ratio when using fluorescent based assays



Specifications

Parameter	Value
Plate length	120.00 ± 0.25 mm
Plate width	80.00 ± 0.25 mm
Plate height	20.70 ± 0.10 mm
Well depth	20.20 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	8.50 ± 0.25 mm
Distance to center of A1 from left edge	10.50 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0710	FrameStar 96 Well Non-Skirted PCR Plate, clear PP wells, purple PC frame, high profile, cut corner H1, 50 plates per case
4ti-0710/B	FrameStar 96 Well Non-Skirted PCR Plate, clear PP wells, blue PC frame, high profile, cut corner H1, 50 plates per case
4ti-0710/C	FrameStar 96 Well Non-Skirted PCR Plate, clear PP wells, clear PC frame, high profile, cut corner H1, 50 plates per case
4ti-0710/G	FrameStar 96 Well Non-Skirted PCR Plate, clear PP wells, green PC frame, high profile, cut corner H1, 50 plates per case
4ti-0710/R	FrameStar 96 Well Non-Skirted PCR Plate, clear PP wells, red PC frame, high profile, cut corner H1, 50 plates per case
4ti-0711	FrameStar 96 Well Non-Skirted PCR Plate, white PP wells, black PC frame, high profile, cut corner H1, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azena Plate Instrument Compatibility Table Pages 162-165.

4titude FrameStar 96 Well Non-Skirted PCR Plate, Low Profile

Low profile, 0.1 ml polypropylene wells, polycarbonate frame, cut corner H1, working volume: <100 µl, total well capacity: 200 µl; universal non-skirted, low profile plate designed for use on all major thermal cyclers

- The low profile wells of this plate are shorter than “standard” profile wells, which decrease the “dead space” between the heated lid of the thermal cycler and the sample in the well
- This eliminates condensation forming on the side wall of the tubes, preventing reduction in PCR volume and increasing the efficiency of the reaction
- This is especially recommended for reaction volumes below 20 µl
- The rigid polycarbonate skirt of this design eliminates warping and distortion during PCR, meaning seal integrity is not compromised and less of your sample is lost through evaporation

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Greatest compatibility with different thermal cyclers

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- <0.1 ml (100 µl) working volume, 0.2 ml (200 µl) total well capacity

Frame


- Rigid polycarbonate frame for added mechanical stability
- Significantly reduced thermal expansion and sample evaporation for improved consistency in PCR results
- Improved seal integrity due to thermal stability of frame
- Alphanumeric grid reference to aid well and sample identification

Use

- Recommended for low volume PCR
- Compatible with standard multichannel pipettes

Options

- Available with the following frame color options for the clear well variety: purple, blue, and clear
- Also available as a black frame with white wells for optimum signal-to-noise ratio when using fluorescent based assays



Specifications

Parameter	Value
Plate length	120.00 ± 0.25 mm
Plate width	80.00 ± 0.25 mm
Plate height	15.60 ± 0.05 mm
Well depth	15.10 ± 0.10 mm
Well diameter	5.50 ± 0.10 mm
Distance to center of A1 from top edge	18.50 ± 0.25 mm
Distance to center of A1 from left edge	10.50 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0720	FrameStar 96 Well Non-Skirted PCR Plate, clear PP wells, purple PC frame, low profile, cut corner H1, 50 plates per case
4ti-0720/B	FrameStar 96 Well Non-Skirted PCR Plate, clear PP wells, blue PC frame, low profile, cut corner H1, 50 plates per case
4ti-0720/C	FrameStar 96 Well Non-Skirted PCR Plate, clear PP wells, clear PC frame, low profile, cut corner H1, 50 plates per case
4ti-0721	FrameStar 96 Well Non-Skirted PCR Plate, clear PP wells, black PC frame, low profile, cut corner H1, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude FrameStar Breakable PCR Plates

FrameStar Breakable PCR Plates

FrameStar Breakable PCR Plates can be easily divided into smaller plate sections, ensuring no tubes are wasted.

The plates combine the advantages of the FrameStar and 8 Well PCR Tube Strip with PC Frame range as well as single tube formats. The result is PCR consumables with thin-walled polypropylene (PP) tubes for optimal PCR results with a rigid polycarbonate frame for easy and reliable handling.

The two-component design minimizes evaporation allowing for the downscaling of reaction volumes and the breakability of the plates offers flexibility to suit your experiment size. Once broken down, strips remain straight and stable for ease of handling and to enable effortless sample tracking if 2D coded.

FrameStar Breakable PCR Plates - Division brings flexibility!

- Dividable vertically, or both vertically and horizontally
- Flexible solutions for every application
- Available as high and low profile plates to suit your reaction volume and instrument format
- Highest instrument compatibility
- Seven frame colors with clear or white tubes available
- Color coding for different workflows
- 2D coded options available
- Error-free sample tracking
- Plate segmentation can be automated
- Also available as pre-cut strips, see 8 Well PCR Tube Strip with PC Frame, page 136


Sealing options

Plates can be sealed with standard heat or adhesive seals and then cut to produce individually sealed strips of wells.


Alternatively, perforated seals or Individual access seals can be used depending on individual tube, 8 or 12 well strips required. Please refer to azenta.com for information on semi-automated, fully automated and individual access heat sealers.




FrameStar Vertically Breakable & Vertically and Horizontally Breakable PCR Plates can be easily divided into smaller plate sections, ensuring no tubes are wasted.



FrameStar Breakable PCR Plates offer you all the flexibility of tube strips, in a plate format




FrameStar Breakable PCR Plates, dividable both horizontally and vertically




AZENTA
LIFE SCIENCES

4titude FrameStar Vertically Breakable PCR Plate




High profile

0.2 ml PP wells, working volume:
<200 μ l, total well capacity: 300 μ lAC



Low profile

0.1 ml PP wells, working volume:
<100 μ l, total well capacity: 200 μ l



96 well semi-skirted plate, vertically scored, snaps easily into strips of 8 tubes or part plates, cut corner A12

- Our FrameStar Vertically Breakable PCR Plates utilize the 2-component design of the FrameStar range, which combines the advantages of thin-walled polypropylene tubes for optimum PCR results and a rigid frame portion for easy and reliable handling

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- 2-component design prevents distortion of tube strips

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- Standard 0.2 ml (200 μ l) working volume, with a 0.3 ml (300 μ l) total well capacity when used with sealing options
- Low profile 0.1ml (100 μ l) working volume, with a 0.2ml (200 μ l) total well capacity when used with sealing options

Frame

- End tabs for easy handling and labeling
- Eliminates strip breakage
- Alphanumeric grid reference to aid well and sample identification

Use

- Rigid PCR plate that can be broken into smaller plate sections
- Fits majority of thermal cyclers
- Compatible with standard multichannel

Options

- Available with the following frame color options for the clear well variety: purple, blue, clear, green, red, black, and white
- Also available as a black frame with white wells for use with optical assays such as qPCR
- 2D coding option
- 2D code reader available

Specifications

Parameter	Value (High Profile)	Value (Low Profile)
Plate length	125.11 \pm 0.25 mm	125.11 \pm 0.25 mm
Plate width	83.22 \pm 0.25 mm	83.22 \pm 0.25 mm
Plate height	20.80 \pm 0.05 mm	15.60 \pm 0.25 mm
Well depth	20.30 \pm 0.10 mm	15.10 \pm 0.10 mm
Well diameter	5.46 \pm 0.10 mm	5.50 \pm 0.10 mm
Distance to center of A1 from top edge	10.11 \pm 0.25 mm	10.11 \pm 0.25 mm
Distance to center of A1 from left edge	13.06 \pm 0.25 mm	13.06 \pm 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm	9.00 mm

2D coding is available for this product. Contact us for more information



AZENTA
LIFE SCIENCES

Ordering Information

FrameStar Vertically Breakable PCR Plate High profile

4ti-1000/P	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, purple PC frame, vertically breakable, high profile, cut corner A12, 50 plates per case
4ti-1000/B	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, blue PC frame, vertically breakable, high profile, cut corner A12, 50 plates per case
4ti-1000/C	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, clear PC frame, vertically breakable, high profile, cut corner A12, 50 plates per case
4ti-1000/G	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, green PC frame, vertically breakable, high profile, cut corner A12, 50 plates per case
4ti-1000/R	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, red PC frame, vertically breakable, high profile, cut corner A12, 50 plates per case
4ti-1000/X	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, black PC frame, vertically breakable, high profile, cut corner A12, 50 plates per case
4ti-1000/W	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, white PC frame, vertically breakable, high profile, cut corner A12, 50 plates per case
4ti-1001	FrameStar 96 Well Semi-Skirted PCR Plate, white PP wells, black PC frame, vertically breakable, high profile, cut corner A12, 50 plates per case

FrameStar Vertically Breakable PCR Plate Low profile


4ti-1200/P	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, purple PC frame, vertically breakable, low profile, cut corner A12, 50 plates per case
4ti-1200/B	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, blue PC frame, vertically breakable, low profile, cut corner A12, 50 plates per case
4ti-1200/C	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, clear PC frame, vertically breakable, low profile, cut corner A12, 50 plates per case
4ti-1200/G	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, green PC frame, vertically breakable, low profile, cut corner A12, 50 plates per case
4ti-1200/R	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, red PC frame, vertically breakable, low profile, cut corner A12, 50 plates per case
4ti-1200/X	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, black PC frame, vertically breakable, low profile, cut corner A12, 50 plates per case
4ti-1200/W	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, white PC frame, vertically breakable, low profile, cut corner A12, 50 plates per case
4ti-1201	FrameStar 96 Well Semi-Skirted PCR Plate, white PP wells, black PC frame, vertically breakable, low profile, cut corner A12, 50 plates per case

FrameStar Breakable PCR plates are not compatible with the ABI 9700 dual block thermal cycler.

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.




4titude FrameStar Vertically and Horizontally Breakable PCR Plate




High profile

0.2 ml PP wells, working volume:
<200 µl, total well capacity: 300 µl



Low profile

0.1 ml PP wells, working volume:
<100 µl, total well capacity: 200 µl



96 well rigid plate, vertically and horizontally scored, snaps easily into part plates, 8 well strips, 12 well strips, part strips or individual tubes

- Azenta's FrameStar Vertically and Horizontally Breakable PCR Plates allow for the most flexible, efficient and cost-effective use of FrameStar PCR plates, ensuring not a single tube is wasted
- Our Vertically and Horizontally Breakable PCR Plates utilize the 2-component design of the FrameStar range, which combines the advantages of thin-walled polypropylene tubes for optimum PCR results and a rigid frame portion for easy and reliable handling

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- Standard 0.2 ml (200 µl) working volume, with a 0.3 ml (300 µl) total well capacity when used with sealing options
- Low profile, 0.1ml wells, <100µl working volume, with a 200 µl total well capacity when used with sealing options

Specifications

Parameter	Value (High Profile)	Value (Low Profile)
Plate length	125.00 ± 0.25 mm	125.00 ± 0.25 mm
Plate width	72.0 ± 0.25 mm	72.0 ± 0.25 mm
Plate height	20.70 ± 0.05 mm	15.60 ± 0.15 mm
Well depth	20.20 ± 0.10 mm	15.10 ± 0.15 mm
Well diameter	5.46 ± 0.10 mm	5.50 ± 0.10 mm
Distance to center of A1 from top edge	4.5 ± 0.25 mm	4.50 ± 0.25 mm
Distance to center of A1 from left edge	8.5 ± 0.25 mm	8.50 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm	9.00 mm

Use

- Plates can be broken up in both 8-strip vertical and 12-strip horizontal directions, resulting in individual strips, part strips, or even individual tubes, giving you the highest range of flexibility
- Breaking a plate is more accurate, more convenient, and safer than cutting it with scissors, as it avoids damaging the sealing rings and contamination of the wells
- Plates can be filled, sealed and separated for storage, processing or distribution
- Easy to break at any temperature
- The plate can be adapted to individual pipetting schemes without the need to waste empty wells
- Separated strips or tubes can be used for the positive control to avoid contamination of the samples
- Compatible with all instruments that fit non-skirted, high profile plates
- NB: In some cases you may have to break off the end tabs for it to fit
- Compatible with standard multichannel pipettes

Options

- 2D coding option
- 2D code reader available

Ordering Information

4ti-1300/X	FrameStar 96 Well Semi-Skirted PCR Plate, clear PP wells, black PC frame, breakable vertically and horizontally, high profile, 50 plates per case
4ti-1400/X	FrameStar 96 Well Semi-skirted PCR Plate, clear PP wells, black PC frame, breakable vertically and horizontally, low profile, 50 plates per case

In some cases it may be required to break off the end tabs to allow for a perfect fit.

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

4titude FrameStar Individual Access 96 Well PCR Plates

FrameStar is our superior technology of making PCR plates with ultra-thin polypropylene wells fitted into a robust polycarbonate frame that provides excellent stability. Our Individual Access plates develop this technology further to supply a novel 96 well plate with individually removable wells combining both flexibility and robustness.

Individual Access plates provide flexibility and robustness for your assay setups. 96 well plate format is suitable for high throughput reagents dispensing on liquid handling robots, while individually sealable and removable wells provide flexibility to end users to accommodate varying throughputs.

The wells are made from medical grade polypropylene which is perfectly suited for use in PCR as well as for long term storage. The exact number of wells required can be used, meaning no wastage of consumables or reagents. Empty frames are also available for tubes to be transferred to, however, frames can be re-used multiple times, depending on their application.


Each well clicks into place within the frame, holding it securely for use with automation and for transport. Additionally, the fit of the tubes is not compromised following a PCR run due to resistance of the rigid frame to thermal expansion.

- **SBS footprint** – suitable for automation
- **Rigid frame** – does not expand during PCR cycles allowing for individual tubes to fit tightly in the frame before and after PCR
- **PP wells** – low binding to nucleic acids and high solvent resistance, ideal for both PCR and storage
- **Thin walled tubes** – optimal heat transfer during PCR
- **Individually sealed using individual seals, specifically designed individual plates** – tubes can be filled, sealed and then single tubes removed
- **Tubes can be removed and inserted again** – ultimate flexibility

Individual Access plates provide the greatest flexibility in PCR plates, affording handling of individual wells in a 96 well plate format.



Individual Access 96 Well Skirted PCR Plate



4titude Individual Access 96 Well Skirted PCR Plates

Low profile, individually removable 0.1 ml polypropylene wells, rigid polycarbonate frame, cut corner H1, working volume: <100 µl, total well capacity: 200 µl

- 96 well, fully skirted, low profile format with each tube inserted separately into the plate frame, allowing for selection and removal of individual tubes from the plate
- Each tube clicks into place in the frame ensuring tubes are secure once inserted
- The rigid polycarbonate skirt of this design eliminates warping and distortion during PCR
- Together with the SBS footprint, this makes the plate suitable for use with robotic systems, ideal for both small to high throughput labs
- Empty frames are available for tubes to be transferred to
- The frames can be re-used multiple times, depending on the application

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Tubes



- Individually removable from the rigid frame, while the rest of the tubes can remain in storage, reducing freeze-thawing of the whole plate
- Thin walled tubes are made from PP for optimal results in (q)PCR whilst also being suitable for storage and incubation
- <100 µl working volume, <200 µl total well capacity

Frame

- Rigid polycarbonate frame provides mechanical stability of the plate
- Alphanumeric grid reference to aid well and sample identification

Use

- Ideal for use with robotic systems
- Fewer transfer steps as the same tube can be used for both (q)PCR and storage
- Each tube can still be separated after sealing
- Each tube can be individually heat sealed using the Individual Access Pierce Heat Seal Strong, such as the contents of each tube is left secure from cross contamination and evaporation for storage and transfer



Options

- Empty frames are available for tubes to be transferred to; frames can be re-used multiple times, depending on the application
- The plate can be sealed in one step using our range of Individual Access Heat Seals on our fully or semi-automated heat sealers (using the 4ti-0613 Individual Access adapter), resulting in individually sealed tubes with pierceable seals to allow for sample access
- Available barcoded upon request
- We also offer customization of the Individual Access Plates such as color coding of the individual wells; please contact us for details

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	15.50 ± 0.25 mm
Well depth	15.10 ± 0.10 mm
Well diameter	5.50 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0960/RA	Individual Access 96 Well Skirted PCR Plate, clear PP wells, white rigid PC frame, individually removable wells, low profile, cut corner H1, 50 plates per case
4ti-0960/RA/F	Individual Access 96 Well Skirted PCR Plate Frame, for individually removable wells, white rigid PC frame, cut corner H1, 10 frames per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude Individual Access 96 Well Skirted Optical Bottom PCR Plates

Low profile, flat optical bottom, 0.1 ml clear polypropylene tubes, black rigid polycarbonate frame, cut corner H1, working volume: <100 µl, total well capacity: 180 µl

- 96 well, fully skirted, low profile format with each tube inserted separately into the plate frame, allowing for selection and removal of individual tubes from the plate
- Each tube clicks into place in the frame ensuring tubes are secure once inserted
- The Individual Access 96 Well Skirted Optical Bottom PCR Plates (next to their counterpart within our FrameStar range, 4ti-0970) are currently unique in the market, being suitable for use in both microscopy and PCR, for example, when single cell sorting is followed by molecular biology applications such as qPCR and sequencing
- Due to their flat bottoms and stackability, these plates are well suited for small sample volume storage (such as compound libraries), with no risk of damaging the seal of the plate below
- Additionally, the small well volume enables excellent sample recovery
- The rigid polycarbonate skirt of this design eliminates warping and distortion during PCR
- Together with the SBS footprint, this makes the plate suitable for use with robotic systems, ideal for both small and high-throughput labs
- Empty frames are available for tubes to be transferred to. The frames can be re-used multiple times, depending on the application

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Tubes

- Individually removable from the rigid frame, while the rest of the tubes can remain in storage, reducing freeze-thawing of the whole plate
- Thin walled tubes are made from PP for optimal results in (q)PCR whilst also being suitable for storage and incubation
- <100 µl working volume, <180 µl total well capacity

Frame

- Rigid polycarbonate frame provides mechanical stability of the plate
- Alphanumeric grid reference to aid well and sample identification

Use

- Suitable for microscopy and small volume sample storage
- Ideal for use with robotic systems
- Fewer transfer steps as the same tube can be used for microscopy, (q)PCR and storage
- Each tube can still be separated after sealing

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



- Each tube can be individually heat sealed using the Individual Access Pierce Heat Seal Strong, such as the contents of each tube is left secure from cross contamination and evaporation for storage and transfer

Options

- Also available as a FrameStar plate (4ti-0970) with fixed wells
- Empty frames are available for tubes to be transferred to; frames can be re-used multiple times, depending on the application
- We also offer customization of the Individual Access Plates such as color coding of the individual wells; please contact us for details
- Available barcoded upon request, please contact us for details
- The plate can be sealed in one step using our range of Individual Access Heat Seals on our Automated Individual Access Heat Sealer (see page 211) or our Semi-Automated Sheet Heat Sealer (see page 213) (using the 4ti-0613 Individual Access adapter), resulting in individually sealed tubes with pierceable seals to allow for sample access

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	15.50 ± 0.25 mm
Well depth	12.00 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0970/RA	Individual Access 96 Well Skirted Optical Bottom PCR Plate, clear PP wells, black PC frame, individually removable wells, low profile, cut corner H1, 50 plates per case
4ti-0970/RA/F	Individual Access 96 Well Skirted PCR Plate Frame, for individually removable wells, black rigid PC frame, cut corner H1, 10 frames per case



AZENTA
LIFE SCIENCES

4titude Individually Removable Well 96 Well Skirted Flat Bottom PCR Plate, 2D Coded

Low profile, flat bottom 2D coded, 0.1 ml clear polypropylene tubes, white rigid polycarbonate frame, linear barcoded, cut corner H1, working volume: <100 µl, total well capacity: 180 µl

- 96 well, fully skirted, low profile format with each tube inserted separately into the plate frame, allowing for selection and removal of individual tubes from the plate
- Each tube clicks into place in the frame ensuring tubes are secure once inserted
- We have now extended our Individual Access range to include PCR plates with 2D coded flat bottoms for applications that require superior sample tracking
- In applications where it is essential for samples to be tracked throughout their processing, e.g. for clinical samples, tube labeling is a far safer and reliable method compared to cap or seal labeling (which can be misplaced, damaged or rendered unreadable after piercing)
- Our plates have a unique 2D code printed on the bottom of the well and a linear barcode on the side of the plate to allow for identification of both the plate and individual tubes within
- The rigid polycarbonate skirt of this design eliminates warping and distortion during PCR
- Together with the SBS footprint, this makes the plate suitable for use with robotic systems, ideal for both small and high throughput labs
- Empty frames are available for tubes to be transferred to
- The frames can be re-used multiple times, depending on the application
- Sample tubes can be identified by the unique 2D code on the base of the wells, removed from the Individual Access plate quickly and easily, and moved to another Individual Access frame for downstream processing or storage
- Quick selection of desired tubes reduces the time needed for the sample plate to be out of the freezer and therefore reduces the chance of other wells defrosting in the meantime
- Flat bottom PCR plates can be stacked, allowing for optimal use of freezer space as multiple plates can be housed in the same space as a storage rack
- The 2D codes on the bottom of the wells can be read by most 2D data-matrix readers including Azenta Rack Readers
- The codes are highly scratch resistant and can withstand cold storage (-80°C), temperatures up to 100°C and solvents such as DMSO
- 2D data-matrix codes utilize data redundancy so even if codes are partly destroyed, the information will still be retained
- Each code is tested for readability and guaranteed to be unique



Key Features


- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Tubes

- Individually removable from the rigid frame, while the rest of the tubes can remain in storage, reducing freeze-thawing of the whole plate
- Thin walled tubes are made from PP for optimal results in (q)PCR whilst also being suitable for storage and incubation
- <100 µl working volume, <180 µl total well capacity

Frame

- Rigid polycarbonate frame provides mechanical stability of the plate
- These plates are labeled with linear Code 128-format barcodes to aid identification and traceability for your samples
- Alphanumeric grid reference to aid well and sample identification



AZENTA
LIFE SCIENCES

4titude Individually Removable Well 96 Well Skirted Flat Bottom PCR Plate, 2D Coded

Use

- Suitable for small volume sample storage, allowing for optimal use of freezer space as multiple plates can be housed in the same space as a storage rack
- Ideal for use with robotic systems
- Fewer transfer steps as the same tube can be used for microscopy, (q)PCR and storage
- Each tube can be individually heat sealed using the Individual Access Pierce Heat Seal Strong, such as the contents of each tube is left secure from cross contamination and evaporation for storage and transfer
- Each tube can still be separated after sealing

Options

- Also available as a Individual Access optical bottom plate (4ti-0970/RA) with black frame, ideal for microscopy and PCR
- Empty frames are available for tubes to be transferred to; frames can be re-used multiple times, depending on the application
- We also offer customization of the Individual Access Plates such as color coding of the individual wells; please contact us for details
- The plate can be sealed in one step using our range of Individual Access Heat Seals on our Automated Individual Access Heat Sealer (see page 211) or on our Semi-Automated Sheet Heat Sealer (see page 213) (using the 4ti-0613 Individual Access adapter), resulting in individually sealed tubes with pierceable seals to allow for sample access

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	15.50 ± 0.25 mm
Well depth	12.00 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm
Code	
Size	12 x 12 px; 2 x 2 mm
Format	white on black
Content	8 digit numeric

Ordering Information

4ti-0975/RA	Individual Access 96 Well Skirted Flat Bottom PCR Plate, clear PP wells, white PC frame, individually removable wells, 2D coded, low profile, cut corner H1, 50 plates per case
4ti-0960/RA/F	Individual Access 96 Well Skirted PCR Plate Frame, for individually removable wells, white rigid PC frame, cut corner H1, 10 frames per case

**Please note, /SBC must be added to 4ti-0960/RA-F to include a standard single barcode.*

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude Individual Access 96 Well Non-Skirted PCR Plate

Low profile, individually removable 0.1 ml polypropylene wells, rigid polycarbonate frame, non-skirted, cut corner H1, working volume: <100 µl, total well capacity: 200 µl

- 96 well, fully skirted, low profile format with each tube inserted separately into the plate frame, allowing for selection and removal of individual tubes from the plate.
- Each tube clicks into place in the frame ensuring tubes are secure once inserted
- The rigid polycarbonate skirt of this design eliminates warping and distortion during PCR
- SBS footprint,
- When used with Non-Skirted PCR Plate Adapter (4ti-0373), it is suitable for use with robotic systems
- Empty frames are available for tubes to be transferred to
- The frames can be re-used multiple times, depending on the application



Specifications

Parameter	Value
Format	96 Well
Plate length	120.00 ± 0.20 mm
Plate width	80.00± 0.20 mm
Plate height	15.60 mm
Well diameter	5.50 mm
Color (frame)	White
Color (tube)	Clear

Ordering Information

4ti-0720/RA	Individual Access 96 Well Non-Skirted PCR Plate, clear PP individually removable wells, white PC frame, low profile, cut corner H1, 50 plates per case
--------------------	---

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

Individual Access Sealing Options


The 96 well Individual Access plate can be sealed in one step using Individual Access seals. These seals result in individually sealed tubes that can be pierceable, allowing for sample access.

Automated Individual Access Heat Sealing

The roll-fed Automated Individual Access Heat Sealer allows for automatic sealing of Individual Access plates using Individual Access Heat Seal rolls.

The roll with the indexed groups of 96 individual sealing discs is automatically fed through the heat sealer. The accurate sealing is controlled by a sensor which gets activated by optical windows in the material feed but can also be adjusted. Sealing temperature, time of sealing and exit delay (for cooling) can be controlled via the instrument's touchscreen.

Custom versions of instrument and sealing material are possible.



The Automated Individual Access Heat Sealer enables researchers to leverage the benefits of the Azenta Individual Access plate and seal range whilst maintaining the gold standard heat sealing.

- Compatible with a wide range of 96 well plate types and tube racks, allowing for the sealing of individually accessible wells and tubes
- Flexibility of seal material choice: Azenta offers a wide range of sealing materials
- Ability to seal custom shapes with custom seals
- Reproducible sealing: fixed pressure and accurate time and temperature controls ensure consistent sealing
- No compressed air required offering flexibility in labs



Individual Access Sealing Options

Semi-Automated Individual Access Heat Sealing

Individual Access Heat Seals are currently available in sheet format for use with the Semi-Automated Sheet Heat Sealer (using the 59-2005 Individual Access adapter).

Individual Access Sealing Procedure using the Semi-Automated Sheet Heat Sealer (59-2000)

1. Place adapter in open drawer of Sealer
2. Place Individual Access Plate on adapter
3. Place Individual Access Seal onto the plate
4. Plate is automatically sealed, remove plate
5. Remove the backing liner from the seal



Cap Mats for PCR Plates

- 96 individual caps in sheet format, blue TPE, pierceable; suitable for sealing all of our 96 well PCR plates
- The caps can be individually applied and removed once detached from the backing liner, making the mats ideally suited for use with our flexible PCR consumables, including Individual Access and divisible plates
- The mats offer an alternative to adhesive and heat sealing, in particular as a temporary solution when samples need to be repeatedly accessed
- They are easily pierceable with pipette tips to access samples, and they are easily removable using 1- and 8-way decappers or, alternatively, using the Azenta Automated Plate Seal Remover if a seal is overlaid on top of the caps

Ordering Information

Semi-Automated Individual Access Heat Sealing

59-2000	Semi-Automated Sheet Heat Sealer , includes adapter (59-2001)
59-2005	Semi-Automated Sheet Heat Sealer Adapter , for Individual Access plates, 1 adapter per case
4ti-05381/RA	Individual Access Pierce Heat Seal Strong , strong heat sealing foil, 96 individual seals in sheet format, 100 sheets (127 x 100mm) per case
4ti-0531/RA	Individual Access Pierce Heat Seal , pierceable heat sealing foil, 96 individual seals in sheet format, 100 sheets (127 x 100mm) per case
4ti-0521/RA-TAB	Individual Access Peel Heat Seal, with tabs , peelable heat sealing foil, 96 individual seals with tabs, sheet format, 100 sheets (127 x 100mm) per case
4ti-0521/RA-8	Individual Access Peel Heat Seal, peelable heat sealing foil , 12 strips, each covering 8 wells, sheet format, 100 sheets (127 x 100mm) per case
4ti-0960/RA	Individual Access 96 Well Skirted PCR Plate , clear PP wells, white rigid PC frame, individually removable wells, low profile, cut corner H1, 50 plates per case
4ti-0753/757	96 Well Skirted PCR Plate with Removable 8 Well Strips , white PP wells, white PC frame, low profile, 50 plates per case
4ti-1200/P	FrameStar 96 Well Semi-Skirted PCR Plate , clear PP wells, purple PC frame, breakable vertically, low profile, cut corner A12, 50 plates per case
4ti-1400/X	FrameStar 96 Well Semi-skirted PCR Plate , clear PP wells, black PC frame, breakable vertically and horizontally, low profile, cut corner A12, 50 plates per case

Automated Individual Access Heat Sealing

59-1000	Automated Individual Access Roll Heat Sealer , for heat sealing of individual tubes or custom shaped consumables, for use with individual seals in roll format, includes adapter A (59-1004)
4ti-0539/RA	Individual Access Pierce Heat Seal Strong , strong heat sealing foil, 96 individual seals in roll format, 1 roll (420m x 100mm)
4ti-0532/RA	Individual Access Pierce Heat Seal , pierceable heat sealing foil, 96 individual seals in roll format, 1 roll (420m x 100mm)
4ti-0522/RA-TAB	Individual Access Peel Heat Seal , with tabs, peelable heat sealing foil, 96 individual seals with tabs, roll format, 1 roll (420m x 100mm)
4ti-0522/RA-8	Individual Access Peel Heat Seal , peelable heat sealing foil, 12 strips of 8 individual seals in roll format, 1 roll (420m x 100mm)
4ti-0960/RA	Individual Access 96 Well Skirted PCR Plate , clear PP wells, white rigid PC frame, individually removable wells, low profile, cut corner H1, 50 plates per case
4ti-0753/757	96 Well Skirted PCR Plate with Removable 8 Well Strips , white PP wells, white PC frame, low profile, 50 plates per case
4ti-1200/P	FrameStar 96 Well Semi-Skirted PCR Plate , clear PP wells, purple PC frame, breakable vertically, low profile, cut corner A12, 50 plates per case
4ti-1400/X	FrameStar 96 Well Semi-skirted PCR Plate , clear PP wells, black PC frame, breakable vertically and horizontally, low profile, cut corner A12, 50 plates per case



PCR Tubes & Strips



4titude PCR Tubes, Strips & Caps

Strips of 8, clear polypropylene PCR tubes (0.2 mL or 0.1 mL) held in a rigid polycarbonate frame.

By molding the frame portion in a more rigid polymer, the mechanical stability is greatly improved compared with traditional single piece products, as seen here.

For color coding of experiments, we offer strips with clear wells in 6 different frame colors. Additionally, for optical assays such as qPCR, we supply strips with white wells in black frames. White wells increase the signal-to-noise ratio by maximizing reflection of light in fluorescent based assays.

PCR Tubes & Strips are available either with cap strips (domed or optically flat) or without cap strips, and are compatible with the majority of thermal cyclers. End tabs allow for easy handling and labelling of the strips and some products are also available with an off-the-shelf 2D code, offering a vast supply of unique code combinations.


8 Well PCR Tube Strips

Features


- 2-component design in an 8 well strip format
- Compatible with the majority of thermal cyclers
- Available with either domed cap strips or with flat, optically clear cap strips
- Available with off-the-shelf 2D code



8 Well PCR Tube Strips,
Plus Strips of Domed Caps



8 Well PCR Tube Strips, Plus
Strips of Flat Optical Caps



8 Well PCR Tube Strips,
Low Profile

8 Well PCR Tube Strips feature a rigid polycarbonate frame for highest mechanical stability




Unlike standard tube strips, the PCR Tubes & Strips will remain straight and stable, even at elevated temperatures and when filled with liquid.

Ordering Information

High profile + Domed Caps

4ti-0785/P	8 Well PCR Tube Strips, clear PP wells, purple PC frame, high profile, plus strips of domed caps, 120 tube strips and cap strips per case
4ti-0785/G	8 Well PCR Tube Strips, clear PP wells, green PC frame, high profile, plus strips of domed caps, 120 tube strips and cap strips per case
4ti-0785/B	8 Well PCR Tube Strips, clear PP wells, blue PC frame, high profile, plus strips of domed caps, 120 tube strips and cap strips per case
4ti-0785/R	8 Well PCR Tube Strips, clear PP wells, red PC frame, high profile, plus strips of domed caps, 120 tube strips and cap strips per case
4ti-0785/X	8 Well PCR Tube Strips, clear PP wells, black PC frame, high profile, plus strips of domed caps, 120 tube strips and cap strips per case
4ti-0785/X/2D	8 Well PCR Tube Strips, clear PP wells, black PC frame, 2D coded, high profile, plus strips of domed caps, 120 tube strips and cap strips per case
4ti-0785/W/2D	8 Well PCR Tube Strips, clear PP wells, white PC frame, 2D coded, high profile, plus strips of domed caps, 120 tube strips and cap strips per case
4ti-0785/XW	8 Well PCR Tube Strips, white PP wells, black PC frame, high profile, plus strips of domed caps, 120 tube strips and cap strips per case
4ti-0785/M	8 Well PCR Tube Strips, clear PP wells, assorted colors PC frame, high profile, plus strips of domed caps, 120 tube strips and cap strips per case
4ti-0785/C	8 Well PCR Tube Strips, clear PP wells, clear PC frame, high profile, plus strips of domed caps, 120 tube strips and cap strips per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

8 Well PCR Tube Strips

Ordering Information

High profile + Flat Caps

4ti-0786/P	8 Well PCR Tube Strips, clear PP wells, purple PC frame, plus strips of flat caps, high profile, 120 tube strips and cap strips per case
4ti-0786/G	8 Well PCR Tube Strips, clear PP wells, green PC frame, plus strips of flat caps, high profile, 120 tube strips and cap strips per case
4ti-0786/B	8 Well PCR Tube Strips, clear PP wells, blue PC frame, high profile, plus strips of flat caps, 120 tube strips and cap strips per case
4ti-0786/R	8 Well PCR Tube Strips, clear PP wells, red PC frame, plus strips of flat caps, high profile, 120 tube strips and cap strips per case
4ti-0786/X	8 Well PCR Tube Strips, clear PP wells, black PC frame, plus strips of flat caps, high profile, 120 tube strips and cap strips per case
4ti-0786/X/2D	8 Well PCR Tube Strips, clear PP wells, black PC frame, plus strips of flat caps, high profile, 2D coded, 120 tube strips and cap strips per case
4ti-0786/W/2D	8 Well PCR Tube Strips, clear PP wells, white PC frame, plus strips of flat caps, high profile, 2D coded, 120 tube strips and cap strips per case
4ti-0786/XW	8 Well PCR Tube Strips, white PP wells, black PC frame, plus strips of flat caps, high profile, 120 tube strips and cap strips per case
4ti-0786/M	8 Well PCR Tube Strips, clear PP wells, assorted colors PC frame, plus strips of flat caps, high profile, 120 tube strips and cap strips per case
4ti-0786/C	8 Well PCR Tube Strips, clear PP wells, clear PC frame, high profile, plus strips of flat caps, 120 tube strips and cap strips per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

Low Profile (no caps)

4ti-0789/B	8 Well PCR Tube Strips, clear PP wells, blue PC frame, low profile, 120 tube strips per case
4ti-0789/XW	8 Well PCR Tube Strips, white PP wells, black PC frame, low profile, 120 tube strips per case
4ti-0789/W	8 Well PCR Tube Strips, clear PP wells, white PC frame, low profile, 120 tube strips per case

High profile (no caps)

4ti-0775/P	8 Well PCR Tube Strips, clear PP wells, purple PC frame, high profile, 120 tube strips per case
4ti-0775/G	8 Well PCR Tube Strips, clear PP wells, green PC frame, high profile, 120 tube strips per case
4ti-0775/B	8 Well PCR Tube Strips, clear PP wells, blue PC frame, high profile, 120 tube strips per case
4ti-0775/R	8 Well PCR Tube Strips, clear PP wells, red PC frame, high profile, 120 tube strips per case
4ti-0775/X	8 Well PCR Tube Strips, clear PP wells, black PC frame, high profile, 120 tube strips per case
4ti-0775/X/2D	8 Well PCR Tube Strips, clear PP wells, black PC frame, 2D coded, high profile, 120 tube strips per case
4ti-0775/W/2D	8 Well PCR Tube Strips, clear PP wells, white PC frame, 2D coded, high profile, 120 tube strips per case
4ti-0775/XW	8 Well PCR Tube Strips, white PP wells, black PC frame, high profile, 120 tube strips per case
4ti-0775/M	8 Well PCR Tube Strips, clear PP wells, assorted color PC frame, high profile, 120 tube strips per case
4ti-0775/C	8 Well PCR Tube Strips, clear PP wells, clear PC frame, high profile, 120 tube strips per case





4titude 8 Well PCR Tube Strip with PC Frame Adapter Lid

A 96 format frame that fits both 8 Well PCR Tube Strips and FrameStar Breakable plates, enabling our flexible PCR consumables to be handled in SBS format; ideal for use in automated or manual manufacturing processes.

Grippers on a robotic deck grip the 8 Well PCR Tube Strip adapter, allowing it to be handled as if it were a plate. The high profile Tube Strips have a push fit with the adapter, ensuring they are held within the plate and do not move around or fall out during handling, thereby ensuring samples are secure throughout the lifecycle.

- SBS footprint for use with robotics – loaded frames can be handled by robotics, increasing efficiency and accuracy during kit assembly
- Dedicated lid with interference fit – samples are securely protected avoiding damage to kits during transport
- Interference fit with high profile tubes – tubes are stable in the frame making them easier to handle and reducing the chance of errors when used with robotics
- Full Breakable plates or individual 8 Well PCR Tube Strips – suitable for use in both high and low throughput environments
- Ability to 2D code end tabs – sample can be tracked from kit manufacture all the way to processing at the end user site
- Locator pins on the deck ensure 8 Well PCR Tube Strips are always loaded in the correct orientation
- Reduces human error
- Adapters can be customized into different colors for custom/OEM opportunities.



Important note: the adapters cannot be placed directly into PCR cyclers due to the height of the skirt. These adapters are designed to provide an automation friendly strip/plate support during processing, after which, the strips/plates can be removed from the adapter and placed into a PCR block.

Specification Adapter

Feature	Information
Format	96 well
Length	127.76mm \pm 0.25mm
Width	85.48mm \pm 0.25mm
Height (without lid)	21.25mm
Height (with lid)	24.00mm
Color (adapter)	White
Color (lid)	Clear

Ordering Information

4ti-0370	8 Well PCR Tube Strip Adapter, 96 well skirted frame, with lid, white PC, cut corner H1, 18 adapters and lids per case
4ti-0371	8 Well PCR Tube Strip Adapter, 96 well skirted frame, white PC, cut corner H1, 18 adapters per case
4ti-0292	8 Well PCR Tube Strip Adapter Lid, without condensation rings, clear PS, low profile, no cut corner, 10 lids per case



AZENTA
LIFE SCIENCES

Non-Skirted PCR Plate Adapter

A 96 format frame that fits non-skirted PCR plates to allow for easy handling and use with automation.

- Cost and space saving design
- Clear polycarbonate
- Universal fit with non-skirted PCR plates
- Compatible with the Universal Microplate Lid
- Non-skirted PCR plates have the widest application range, fitting most cyclers
- However, they don't offer the same level of compatibility with robotic platforms and are not as multichannel pipetting friendly as fully- or semi-skirted plates are
- The Azenta Non-Skirted PCR Plate Adapter was developed to address all that, improving the overall ease of use of non-skirted plates
- Azenta adapters are manufactured in class 7 ISO certified cleanroom, and all lids are certified free of DNase, RNase, human DNA, bacterial and eukaryotic cells, dust and endotoxins/pyrogens



Key Features

- Accommodates both low and high profile non-skirted plates
- No locator pins, ensuring compatibility with all non-skirted plates and tubestrips
- Compatible with the Universal Microplate Lid, for a quick and easy sealing solution to protect samples from contamination and evaporation

Specification Non-Skirted PCR Plate Adapter

Feature	Information
Format	96 well
Length	127.76mm \pm 0.25mm
Width	85.48mm \pm 0.25mm
Height (without lid)	21.2mm
Height (with lid)	24.00mm
Color	Clear

Ordering Information

4ti-0373	Non-Skirted PCR Plate Adapter , 96 well skirted frame, clear PC, cut corner H1, 18 adapters per case
Combi Pack	
4ti-0372	Non-Skirted PCR Plate Adapter , 96 well skirted frame, with lid , clear PC, cut corner H1, 18 adapters and lids per case
Compatible lids	
4ti-0290	Universal Microplate Lid , without condensation rings, clear, low profile, no cut corner, 50 lids per case



4titude 96 Well PCR Plates with 8 Well Removeable Tube Strips

PCR tube strips in frames - excellent flexibility

The 96 Well PCR Plates with 8 Well Removeable Tube Strip system offers total flexibility in plate usage. It allows the user to insert or remove strips of 8 tubes from a 96 well plate frame.

- 8 Well Removeable Tube Strips can be used as a stand-alone product or inserted in one of two PCR plate frames for ease of handling
- 96 Well PCR Plate frames are available for use on Roche LightCycler® 480 and for universal use
- Pre-loaded frames are available
- Sealable using cap strips, adhesive seals or heat seals



8 Well Removeable Tube Strips can be used as a stand-alone product or with a frame for ease of handling.

8 Well PCR Tube Strips

Low profile, 0.1 ml polypropylene wells, working volume: <100 µl, total well capacity: 200 µl; can be used as a stand-alone product or inserted in any of the 96 Well PCR Plate frames available for ease of handling

- These low profile PCR tube strips are molded from virgin polypropylene under cleanroom conditions
- They are available in either clear polypropylene for standard PCR techniques, or in white polypropylene for use in fluorescent detection applications like qPCR, as they give the highest sensitivity and consistency as all the fluorescence signal is reflected back to the detector

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Tubes

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- Low profile 0.1 ml (100 µl) working volume, with a 0.2 ml (200 µl) total well capacity when used with sealing options
- Raised well rims prevent cross contamination and facilitate effective sealing for reduced evaporation

Specifications

Parameter	Value
Strip length	82.00 ± 0.10 mm
Strip width	8.20 mm
Strip height	15.60
Well depth	14.20 mm
Well diameter	5.50 ± 0.10 mm
Pitch (distance between A1 and A2)	9.00 mm



Use

- Can be used as individual strips, or in a 96 Well PCR Plate frame for easier handling
- Compatible with standard multichannel pipettes

Use

- 8 Well Removeable Tube Strips are available with clear wells
- 8 Well Removeable Tube Strips also available with white wells, giving the highest sensitivity and consistency for fluorescent detection during qPCR
- Combi packs available with strips of 8 flat optical caps
- 8 Well Removeable Tube Strips can be ordered separately or pre-loaded on 96 Well PCR Plate frames

Ordering Information

4ti-0753	8 Well PCR Tube Strips, white PP wells, low profile, 120 tube strips per case
4ti-0753/C	8 Well PCR Tube Strips, clear PP, low profile, 120 tube strips per case
Combi Pack	
4ti-0754/C	8 Well PCR Tube Strips, clear PP, low profile, plus cap strips, 4ti-0753 plus 4ti-0751 combi pack, 120 tube strips and cap strips per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude 96 Well Skirted PCR Plate Frame for Removable 8 Well Tube Strips

Rigid polycarbonate frame, for use with Removable 8 Well Tube Strips available for use on Roche LightCycler® 480 and for universal use

- These 96 Well Skirted PCR Plate frames are molded from rigid polycarbonate for use with our polypropylene Removable 8 Well Tube Strips
- The Removable 8 Well Tube Strips are inserted into the 96 Well Skirted PCR Plate frames to form a complete plate, part plate, or individual strip with an easy-to-handle frame
- As the Removable 8 Well Tube Strips are interchangeable with the 96 Well Skirted PCR Plate frames, the frames themselves are reusable; simply purchase more Removable 8 Well Tube Strips to use
- 96 Well Skirted PCR Plate frames can be purchased separately, or pre-loaded with Removable 8 Well Tube Strips

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Frames

- Rigid polycarbonate frames, which reduce thermal expansion and sample evaporation during PCR, leading to improved consistency in PCR results
- Alphanumeric grid reference to aid well and sample identification


Use

- Compatible with standard multichannel pipettes


Specifications

Parameter	Value	
	96 Well Skirted PCR Plate Frame	96 Well Semi-Skirted PCR Plate Frame, Roche Style
Plate length	127.76 ± 0.25 mm	127.70 ± 0.25 mm
Plate width	85.48 ± 0.25 mm	85.48 ± 0.25 mm
Plate height	16.10 ± 0.25 mm	15.60 ± 0.25 mm
Well depth	15.10 ± 0.10 mm	15.10 ± 0.10 mm
Well diameter	5.50 ± 0.10 mm	5.50 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm	9.00 mm

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



4ti-0757-F



4ti-0753/757 &
4ti-0753/C/757

Options

- Available as a 96 Well Skirted Frame for universal use
- Also available as a 96 Well Semi-Skirted Frame, Roche Style for use with the Roche LightCycler® 480
- Frames can be purchased on their own or pre-loaded with Removable 8 Well Tube Strips
- Removable 8 Well Tube Strips are available with clear wells
- Removable 8 Well Tube Strips also available with white wells, giving the highest sensitivity and consistency for fluorescent detection during qPCR
- Available barcoded upon request

Ordering Information

96 Well Skirted PCR Plate Frame

4ti-0757-F	96 Well Skirted PCR Plate for Removable 8 Well Strips, white PC frame, 50 plates per case
4ti-0753/757	96 Well Skirted PCR Plate with Removable 8 Well Strips, white PP wells, white PC frame, low profile, 50 plates per case
4ti-0753/C/757	96 Well Skirted PCR Plate with Removable 8 Well Strips, clear PP wells, white PC frame, low profile, 50 plates per case

96 Well Semi-Skirted PCR Plate Frame, Roche Style

4ti-0950W-F	96 Well Semi-Skirted PCR Plate for Removable 8 Well Strips, Roche style, white PC frame, low profile, cut corner H12, 10 frames per case
4ti-0753/950W	96 Well Skirted PCR Plate with Removable 8 Well Strips, Roche style, white PP wells, white PC frame, low profile, 50 plates per case



AZENTA
LIFE SCIENCES

4titude Standard PCR Plates, Strips & Tubes

Azenta offers a wide range of PCR consumables for low to medium throughput applications. Our standard one-piece PCR consumables are manufactured from virgin polypropylene in our Class 7 ISO certified clean-room production facility, and comply to the same stringent QC requirements as our FrameStar range.

The ultra-thin walled tubes of our standard PCR plates maximize heat transfer and the raised rims facilitate sealing. Our range of plates includes fully skirted, semi-skirted, and non-skirted plates, available in clear or white (for qPCR), with additional colors for non-skirted plates available.

All our PCR consumables are certified free from RNase, DNase, and human genomic DNA.

These plates fit most thermal cyclers; for a complete list please see the compatibility table.

- **Clean-room Injection Molding - Class 7 ISO Certification**

No contamination and 10 fold lower amount of air particles compared to most PCR plate manufacturers

- **Virgin, Medical Grade Polymers**

No leakage of substances which may have a detrimental effect on product purity

- **Certified RNase-, DNase-, DNA, and Pyrogen-free**

Inhibitor free consumables

- **Ultra-thin and consistent wall thickness**

Fast and precise thermal transfer

Azenta's standard PCR plates, strips and tubes are manufactured from virgin polypropylene under ISO certified cleanroom conditions.



4titude 384 Well Skirted PCR Plate

Polypropylene, cut corner A24, working volume: <30 µl, total well capacity: 55 µl; designed for use on standard 384 well thermal cyclers

- All of our PCR plates are molded from virgin polypropylene under ISO certified cleanroom conditions in our UK-based production facility, and as such comply with the same stringent requirements as our FrameStar range

Key Features

- Subject to a strict QC procedure which includes a visual check and leak testing of every tube
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time quantitative PCR (RT-qPCR) results
- <30 µl working volume, 55 µl total well capacity when used with sealing options
- Maximum thermal conductivity for efficient heat transfer and precise thermal cycling
- Raised well rims prevent cross contamination and facilitate effective sealing for reduced evaporation

Frame

- Alphanumeric grid reference to aid well and sample identification

Use

- Recommended for low volume PCR
- Compatible with blocks designed for 384 well PCR plates
- Compatible with standard multichannel pipettes

Options

- Super clear well option available to maximize sample visibility
- Also available as a frosted plate for increased qPCR signal intensities and improved detection sensitivity
- Available barcoded upon request
- Available as a FrameStar 2-component PCR plate with the same dimensions, but with the added benefits of a polycarbonate frame that eliminates plate distortion and warping during the PCR process, leading to less evaporation and sample loss



Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	10.30 ± 0.05 mm
Well depth	9.20 ± 0.10 mm
Well diameter	3.00 ± 0.10 mm
Distance to center of A1 from top edge	8.99 ± 0.25 mm
Distance to center of A1 from left edge	12.13 ± 0.25 mm
Pitch (distance between A1 and A2)	4.50 mm

Ordering Information

4ti-1384	384 Well Skirted PCR Plate, clear PP, cut corner A24, 50 plates per case
4ti-1387	384 Well Skirted PCR Plate, frosted PP, cut corner A24, 50 plates per case
4ti-1385	384 Well Skirted PCR Plate, white PP, cut corner A24, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude 384 Well Skirted PCR Plate, Roche Style

White polypropylene, cut corners A24 and P24; working volume: <30 µl, total well capacity: 55 µl; designed for use on the Roche LightCycler® 480 with 384 well block

- The dimensions of these plates are designed for optimum compatibility with the Roche LightCycler® 480, and are in a 384 well format for reaction volumes of up to 30 µl
- All of our PCR plates are molded from virgin polypropylene under ISO certified cleanroom conditions, and as such comply with the same stringent requirements as our FrameStar range

Key Features

- Subject to a strict QC procedure which includes a visual check and leak testing of every tube
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene wells for optimum PCR and real-time (RT-qPCR) results
- <30 µl working volume, 55 µl total well capacity when used with sealing options
- Raised well rims prevent cross contamination and facilitate effective sealing for reduced evaporation
- Maximum thermal conductivity for efficient heat transfer and precise thermal cycling

Frame


- Alphanumeric grid reference to aid well and sample identification

Use

- Recommended for low volume PCR
- Compatible with Roche LightCycler® 480 with 384 well block
- Compatible with standard multichannel pipettes

Options

- Available as a plate made from white polypropylene for optimum signal-to-noise ratio when using fluorescent based assays
- Available barcoded upon request
- Also available as a FrameStar 2-component PCR plate with the same dimensions, but with the added benefits of a polycarbonate frame that eliminates plate distortion and warping during the PCR process, leading to less evaporation and sample loss



Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	10.30 ± 0.05 mm
Well depth	9.20 ± 0.10 mm
Well diameter	3.00 ± 0.10 mm
Distance to center of A1 from top edge	8.99 ± 0.25 mm
Distance to center of A1 from left edge	12.13 ± 0.25 mm
Pitch (distance between A1 and A2)	4.50 mm

Ordering Information

4ti-1381	384 Well Skirted PCR Plate, Roche style, white PP, cut corner A24/P24, 50 plates per case
-----------------	--

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude 96 Well Skirted PCR Plate

Low profile, 0.1 ml wells, polypropylene, cut corner H1, working volume: <100 µl, total well capacity: 200 µl; universal 96 well skirted plate, designed for use on standard thermal cyclers

- All of our PCR plates are molded from virgin polypropylene under ISO certified cleanroom conditions in our UK-based production facility, and as such comply with the same stringent requirements as our FrameStar range.

Key Features

- Subject to a strict QC procedure which includes a visual check and leak testing of every tube
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- Low profile 0.1 ml (100 µl) working volume, with a 0.2 ml (200 µl) total well capacity when used with sealing options
- Maximum thermal conductivity for efficient heat transfer and precise thermal cycling
- Raised well rims prevent cross contamination and facilitate effective sealing for reduced evaporation

Frame


- Alphanumeric grid reference to aid well and sample identification

Use

- Compatible with blocks designed for skirted PCR plates
- Compatible with standard multichannel pipettes

Options

- Super clear wells maximize sample visibility
- Also available as a white plate, ideal for qPCR, giving optimal signal-to-noise ratio for fluorescence-based assays
- Available barcoded upon request
- Available as a FrameStar 2-component PCR plate with the same dimensions, but with the added benefits of a polycarbonate frame that eliminates plate distortion and warping during the PCR process, leading to less evaporation and sample loss



Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	15.50 ± 0.05 mm
Well depth	15.00 ± 0.10 mm
Well diameter	5.50 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0740	96 Well Skirted PCR Plate, clear PP, low profile, cut corner H1, 50 plates per case
4ti-0741	96 Well Skirted PCR Plate, white PP, low profile, cut corner H1, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude 96 Well Semi-Skirted PCR Plate

High profile, 0.2 ml wells, polypropylene, cut corner A12, working volume: <200 µl, total well capacity: 300 µl; universal semi-skirted plate, designed for use on standard thermal cyclers

- All of our PCR plates are molded from virgin polypropylene under ISO certified cleanroom conditions in our UK-based production facility, and as such comply with the same stringent requirements as our FrameStar range.

Key Features

- Subject to a strict QC procedure which includes a visual check and leak testing of every tube
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- Standard 0.2 ml (200 µl) working volume, with a 0.3 ml (300 µl) total well capacity when used with sealing options
- Maximum thermal conductivity for efficient heat transfer and precise thermal cycling
- Raised well rims prevent cross contamination and facilitate effective sealing for reduced evaporation

Frame



- Alphanumeric grid reference to aid well and sample identification

Use

- Compatible with universal standard block thermal cyclers and sequencers
- Compatible with standard multichannel pipettes

Options

- Super clear wells maximize sample visibility
- Also available as a white plate, ideal for qPCR, giving optimal signal-to-noise ratio for fluorescence-based assays
- Available barcoded upon request
- Available as a FrameStar 2-component PCR plate with the same dimensions, but with the added benefits of a polycarbonate frame that eliminates plate distortion and warping during the PCR process, leading to less evaporation and sample loss



Specifications

Parameter	Value
Plate length	124.42 ± 0.25 mm
Plate width	84.02 ± 0.25 mm
Plate height	20.70 ± 0.25 mm
Well depth	20.20 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	10.51 ± 0.10 mm
Distance to center of A1 from left edge	12.71 ± 0.10 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0760	96 Well Semi-Skirted PCR Plate, clear PP, high profile, cut corner A12, 50 plates per case
4ti-0761	96 Well Semi-Skirted PCR Plate, white PP, high profile, cut corner A12, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude 96 Well Semi-Skirted PCR Plate, Roche Style

Low profile, 0.1 ml wells, polypropylene, cut corner H12, working volume: <100 µl, total well capacity: 200 µl; designed for use on Roche LightCycler 480®

- This Roche Style plate is designed to achieve optimized assay conditions on the Roche LightCycler® 480
- This particular style of plate is in a low profile 96-well format, perfect for reaction volumes of 10-100 µl
- The wells are shorter than “standard” profile wells, which decrease the “dead space” between the heated lid of the thermal cycler and the sample in the well
- This eliminates condensation forming on the side wall of the tubes, preventing reduction in PCR volume and increasing the efficiency of the reaction
- This is especially recommended for reaction volumes below 20 µl

Key Features

- Subject to a strict QC procedure which includes a visual check and leak testing of every tube
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- Low profile 0.1 ml (100 µl) working volume, with a 0.2 ml (200 µl) total well capacity when used with sealing options
- Maximum thermal conductivity for efficient heat transfer and precise thermal cycling

Frame


- Alphanumeric grid reference to aid well and sample identification

Use

- Optimized for use with the Roche LightCycler® 480
- Recommended for low volume PCR
- Compatible with standard multichannel pipettes

Options

- Available with white wells give optimum signal-to-noise ratio when using fluorescent-based assays
- Combi packs available with qPCR Seal (0560)
- Available barcoded upon request
- Available as a FrameStar 2-component PCR plate with the same dimensions, but with the added benefits of a polycarbonate frame that eliminates plate distortion and warping during the PCR process, leading to less evaporation and sample loss



Specifications

Parameter	Value
Plate length	127.70 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	15.60 ± 0.10 mm
Well depth	15.10 ± 0.10 mm
Well diameter	5.50 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0955	96 Well Semi-Skirted PCR Plate, Roche style, white PP, low profile, cut corner H12, 50 plates per case
Combi packs	
4ti-0955/0560	96 Well Semi-Skirted PCR Plate, Roche style, plus qPCR Seal, 4TI-0955 plus 4TI-0560, combi pack, 50 plates and seals per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

4titude 96 Well Semi-Skirted PCR Plate with Upstand ABI® Style

High profile, 0.2 ml wells, polypropylene, cut corner A12, working volume: <200 µl, total well capacity: 300 µl; designed for use on ABI® thermal cyclers

- We recommend this semi-skirted plate for use with ABI® thermal cyclers and sequencers
- It can be used directly with ABI® instruments with no adapters and no re-calibration of the instruments necessary
- The only case where this is not true is with the ABI Fast Block thermal cyclers, in which case using our FrameStar FastPlate is recommended instead
- All of our PCR plates are molded from virgin polypropylene under ISO certified cleanroom conditions in our UK-based production facility, and as such comply with the same stringent requirements as our FrameStar range

Key Features

- Subject to a strict QC procedure which includes a visual check and leak testing of every tube
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells


- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time qPCR results
- Standard 0.2 ml (200 µl) working volume, with a 0.3 ml (300 µl) total well capacity when used with sealing options
- Maximum thermal conductivity for efficient heat transfer and precise thermal cycling
- Raised well rims prevent cross contamination and facilitate effective sealing for reduced evaporation

Frame

- Alphanumeric grid reference to aid well and sample identification

Use

- Designed for and compatible with ABI® thermal cyclers and Real Time PCR instruments
- Compatible with standard multichannel pipettes



Options

- Clear version has super clear wells for maximum sample visibility
- Also available with frosted wells for increased qPCR signal intensities and improved detection sensitivity
- Available barcoded upon request
- Available as a FrameStar 2-component PCR plate with the same dimensions, but with the added benefits of a polycarbonate frame that eliminates plate distortion and warping during the PCR process, leading to less evaporation and sample loss

Specifications

Parameter	Value
Plate length	126.00 ± 0.25 mm
Plate width	86.00 ± 0.25 mm
Plate height	22.10 ± 0.10 mm
Well depth	20.20 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	11.50 ± 0.10 mm
Distance to center of A1 from left edge	13.50 ± 0.10 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0735	96 Well Semi-Skirted PCR Plate, ABI style, clear PP, with upstand, high profile, cut corner H12, 50 plates per case
4ti-0736	96 Well Semi-Skirted PCR Plate, ABI style, frosted PP, with upstand, high profile, cut corner H12, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude 96 Well Non-Skirted PCR Plate

High profile, 0.2 ml wells, polypropylene, cut corner H12, working volume: <200 µl, total well capacity: 300 µl; universal non-skirted plate, designed for use on standard thermal cyclers

- All of our PCR plates are molded from virgin polypropylene under ISO certified cleanroom conditions in our UK-based production facility, and as such comply with the same stringent requirements as our FrameStar range

Key Features

- Subject to a strict QC procedure which includes a visual check and leak testing of every tube
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- 0.2 ml (200 µl) working volume, with a 0.3 ml (300 µl) total well capacity when used with sealing options
- Maximum thermal conductivity for efficient heat transfer and precise thermal cycling
- Raised well rims prevent cross contamination and facilitate effective sealing for reduced evaporation

Frame



- Alphanumeric grid reference to aid well and sample identification

Use

- Compatible with universal standard block thermal cyclers and sequencers
- Compatible with standard multichannel pipettes

Specifications

Parameter	Value
Plate length	120.00 ± 0.25 mm
Plate width	80.00 ± 0.25 mm
Plate height	20.20 ± 0.10 mm
Well depth	20.20 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	8.75 ± 0.25 mm
Distance to center of A1 from left edge	10.75 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm



Options

- Available with the following color options: clear, blue, red, green, yellow, purple, and white
- Super clear well version maximizes sample visibility
- White plate ideal for qPCR, giving optimal signal-to-noise ratio for fluorescence-based assays

Ordering Information

4ti-0750	96 Well Non-Skirted PCR Plate, clear PP, high profile, cut corner H12, 50 plates per case
4ti-0750-25	96 Well Non-Skirted PCR Plate, clear PP, high profile, cut corner H12, 25 plates per case
4ti-0750/P	96 Well Non-Skirted PCR Plate, purple PP, high profile, cut corner H12, 50 plates per case
4ti-0750/B	96 Well Non-Skirted PCR Plate, blue PP, high profile, cut corner H12, 50 plates per case
4ti-0750/G	96 Well Non-Skirted PCR Plate, green PP, high profile, cut corner H12, 50 plates per case
4ti-0750/R	96 Well Non-Skirted PCR Plate, red PP, high profile, cut corner H12, 50 plates per case
4ti-0750/W	96 Well Non-Skirted PCR Plate, white PP, high profile, cut corner H12, 50 plates per case
4ti-0750/Y	96 Well Non-Skirted PCR Plate, yellow PP, high profile, cut corner H12, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

4titude Non-Skirted PCR Plate Segments

High profile, 0.2 ml wells, polypropylene, cut corner H12, working volume: <200 µl, total well capacity: 300 µl; 96 Well Non-Skirted PCR Plate divided into 8 well, 16 well, 24 well, 32 well or 48 well segments

- These plates are produced by dividing our 96 Well Non-Skirted PCR Plate into smaller segments, a versatile solution for when a whole PCR plate may not be needed
- All of our PCR plates are molded from virgin polypropylene under ISO certified cleanroom conditions in our UK-based production facility, and as such comply with the same stringent requirements as our FrameStar range

Key Features



- Subject to a strict QC procedure which includes a visual check and leak testing of every tube
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- Standard 0.2 ml (200 µl) working volume, with a 0.3 ml (300 µl) total well capacity when used with sealing options
- Maximum thermal conductivity for efficient heat transfer and precise thermal cycling
- Raised well rims prevent cross contamination and facilitate effective sealing for reduced evaporation

Frame

- Alphanumeric grid reference to aid well and sample identification



Use

- Compatible with universal standard block thermal cyclers and sequencers
- Compatible with standard multichannel pipettes

Options

- Available with the following color options: clear, blue, red, green, yellow, and purple
- Super clear well version maximizes sample visibility

Specifications

Parameter	Value
Plate length	120.00 ± 0.25 mm
Plate width	80.00 ± 0.25 mm
Plate height	20.20 ± 0.10 mm
Well depth	20.20 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	8.75 ± 0.25 mm
Distance to center of A1 from left edge	10.75 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

4titude Non-Skirted PCR Plate Segments

Ordering Information

Clear 8 Well Non-Skirted PCR Plate Segment

4ti-0750/8	8 Well Non-Skirted PCR Plate Segment, clear PP, high profile, 600 segments per case
4ti-0750/16	16 Well Non-Skirted PCR Plate Segment, clear PP, high profile, 300 segments per case
4ti-0750/24	24 Well Non-Skirted PCR Plate Segment, clear PP, high profile, 200 segments per case
4ti-0750/32	32 Well Non-Skirted PCR Plate Segment, clear PP, high profile, 150 segments per case
4ti-0750/48	48 Well Non-Skirted PCR Plate Segment, clear PP, high profile, 100 segments per case

Blue 8 Well Non-Skirted PCR Plate Segment

4ti-0750/8/B	8 Well Non-Skirted PCR Plate Segment, blue PP, high profile, 600 segments per case
4ti-0750/16/B	16 Well Non-Skirted PCR Plate Segment, blue PP, high profile, 300 segments per case
4ti-0750/24/B	24 Well Non-Skirted PCR Plate Segment, blue PP, high profile, 200 segments per case
4ti-0750/32/B	32 Well Non-Skirted PCR Plate Segment, blue PP, high profile, 150 segments per case
4ti-0750/48/B	48 Well Non-Skirted PCR Plate Segment, blue PP, high profile, 100 segments per case

Purple 8 Well Non-Skirted PCR Plate Segment

4ti-0750/8/P	8 Well Non-Skirted PCR Plate Segment, purple PP, high profile, 600 segments per case
4ti-0750/16/P	16 Well Non-Skirted PCR Plate Segment, purple PP, high profile, 300 segments per case
4ti-0750/24/P	24 Well Non-Skirted PCR Plate Segment, purple PP, high profile, 200 segments per case
4ti-0750/32/P	32 Well Non-Skirted PCR Plate Segment, purple PP, high profile, 150 segments per case
4ti-0750/48/P	48 Well Non-Skirted PCR Plate Segment, purple PP, high profile, 100 segments per case

Yellow 8 Well Non-Skirted PCR Plate Segment

4ti-0750/8/Y	8 Well Non-Skirted PCR Plate Segment, yellow PP, high profile, 600 segments per case
4ti-0750/16/Y	16 Well Non-Skirted PCR Plate Segment, yellow PP, high profile, 300 segments per case
4ti-0750/24/Y	24 Well Non-Skirted PCR Plate Segment, yellow PP, high profile, 200 segments per case
4ti-0750/32/Y	32 Well Non-Skirted PCR Plate Segment, yellow PP, high profile, 150 segments per case
4ti-0750/48/Y	48 Well Non-Skirted PCR Plate Segment, yellow PP, high profile, 100 segments per case

Red 8 Well Non-Skirted PCR Plate Segment

4ti-0750/8/R	8 Well Non-Skirted PCR Plate Segment, red PP, high profile, 600 segments per case
4ti-0750/16/R	16 Well Non-Skirted PCR Plate Segment, red PP, high profile, 300 segments per case
4ti-0750/24/R	24 Well Non-Skirted PCR Plate Segment, red PP, high profile, 200 segments per case
4ti-0750/32/R	32 Well Non-Skirted PCR Plate Segment, red PP, high profile, 150 segments per case
4ti-0750/48/R	48 Well Non-Skirted PCR Plate Segment, red PP, high profile, 100 segments per case

Green 8 Well Non-Skirted PCR Plate Segment

4ti-0750/8/G	8 Well Non-Skirted PCR Plate Segment, green PP, high profile, 600 segments per case
4ti-0750/16/G	16 Well Non-Skirted PCR Plate Segment, green PP, high profile, 300 segments per case
4ti-0750/24/G	24 Well Non-Skirted PCR Plate Segment, green PP, high profile, 200 segments per case
4ti-0750/32/G	32 Well Non-Skirted PCR Plate Segment, green PP, high profile, 150 segments per case
4ti-0750/48/G	48 Well Non-Skirted PCR Plate Segment, green PP, high profile, 100 segments per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

4titude 8 Well PCR Tube Strips

High profile, 0.2 ml wells, clear polypropylene, working volume: <200 µl, total well capacity: 300 µl; suitable for all standard 0.2 ml block thermal cyclers

- These PCR tubes are molded from virgin polypropylene in our UK-based Class 7 ISO certified cleanroom production facility, and comply to the same stringent requirements as our FrameStar range
- Recommended for low to medium throughput applications

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Tubes

- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results
- Standard 0.2 ml (200 µl) working volume, with a 0.3 ml (300 µl) total well capacity when used with sealing options
- Individually numbered tubes

Use

- Suitable for all standard 0.2 ml block thermal cyclers
- Can be cut into sections

Options

- Available with strips of domed (4ti-0780) or flat optical (4ti-0784) sealing caps*

* See website for the latest Important Product Information



Ordering Information

4ti-0781	8 Well PCR Tube Strips, clear PP, high profile, 125 strips per case
Code Combi Pack	
4ti-0780	8 Well PCR Tube Strips, plus Strips of Domed Caps, 4ti-0781 plus 4ti-0782, combi pack, 125 tube strips and cap strips per case
4ti-0784	8 Well PCR Tube Strips, plus Strips of Flat Caps, 4ti-0781 plus 4ti-0783, combi pack, 125 tube strips and cap strips per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude 8 Well PCR Tube Strips, With Attached Caps

Low profile 0.1 ml (working volume: <100 µl, total well capacity: 200 µl), clear polypropylene, with attached flat optically clear caps; and standard 0.2 ml (working volume: <200 µl, total well capacity: 300 µl), clear polypropylene, with attached domed or flat optically clear caps

- These PCR tubes with attached caps are molded from virgin polypropylene, comply to the same stringent requirements as our FrameStar range and are free from RNase, DNase, and human genomic DNA
- Recommended for low to medium throughput applications

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Tubes

- The tubes on each strip are joined together either by three links (flat cap variety - 4ti-0792 and 4ti-0793) or by one strong link (domed cap variety - 4ti-0794) to make the strip more rigid and help reduce any chance of spillage
- Each cap is separately joined to a tube, making it impossible to either cross contaminate another tube with the wrong cap, or to lose a cap altogether

Use

- The tube strips can be easily separated by cutting the links, to make smaller sections or individual tubes
- The flat optically clear caps enable light signals, such as fluorescence, to pass through without affecting the signal, and are suitable for imaging techniques including RT-qPCR

Options


- The strips are available with either tethered flat optically clear caps or domed caps



Ordering Information

4ti-0792	8 Well PCR Tube Strips, with Attached Flat Caps, 0.2ml wells, clear PP, high profile, 120 tube strips per case
4ti-0793	8 Well PCR Tube Strips, with Attached Flat Caps, 0.1ml wells, clear PP, low profile, 120 tube strips per case
4ti-0794	8 Well PCR Tube Strips, with Attached Domed Caps, 0.2ml wells, clear PP, high profile, 120 tube strips per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

4titude 4 Well PCR Tube Strips, Rotor-Gene® Style, With Caps

0.1 ml wells, clear polypropylene, with strips of 4 caps, working volume: <100 µl, total well capacity: 200 µl; designed for Qiagen/Corbett Rotor-Gene® instruments

- Azenta offers a range of PCR consumables for low to medium throughput applications
- These PCR tubes and caps are molded from virgin polypropylene to prevent any PCR background signal
- They are especially suited for use with Qiagen/Corbett Rotor-Gene® instruments

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Tubes

- Optically clear tubes ideally suited for qPCR
- Low profile 0.1 ml (100 µl) working volume, 0.2 ml (200 µl) total well capacity
- Ultra-smooth, uniform, thin-walled polypropylene tubes for optimum PCR and real-time (RT-qPCR) results

Caps

- Frosted cap extensions allow for easy handling and labeling

Use

- Designed for Qiagen/Corbett Rotor-Gene® instruments
- Tube and cap strips can be separated for individual use
- Pack of 250 strips of tubes and caps sufficient for 1,000 reactions



Ordering Information

4ti-0796

4 Well PCR Tube Strips, Rotor-Gene style, plus strips of caps, clear PP, 250 tube strips and cap strips per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

Cap² 0.2ml Dual-Cap Sample Collection PCR Tube

PCR Tube Innovation for Workflow Versatility

A unique PCR tube featuring a patent pending dual-cap mechanism, 2D datamatrix code as well as a human readable identification, allowing for wide-ranging workflow versatility. Innovative features allow for utilization as a consumable suitable for sample collection/transport that can be directly processed for analysis, ideal for genetic testing workflows.

Key Features

Application and Workflow Flexibility

- PCR tube with patent pending cap-in-cap design, sample identifiers and automation friendly capabilities offer a flexible, scalable solution
- Hinged lid allows for flexibility for use in manual workflows, alongside mitigating against contamination
- Orange Screw cap enables automated opening and closing using an Azenta fully automated or semi-automated handheld capper
- Tube-specific PCR Rack in SBS footprint for PCR cycling as well as liquid handler use

Sample Identification

- Tab on the side of the tube with 2D datamatrix code and human readable code
- Full rack decoding possible with the Verification Camera of Azenta IntelliXcap Acoustic (part 46-8014)

High Quality Tubes

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen
- Leak tested
- 0.2 ml (200 µl) total well capacity
- Polypropylene, ultra-thin walled, low DNA binding tubes make for optimum PCR results and mitigation against
- DNA loss from samples



Ordering Information

4ti-0791	Cap² 0.2ml Dual-Cap, Sample Collection, PCR Tube. Clear polypropylene, attached flip cap with orange screw cap, 2D code and human readable number on tab, 500 tubes per case
4ti-0374	Cap² 0.2ml Dual-Cap Sample Collection PCR Tube Rack. 48 well skirted frame, SBS format, white PC, cut corner H1/H12, suitable for 0.2ml Dual-Cap PCR Tubes (4TI-0791), 10 racks per case

4titude Individual PCR Tubes, With Attached Caps, With Or Without 2D Barcode

High profile, 0.2 ml wells, polypropylene, working volume <200 µl, total well capacity: 300 µl; with either flat or domed attached caps; also available with unique 2D coded flat caps

- These PCR tubes are molded from virgin polypropylene in our UK-based Class 7 ISO certified cleanroom production facility, and comply to the same stringent requirements as our FrameStar range
- Recommended for low to medium throughput applications
- Our individual tubes with flat caps are also available with a unique 2D code applied to the top of each tube
- This allows for clear labeling of tubes, despite the limited space available
- Printed 2D codes allow for easy sample tracking and are more reliable than the use of adhesive stickers
- Our compact hand held barcode scanner, 4ti-4060, delivers the speed typical of laser scanners on any barcode, including both 1D and 2D codes



Key Features

- Leak tested
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Tubes

- 0.2 ml (200 µl) working volume, 0.3 ml (300 µl) total well capacity
- Snap-shut cap design

Uses

- Suitable for all standard 0.2 ml block thermal cyclers

Options

- Flat and domed cap designs
- Our individual tubes with flat caps (4ti-0790) are also available with unique 2D code applied to the top of each tube (use product code 4ti-0790/2D when ordering); the barcode is suitable for low temperature storage and high temperature thermal cycling
- Our hand held barcode scanner, 4ti-4060 together with the user-friendly software allows for reliable and convenient sample management

Ordering Information

4ti-0790	Individual PCR Tubes with Attached Flat Caps, 0.2ml, clear PP wells, 1000 tubes per case
4ti-0790/2D	Individual PCR Tubes with Attached Flat Caps, 0.2ml, clear PP, 2D coded, high profile, 960 tubes per case
4ti-0795	Individual PCR Tubes with Attached Domed Caps, 0.2ml, clear PP, high profile, 1000 tubes per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.

4titude Non-Skirted PCR Plates Breakable Horizontally or Vertically

Non-Skirted PCR Plates Breakable Horizontally or Vertically allow for the most flexible, efficient and cost-effective use of a PCR plate. Avoid the costly use of half-empty plates or the fiddly separation of plates with scissors. Cutting plates can damage wells and sealing rings, risking evaporation and sample contamination.

Based on our standard non-skirted PCR plate (4ti-0750), Non-Skirted PCR Plates Breakable Horizontally or Vertically can be quickly and easily divided along the perforations between the rows. The correct number of wells can be separated off for each experiment, saving time and costs.

The Non-Skirted PCR Plates Breakable Horizontally or Vertically are available perforated either in the vertical direction, tearing into 8 well strips, or in the horizontal direction, tearing into 12 well strips. Both versions maintain all the benefits of our standard non-skirted PCR plate, but with increased flexibility.

- **Allows for the most flexible and efficient use of a PCR plate** — no need to run half-empty plates, so reducing costs
- **Plate is perforated to enable accurate tearing into either 8 well or 12 well strips** — no tricky cutting of plates with scissors risking perforating wells, damaging sealing rings and contamination
- **Black grid reference on all strips** — No sample identification errors
- **Non-skirted plates** — Universal cyclers and sequencer compatibility
- **8 well version is easily divided into 24 and 48 well plates** to fit a 24 or 48 well thermal cycler block
- **12 well version perfectly suited for gradient cyclers**
- **White version available for superior qPCR performance**



96/12 Non-Skirted PCR Plates Breakable Horizontally or Vertically allow you to make full use of your gradient PCR instruments. The temperature gradient is typically created along the horizontal direction of the block, thus 12 well strips or sections are ideal.



How trustworthy are your scissors?

Scissors are widely used by everyone in the lab and are typically highly contaminated with substances including bacteria and DNA. Cutting PCR plates with scissors should be avoided as it may lead to contamination of the wells.



AZENTA
LIFE SCIENCES

4titude 96 Well Non-Skirted PCR Plate

Breakable Horizontally or Vertically

96 well non-skirted PCR plate, vertically or horizontally perforated, high profile, 0.2 ml wells, polypropylene, cut corner H12; working volume: <200 µl, total well capacity: 300 µl; tears easily into strips or part plates; universal cycler and sequencer compatibility



- Non-Skirted PCR Plates Breakable Horizontally or Vertically allow for the most flexible, efficient and cost-effective use of a PCR plate
- Avoid the costly use of half-empty plates or the fiddly separation of plates with scissors
- Scissors are widely used by everyone in the lab for cutting diverse materials, so are typically highly contaminated with substances including bacteria and DNA
- Cutting plates with scissors should be avoided as it can perforate wells and damage sealing rings, risking evaporation and sample contamination
- Based on our 96 Well Non-Skirted PCR Plate, plates can be quickly and easily divided along the perforations between the rows
- The correct number of wells can be separated off for each experiment, saving time and costs
- Available perforated either vertically, tearing into 8-well strips (96/8), or horizontally, tearing into 12-well strips (96/12)
- The new 12-strip 96/12 plates allow you to make full use of your gradient PCR instruments
- The temperature gradient is typically created along the horizontal direction of the block, thus 12-well strips or sections are ideal
- All of our PCR plates are molded from virgin polypropylene under ISO-certified cleanroom conditions in our UK-based production facility, and as such comply with the same stringent requirements as our FrameStar range
- This plate is perforated under the same cleanroom conditions post-manufacture so it can be separated into part plates or individual strips of tubes
- The Breakable Horizontally or Vertically Plate maintains all benefits of the 96 Well Non-Skirted PCR Plate, but adds extra flexibility

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Frame

- Non-skirted plates: universal cycler and sequencer compatibility
- Alphanumeric grid reference on all strips: no sample identification errors
- Plate is perforated to enable accurate tearing into 8-well or 12-well strips: no tricky cutting of plates with scissors, risking perforated wells, damaging sealing rings, and contamination



Use

- Allows for the most flexible and efficient use of a PCR plate: no need to run half-empty plates, so reducing costs
- Snaps into strips for lower throughput: cost effective

Options

- The 96/8 8-well version is easily divided into 24 and 48-well plates to fit a 24 or 48-well thermal cycler block
- The 96/12 12-well version perfectly suited for gradient cyclers
- White version available for superior qPCR performance
- Sealable with Azenta PCR cap strips, adhesive seals and heat seals

Specifications

Parameter	Value
Plate length	120.00 ± 0.25 mm
Plate width	80.00 ± 0.25 mm
Plate height	20.20 ± 0.10 mm
Well depth	20.20 ± 0.10 mm
Well diameter	5.46 ± 0.10 mm
Distance to center of A1 from top edge	8.75 ± 0.25 mm
Distance to center of A1 from left edge	10.75 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0750/TA	96 Well Non-Skirted PCR Plate, clear PP, breakable vertically, high profile, cut corner H12, 50 plates per case
4ti-0750/W/TA	96 Well Non-Skirted PCR Plate, white PP, breakable vertically, high profile, cut corner H12, 50 plates per case
4ti-0750/TA/12	96 Well Non-Skirted PCR Plate, clear PP, breakable horizontally, high profile, cut corner H12, 50 plates per case

For recommended plate options depending on manufacturer, block type, and instrument please refer to the Azenta Plate Instrument Compatibility Table Pages 162-165.



AZENTA
LIFE SCIENCES

Thermal Cycler and PCR Product Compatibility Table



Instrument Compatibility Table

FrameStar Polypropylene & Polycarbonate PCR Plates

Skirt			skirted			semi-skirted				
Number of wells			384	384	96	96	96	96	96	96
High/Low Profile			/	/	L	L	L	H	H	H
Product code/Short description			4ti-0384 FrameStar 384	4ti-0380 FrameStar 384 Roche	4ti-0960 FrameStar 96	4ti-0950 & 4ti-0954 FrameStar 96 Roche	4ti-0910 FrameStar 96 FastPlate	4ti-0730 FrameStar 96 ABI®	4ti-0770 FrameStar 96 ABI®	4ti-0900 FrameStar 96
ABI® / LIFE TECHNOLOGIES / THERMO FISHER SCIENTIFIC										
Thermal Cyclers	96 well standard block	Veriti, Proflex, Simplicamp						✓	△	
		GeneAmp® 2700 / 2720 / 9600 / 9700					✓	△		
	96 well FAST block	GeneAmp® 9800 FAST, Veriti FAST					△			
	384 well block	GeneAmp® 9700, Veriti, Proflex, Multiblock system	△							
qPCR Cyclers	96 well standard block	7000, 7300, 7500, 7700, 7900 HT						✓	△	
		QuantStudio™ 3 / 5 / 6 / 7 / 12K, ViiA7™					✓	△		
	96 well FAST block	StepOne								
		StepOne Plus™					△			
		7500 FAST, 7900 HT FAST					△			
		QuantStudio™ 3 / 5 / 6 / 7 / 12K, ViiA7™					△			
384 well block	QuantStudio™ 5 / 6 / 7 / 12K, ViiA7™ 7900 HT FAST	△								
Sequencers	96 well block	3100, 3130XL, 3500, 3500XL, 3730, 3730XL						△	✓	
	384 well block	3100, 3130XL, 3500, 3500XL, 3730, 3730XL	△							
AGILENT / STRATAGENE										
Thermal Cyclers	96 well block	Surecycler 8800								
	384 well block	SureCycler 8800	△							
	96 well block	Robocycler Gradient			△					✓
qPCR Cyclers	96 well block	AriaMx			△					
		Mx3000P™ Mx3005P™							△	✓
		Mx4000™						✓	△	
ANALYTIK JENA / BIOMETRA										
Thermal Cyclers	Strips only	TRIO, Tpersonal, T3 Thermocycler								
	96 well block	Flexcycler2, T1 Thermocycler, Tgradient, Tone, Tadvanced, TProfessional (Standard/Basic) Gradient/XL			△				✓	✓
		Trobot 96, SpeedCycler2 (SP, SPR)			△				✓	✓
	384 well block	Flexcycler2, T1 Thermocycler, Tadvanced, TProfessional, Trobot 96	△							
qPCR Cyclers	96 well block	qTOWER³ / G / touch, Toptical							△	✓
	384 well block	qTOWER³ 84 / 84G	△							
BIOER TECHNOLOGIES										
Thermal Cyclers	Strips only	GeneQ								
	96 well block	Gene Touch 96							△	✓
	384 well block	Gene Touch 384	△							
BIO-RAD										
Thermal Cyclers	Strips only	Genecycler								
	96 well block	C1000 Touch, S1000			△				✓	✓
		iCycler™ MyCycler™, T100							△	✓
384 well block	C1000 Touch, S1000	△								
qPCR Cyclers	96 well block	CFX96 Touch, CFX96 Touch Deep Well, CFX connect			△					
		MyiQ™ iCycler™ IQ / IQ 4 / IQ 5							△	✓
	384 well block	CFX384 Touch	△							

¹ Short product code shown only (without details on frame and well color), please refer to the corresponding product page for ordering details of all variations available.

² For compatibility information of the PCR Plate for Removable 8 Well Tube Strips please refer to the table entry for the respective PCR Plate for Removable 8 Well Tube Strips frame.

△ Recommended option ✓ Compatible • Should be compatible, please check with your specific instrument/block ✱ Compatible with PCR Plates Breakable Horizontally or Vertically only

Azenta US, Inc. recognizes that designated trademarks and brands are the property of their respective owners.



AZENTA
LIFE SCIENCES

FrameStar PCR Plates, Individually Accessible Plate, PCR Tube Strips								Standard PCR Plates and Strips												
non skirted		Breakable Vertically		Breakable Vertically & Horizontally		skirted	strips	skirted			semi-skirted			non-skirted	strips					strips
96	96	96	96	96	96	96	8	384	384	96	96	96	96	96	8	8	8	8	8	4
H	L	H	L	H	L	L	H	/	/	L	H	L	H	H	H	L	H	L	H	/
4ti-0710 FrameStar 96	4ti-0720 FrameStar 96	4ti-1000 FrameStar Breakable Vertically	4ti-1200 FrameStar Breakable Vertically	4ti-1300 FrameStar Breakable Vertically & Horizontally	4ti-1400 FrameStar Breakable Vertically & Horizontally	4ti-0960/RA Individual Access 96	4ti-0785 & 4ti-0786 PCR Tube Strips 8 Well	4ti-1384 384 Well	4ti-1381 384 Well Roche	4ti-0740 96 Well	4ti-0760 96 Well	4ti-0955 96 Well Roche	4ti-0735 96 Well ABI®	4ti-0750 - 96 Well 4ti-0750/7A- dividable	4ti-0781 8 Well	4ti-0753 PCR Tube Strip	4ti-0792 8 Well, Attached Flat Caps	4ti-0793 8 Well, Attached Flat Caps	4ti-0794 8 Well, Attached Domed C.	4ti-0796 4 Well Rotor-Gene®
		✓		✓			✓				✓		✓	✓	✓		✓		✓	
✓		✓		✓			✓				✓		✓	✓	✓		✓		✓	
			✓		✓								✓	✓	✓		✓		✓	
								✓												
✓		✓		✓			✓				✓		✓	✓	✓		✓			
		✓		✓			✓				✓		✓	✓	✓		✓			
			✓		✓								✓	✓						
					✓											△				
					✓											✓			✓	
					✓														✓	
					✓														✓	
					✓														✓	
✓		✓		✓							✓		✓							
								✓												
△														✓	✓					
								✓												
✓			✓		✓	✓				✓				✓		✓		✓		
	✓		✓		✓	✓				✓						✓		✓		
✓		✓		✓			✓				✓			✓	✓	✓		✓		
✓		✓		✓			✓				✓			✓	✓	✓		✓		
		△		✓			✓							*	✓		✓		✓	
✓		✓		✓			✓				✓			✓	✓		✓		✓	
								✓												
		△		✓			✓							*	✓		✓		✓	
✓	✓	✓	✓	✓	✓	✓	✓			✓	✓			✓	✓	✓	✓	✓	✓	
✓		✓		✓			✓				✓			✓	✓		✓		✓	
								✓												
	✓		✓		✓	✓				✓						✓		✓		
✓		✓		✓			✓				✓					✓		✓		
								✓												
								✓												

The individual instrument-plate compatibility information is for guidance only. Samples of all plates are available before purchase to ensure compatibility. Please check installed heat block, refer to instruments manual for details or contact our technical support.

Azenta recognises that designated trademarks and brands of the Instrument Compatibility Table are the property of their respective owners.



Instrument Compatibility Table

FrameStar Polypropylene & Polycarbonate PCR Plates

Skirt			skirted			semi-skirted				
Number of wells			384	384	96	96	96	96	96	96
High/Low Profile			/	/	L	L	L	H	H	H
Product code/Short description			4ti-0384 FrameStar 384	4ti-0380 FrameStar 384 Roche	4ti-0960 FrameStar 96	4ti-0950 & 4ti-0954 FrameStar 96 Roche	4ti-0910 FrameStar 96 FastPlate	4ti-0730 FrameStar 96 ABI®	4ti-0770 FrameStar 96 ABI®	4ti-0900 FrameStar 96
BIO-RAD MJ RESEARCH										
Thermal Cyclers	Strips only	Mini Gradient								
	96 well block	Personal								
		PTC100™ / 200™ / 220™ / 221™ / 225™ / 240™							✓	✓
	384 well block	PTC200™ / 220™ / 221™ / 225™ / 240™								
qPCR Cyclers	Strips only	MiniOpticon								
	96 well block	Chromo4™							✓	✓
		Opticon2™								
CORBETT RESEARCH										
Thermal Cyclers	96 well block	(Qiagen) Palm Cycler								✓
	384 well block	(Qiagen) Palm Cycler 384								
qPCR Cycl.	Strips only	Rotor-Gene series								
EPPENDORF										
Thermal Cyclers	96 well block	MasterCycler® ep / ep gradient / Pro / Pro S / nexus / nexus gradient / nexus SX1 / nexus GSX1							✓	✓
		MasterCycler® nexus X2 / GX2 / GX2e / X2e								
	384 well block	MasterCycler® ep 384 / Pro 384								
qPCR Cycl.	96 well block	Mastercycler™ ep realplex							✓	✓
GE HEALTHCARE / AMERSHAM										
Sequencers	96 well block	MegaBACE™ 500, MegaBACE™ 1000 mark 2								
	384 well block	MegaBACE™ 4000								
PEQLAB / VWR										
Thermal Cyclers	Strips only	peqSTAR XS, peqSTAR 2X								
	96 well block	peqSTAR 96X							✓	✓
	384 well block	peqSTAR 384X								
ROCHE										
qPCR Cyclers	96 well block	LC96, LC480								
	384 well block	LC480								
	Strips only	Nano								
SENSOQUEST										
Thermal Cyclers	96 well block	Labcycler			✓		✓	✓	✓	✓
	384 well block	Labcycler								
TAKARA										
Thermal Cycl.	96 well block	Dice touch, Gradient								✓
TECHNE										
Thermal Cyclers	Strips only	3Prime, 3PrimeG, 3PrimeX								
	96 well block	Prime, PrimeG, Prime Elite, Prime Elite Satellite								✓
		PCRmax Alpha cycler 1 / 2 / 4								✓
		TC412, TC512, Genius, Genius Quad, Touchgene, Touchgene Gradient, Flexigene							✓	✓
	384 well block	Prime, PrimeG, Prime Elite, Prime Elite Satellite								
		PCRmax Alpha cycler 1 / 2 / 4								
	TC412, TC512, Genius, Genius Quad, Flexigene									
qPCR Cycl.	96 well block	Quantica								

¹ Short product code shown only (without details on frame and well color), please refer to the corresponding product page for ordering details of all variations available.

² For compatibility information of the PCR Plate for Removable 8 Well Tube Strips please refer to the table entry for the respective PCR Plate for Removable 8 Well Tube Strips frame.

△ Recommended option ✓ Compatible • Should be compatible, please check with your specific instrument/block ✱ Compatible with PCR Plates Breakable Horizontally or Vertically only

Azenta US, Inc. recognizes that designated trademarks and brands are the property of their respective owners.



AZENTA
LIFE SCIENCES

FrameStar PCR Plates, Individually Accessible Plate, PCR Tube Strips								Standard PCR Plates and Strips												
non skirted		Breakable Vertically		Breakable Vertically & Horizontally		skirted	strips	skirted			semi-skirted			non-skirted	strips					strips
96	96	96	96	96	96	96	8	384	384	96	96	96	96	96	8	8	8	8	8	4
H	L	H	L	H	L	L	H	/	/	L	H	L	H	H	H	L	H	L	H	/
41i-0710 FrameStar 96	41i-0720 FrameStar 96	41i-1000 FrameStar Breakable Vertically	41i-1200 FrameStar Breakable Vertically	41i-1300 FrameStar Breakable Vertically & Horizontally	41i-1400 FrameStar Breakable Vertically & Horizontally	41i-0960/RA Individual Access 96	41i-0785 & 41i-0786 PCR Tube Strips 8 Well	41i-1384 384 Well	41i-1381 384 Well Roche	41i-0740 96 Well	41i-0760 96 Well	41i-0955 96 Well Roche	41i-0735 96 Well ABI®	41i-0750 - 96 Well 41i-0750/TA - dividable	41i-0781 8 Well	41i-0753 PCR Tube Strip	41i-0792 8 Well, Attached Flat Caps	41i-0793 8 Well, Attached Flat Caps	41i-0794 8 Well, Attached Domed C.	41i-0796 4 Well Rotor-Gene®
		△		✓			✓							*	✓		✓		✓	
✓																				
✓	✓	✓	✓	✓	✓	✓	✓			✓	✓			✓	✓	✓	✓	✓	✓	
								✓												
			✓		✓											△		✓		
✓	✓	✓		✓		✓	✓			✓	✓			✓	✓	✓	✓	✓		
	✓		✓		✓	✓				✓						✓		✓		
✓		✓		✓			✓				✓			✓	✓		✓		✓	
								✓												△
✓	✓	✓	✓	✓	✓	✓	✓			✓	✓			✓	✓	✓	✓	✓	✓	
		△	✓	✓	✓		✓							✓	✓	✓	✓	✓	✓	
								✓												
✓	✓	✓	✓	✓	✓	✓	✓			✓	✓			✓	✓	✓	✓	✓		
	✓		✓		✓					✓						✓				
								✓												
									✓							△				
✓	✓	✓	✓	✓	✓	✓	✓			✓	✓		✓	✓	✓	✓	✓	✓	✓	
								✓												
		✓		✓			✓				✓			✓	✓				✓	
✓		✓		✓			✓				✓			✓	✓		✓		✓	
		△		✓			✓							*	✓		✓		✓	
✓		✓		✓			✓				✓			✓	✓		✓		✓	
✓		✓		✓			✓				✓			✓	✓		✓		✓	
✓	✓	✓	✓	✓	✓	✓	✓			✓	✓			✓	✓	✓	✓	✓	✓	
								✓												
								✓												
								✓												
	✓		✓		✓	✓				✓						✓		✓		

The individual instrument-plate compatibility information is for guidance only. Samples of all plates are available before purchase to ensure compatibility. Please check installed heat block, refer to instruments manual for details or contact our technical support.

Azenta recognises that designated trademarks and brands of the Instrument Compatibility Table are the property of their respective owners.



AZENTA
LIFE SCIENCES

Storage Plates



4titude 384 Square Deep Well Storage Microplate

190 µl square wells, V shaped bases, clear polypropylene

- This plate is designed for high density sample collection and storage, for a wide array of applications within cell biology, molecular biology and drug discovery
- The V shaped wells allow for complete sample retention with pipettes, and are ultra-flat for a completely uniform seal, either with adhesive or heat seals, or (4ti-0139)
- Our plates are manufactured in cleanroom facilities which are certified free of RNase, DNase, human genomic DNA and endotoxins
- We use only the highest medical grade virgin polypropylene with high resistance against chemicals such as DMSO, phenol, and chloroform

Key Features

- Made from polypropylene, a material with very low biomolecular binding characteristics, a high temperature tolerance and resistance to many standard laboratory chemicals, including DMSO
- Autoclavable
- Free from DNase, RNase, human genomic DNA, and endotoxin

Wells


- Square wells, to make the best use of space, and to improve mixing
- Conical V-bottoms

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint

Use

- Clear plate, perfect for storage applications
- Temperature range for use: -80°C to 120°C
- Best for precipitation, centrifugation and small volume recovery due to the conical V bottom, as the liquid gathers at the lowest point of the V for easy pipetting
- Stackable
- Suitable for adhesive and heat sealing



Options

- Non-gamma treated as standard, gamma treatment available upon request; contact us for details and pricing
- Available barcoded upon request
- A clear silicone sealing mat (4ti-0139) is also available to fit this product, with 384 sections to correspond with the wells of this plate
- Ultra-low DNA binding option available (4ti-LB0147) for sensitive applications with ultra-low DNA input and for maximum DNA recovery after low temperature storage and high temperature incubation

Ordering Information

4ti-0147	384 Square Deep Well Storage Microplate, 190ul square wells, V-shaped bottom, clear PP, 100 plates per case
4ti-LB0147	384 Square Deep Well Storage Microplate, low binding, 190µl square wells, V-shaped bottom, clear PP, 100 plates per case
Related Products	
4ti-0139	384 Square Well Sealing Cap Mat, clear silicone, for use with square 384 well microplates and deep well storage microplates, 50 mats per case



4titude 96 Square Deep Well Storage Microplate

1.2 ml and 2.2 ml square wells, U and V shaped bottom, clear polypropylene

- These storage plates are designed for high density sample collection and storage for a wide array of applications within cell biology, molecular biology and drug discovery
- They come in three formats: 1.2 ml volume U shaped well plates (4ti-0126), 2.2 ml volume V shaped (4ti-0132), and 2.2 ml volume U shaped (4ti-0136) plates
- All of these plates are ultra-flat for a completely uniform seal, either with adhesive or heat seals, or with our accompanying silicone sealing mat (4ti-0137)
- These plates are manufactured in cleanroom facilities which are certified free of RNase, DNase, human genomic DNA and endotoxins
- We use only the highest medical grade virgin polypropylene with high resistance against chemicals such as DMSO, phenol, and chloroform

Key Features

- Made from polypropylene, a material with very low biomolecular binding characteristics, a high temperature tolerance and resistance to many standard laboratory chemicals, including DMSO
- Autoclavable
- Free from RNase, DNase, human genomic DNA and endotoxins

Wells


- Conical V-bottoms (4ti-0132) and U-shaped bottoms (4ti-0126 and 4ti-0136)
- Square wells, to make the best use of space, and to improve mixing

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint

Use

- Clear plate, perfect for storage applications
- Temperature range for use: -80°C to 120°C
- Best for precipitation, centrifugation and small volume recovery due to the conical V bottom, as the liquid gathers at the lowest point of the V for easy pipetting
- U-shaped bottom plates are best for washing, mixing and pelleting and give high surface area
- Stackable
- Suitable for adhesive and heat sealing




Options

- Available with a U-shaped bottom in 1.2 ml and 2.2 ml
- Also available with a V-shaped bottom in 2.2 ml
- Non-gamma treated as standard, gamma treatment available upon request; contact us for details and pricing
- Available barcoded upon request
- A clear silicone sealing mat (4ti-0137) is also available to fit these plates

Ordering Information

4ti-0126	96 Square Deep Well Storage Microplate, 1.2ml square wells, U-shaped bottom, clear PP, 100 plates per case
4ti-0132	96 Square Deep Well Storage Microplate, 2.2ml square wells, V-shaped bottom, clear PP, 50 plates per case
4ti-0136	96 Square Deep Well Storage Microplate, 2.2ml square wells, U-shaped bottom, clear PP, 50 plates per case
Related Products	
4ti-0137	96 Square Well Sealing Cap Mat, clear silicone, for use with square 96 well microplates and deep well storage microplates, 50 mats per case



AZENTA
LIFE SCIENCES

4titude 96 Round Deep Well Storage Microplate

1.2 ml or 2.0 ml round wells, U-shaped bottom, clear polypropylene

- These 96 deep well storage microplates are suitable for many manual and automated protocols, such as the Illumina® library and sample preparations, due to the plates' SBS footprint
- Additionally, they are suitable for use on magnetic plates for bead separation protocols
- Our round deep well storage microplates are available in two types of well volumes: 1.2 ml wells (4ti-0120) and 2.0 ml wells (4ti-0130)
- Both of these plates are ultra-flat for a completely uniform seal, either with our adhesive or heat seals, or with our accompanying silicone sealing mats (4ti-0135 mat for the 1.2 ml plate ; 4ti-0138 mat for the 2.0 ml plate)
- Our plates are manufactured in cleanroom facilities which are certified free of RNase, DNase, human genomic DNA and endotoxins
- We use only the highest medical grade virgin polypropylene with high resistance against chemicals such as DMSO, phenol, and chloroform

Key Features

- Made from polypropylene, a material with very low biomolecular binding characteristics, a high temperature tolerance and resistance to many standard laboratory chemicals, including DMSO
- Autoclavable
- Free from RNase, DNase, human genomic DNA and endotoxins

Wells

- Round wells, suitable for most applications as they reduce droplet effects and wicking
- U-shaped bottoms

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint

Use

- Clear plate, perfect for storage applications
- Suitable as a collection plate from filter systems
- Temperature range for use: -80°C to 120°C
- U-shaped bottom plates are best for washing, mixing and pelleting and give high surface area
- Stackable
- Suitable for adhesive and heat sealing



Options

- Available with 1.2 ml and 2.0 ml wells
- Non-gamma treated as standard, gamma treatment available upon request; contact us for details and pricing
- Available barcoded upon request
- 2 clear silicone sealing mats (4ti-0135 and 4ti-0138) are also available to fit these plates

Ordering Information

4ti-0120	96 Round Deep Well Storage Microplate, 1.2ml round wells, U-shaped bottom, clear PP, 50 plates per case
4ti-0130	96 Round Deep Well Storage Microplate, 2.0ml round wells, U-shaped bottom, clear PP, 50 plates per case
Related Products	
4ti-0135	96 Round Well Sealing Cap Mat, white silicone, for use with 4TI-0120 only, 100 mats per case
4ti-0138	96 Round Well Sealing Cap Mat, clear silicone, for use with round 96 well microplates and deep well storage microplates (not for use with 4ti-0120 and 4ti-0110), 50 mats per case



AZENTA
LIFE SCIENCES

4titude 96 Round Deep Well Storage Microplate, For Magnetic Separators

1 ml round wells, V-shaped bottom, clear polypropylene

- These 96 deep well storage microplates are especially designed for use with magnetic separators for bead separation protocols
- Due to the special shape of the stacking ribs, the plate sits much lower on the magnetic separator than standard round well plates thus facilitating the speed and efficiency of the separation process
- The plates are ultra-flat for a completely uniform seal; sealing possible with either adhesive or heat seals

Key Features

- Made from polypropylene, a material with very low biomolecular binding characteristics, a high temperature tolerance and resistance to many standard laboratory chemicals, including DMSO
- Autoclavable
- Free from RNase, DNase, human genomic DNA and endotoxins

Wells


- Round wells: suitable for most applications, as they reduce droplet effects and wicking
- V-shaped bottom

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint

Use

- Perfect fit to magnetic separators: the magnetic beads separate perfectly into e.g. a ring (if a ring magnet is used) allowing easy removal of supernatant from the center of the well by manual or automated pipetting
- Replaces the 96-well storage plates, round well, 0.8 ml (MIDI plate, Fisher Scientific® part number AB-0859), in e.g. Illumina® protocols
- Best for small volume recovery due to the V-shaped bottom, as the liquid gathers at the lowest point of the V for easy pipetting
- Compatible with plate readers and ideal for use with automation
- Temperature range for use: -80°C to 120°C



Options

- Ultra-low DNA binding option available (4ti-LB0125) for sensitive applications with ultra-low DNA input and for maximum DNA recovery after low temperature storage and high temperature incubation
- Available barcoded upon request
- Non-gamma treated as standard, gamma treatment available upon request; contact us for details and pricing
- Silicone sealing cap mat (4ti-0124) available separately: efficiently prevents cross contamination and sample evaporation to ensure a high degree of sample security

Ordering Information

4ti-0125	96 Round Deep Well Storage Microplate, for magnetic separators, 1.0ml round wells, V-shaped bottom, clear PP, 50 plates per case
4ti-LB0125	96 Round Deep Well Storage Microplate, for magnetic separators, low binding, 1.0ml round wells, V-shaped bottom, clear PP, 50 plates per case
Related Products	
4ti-0124	96 Round well Sealing Cap Mat, clear silicone, for use with 4ti-0125, 50 mats per case

4titude 96 Round Well Storage Microplate

200 µl, 300 µl, 330 µl or 350 µl round wells, U or V shaped bottom, clear polypropylene; also available as a low binding plate with ultra-low DNA binding properties

- These shallow 96 well storage microplates are particularly suitable for collection and preservation of samples widely used in cell biology research, molecular biology research, and drug discovery
- Four formats of this plate are available: a U-bottom plate with 300 µl or 330 µl well volumes (4ti-0110 and 4ti-0116 respectively), and a V-bottom plate with 200 µl or 330 µl well volume (4ti-LB0109 and 4ti-0117 respectively)
- Also available as a low binding plate with ultra-low DNA binding properties (4ti-LB0109) for sensitive applications such as Next Generation Sequencing sample preparation
- All of these plates are ultra-flat for a completely uniform seal, either with adhesive or heat seals, or with our accompanying silicone sealing mat (4ti-0138)
- These plates are manufactured in cleanroom facilities which are certified free of RNase, DNase, human genomic DNA and endotoxins
- We use only the highest medical grade virgin polypropylene with high resistance against chemicals such as DMSO, phenol, and chloroform

Key Features

- Made from polypropylene, a material with very low biomolecular binding characteristics, a high temperature tolerance and resistance to many standard laboratory chemicals, including DMSO
- Autoclavable
- Free from RNase, DNase, human genomic DNA and endotoxins
- Low binding option made of selected low bind polymers, no coating is used to achieve the binding characteristics; they feature stacking ribs under the deck for improved stability, strengthening the plate for use in robotic automation applications

Wells

- Conical V-bottoms (4ti-LB0109 and 4ti-0117) and U-shaped bottoms (4ti-0110 and 4ti-0116)
- Round wells, suitable for most applications as they reduce droplet effects and wicking

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint

Use

- Clear plate, perfect for storage applications
- Suitable as a collection plate from filter systems
- Temperature range for use: -80°C to 120°C



- Best for precipitation, centrifugation and small volume recovery due to the V-shaped bottom, as the liquid gathers at the lowest point of the V for easy pipetting
- U-shaped bottom plates are best for washing, mixing and pelleting and give high surface area
- Stackable
- Suitable for adhesive and heat sealing

Options

- Available with a U-shaped bottom in 300 µl and 330 µl, and with V-shaped bottom in 200 µl and 330 µl
- Ultra-low DNA binding option available (4ti-LB0109) for sensitive applications with ultra-low DNA input and for maximum DNA recovery after low temperature storage and high temperature incubation; learn more about our low binding range (not suitable for use on PCR blocks)
- Non-gamma treated as standard, gamma treatment available upon request; contact us for details and pricing
- Available barcoded upon request
- A clear silicone sealing mat (4ti-0138) is also available to fit these plates

Ordering Information

4ti-LB0109	96 Round Well Microplate, low binding, 200ul round wells, V-shaped bottom, clear PP, 50 plates per case
4ti-0110	96 Round Well Microplate, 300ul round wells, U-shaped bottom, clear PP, 100 plates per case
4ti-0116	96 Round Well Microplate, 350ul round wells, U-shaped bottom, clear PP, 100 plates per case
4ti-0117	96 Round Well Microplate, 330ul round wells, V-shaped bottom, clear PP, 100 plates per case
Related Products	
4ti-0138	96 Round Well Sealing Cap Mat, clear silicone, for use with round 96 well microplates and deep well storage microplates (not for use with 4ti-0120 and 4ti-0110), 50 mats per case



AZENTA
LIFE SCIENCES

4titude 12 Channel Reservoir Plate

Pyramid channel bottom, clear polypropylene

- Our Reservoir Plates are perfect for storing volumes of samples to be pipetted into other microplates for further applications
- Both plate formats - open format and 12 column format - are compatible with standard 12 and 96 well channel pipettes, and are suitable for automated systems, for instance the PerkinElmer® Sciclone® NGS Workstation
- These plates are manufactured in cleanroom facilities which are certified free of RNase, DNase, human genomic DNA and endotoxins
- We use only the highest medical grade virgin polypropylene with high resistance against chemicals such as DMSO, phenol, and chloroform
- Pyramid channel bases are best for small volume recovery
- The standard SBS footprint ensures its compatibility with automation

Key Features

- Made from polypropylene, a material with very low biomolecular binding characteristics, a high temperature tolerance and resistance to many standard laboratory chemicals, including DMSO
- Stackable
- Free from RNase, DNase, human genomic DNA and endotoxins

Channels

- Pyramid channel bottom, perfect for maximum sample retrieval; both the open format plate (4ti-0133) and the 12 channel plate (4ti-0131) have pyramid bases
- Rectangular channels (4ti-0131), to ensure the best use of space and allow the use of multichannel pipettes

Frame

- SBS footprint

Use

- Perfect for maximum sample retrieval; even the smallest volume can be retrieved as it gathers in the base of the pyramid, and can be easily taken in by a pipette
- Compatible with plate readers and ideal for use with automation
- Suitable as a collection plate from filter systems
- Temperature range for use -80°C to 120°C




Options

- Available barcoded upon request
- Non-gamma treated as standard, gamma treatment available upon request; contact us for details and pricing

Ordering Information

4ti-0131	12 Channel Reservoir Plate, 21ml channels, pyramid bottom, clear PP, 25 plates per case
4ti-0133	12 Channel Reservoir Plate, 290ml open format, 12 channel pyramid bottom, clear PP, 25 plates per case



4titude 96 Square Well Microplate, KingFisher™ Style

200 µl square wells, V-shaped bottom, clear polypropylene

- Suitable for use with KingFisher Flex, Apex and Presto
- Ultra-flat deck and standard SBS footprint ensures its compatibility with automation

Key Features

- Made from polypropylene, a material with very low biomolecular binding characteristics, a high temperature tolerance and resistance to many standard laboratory chemicals, including DMSO

Wells

- Conical V-bottoms
- Square wells, to make the best use of space, and to improve mixing

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint

Use

- Stackable
- Compatible with KingFisher system

Options

- Non-gamma treated as standard, gamma treatment available upon request; contact us for details and pricing
- Available barcoded upon request



Ordering Information

4ti-0151

Microplate KingFisher Style, 96 square
200ul wells, V-shaped bottom, 50 plates
per case

4titude 96 Square Deep Well Microplate, KingFisher™ Style

2.0 ml square wells, V-shaped bottom, clear polypropylene

- Suitable for use with KingFisher Duo Prime, Flex, Apex and Presto
- Designed for use with magnetic bead separators
- Ultra-flat deck and standard SBS footprint ensures its compatibility with automation

Key Features

- Made from polypropylene, a material with very low biomolecular binding characteristics, a high temperature tolerance and resistance to many standard laboratory chemicals, including DMSO

Wells

- Conical V-bottoms
- Square wells, to make the best use of space, and to improve mixing

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint

Use

- Stackable
- Compatible with KingFisher system

Options

- Non-gamma treated as standard, gamma treatment available upon request; contact us for details and pricing
- Available barcoded upon request



Ordering Information

4ti-0150

Deep Well Microplate KingFisher Style,
96 square 200ul wells, V-shaped bottom,
50 plates per case



AZENTA
LIFE SCIENCES

4titude 96 Tip Comb for Deep Well Magnets, KingFisher™ Style

Clear polypropylene

- Suitable for use with KingFisher Flex, Apex and Presto

Key Features

- Made from polypropylene, a material with very low biomolecular binding characteristics, a high temperature tolerance and resistance to many standard laboratory chemicals, including DMSO

Use

- Suitable for use with 96 deep well magnets
- Compatible with KingFisher system

Options

- Non-gamma treated as standard, gamma treatment available upon request; contact us for details and pricing
- Available barcoded upon request



Ordering Information

4ti-0152

Comb for Deep Well Magnets KingFisher
Style, 96 tips, 50 plates per case



AZENTA
LIFE SCIENCES

Assay Plates



4titude 384 Well Optically Clear Tissue Culture Plate

190 µm clear base imaging microplate, black frame

- Optically Clear Tissue Culture Plates have been designed for high content screening (HCS) assays in drug development and related areas
- They are also suitable for homogeneous assays, employing fluorescence intensity, FRET and TR-FRET where measurements are bottom-read
- The high-quality optical base plate assures the necessary accuracy and consistency for automated high throughput systems, generating optimum signal-to-noise ratios
- Using state-of-the-art manufacturing technology, we have developed a product which offers several key advantages to the end user
- Optically Clear Tissue Culture Plates are assembled using unique patented laser welding technology which reduces autofluorescence and does not inhibit cell growth
- Available with 190 µm polystyrene base
- Available TC treated: our advanced tissue-culture (TC) treatment method evenly coats each well for optimal cell adhesion properties for highest reproducibility between wells, plates and batches

Key Features

- Leak free
- Free from DNase, RNase, human genomic DNA, and cyto-toxic

Wells

- Optimum signal-to-noise ratios
- Good cell adhesion
- Wicking and bubble formation eliminated

Frame

- SBS footprint
- Alphanumeric grid reference

Use

- Recommended for confocal microscopy due to reduced cross-talk and superior base flatness
- Reduced autofluorescence ensures their suitability for fluorescent assays
- Suitable for assays that measure absorbance in the visible light range (400-900 nm wavelengths)
- Suitable for colorimetric assays
- Suitable for homogeneous assays
- Suitable for use in BMG lab tech, Molecular Devices, Promega Glomax and other plate readers
- Suitable for adhesive and heat sealing: our Moisture Barrier Seal 384 (4ti-0516/384) is recommended as it is optically clear and allows repeated imaging without removal, reducing contamination risks



Options

- Available with 190 µm polystyrene base
- All plates come gamma irradiated, apart from 4ti-0204
- TC treated option available (4ti-0201)
- Collagen 1 treated option available (4ti-0205)
- Poly D-Lysine treated option available (4ti-0206)
- Available barcoded upon request

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	14.35 ± 0.25 mm
Well depth	11.35 ± 0.25 mm
Well diameter	3.70 ± 0.10 mm
Distance to center of A1 from top edge	8.99 ± 0.25 mm
Distance to center of A1 from left edge	12.13 ± 0.25 mm
Pitch (distance between A1 and A2)	4.50 mm

Ordering Information

4ti-0201	Optically Clear Tissue Culture Plate, 384 well, tissue culture treated, gamma treated, with lid, 190µm clear base, black frame, 24 plates and lids per case
4ti-0203	Optically Clear Tissue Culture Plate, 384 well, gamma treated, 190µm clear base, black frame, 30 plates per case
4ti-0204	Optically Clear Tissue Culture Plate, 384 well, 190µm clear base, black frame, 30 plates per case
Lids	
4ti-0280	Microplate 384 Lid, without condensation rings, clear, low profile, cut corner A1/P1, 100 lids per case
4ti-0281	Microplate 384 Lid, gamma treated, without condensation rings, clear, low profile, cut corner A1/P1, 100 lids per case



AZENTA
LIFE SCIENCES

4titude 96 Well Optically Clear Tissue Culture Plate

190 µm clear base imaging microplate, black frame

- Optically Clear Tissue Culture Plates have been designed for high content screening (HCS) assays in drug development and related areas
- They are also suitable for homogeneous assays, employing fluorescence intensity, FRET and TR-FRET where measurements are bottom-read
- The high-quality optical base plate assures the necessary accuracy and consistency for automated high throughput systems, generating optimum signal-to-noise ratios
- Optically Clear Tissue Culture Plates are assembled using unique patented laser welding technology which reduces autofluorescence and does not inhibit cell growth
- Available with 190 µm polystyrene base
- Available TC treated: our advanced tissue-culture (TC) treatment method evenly coats each well for optimal cell adhesion properties for highest reproducibility between wells, plates and batches

Key Features

- Leak free
- Free from DNase, RNase, human genomic DNA, and cyto-toxic

Wells

- Optimum signal-to-noise ratios
- Good cell adhesion
- Wicking and bubble formation eliminated

Frame

- SBS footprint
- Alphanumeric grid reference

Use

- Recommended for confocal microscopy due to reduced cross-talk and superior base flatness
- Reduced autofluorescence ensures their suitability for fluorescent assays
- Suitable for assays that measure absorbance in the visible light range (400-900 nm wavelengths)
- Suitable for colorimetric assays
- Suitable for homogeneous assays
- Suitable for use in BMG lab tech, Molecular Devices, Promega Glomax and other plate readers
- Suitable for adhesive and heat sealing: our Moisture Barrier Seal 96 (4ti-0516/96) is recommended as it is optically clear and allows repeated imaging without removal, reducing contamination risks



Options

- Available with 190 µm polystyrene base
- All plates come gamma irradiated, apart from 4ti-0224
- TC treated option available (4ti-0221)
- Collagen 1 treated option available (4ti-0225)
- Poly D-Lysine treated option available (4ti-0226)
- Available barcoded upon request

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	14.35 ± 0.25 mm
Well depth	10.8 ± 0.25 mm
Well diameter	6.3 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0221	Optically Clear Tissue Culture Plate, 96 well, tissue culture treated, gamma treated, with lid, 190µm clear base, black frame, 24 plates and lids per case
4ti-0223	Optically Clear Tissue Culture Plate, 96 well, gamma treated, 190µm clear base, black frame, 30 plates per case
4ti-0224	Optically Clear Tissue Culture Plate, 96 well, 190µm clear base, black frame, 30 plates per case
Lids	
4ti-0282	Microplate 96 Lid, with condensation rings, clear, low profile, cut corner A1/H1, 80 lids per case
4ti-0283	Microplate 96 Lid, gamma treated, with condensation rings, clear, low profile, cut corner A1/H1, 80 lids per case



AZENTA
LIFE SCIENCES

4titude 24 Well Optically Clear Tissue Culture Plate

190 µm clear base imaging microplate, black frame

- Optically Clear Tissue Culture Plates have been designed for high content screening (HCS) assays in drug development and related areas
- They are also suitable for homogeneous assays, employing fluorescence intensity, FRET and TR-FRET where measurements are bottom-read
- The high-quality optical base plate assures the necessary accuracy and consistency for automated high throughput systems, generating optimum signal-to-noise ratios
- Using state-of-the-art manufacturing technology, we have developed a product which offers several key advantages to the end user
- Optically Clear Tissue Culture Plates are assembled using unique patented laser welding technology which reduces autofluorescence and does not inhibit cell growth
- Available with 190 µm polystyrene base
- Available TC treated: our advanced tissue-culture (TC) treatment method evenly coats each well for optimal cell adhesion properties for highest reproducibility between wells, plates and batches

Key Features

- Leak free
- Free from DNase, RNase, human genomic DNA, and cyto-toxic

Wells

- Optimum signal-to-noise ratios
- Good cell adhesion
- Wicking and bubble formation eliminated

Frame

- SBS footprint
- Alphanumeric grid reference

Use

- Recommended for confocal microscopy due to reduced cross-talk and superior base flatness
- Reduced autofluorescence ensures their suitability for fluorescent assays
- Suitable for assays that measure absorbance in the visible light range (400-900 nm wavelengths)
- Suitable for colorimetric assays
- Suitable for homogeneous assays
- Suitable for use in BMG lab tech, Molecular Devices, Promega Glomax and other plate readers
- Suitable for adhesive and heat sealing: our Moisture Barrier Seal 24 (4ti-0516/24) is recommended as it is optically clear and allows repeated imaging without removal, reducing contamination risks



Options

- Available with 190 µm polystyrene base
- All plates come gamma irradiated, apart from 4ti-0244
- TC treated option available (4ti-0241)
- Collagen 1 treated option available (4ti-0245)
- Poly D-Lysine treated option available (4ti-0246)
- Available barcoded upon request

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	15.00 ± 0.25 mm
Well depth	12.50 ± 0.25 mm
Well diameter	14.50 ± 0.10 mm
Distance to center of A1 from top edge	15.74 ± 0.25 mm
Distance to center of A1 from left edge	18.88 ± 0.25 mm
Pitch (distance between A1 and A2)	18.00 mm

Ordering Information

4ti-0241	Optically Clear Tissue Culture Plate, 24 well, tissue culture treated, gamma treated, with lid, 190µm clear base, black frame, 24 plates and lids per case
4ti-0243	Optically Clear Tissue Culture Plate, 24 well, gamma treated, 190µm clear base, black frame, 30 plates per case
4ti-0244	Optically Clear Tissue Culture Plate, 24 well, 190µm clear base, black frame, 30 plates per case
Lids	
4ti-0284	Microplate 24 Lid, with condensation rings, clear, low profile, cut corner A1, 80 lids per case
4ti-0286	Microplate 24 Lid, gamma treated, with condensation rings, clear, low profile, cut corner A1, 80 lids per case



AZENTA
LIFE SCIENCES

4titude 384 Well Ultra Optically Clear Plates

Ultra-clear based imaging microplate, black frame

- The main advantage of using Ultra Optically Clear Plates is in the optically ultra-clear base of the plate, which gives superior results by delivering low absorbance and high transmission, together with low background signals
- Not only this, but due to its incredible optical clarity, the ultra-clear base of the Ultra Optically Clear Plate delivers improved transmission of signals at low wavelengths compared to standard optical films
- It allows DNA measurements at 260/280 nm wavelengths in a medium or high throughput

Key Features

- Free from DNase, RNase and human genomic DNA

Base

- Ultra-clear base improves transmission for low wavelengths
- Peel-back film on the base for scratch free surface
- Optimum signal-to-noise ratios

Frame

- Alphanumeric grid reference to aid well and sample identification

Use

- Suitable for adhesive and heat sealing

Options

- Available barcoded upon request
- Low profile lid available (4ti-0280)



Ordering Information

4ti-0214	Ultra Optically Clear Plate, 384 well, clear base, black frame, 30 plates per case
Lids	
4ti-0280	Microplate 384 Lid, without condensation rings, clear, low profile, cut corner A1/P1, 100 lids per case



AZENTA
LIFE SCIENCES

4titude 96 Well UV Ultra Optically Clear Cap Plates

Ultra-clear based imaging microplate

- The main advantage of using Ultra Optically Clear Plates is in the optically ultra-clear base of the plate, which gives superior results by delivering low absorbance and high transmission, together with low background signals
- Not only this, but due to its incredible optical clarity, the ultra-clear base of the Ultra Optically Clear Plate delivers improved transmission of signals at low wavelengths compared to standard optical films
- It allows DNA measurements at 260/280 nm wavelengths in a medium or high throughput

Key Features

- Free from RNase, DNase and human genomic DNA

Base

- Ultra-clear base improves transmission for low wavelengths
- Optimum signal-to-noise ratios

Frame


- Alphanumeric grid reference to aid well and sample identification

Use

- Suitable for adhesive and heat sealing

Options

- Available barcoded upon request
- Lid available (4ti-0290)



Ordering Information

4ti-0234	Ultra Optically Clear Plate, 96 well, clear base, clear frame, 30 plates per case
Lids	
4ti-0290	Universal Microplate Lid, without condensation rings, clear, low profile, no cut corner, 50 lids per case



4titude 24 Well Assay Plate

1.88 ml round wells, flat base, polystyrene, cut corner A1

- The Azenta 24 well polystyrene assay plates have flat bottom wells and have been designed for fluorescence applications
- The plates give optimum results for most top-reading instruments and conform to standard SBS footprint
- Black plates have a low background fluorescence and minimize light scattering
- This high quality plate assures the necessary accuracy and consistency for automated high throughput systems



Key Features

- Made from polystyrene, a hard material with smooth surfaces
- Non-gamma treated as standard
- Leak-free
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Flat bottom wells, designed for fluorescence and cell culture application
- Round wells, for reduced droplet effects and wicking

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint

Use

- Best used for cell culture, fluorescence
- Compatible with top reading plate readers and ideal for use with automation

Options

- Available as a black microplate, ideal for low background fluorescence; it minimizes light scattering, suitable for fluorescent assays, recommended for top reading fluorescence instrumentation
- Gamma treatment available upon request
- Available barcoded upon request
- Our Moisture Barrier Seal 24 (4ti-0516/24) is recommended as sealing option, as it is optically clear and allows repeated imaging without removal, reducing contamination risks
- Lids available

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	15.00 ± 0.25 mm
Well depth	12.50 ± 0.25 mm
Well diameter	14.50 ± 0.10 mm
Distance to center of A1 from top edge	15.74 ± 0.25 mm
Distance to center of A1 from left edge	18.88 ± 0.25 mm
Pitch (distance between A1 and A2)	18.00 mm

Ordering Information

4ti-0262	24 Well Assay Plate, 1.88ml round wells, flat base, black PS, cut corner A1, 100 plates per case
Lids	
4ti-0284	Microplate 24 Lid, with condensation rings, clear, low profile, cut corner A1, 80 lids per case
4ti-0286	Microplate 24 Lid, gamma treated, with condensation rings, clear, low profile, cut corner A1, 80 lids per case



4titude 96 Well Assay Plate

0.35 ml round wells, flat base, polystyrene, cut corner A1/H1

- The Azenta 96 well black assay plate has been specifically designed for fluorescence
- It is also suitable for homogeneous assays employing fluorescence intensity, FRET and TR-FRET where measurements are top-read
- This high quality plate assures the necessary accuracy and consistency for automated high throughput systems, generating optimum signal to noise ratios
- The white assay microplates have been designed for luminescence applications, such as Luciferase Reporter Assays and scintillation applications
- The white plate maximizes signal intensity in cases of low signal from some or all the wells, and it is designed to give optimum results for most top reading instruments, and to conform to standard SBS footprint

Key Features

- Made from polystyrene, a hard polished material
- Leak-free
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Flat bottom wells, suitable for top reading and cell culture application
- Round wells, for reduced droplet effects and wicking

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint

Use

- Best used for cell culture, fluorescence, luminescence, and light assays, ELISA, and homogeneous assays
- Recommended for top reading fluorescence instrumentation
- Compatible with plate readers and ideal for use with automation

Options

- Available as a black plate, for low background fluorescence and minimum light scattering
- White microplate, suitable for visible light assays and imaging, and recommended for top reading instrumentation
- Non-gamma treated as standard; gamma treatment available upon request
- Available barcoded upon request



- Our Moisture Barrier Seal 96 (4ti-0516/96) is recommended as a sealing option, as it is optically clear and allows repeated imaging without removal, reducing contamination risks
- Lids available

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	14.35 ± 0.25 mm
Well depth	10.8 ± 0.25 mm
Well diameter	6.3 ± 0.10 mm
Distance to center of A1 from top edge	11.24 ± 0.25 mm
Distance to center of A1 from left edge	14.38 ± 0.25 mm
Pitch (distance between A1 and A2)	9.00 mm

Ordering Information

4ti-0263	96 Well Assay Plate, 0.35ml round wells, flat base, black PS, cut corner A1/H1, 100 plates per case
4ti-0273	96 Well Assay Plate, 0.35ml round wells, flat base, white PS, cut corner A1/H1, 100 plates per case
Lids	
4ti-0282	Microplate 96 Lid, with condensation rings, clear, low profile, cut corner A1/H1, 80 lids per case
4ti-0283	Microplate 96 Lid, gamma treated, with condensation rings, clear, low profile, cut corner A1/H1, 80 lids per case



AZENTA
LIFE SCIENCES

4titude 384 Well Assay Plate, Solid Wells

0.12 ml rounded square wells, flat base, polystyrene, cut corner A1/P1

- The Azenia 384 well black assay plate has been specifically designed for fluorescence
- It is also suitable for homogeneous assays employing fluorescence intensity, FRET and TR-FRET where measurements are top-read
- This high quality plate assures the necessary accuracy and consistency for automated high throughput systems, generating optimum signal to noise ratios
- The 384 well white solid bottom assay microplate has been specifically designed for luminescence applications, such as Luciferase Reporter Assays and scintillation applications
- It reduces well-to-well crosstalk, and the solid white color boosts signal in cases of low signal from some or all the wells
- The non-treated clear microplate is ideal for colorimetric assays and sample storage
- The rounded square wells eliminate wicking (capillary action)
- The flat bottom is ideal for optical reading

Key Features

- Made from polystyrene, a hard material with optical clarity
- Leak-free
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Wells

- Flat bottom wells, suitable for top reading and cell culture application
- Rounded square wells, to ensure the best use of space and improve sample mixing, suitable for small volumes
- Wicking and bubble formation eliminated

Frame

- Alphanumeric grid reference to aid well and sample identification
- SBS footprint



Use

- Best used for cell culture, fluorescence, luminescence, imaging and light assays, ELISA, and homogeneous assays
- Recommended for top reading fluorescence instrumentation
- Compatible with plate readers and ideal for use with automation
- Suitable for use in BMG labtech, Molecular Devices, Promega Glomax and other plate readers; for full instrument compatibility please contact us

Options

- Black assay microplate, for low background fluorescence and minimum light scattering
- White microplate, suitable for visible light assays and imaging, and recommended for top reading instrumentation
- The clear microplate offers the best solution for absorption, ELISA, spectrophotometric and colorimetric assays, and storage applications
- Non-gamma treated as standard; gamma treatment available upon request
- Available barcoded upon request
- Our Moisture Barrier Seal 384 (4ti-0516/384) is recommended as sealing option, as it is optically clear and allows repeated imaging without removal, reducing contamination risks
- Lids available

4titude 384 Well Assay Plate

Specifications

Parameter	Value
Plate length	127.76 ± 0.25 mm
Plate width	85.48 ± 0.25 mm
Plate height	14.35 ± 0.25 mm
Well depth	11.35 ± 0.25 mm
Well diameter	3.70 ± 0.10 mm
Distance to center of A1 from top edge	8.99 ± 0.25 mm
Distance to center of A1 from left edge	12.13 ± 0.25 mm
Pitch (distance between A1 and A2)	4.50 mm

Ordering Information

4ti-0254	384 Well Assay Plate, 0.12ml rounded square wells, flat base, clear PS, cut corner A1/P1, 100 plates per case
4ti-0264	384 Well Assay Plate, 0.12ml rounded square wells, flat base, black PS, cut corner A1/P1, 100 plates per case
4ti-0274	384 Well Assay Plate, 0.12ml rounded square wells, flat base, white PS, cut corner A1/P1, 100 plates per case
Lids	
4ti-0280	Microplate 384 Lid, without condensation rings, clear, low profile, cut corner A1/P1, 100 lids per case
4ti-0281	Microplate 384 Lid, gamma treated, without condensation rings, clear, low profile, cut corner A1/P1, 100 lids per case

4titude Microplate Sealing Consumables

Azenta offers a wide range of plate sealing solutions for whatever your requirements are. You can choose between sealing with strip caps, mats, lids, adhesive seals in strip or plate format, and heat seals in flexible formats up to plate size. The choice of an optimized sealing solution is particularly important for (q)PCR because thermal cycling can be associated with evaporation of reaction reagents.



Heat Sealing & Adhesive Sealing

Dependent on your application requirements we offer a wide selection of materials to choose from within both our adhesive seal and heat seal ranges. You have the option to choose your seal based on a wide variety of properties offered, including peelability, pierceability, gas permeability, optical clarity, temperature stability and solvent resistance.

The best sealing results can be obtained by using flat rigid plates like FrameStar plates and heat seals, using reliable, high quality sealing instruments for seal application, like the Semi-Automated Sheet Heat Sealer or the Automated Roll Heat Sealer.

Alternatively, popular adhesive seals need to be applied well using a seal roller or seal applicator.

Azenta seals are produced and processed under strictly controlled environmental conditions and according to our ISO standard manufacturing. All of our seals are DNase, RNase, human genomic DNA, dust, endotoxin/pyrogen free. Dimensional and functional tests are performed on all production lots. If your seal of choice is not offered gamma treated as a stock product, then please contact us; we can offer custom gamma treatment of any seals if required.

Caps, Lids & Mats

As an alternative to sealing films, Azenta offers multiple types of cap strips and individual well cap mats for sealing of both plates and tubes - domed and flat, strips of 8, strips of 12, and our new optically superior Strips of 8 Flat optically superior Crystal Clear caps.

A variety of rigid polystyrene lids are available for PCR plates and microplates, as well as silicone sealing mats for our storage plate range. We offer lids that are compatible with our FrameStar PCR plate range, and all assay plate ranges including the Optically Clear Plates range and Ultra Optically Clear Plates. The silicone cap mats are for use with our storage plates, and come in a variety of formats depending on the plate.

Heat Sealing Consumables & Instrumentation

Heat sealing is the gold standard method of plate and tube sealing. It prevents sample loss and maximizes sample security, by ensuring a complete seal and preventing evaporation, leakage and contamination.

The sealing performance of heat seals is superior to all other methods including cap, mat and adhesive sealing. The variability of sealing integrity seen when using adhesive seals, caps or mats is reduced. The optimized sealing performance of a heat seal allows the use of smaller reagent volumes, leading to reagent cost savings and thus making heat sealing the most cost efficient sealing method for a wide range of applications.

Heat seals are available as sheets, for manual or semi-automated heat sealers, as well as in different roll formats for automated sealers. Azenta sealing consumables are compatible with a wide range of heat sealers, please refer to the Instrument Compatibility Table on pages 162-165.

Depending on throughput, we recommend our Semi-Automated Sheet Heat Sealer (page 211) or Automated Roll Heat Sealer (page 213) for applying your heat seal.

In addition to instruments and consumables we also offer a Thermal Test Film for the optimization and troubleshooting of heat sealing applications.

Economic efficiency

Investing in heat sealing solutions leads to a per plate reduction of sealing costs for all applications including storage, PCR and qPCR as shown in the table below for (q)PCR. The cost is further reduced by changing to the use of rolls seals rather than sheet seals.

Benefit of Heat Sealing	Why this is important
Enables a tight seal around each sample through melding of seal with sealing rings	Maintains sample integrity and minimizes evaporation
Achieves a consistent seal across the plate through unrivalled optimized sealing performance of Azenta heat sealers	Provides the ability to use the whole plate; Eliminates any potential for variation between wells as a result of evaporation
Fast and convenient application through touch of a button	Ability to seal a number of plates quickly and efficiently
No sample contamination through sealing materials, as seal surface is identical to well material	Sample contamination may affect the sample and resulting experiment

Save your time

The application of heat seals is also easier and faster than when using caps or adhesive seals. Typical sealing times of a semi-automatic heat sealer, such as the Azenta Semi-Automated Sheet Heat Sealer, are around 2.5 seconds. A fully automatic roll heat sealer, such as the Azenta Automated Roll Heat Sealer, allows for sealing cycle times of less than 15 seconds.

Set your standard

Reproducible sealing quality can be guaranteed by standardizing the sealing parameters: time, temperature and pressure. Azenta offers the widest choice of heat seal materials available, with sheet formats for manual or semi-automatic heat sealers and roll formats for automated heat sealers.

Choose your application

Depending on the material of the plate (PP, PE, PS, COC, PC), the presence of solvents like DMSO in your sample and the storage or application temperatures required, we can offer a wide range of seals covering 100% DMSO storage and sealing integrity temperatures between -200°C to +120°C.

Highest quality standards

Our seals are produced according to our ISO standard manufacturing. All of our seals are DNase, RNase, human genomic DNA, dust, endotoxin and pyrogen free. If your seal of choice is not offered gamma treated as a stock product then please contact us as we can offer custom gamma treatment of any seals if required.



Optically Clear Peelable Heat Seal

Peelable heat sealing film, optically clear; suitable for qPCR and optical applications

- Heat sealing offers a 100% effective method for plate sealing for a complete seal integrity, as well as being quick and cost effective
- Our Clear Heat Seal is an optically clear laminate film forming a peelable seal to polypropylene, polyethylene, polystyrene and polycarbonate plates
- Samples can also be accessed by pre-piercing with a blade, needle or our Pierce Plate (4ti-0398)
- The optical clarity of this seal enables its use for sealing plates required for imaging use, including fluorescent detection methods such as qPCR and colorimetric assays
- The Clear Heat Seal forms a complete seal to a plate enabling both low temperature uses, including low temperature storage, and high temperature uses, such as PCR (when used with a pressurized heated lid)
- This seal demonstrates moderate solvent resistance and can be utilized for short term compound storage at room temperature
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers such as our Automated Roll Heat Sealer
- Our sheet seals are inter-leaved with paper sheets to help denote which side is the sealing side and aid removal of one sheet at a time from the pack
- For applications requiring high tensile strength (e.g. bead mill applications) please see our Clear Heat Seal Plus

Key Features

- Peelable
- Seal integrity range: -80°C to 80 °C (to 110°C with a pressurized heated PCR lid)
- Optically clear
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Applications: imaging, fluorescent detection, and colorimetric assays
- Suitable for PCR and qPCR
- Seals polypropylene, polyethylene, polystyrene and polycarbonate plates



*110°C with pressurized heated lid



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Non-gamma treated as standard; gamma treatment available upon request
- Also available with our FrameSeal™ technology, a disposable, rigid, plastic frame perfect for use with a robotic gripper to be used within an automation cell

Ordering Information

4ti-0540	Clear Heat Seal, peelable heat sealing film, 1 roll (500m x 78mm) ¹
4ti-0540/80	Clear Heat Seal, peelable heat sealing film, 1 roll (80m x 78mm) ¹
4ti-0540/REMP	Clear Heat Seal, peelable heat sealing film, 1 roll (500m x 78mm) ³
4ti-0540S	Clear Heat Seal, peelable heat sealing film, 1 sample roll (5m x 78mm) ¹
4ti-0542	Clear Heat Seal, peelable heat sealing film, 1 roll (350m x 115mm) ²
4ti-0542/REMP	Clear Heat Seal, peelable heat sealing film, 1 roll (350m x 115mm) ⁴
4ti-0542S	Clear Heat Seal, peelable heat sealing film, 1 sample roll (5m x 115mm) ²
4ti-0541	Clear Heat Seal, peelable heat sealing film, 100 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®.

³ Compatible with Azenta Portrait Heat Sealer.

⁴ Compatible with Azenta Landscape / Stacking Heat Sealers



AZENTA
LIFE SCIENCES

Optically Clear Weld Heat Seal

Optically clear heat sealing film, non-peelable, difficult to pierce; suitable for qPCR, optical applications and storage


- Heat sealing offers a 100% effective method for plate sealing for a complete seal integrity, as well as being quick and cost effective
- Our Clear Weld Heat Seal is an optically clear polymer film forming a permanent seal to polypropylene plates
- The optical clarity of this seal enables its use for sealing plates required for imaging use, including fluorescent detection methods such as qPCR and colorimetric assays
- The Clear Weld Heat Seal forms a complete seal to a plate, enabling both low and very high temperature uses, including low temperature storage and high temperature incubations
- This seal is suitable for PCR/qPCR, even without the use of a pressurized heated lid, and is 100% effective when used in water bath thermal cyclers
- The permanent nature of this 100% effective seal renders it suitable for the storage and disposal of hazardous material
- Clear Weld Heat Seal demonstrates a good solvent resistance and can be utilized for long term compound storage
- Samples can be accessed by pre-piercing with a blade, needle or our Pierce Plate (4ti-0398)
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer

Key Features

- Permanent seal
- Difficult to pierce
- Non-peelable
- Seal integrity range: -80°C to 110°C
- DMSO and solvent resistant
- Autoclavable (121°C)* Once autoclaved, not recommended for PCR applications
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use


- Applications: suitable for PCR and qPCR, long term storage, and disposal of hazardous materials
- Seals polypropylene plates
- Compatible with water bath thermal cyclers



Seal Integrity
Temperature Range

-80°C
110°C

*Once autoclaved, not recommended for PCR



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Non-gamma treated as standard; gamma treatment available upon request
- Also available with our FrameSeal™ technology, a disposable, rigid, plastic frame perfect for use with a robotic gripper to be used within an automation cell

Ordering Information

4ti-0573	Clear Weld Heat Seal, optically clear heat sealing film, 1 roll (610m x 78mm) ¹
4ti-0573/122	Clear Weld Heat Seal, optically clear heat sealing film, 1 roll (122m x 78mm) ¹
4ti-0573S	Clear Weld Heat Seal, optically clear heat sealing film, 1 sample roll (5m x 78mm) ¹
4ti-0574	Clear Weld Heat Seal, optically clear heat sealing film, 1 roll (500m x 115mm) ²
4ti-0574S	Clear Weld Heat Seal, optically clear heat sealing film, 1 sample roll (5m x 115mm) ²
4ti-0575**	Clear Weld Heat Seal, optically clear heat sealing film, 100 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®

** Seals are inter-leaved with paper sheets. This helps to denote which side is the sealing side, plus aids removal of one sheet at a time from the pack



AZENTA
LIFE SCIENCES

Strong Clear Heat Seal Plus

High tensile strength heat sealing film, optically clear, peelable; suitable for bead mill applications

- Heat sealing offers a 100% effective method for plate sealing for a complete seal integrity, as well as being quick and cost effective
- Our Clear Heat Seal Plus is a clear polymer film forming a peelable seal to polypropylene, polystyrene and cyclic olefin copolymer (COC) plates
- The excellent tensile strength of this seal enables its use for sealing of microplates during homogenisation or disruption of seeds or other material such as bead mill applications
- Clear Heat Seal Plus forms a complete seal to a plate also enabling short term sample storage
- The Clear Heat Seal Plus is thicker than our standard Clear Heat Seal for application requiring high tensile strength, but is not as suitable for optical applications due to decreased transmission through the thicker sealing layers
- For optical applications please refer to our Clear Heat Seal
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer

Key Features

- Peelable
- Seal integrity range: -80°C to 80°C (110°C with a pressurized heated lid)
- High tensile strength for bead mill applications
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Applications: suitable for bead mill applications and PCR



Seal Integrity Temperature Range	-80°C 80°C
----------------------------------	---------------

*110°C with pressurized heated lid



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in two roll dimensions, to suit your choice of automatic heat sealing equipment
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0549	Clear Heat Seal Plus, high tensile strength heat sealing film, 1 roll (250m x 78mm) ¹
4ti-0549/S	Clear Heat Seal Plus, high tensile strength heat sealing film, 1 sample roll (5m x 78mm) ¹
4ti-0548	Clear Heat Seal Plus, high tensile strength heat sealing film, 1 roll (250m x 115mm) ²
4ti-0548/S	Clear Heat Seal Plus, high tensile strength heat sealing film, 1 sample roll (5m x 115mm) ²
4ti-05481	Clear Heat Seal Plus, high tensile strength heat sealing film, 100 sheets (125 x 78mm) per case

¹ Compatible with Agilent (Velocity 11) PlateLoc®

² Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube



AZENTA
LIFE SCIENCES

Clear Heat Seal Easily Pierceable

Thin polyester heat sealing film, easily pierceable with autosampler needles/ABI® 3730; suitable for PCR, qPCR and optical applications

- Heat sealing offers a 100% effective method for plate sealing for a complete seal integrity, as well as being quick and cost effective
- Our Clear Heat Seal Easily Pierceable is an optically clear polyester backed film, forming a pierceable seal to polypropylene, polyethylene, polystyrene and cyclic olefin copolymer (COC) plates
- The optical clarity of this seal enables its use for sealing plates required for imaging use, including fluorescent detection methods such as qPCR and colorimetric assays
- Its pierceability renders it useful for automation and for use on needle, capillary and tip based liquid handling systems
- Effective on the ABI® 3730 capillary sequencer, removing the need for the use of expensive septa mats
- The Clear Heat Seal Easily Pierceable forms a complete seal to a plate enabling moderately low and high temperature uses, including PCR when using a pressurized heated lid thermal cycler
- Demonstrates a moderate solvent resistance and can be utilized for short term compound storage
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer
- A perforated roll is available, for easy removal of sheets, for use with manual and semi-automated sealers

Key Features

- Permanent seal
- Easily pierceable with autosampler needles/ABI® 3730
- Seal integrity range: -20°C to 80°C (or 110°C when used with pressurized heated PCR lid)
- Moderate solvent resistance
- Optically clear
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Seals polypropylene, polyethylene, polystyrene and cyclic olefin copolymer (COC) plates
- Suitable for capillary sequencers and automated liquid handlers, e.g. ABI® 3730
- Suitable for PCR and qPCR, and short term storage



Seal Integrity Temperature Range	-80°C
	80°C

*110°C with pressurized heated lid



Options

- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Sheet format: 125 x 78 mm, available from a roll with 1,000 perforated sheets; seals all SBS plate formats, from 12 well to 1536 well
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0580	Clear Heat Seal, easily pierceable with autosampler needles, thin polyester heat sealing film, 1 roll (610m x 78mm)¹
4ti-0580/122	Clear Heat Seal, easily pierceable with autosampler needles, thin polyester heat sealing film, 1 roll (122m x 78mm)¹
4ti-0580S	Clear Heat Seal, easily pierceable with autosampler needles, thin polyester heat sealing film, 1 sample roll (5m x 78mm)¹
4ti-0582	Clear Heat Seal, easily pierceable with autosampler needles, thin polyester heat sealing film, 1 roll (500m x 115mm)²
4ti-0582S	Clear Heat Seal, easily pierceable with autosampler needles, thin polyester heat sealing film, 1 sample roll (5m x 115mm)²
4ti-0581	Clear Heat Seal, easily pierceable with autosampler needles, thin polyester heat sealing film, 1 perforated roll with sheets (125 x 78mm)

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®



AZENTA
LIFE SCIENCES

4titude Peel Foil Heat Seal

Peelable heat sealing foil; suitable for low temperature storage, high temperature uses and PCR

- Heat sealing offers a 100% effective method of plate sealing, for complete seal integrity, as well as being quick and cost effective
- Our Peel Heat Seal is a laminate seal compatible with polypropylene plates
- It can be removed from polypropylene plates by peeling, even with a plate which has been removed directly from -80°C storage
- Peel Heat Seal forms a complete seal to a plate enabling very low temperature uses, including very low temperature storage, and high temperature uses, such as PCR (when used with a pressurized heated lid)
- The seal demonstrates moderate solvent resistance and can be utilized for short term compound storage at room temperature
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer

Key Features

- Peelable from polypropylene and COC plates
- Seal integrity range: -80°C to 90°C (110°C when used with pressurized heated lid)
- Good solvent resistance including DMSO
- Autoclavable (121°C)
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Suitable for short term room temperature storage
- Suitable for very low temperature storage



Seal Integrity Temperature Range	-80°C
	90°C

*110°C with pressurized heated lid

Peelable	Autoclavable

Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0520	Peel Heat Seal, peelable heat sealing foil, 1 roll (610m x 78mm)¹
4ti-0520/122	Peel Heat Seal, peelable heat sealing foil, 1 roll (122m x 78mm)¹
4ti-0520S	Peel Heat Seal, peelable heat sealing foil, 1 sample roll (5m x 78mm)¹
4ti-0522	Peel Heat Seal, peelable heat sealing foil, 1 roll (500m x 115mm)²
4ti-0522S	Peel Heat Seal, peelable heat sealing foil, 1 sample roll (5m x 115mm)²
4ti-0521	Peel Heat Seal, peelable heat sealing foil, 100 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®

4titude Universal Foil Peel Heat Seal

Peelable heat sealing foil with wide material compatibility; suitable for low temperature storage, high temperature uses and PCR

- Heat sealing offers a 100% effective method of plate sealing, for complete seal integrity, as well as being quick and cost effective
- Our Universal Peel Heat Seal is a laminate seal compatible with PP, PE, PS, COC and PC plates, providing the highest flexibility in plate material choice
- Sample access is possible by peeling from the compatible materials and also by piercing the seal with needles, but not with plastic tips
- It is resealable by applying another Universal Peel Heat Seal directly on top of a previously pierced seal
- Universal Peel Heat Seal forms a complete seal to a plate enabling low temperature uses as well as high temperature uses, such as PCR (when used with a pressurized heated lid)
- The seal demonstrates moderate solvent resistance
- Use of Universal Peel Heat Seal makes roll changes unnecessary even for customers using different microplate materials, because with minor adjustments of sealing parameters all common microplate materials can be sealed to
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer

Key Features

- Pierceable with needle, but not with standard pipette tips
- Peelable
- Seal integrity temperature range: -80°C to 90°C (110°C when used with a pressurized heated lid)
- Autoclavable (121°C)
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Suitable for low temperature sample storage and high temperature uses, such as PCR
- Compatible with PP, PE, PS, COC and PC plates
- Wide material compatibility allows for high throughput sealing of different plates without the need for roll changes
- Resealable by applying another Universal Peel Heat Seal directly on top of a previously pierced seal



Seal Integrity Temperature Range	-80°C 90°C
----------------------------------	---------------

*110°C with pressurized heated lid

Peelable	Autoclavable

Options


- Sheet format: 125 mm x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in two roll dimensions, to suit your choice of automatic heat sealing equipment

Ordering Information

4ti-0523	Universal Peel Heat Seal, peelable heat sealing foil with wide material compatibility, 1 roll (610m x 78mm) ¹
4ti-0523S	Universal Peel Heat Seal, peelable heat sealing foil with wide material compatibility, 1 sample roll (5m x 78mm) ¹
4ti-0524	Universal Peel Heat Seal, peelable heat sealing foil with wide material compatibility, 1 roll (500m x 115mm) ²
4ti-0524S	Universal Peel Heat Seal, peelable heat sealing foil with wide material compatibility, 1 sample roll (5m x 115mm) ²
4ti-05231	Universal Peel Heat Seal, peelable heat sealing foil with wide material compatibility, 100 sheets (125 x 78mm) per case

¹ Compatible with Azena Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®



AZENTA
LIFE SCIENCES

4titude Individual Access Peel Heat Seal

Peelable heat sealing foil, 96 individual seals with tabs, or 12 strips each covering 8 wells; suitable for very low temperature storage/high temperature uses/PCR


- Heat sealing offers a 100% effective method of plate sealing, for complete seal integrity, as well as being quick and cost effective
- Our Individual Access Peel Heat Seal is a laminate seal compatible with polypropylene plates, featuring 96 individual foil seal spots or 12 strips of individual spots on a removable backing
- These seals result in individually sealed tubes/strips, and they can be removed from polypropylene plates by peeling, even with a plate which has been removed directly from -80°C storage
- Individual Access Peel Heat Seal forms a complete seal to a plate enabling very low temperature uses, including very low temperature storage, and high temperature uses, such as PCR (when used with a pressurized heated lid)
- The seal demonstrates moderate solvent resistance and can be utilized for short term compound storage at room temperature
- This seal is available as sheets, for use with manual and semi-automated sealers, such as our Semi-Automated Sheet Heat Sealer (using the 4ti-0613 Individual Access adapter)

Key Features

- 96 individual foil seal spots or 12 strips of 8 spots on a removable backing
- 4 pin holes for exact positioning in special adapters of the Semi-Automated Sheet Heat Sealer
- Peelable from polypropylene and COC plates
- Seal integrity range: -80°C to 90°C (110°C when used with pressurized heated lid)
- Good solvent resistance including DMSO
- Autoclavable (121°C)
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Suitable for short term room temperature storage
- Suitable for very low temperature storage
- Best used in combination with our Individual Access plates and Breakable Vertically PCR Plates




Seal Integrity Temperature Range	-80°C 90°C
*110°C with pressurized heated lid	
Peelable	Autoclavable

Options

- Available as 96 individual seals with tabs per sheet
- Also available as 12 strips of 8 individual seals per sheet
- Sheet format: 127 x 100 mm
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0521/RA-TAB	Individual Access Peel Heat Seal, with tabs, peelable heat sealing foil, 96 individual seals with tabs, sheet format, 100 sheets (127 x 100mm) per case
4ti-0521/RA-8	Individual Access Peel Heat Seal, peelable heat sealing foil, 12 strips, each covering 8 wells, sheet format, 100 sheets (127 x 100mm) per case



AZENTA
LIFE SCIENCES

DMSO Resistant Peel Heat Seal

Solvent resistant heat sealing foil, peelable; suitable for low and room temperature compound storage

- Our DMSO Resistant Peel Heat Seal is a foil seal compatible with polypropylene and forming an excellent seal to cyclic olefin copolymer (COC) plates
- The solvent resistance of this seal enables its use for low and room temperature compound storage in Dimethyl Sulfoxide (DMSO) and organic solvents
- 100% DMSO can be stored at room temperature for 12 months without deterioration of the seal
- It forms a weld seal to polyethylene plates and cannot be peeled off
- Access is by piercing using a blade, needle, or an Azenta Pierce Plate (4ti-0398)
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer
- The sheet seals are inter-leaved with paper sheets to help denote which side is the sealing side and to aid removal of one sheet at a time from the pack

Key Features

- Permanent seal to polyethylene
- Peelable seal to polypropylene and COC
- Seal integrity range: -80°C to 40°C
- High solvent resistance, including 100% DMSO
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Suitable for long term storage



Seal Integrity
Temperature Range

-80°C
40°C



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0585	DMSO Resistant Peel Heat Seal, solvent resistant heat sealing foil, 1 roll (500m x 78mm) ¹
4ti-0585/100	DMSO Resistant Peel Heat Seal, solvent resistant heat sealing foil, 1 roll (100m x 78mm) ¹
4ti-0585/REMP	DMSO Resistant Peel Heat Seal, solvent resistant heat sealing foil, 1 roll (500m x 78mm) ³
4ti-0585S	DMSO Resistant Peel Heat Seal, solvent resistant heat sealing foil, 1 sample roll (5m x 78mm) ¹
4ti-0586	DMSO Resistant Peel Heat Seal, solvent resistant heat sealing foil, 1 roll (500m x 115mm) ²
4ti-0586/REMP	DMSO Resistant Peel Heat Seal, solvent resistant heat sealing foil, 1 roll (500m x 115mm) ⁴
4ti-0586S	DMSO Resistant Peel Heat Seal, solvent resistant heat sealing foil, 1 sample roll (5m x 115mm) ²
4ti-0587	DMSO Resistant Peel Heat Seal, solvent resistant heat sealing foil, 1 perforated roll with 100 sheets (125 x 78mm)

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®

³ Compatible with Azenta Portrait Heat Sealer

⁴ Compatible with Azenta Landscape / Stacking Heat Sealers



AZENTA
LIFE SCIENCES

Pierceable Heat Seal

Pierceable heat sealing foil, high solvent resistance, resealable, suitable for PCR/storage/shipping



- Heat sealing offers a 100% effective method of plate sealing, for complete seal integrity, as well as being quick and cost effective
- Pierce Heat Seal is compatible with polypropylene and polystyrene plates
- This seal demonstrates good solvent resistance and can be used for low temperature and room temperature compound storage in DMSO and organic solvents
- Pierce Heat Seal can be pierced with a pipette tip manually, or by a liquid handling robot
- This seal can be resealed by applying another Pierce Seal straight on top of a previously pierced seal
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer
- A blue stripe on the sheet foils clearly indicates the non-sealing surface, for ease of seal orientation and handling
- Pierce Heat Seal sheets are also available with a printed grid reference on the non-sealing surface

Key Features

- Pierceable
- Resealable
- Seal integrity range: -20°C to 120°C
- Good solvent resistance
- Autoclavable (121°C)
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Seals polypropylene and polystyrene plates
- Suitable for long term storage



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Available with printed grid reference (in sheet format, 4ti-0531/GR)
- Non-gamma treated as standard; gamma treatment available upon request
- Custom printing available on request

Ordering Information

4ti-0530	Pierce Heat Seal, pierceable heat sealing foil, 1 roll (610m x 78mm) ¹
4ti-0530/122	Pierce Heat Seal, pierceable heat sealing foil, 1 roll (122m x 78mm) ¹
4ti-0530/REMP	Pierce Heat Seal, pierceable heat sealing foil, 1 roll (610m x 78mm) ³
4ti-0530S	Pierce Heat Seal, pierceable heat sealing foil, 1 sample roll (5m x 78mm) ¹
4ti-0532	Pierce Heat Seal, pierceable heat sealing foil, 1 roll (500m x 115mm) ²
4ti-0532/REMP	Pierce Heat Seal, pierceable heat sealing foil, 1 roll (500m x 115mm) ⁴
4ti-0532S	Pierce Heat Seal; pierceable heat sealing foil, 1 sample roll (5m x 115mm) ²
4ti-0531	Pierce Heat Seal, pierceable heat sealing foil, 100 sheets (125 x 78mm) per case
4ti-0531/GR	Pierce Heat Seal, with Grid Reference, pierceable heat sealing foil, 100 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®

³ Compatible with Azenta Portrait Heat Sealer

⁴ Compatible with Azenta Landscape / Stacking Heat Sealers

Individual Access Pierceable Heat Seal

Pierceable heat sealing foil, 96 individual seals in sheet format, high solvent resistance, resealable; suitable for PCR/storage/shipping


- Heat sealing offers a 100% effective method of plate sealing, for complete seal integrity, as well as being quick and cost effective
- Individual Access Pierce Heat Seal is compatible with polypropylene and polystyrene plates, featuring 96 individual foil seal spots on a removable backing
- This seal demonstrates good solvent resistance and can be used for low temperature and room temperature compound storage in DMSO and organic solvents
- These seals result in individually sealed tubes, and they can be pierced with a pipette tip manually, or by a liquid handling robot
- Individual Access Pierce Heat Seal can be resealed by applying another Individual Access Pierce Seal straight on top of a previously pierced seal
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer

Key Features

- 96 individual foil seal spots on a removable backing
- Pierceable
- 4 pin holes for exact positioning in special adapters of the Semi-Automated Sheet Heat Sealer
- Resealable
- Seal integrity range: -20°C to 120°C
- Good solvent resistance
- Autoclavable (121°C)
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Seals polypropylene and polystyrene plates
- Suitable for long term storage
- Best used in combination with our Individual Access plates




Options

- Sheet format: 127 x 100 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0531/RA

Individual Access Pierce Heat Seal, pierceable heat sealing foil, 96 individual seals in sheet format, 100 sheets (127 x 100mm) per case



AZENTA
LIFE SCIENCES

Pierce Heat Seal Strong

Strong heat sealing foil, peelable from COC plates, pierceable, suitable for PCR, sample shipping, compound storage

- Heat sealing offers an 100% effective method of plate sealing, for complete seal integrity, as well as being quick and cost effective
- Pierce Heat Seal Strong is compatible with polypropylene and COC plates
- This seal is peelable from COC plates and gives a weld seal to polypropylene plates
- Demonstrating good solvent resistance it can be used for low temperature and room temperature compound storage in DMSO and organic solvents
- The seal can be pierced with a pipette tip manually, or by a liquid handling robot
- Applications include PCR, sample shipping, low and room temperature compound storage with DMSO and other organic solvents
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer

Key Features

- Pierceable
- Peelable from COC plates
- Seal integrity range: -20°C to 120°C
- Good solvent resistance
- Gives a weld seal to polypropylene plates
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Seals polypropylene and COC plates
- Recommended for PCR, sample shipping, low and room temperature compound storage with DMSO and other organic solvents
- Suitable for long term storage



Seal Integrity
Temperature Range

-20°C
120°C



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in two roll dimensions: 610 m x 78 mm, and 500 m x 115 mm
- Custom printing available upon request
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0538	Pierce Heat Seal Strong, strong heat sealing foil, 1 roll (610m x 78mm) ¹
4ti-0538S	Pierce Heat Seal Strong, strong heat sealing foil, 1 sample roll (5m x 78mm) ¹
4ti-0539	Pierce Heat Seal Strong, strong heat sealing foil, 1 roll (500m x 115mm) ²
4ti-0539S	Pierce Heat Seal Strong, strong heat sealing foil, 1 sample roll (5m x 115mm) ²
4ti-05381	Pierce Heat Seal Strong, strong heat sealing foil, 100 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®



AZENTA
LIFE SCIENCES

Individual Access Pierce Heat Seal Strong

Sheets of 96 foil seal spots for sealing of individual wells, suitable for storage and PCR, pierceable

- Individual Access Pierce Heat Seal Strong features 96 individual foil seal spots on a removable backing
- These seals result in individually sealed tubes that are pierceable, allowing for sample addition straight into pre-dispensed reagents, without fiddly removal of the seal
- Individual Access Pierce Heat Seal Strong is compatible with polypropylene and COC plates
- This seal is peelable from COC plates and gives a weld seal to polypropylene plates
- Demonstrating good solvent resistance, it can be used for low temperature and room temperature compound storage in DMSO and organic solvents
- The seal can be pierced with a pipette tip manually, or by a liquid handling robot
- Applications include sample shipping, storage and PCR
- This seal is available as sheets, for use with manual and semi-automated sealers, such as our Semi-Automated Sheet Heat Sealer (using the 4ti-0613 Individual Access adapter)
- The 96 well Individual Access plate can be sealed in one step resulting in individually sealed tubes that are pierceable, allowing for sample access

Key Features

- 96 individual foil seal spots on a removable backing
- 4 pin holes for exact positioning in special adapters of the Semi-Automated Sheet Heat Sealer
- Seals polypropylene (weld seal) and COC (peelable seal) plates
- Pierceable
- Seal integrity range: -20°C to 120°C
- Good solvent resistance
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Recommended for PCR, sample shipping, low and room temperature compound storage with DMSO and other organic solvents
- Best used in combination with our Individual Access 96 Well Skirted PCR Plate (4ti-0960/RA)



Seal Integrity
Temperature Range

-20°C
120°C



Options

- Sheet format: 127 mm x 100 mm

Ordering Information

4ti-05381/RA	Individual Access Pierce Heat Seal Strong, strong heat sealing foil, 96 individual seals in sheet format, 100 sheets (127 x 100mm) per case
4ti-0539/RA	Individual Access Pierce Heat Seal Strong, strong heat sealing foil, 96 individual seals in roll format, 1 roll (420m x 100mm)



AZENTA
LIFE SCIENCES

Foil Heat Seal

Aluminium heat sealing foil, resealable, peelable, pierceable; suitable for compound storage, PCR


- Our Foil Heat Seal is compatible with polypropylene and polystyrene plates
- This seal demonstrates moderate solvent resistance and can be used for low temperature compound storage in DMSO and organic solvents and short term room temperature storage
- The Foil Heat Seal can be pierced with a pipette tip, manually or by liquid handling robots, or it can be removed by peeling
- It can be resealed by applying another Foil Heat Seal directly on top of a previously pierced seal
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer
- A red stripe on the sheet foils clearly indicates the non-sealing surface, for ease of seal orientation and handling

Key Features

- Pierceable
- Peelable
- Resealable
- Seal integrity range: -20°C to 110°C
- Moderate solvent resistance
- Autoclavable (121°C)
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Seals polypropylene and polystyrene plates
- Suitable for short term storage



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0535	Foil Heat Seal, aluminium heat sealing foil, 1 roll (610m x 78mm) ¹
4ti-0535/122	Foil Heat Seal, aluminium heat sealing foil, 1 roll (122m x 78mm) ¹
4ti-0535/REMP	Foil Heat Seal, aluminium heat sealing foil, 1 roll (610m x 78mm) ³
4ti-0535S	Foil Heat Seal, aluminium heat sealing foil, 1 sample roll (5m x 78mm) ¹
4ti-0537	Foil Heat Seal, aluminium heat sealing foil, 1 roll (500m x 115mm) ²
4ti-0537/REMP	Foil Heat Seal, aluminium heat sealing foil, 1 roll (500m x 115mm) ⁴
4ti-0537S	Foil Heat Seal, aluminium heat sealing foil, 1 sample roll (5m x 115mm) ²
4ti-0536	Foil Heat Seal, aluminium heat sealing foil, 100 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®

³ Compatible with Azenta Portrait Heat Sealer

⁴ Compatible with Azenta Landscape | Stacking Heat Sealers

Polystyrene Plate Foil Heat Seal

Peelable heat sealing foil, seals to polystyrene plates, resealable, pierceable; suitable for compound storage

- Heat sealing offers a 100% effective method of plate sealing, for complete seal integrity, as well as being quick and cost effective
- Our Polystyrene Foil Heat Seal produces a stronger seal to polystyrene plates than our standard Foil Heat Seal
- Compatible with polypropylene, polystyrene and polycarbonate plates
- This seal demonstrates moderate solvent resistance and can be used for low temperature compound storage, in DMSO and organic solvents, and short term room temperature storage
- Polystyrene Foil Heat Seal can be pierced with a pipette tip manually, by a liquid handling robot, using the Azenta Pierce Plate (4ti-0398), or it can be removed by peeling. It can be resealed by applying another Polystyrene Foil Heat Seal directly on top of a previously pierced seal
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer

Key Features

- Pierceable
- Peelable
- Resealable foil on foil
- Seal integrity range: -20°C to 110°C
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Seals polypropylene, polystyrene and polycarbonate plates
- Suitable for short term storage
- Suitable for low temperature compound storage



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0545	Polystyrene Foil Heat Seal, peelable heat sealing foil, 1 roll (610m x 78mm) ¹
4ti-0545/122	Polystyrene Foil Heat Seal, peelable heat sealing foil, 1 roll (122m x 78mm) ¹
4ti-0545S	Polystyrene Foil Heat Seal, peelable heat sealing foil, 1 sample roll (5m x 78mm)
4ti-0546	Polystyrene Foil Heat Seal, peelable heat sealing foil, 1 roll (500m x 115mm) ²
4ti-0547	Polystyrene Foil Heat Seal, peelable heat sealing foil, 100 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®



AZENTA
LIFE SCIENCES

Thermal Bond Heat Seal

Heavy duty heat sealing foil, peelable; suitable for long term storage, transportation

- Heat sealing offers a 100% effective method of plate sealing, for complete seal integrity, as well as being quick and cost effective
- Our Thermal Bond Heat Seal is a heavy duty laminate foil seal suitable for providing a very strong, but peelable seal
- Compatible with polypropylene plates to provide a high degree of sample protection
- Demonstrates very good solvent resistance and can be used for very low temperature compound storage, in DMSO and organic solvents, and long term room temperature storage such that it is recommended as suitable for sample transportation
- The seal can be pierced only by using a blade or using the Azenta Pierce Plate (4ti-0398)
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer

Key Features

- Peelable
- Seal integrity range: -200°C to 110°C
- High solvent resistance
- Autoclavable (121°C)
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen


Use

- Recommended for PCR, including with water bath thermal cyclers
- Seals polypropylene plates
- Suitable for long term storage and transportation
- Suitable for very low temperature storage



Seal Integrity
Temperature Range

-200°C
110°C



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0590	Thermal Bond Heat Seal, heavy duty heat sealing foil, 1 roll (500m x 78mm) ¹
4ti-0590/100	Thermal Bond Heat Seal, heavy duty heat sealing foil, 1 roll (100m x 78mm) ¹
4ti-0590S	Thermal Bond Heat Seal, heavy duty heat sealing foil, 1 sample roll (5m x 78mm) ¹
4ti-0592	Thermal Bond Heat Seal, heavy duty heat sealing foil, 1 roll (300m x 115mm) ²
4ti-0592S	Thermal Bond Heat Seal, heavy duty heat sealing foil, 1 sample roll (5m x 115mm) ²
4ti-0591	Thermal Bond Heat Seal, heavy duty heat sealing foil, 100 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Not compatible with Agilent (Velocity 11) PlateLoc®

Gas Permeable Heat Seal

Heat sealing membrane, limits evaporation, peelable, pierceable; suitable for cell culture, seed and insect storage

- Heat sealing offers a 100% effective method for plate sealing for a complete seal integrity, as well as being quick and cost effective
- Our Gas Permeable Heat Seal is made from a woven material and is designed for use in cell culture, due to its porous nature
- The small pore size (<20 µm) of this material enables gas exchange, per 24 hours of >20 g/m², whilst evaporation is reduced to a minimum
- Compatible with polypropylene, polystyrene and cyclic olefin copolymer (COC) plates
- It can be removed by peeling, or it can be pierced with a pipette tip manually, using a liquid handling robot or with our Pierce Plate (4ti-0398)
- Gas Permeable Heat Seal can be utilized for effective overnight incubations, during which it demonstrates significant reductions in evaporation compared to lids
- It can also be used for insect and seed storage as it enables gas exchange, whilst providing an inert surface with no adhesive to interfere with the well contents
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer

Key Features


- Peelable
- Pierceable
- Seal integrity range: -20°C to 80°C
- Gas permeability rate: 180 m³/m²/day
- Moisture vapor transmission rate: 20 g/m²/day
- Autoclavable (121°C)
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Seals polypropylene, polystyrene and cyclic olefin copolymer (COC) plates
- Suitable for cell culture, overnight incubations, as well as insect and seed storage



Seal Integrity
Temperature Range
-20°C
80°C



Options

- Sheet format: 125 x 78 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: available in a variety of roll dimensions, to suit your choice of automatic heat sealing equipment
- Available gamma treated from stock: gamma treated variations have /ST added to the end of their product code

Ordering Information

4ti-0598	Gas Permeable Heat Seal, heat sealing membrane, 1 roll (610m x 78mm) ¹
4ti-0598/122	Gas Permeable Heat Seal, heat sealing membrane, 1 roll (122m x 78mm) ¹
4ti-0598S	Gas Permeable Heat Seal, heat sealing membrane, 1 sample roll (5m x 78mm) ¹
4ti-0599	Gas Permeable Heat Seal, heat sealing membrane, 1 roll (500m x 115mm) ²
4ti-0599S	Gas Permeable Heat Seal, heat sealing membrane, 1 sample roll (5m x 115mm) ²
4ti-0597	Gas Permeable Heat Seal, heat sealing membrane, 100 sheets (125 x 78mm) per case
4ti-0597/ST	Gas Permeable Heat Seal, Gamma treated, heat sealing membrane, 10 x 10 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube

² Compatible with Agilent (Velocity 11) PlateLoc®



AZENTA
LIFE SCIENCES

Clear Heat Seal, Peelable Film

Peelable heat sealing film, optically clear, with 3mm slits for gas transfer; suitable for insect studies, seed storage


- Heat sealing is a quick and cost effective method of plate sealing
- Our Clear Heat Seal, Peelable Film is based on our Clear Heat Seal, with the addition of 3mm slits across the width of the seal
- These slits render the seal gas permeable, whilst retaining evaporation to a minimum, compared to the use of lids
- The Clear Heat Seal, Peelable Film is compatible with polypropylene, polyethylene and polystyrene plates
- The seal can be removed by peeling, or it can be pierced with a pipette tip manually, using a liquid handling robot or with our Pierce Plate (4ti-0398)
- The Clear Heat Seal, Peelable Film has a wider seal integrity temperature range, from -80°C to 110°C, than our Gas Permeable Heat Seal
- It can be used for insect and seed storage, as it enables gas exchange, whilst providing an inert surface with no adhesive to interfere with the well contents
- This seal is available as sheets, for use with manual and semi-automated sealers such as our Semi-Automated Sheet Heat Sealer
- Also available in multiple roll formats compatible with specified automated heat sealers, such as our Automated Roll Heat Sealer

Key Features

- Peelable
- Seal integrity range: -80°C to 110°C
- 3mm slits for gas transfer
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Seals polypropylene, polyethylene and polystyrene plates
- Suitable for insect and seed storage



Options

- Sheet format: 125 x 78mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates, from 12 well to 1536 well
- Roll format: 500m x 78mm; approx. 4,200 seals
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0540/SLIT	Clear Heat Seal, gas permeable, clear heat sealing film, with 3mm slits for gas transfer, 1 roll (450m x 78mm) ¹
4ti-0540/SLIT/S	Clear Heat Seal, gas permeable, clear heat sealing film, with 3mm slits for gas transfer, 1 sample roll (5m x 78mm)
4ti-0541/SLIT	Clear Heat Seal, gas permeable, clear heat sealing film, with 3mm slits for gas transfer, 100 sheets (125 x 78mm) per case

¹ Compatible with Azenta Semi-Automated Sheet Heat Sealer / Thermo Fisher ALPS 300™ and ALPS 3000™ / KBiosystems Wasp™ / KBioscience FlexiSeal and Cube



Thermal Individual Tube Seal, 96-Format

Thermal Individual Tube Seal was designed for the heat sealing of 96 well microplates or 96 format cluster tubes, preserving sample integrity and audit trail in compound library and other high-throughput applications

- 96 individual sealing spots on a backing liner
- Creates an air-tight seal
- Peelable and pierceable
- “Embossed” seal construction
- Available with or without tabs for easy peeling
- For use with heat sealers
- Heat sealing offers a 100% effective method of microplate and tube sealing, for complete seal integrity, as well as being quick and cost effective
- Thermal Individual Tube Seal, 96-Format was designed for the heat sealing of 96 well microplates or 96 format cluster tubes, preserving sample integrity and audit trail in compound library and other high-throughput applications
- This is a foil-based heat seal which consists of 96 individual round seals held on a convenient to handle sealing sheet
- The seal is applied using a manual or semi-automated heat sealer, such as our Semi-Automated Sheet Heat Sealer
- After sealing, the microplate frame or the rack can be removed to leave 96 individually sealed wells or tubes
- The individual seals can be removed as required, by hand or using forceps, via the seal removal tab; no polymer residue is left on the well/tube following removal
- The individual seals can also be pierced with a pipette tip manually, using a liquid handling robot

Key Features

- Peelable
- Pierceable
- Seals polypropylene
- 96 individual foil-based seals held on an easy to handle sealing sheet
- Good solvent resistance, including DMSO
- Free from DNase, RNase, and human genomic DNA, endotoxin/pyrogen free



Use

- Air-tight sealing that works as an impenetrable barrier for added sample security
- Compatible with 96 well microplates and 96 format racked tubes
- For use with manual or semi-automated heat sealers, such as our Semi-Automated Sheet Heat Sealer
- Tabs for easy removal: hold the tube and peel the tab upwards with forceps or fingers
- Engineered to come away cleanly, leaving no residue, for easy resealing
- Suitable for long term storage and sample shipping
- Suitable for high temperature incubations
- Suitable for cryogenic storage*

Options

- Sheet format: individual spots with tabs for easy peeling
- Sheet format: individual spots without tabs

Ordering Information

66-1001	Thermal Seal, Individual, with tab for manual peeling, 96 Tube Seals per sheet, 50 sheets per case
66-1021	Thermal Seal, Individual, without tab, 96 Tube Seals per sheet, 50 sheets per case

**Not for use in liquid phase nitrogen*

Heat Sealing Consumables Comparison & Instrument Compatibility Table

	Clear Heat Seals				Peelable Heat Seals		
Name	Clear Heat Seal	Clear Weld Heat Seal	Clear Heat Seal Easily Pierceable	Clear Heat Seal Plus	Peel Heat Seal	Universal Peel Heat Seal	DMSO Resistant Peel Heat Seal
Specifications							
Application	qPCR Short term compound storage	PCR, esp. water bath cyclers qPCR Storage & disposal of hazardous materials	qPCR and for use with ABI 3730 Sequencer	Homogenisation or disruption of seeds or other material, e.g. bead mill applications	Low temperature compound storage Short term room temperature compound storage (<5 days) PCR	Low temperature compound storage High temperature applications PCR	Low/room temperature compound storage with DMSO & other organic solvents
Special Properties	Good optical clarity Moderate solvent resistance	Good optical clarity Resistance to DMSO	Good optical clarity Some solvent resistance	High tensile strength	Can be peeled directly from -80°C freezer Moderate resistance to solvents at room temperature	Moderate solvent resistance Re-sealable with another Universal Peel Heat Seal	Can be peeled directly from -80°C freezer High resistance to solvents even at elevated temperatures
Seal Integrity Min Temperature	-80°C	-80°C	-20°C	-80°C	-80°C	-80°C	-80°C
Seal Integrity Max Temperature	80°C (or 110°C with pressurized heated PCR lid)	110°C	80°C (or 110°C with pressurized heated PCR lid)	80°C (or 110°C with pressurized heated PCR lid)	90°C (or 110°C with pressurized heated PCR lid)	90°C (or 110°C with pressurized heated PCR lid)	40°C
Pierceable			✓			(✓)	
Peelable	✓			✓	✓	✓	✓
RNase/DNase free	✓	✓	✓	✓	✓	✓	✓
Material	Laminate	Polymer	Polymer	Polymer	Laminate	Laminate	Laminate
Seals to	PP, PE, PS, PC, COC	PP	PP, PE, PS, COC	PP, PS, COC	PP, COC	PP, PE, PS, PC, COC	PP, PE, COC
Sealing parameters with 96 Well PP Plates	175-185°C 2-3 s	175-185°C 2-3 s	165-175°C 3 s	175-185°C 2-3 s	175-185°C 3 s	175°C 2 s	175-185°C 3 s
Sealing parameters with 384 Well PP Plates	165-180°C 3 s	170-175°C 2-3 s	165-175°C 2 s	165-180°C 3 s	170-175°C 2-3 s	175°C 2 s	170-175°C 2-3 s
Sealing parameters with Optically Clear Tissue Culture Plates	185-200°C 3 s	N/A	175-185°C 2-3 s	185-200°C 3 s	N/A	180°C 2 s	N/A
Product Codes/Instrument Compatibility							
Compatible with Azenta Automated Roll Heat Sealer Thermo Fisher ALPS 300™ and ALPS 3000™ KBiosystems Wasp™ and Chameleon™ KBioscience FlexiSeal and Cube							
Roll, 78 mm width	4ti-0540	4ti-0573	4ti-0580	4ti-0549	4ti-0520	4ti-0523	4ti-0585
Roll, 78 mm width, short roll*	4ti-0540/80	4ti-0573/122	4ti-0580/122		4ti-0520/122		4ti-0585/100
Sample roll, 78 mm width	4ti-0540S	4ti-0573S	4ti-0580S	4ti-0549/S	4ti-0520S	4ti-0523S	4ti-0585S
Compatible with Agilent (Velocity 11) PlateLoc®							
Roll, 115 mm width	4ti-0542	4ti-0574	4ti-0582	4ti-0548	4ti-0522	4ti-0524	4ti-0586
Sample roll, 115 mm width	4ti-0542S	4ti-0574S	4ti-0582S	4ti-0548/S	4ti-0522S	4ti-0524S	4ti-0586S
Compatible with Azenta Portrait Heat Sealer							
Roll, 78 mm width, large core	4ti-0540/REMP						4ti-0585/REMP
Compatible with Azenta Landscape Stacking Heat Sealers							
Roll, 115 mm width, large core	4ti-0542/REMP						4ti-0586/REMP
Compatible with Azenta Semi-Automated Sheet Heat Sealer Thermo Fisher ALPS™ 25 and ALPS™ 50 KBiosystems E-Fly 2 Azenta Easy Sealer							
Sheets	4ti-0541	4ti-0575	4ti-0581	4ti-05481	4ti-0521	4ti-05231	4ti-0587
Compatible with Azenta Semi-Automated Sheet Heat Sealer							
Individual Access, sheets					4ti-0521/RA-TAB 4ti-0521/RA-8		
Compatible with Azenta Automated Individual Access Heat Sealer							
Individual Access, roll, 100 mm width					4ti-0522/RA-TAB		

* For use with the Azenta Automated Roll Heat Sealer when using lower roll position and the optional dust cover for protection of the roll - ** NOT compatible with Agilent (Velocity 11) PlateLoc®

Azenta US, Inc. recognizes that designated trademarks and brands are the property of their respective owners.



AZENTA
LIFE SCIENCES

Pierceable Heat Seals		Foil Heat Seals			Gas Permeable Heat Seals	
Pierce Heat Seal	Pierce Heat Seal Strong	Foil Heat Seal	Polystyrene Foil Heat Seal	Thermal Bond Heat Seal	Gas Permeable Heat Seal	Gas Permeable Clear Heat Seal
PCR Compound storage Sample shipping	PCR Compound storage Sample shipping	Low temperature compound storage Short-term room temperature compound storage PCR	Low temperature compound storage Short-term room temperature compound storage PCR	Low temperature transportation & storage PCR, esp. water bath cyclers Storage of organic solvents, acids & alkalines	Cell culture Over night incubation Seed and insect storage	Storage e.g. for seeds or insects
Easily pierceable Resistant to DMSO Re-sealable with another Pierce Heat Seal Color print identifies non-sealing surface	Easily pierceable Resistant to DMSO Re-sealable with another Pierce Heat Seal Color print identifies non-sealing surface	Re-sealable with another Foil Seal Resistant to DMSO Color print identifies non-sealing surface	Re-sealable with another Foil Seal Resistant to DMSO	Very strong seal with PP Resistant to DMSO and other solvents	Small pore size of 20 µm allows gaseous exchange & limits evaporation Gas permeability: 180 m³/m²/day Moisture vapor transmission: 20g/m²/day	3 mm slits across entire surface of seal makes this permeable to gases
-20°C	-20°C	-20°C	-20°C	-200°C	-20°C	-80°C
120°C	120°C	110°C	110°C	110°C	80°C	100°C material integrity (not seal)
✓	✓	✓	✓		✓	✓
		✓	✓	✓	✓	✓
✓	✓	✓	✓	✓	✓	✓
Foil	Foil	Foil	Foil	Laminate	Woven material	Laminate
PP, PS	PP, COC	PP, PS	PP, PS, PC	PP	PP, PS, COC	PP, PE, PS
160-175°C 2 s	170-180°C 2 s	165-180°C 2 s	165-180°C 2 s	170-180°C 2-3 s	170°C 2 s	175-185°C 2-3 s
160-175°C 2 s	170-180°C 2 s	165-175°C 2-3 s	165-175°C 2-3 s	160-170°C 2 s	170°C 2 s	165-180°C 3 s
185-200°C 3 s	180-200°C 3 s	185-200°C 3 s	185-200°C 3 s	N/A	170°C 2 s	185-200°C 3 s
4ti-0530	4ti-0538	4ti-0535	4ti-0545	4ti-0590	4ti-0598	4ti-0540/SLIT
4ti-0530/122		4ti-0535/122	4ti-0545/122	4ti-0590/100	4ti-0598/122	
4ti-0530S	4ti-0538S	4ti-0535S	4ti-0545S	4ti-0590S	4ti-0598S	4ti-0540/SLIT/S
4ti-0532	4ti-0539	4ti-0537	4ti-0546	4ti-0592**	4ti-0599**	
4ti-0532S	4ti-0539S	4ti-0537S	4ti-0546/S	4ti-0592S**	4ti-0599S**	
4ti-0530/REMP		4ti-0535/REMP				
4ti-0532/REMP		4ti-0537/REMP				
4ti-0531	4ti-05381	4ti-0536	4ti-0547	4ti-0591	4ti-0597	4ti-0541/SLIT
4ti-0531/RA	4ti-05381/RA					
4ti-0532/RA	4ti-0539/RA					

Azenta recognizes that designated trademarks and brands of the Instrument Compatibility Table are the property of their respective owners.




Thermosensitive Color Forming Film

Thermosensitive color-forming film; for evaluation of consistent temperature across a heating block; to be used in conjunction with roll fed or sheet fed heat sealing instruments

- Azenta Thermal Test Film (TTF) can be used in conjunction with a roll fed (e.g. Azenta Automated Roll Heat Sealer) or sheet fed (e.g. Azenta Semi-Automated Sheet Heat Sealer)
- It checks the uniformity and reproducibility of the heat sealing block of the instrument
- The film can be used to effectively test the temperature of the heating block between 160°C and 200°C

Key Features

- The Thermal Test Film has a thermosensitive color-forming layer plus a protective layer, both attached to the base material
- Depending on the temperature applied to the film, a color is produced in varying density and hue, giving a perfect image of heat distribution across the heating block of your heat sealer
- The color varies according to dwell time and temperature
- The shorter the duration, the paler and more blueish the color
- The longer the duration, the more saturated and reddish the color



Specifications

Parameter	Value
Sealing temperature range	160 °C to 200 °C
Recommended ambient temperature	15 °C to 30 °C
Recommended ambient humidity	35% RH to 80% RH


Ordering Information

4ti-0640	Thermosensitive Color Forming Film, used to check uniformity and reproducibility of a heat sealing block, 25 sheets (125 x 80 mm)
4ti-0641	Thermosensitive Color Forming Film, starter kit, 10 sheets (125 x 80 mm) and 1 silicone pad
4ti-0642	Thermosensitive Color Forming Film, roll kit, 1 roll (1m x 80 mm) and 1 silicone pad

Thermal Test Film Color Chart

Temperature/Duration	150°C	160°C	170°C	180°C	190°C	200°C	210°C
1 second							
10 seconds							
60 seconds							

Note: This color chart is just an example based on the results of the tests performed in our laboratories. Before using the Thermal Test Film, a similar chart needs to be created that is based on your actual measurement conditions



Automated Roll Heat Sealer

The Automated Roll Heat Sealer enables up to 5,000 perfect seals without manual user intervention, providing a true walk-away system. The main advantage of the Automated Roll Heat Sealer is that it is powered by electric motors, with no requirement for any air supply. This enables the Automated Roll Heat Sealer to generate a reliable and consistent sealing pressure, resulting in superior seal uniformity. In addition, it provides the user the flexibility to use this instrument without the need for an external air supply, enabling ease-of-use as a stand-alone unit on a lab bench or within integrated robotic set-ups. Consistent sealing is achieved through a fixed high sealing pressure and accurate time and temperature controls, ensuring reproducible seal uniformity on the widest range of plates. The high-performance heating block design enables rapid heating with no delay between seals, and provides a uniform temperature across the entire heating block. The Automated Roll Heat Sealer is compatible with a wide range of SBS footprint plates, including all PCR plate formats from 96 to 1536 wells, assay plates, deep well storage plates and microplates. The Automated Roll Heat Sealer can be used with the wide range of sealing materials that Azenta offers, including gas permeable seals with no further instrument modification, enabling a wide range of applications.

- Unrivalled sealing performance and consistency; 5,000 perfect seals without intervention
- SiLA compatible; Easy integration into robotic systems
- Powerful electric motor generates consistent sealing pressure for superior sealing uniformity
- High-performance heating block design generates uniform temperature across entire heating block

Features

- Powerful electric motor drives plate sealing mechanism
- SiLA compatible
- Compatible with a wide range of plates and seals; 2 positions for different roll sizes
- Variable time and temperature controls
- Seal cycle time < 15 seconds
- High-performance heating block design
- Color touch screen with intuitive user interface
- Unlimited password protected protocols
- Auto standby mode
- Optional roll cover
- 2 year warranty



Color touch screen with intuitive user interface – ease of use



Benefits

- Reliable and consistent sealing pressure; Superior sealing uniformity; No requirement for any air supply
- Easy integration into robotic systems through plug & play
- Flexibility with a wide range of consumables including gas permeable seals without instrument modification
- Enables optimization of any seal / plate type
- Time-saving
- Rapid heating, no delay between seals; Uniform temperature across entire heating block, +/- 1 degree edge-to-edge, corner-to-corner
- Ease of use, saving time
- Save personalized and SOP-set temperatures and times
- Energy-saving; prolongs component life
- Seal protection for sensitive applications
- Instrument reliability



AZENTA
LIFE SCIENCES

Automated Roll Heat Sealer

Providing The Perfect Solution For Automated Heat Sealing

Plug and Play Robotic Integration

- The automated Roll Heat Sealer is compliant with SiLA standards for rapid integration of automated systems (www.sila-standard.org). This means it can be “plug and play” connected with other instruments, such as readers, robotic arms and liquid handlers, to give a custom automated system without the need to write expensive custom drivers.
- The Automated Roll Heat Sealer is extremely versatile leaving you the freedom to expand and reconfigure your systems such as adding robotic plate handling.
- Full communication protocol available — operate the instrument and record sealing conditions for each seal using the RS232 communication port
- SiLA compatible — quick and easy integration with other SiLA compliant devices using the SiLA driver




Technical Specifications

Parameter	Value
Dimensions (W x L x H)	230 x 507 x 276 mm
<i>Please note: Additional space is required if large seal rolls are used</i>	
Sealing Temperature Range	100-195°C
Sealing Time Range	0.1-10 sec
Weight (without roll)	27 kg
Power Supply	V in: AC 100-240 V V out: DC 24 V 320 W
Power Consumption	700 W (max)
Working Temperature Range	10-30°C
Operating Humidity (RH)	0-85%
Connection	RS-232 serial port, USB port



*Flexibility of seal material
choice - Azenta offers a wide
range of sealing films and foils*



*Compatibility with a wide
range of plate types - two
adapters supplied for optimal
sealing results, with custom
adapters available on request*

Ordering Information

Automated Roll Heat Sealer

Includes: Power cord, manual, plate support adapters A (4ti-0665-2) and B (4ti-0665-3), 24 months parts and labor warranty

4ti-0665	Automated Roll Heat Sealer, for use with adapters A and B
Accessories	
4ti-0665-1	Seal Loading Tool; 1 tool
4ti-0665-2	Plate Support Adapter A, for mid height plates; 1 adapter
4ti-0665-3	Plate Support Adapter B, for 96/384 well PCR plates, 1536 well plates and standard height plates; 1 adapter
4ti-0665-56	Elevator Adapter, for very low skirted microplate; 1 adapter
4ti-0665-57	Lightweight Adapter, for low skirted microplates; 1 adapter
4ti-0665-58	Lightweight Adapter, for high skirted microplates; 1 adapter
4ti-0665-59	Lightweight Adapter, for 384 well microplates; 1 adapter
4ti-0665-60	Lightweight Adapter, for 96 well microplates; 1 adapter
4ti-0665-4	Roll Holder Set; includes spindle, clamping wheels, locking nut; 1 set
4ti-0665-5	Vacuum Cups, Front; clear, set of 2; 1 set
4ti-0665-6	Vacuum Cups, Back; black, set of 2; 1 set
4ti-0665-7	Plastic Tweezers; 1 tweezers
4ti-0665-8	Clear Plastic Roll Dust Cover; 1 cover
4ti-0665-41	SiLA driver; 1 license



AZENTA
LIFE SCIENCES

Automated Individual Access Roll Heat Sealer

The Automated Individual Access Heat Sealer is an automated roll heat sealer for higher throughput capable of sealing individual wells or tubes, enabling researchers to leverage the benefits of the Azenta Individual Access range whilst maintaining the gold standard heat sealing provided by the Azenta Automated Roll Heat Sealer.

Concept

Individual access utilizes a plate with individually removable wells, together with seals consisting of individual foil seal spots. This enables sealing of individually accessible tubes and thereby provides flexibility for single access or placement of tubes within a rack. In addition, the Automated Individual Access Heat Sealer also has the ability to seal custom shaped consumables with custom shaped seals to accommodate tailor-made needs. The resulting individual access of tubes and consumables through sealing of individual tubes and custom shaped consumables enables high throughput manufacturing but shipment of individual product. End-users can then utilize tubes individually as needed by taking one well at a time.


Features

- Individual sealing of plate wells or tubes
- Color touch screen with intuitive user interface
- Variable time and temperature controls
- Compatibility with wide range of plates and seals
- Ability to seal custom shapes with custom seals
- Unlimited password protected protocols
- Rapid heating
- Auto standby mode
- 2 year warranty



A solution for:

- Diagnostic Kit Manufacturers
- Oligonucleotide Production
- Synthetic Biology
- Sample Storage
- Individual Access Users
- Antibody Extraction



Benefits

- Enables individual access to plate wells or tubes and allows high throughput manufacturing but shipment of individual product
- Ease of use, saving time
- Enables optimization of any seal and plate type
- Enables specific requirements and unique applications
- Custom potential for tailor-made consumables
- Save personalized and SOP-set temperatures and times
- Easy integration into robotic systems through plug & play
- Fast start-up time; block uniformity maintained to $\pm 1^{\circ}\text{C}$
- Energy saving; prolongs component life
- Instrument reliability

Ordering Information

Automated Individual Access Heat Sealer Includes: Power cord, manual, plate support adapters, 24 months parts and warranty.	
59-1000	Automated Individual Access Heat Sealer
Accessories	
59-1001	Automated Individual Access Heat Sealer Seal Loading Tool
59-1004	Automated Individual Access Heat Sealer Plate Support Adapter, A
59-1003	Automated Individual Access Heat Sealer Waste Collection Core
59-1002	Automated Individual Access Heat Sealer Spindle Support

Semi-Automated Sheet Heat Sealer



Semi-Automated Sheet Heat Sealer

The Semi-Automated Sheet Heat Sealer is compatible with a wide range of seals and plates of differing designs and heights.

With variable temperature and time settings, sealing conditions are easily optimized to produce a tight seal, eliminating sample loss.

Plate and seal are placed on the holder, the “*Operate*” button pressed and the drawer automatically closes. The sealing process is controlled by an electric mechanism.

Heat sealers provide a mechanism for controlled plate sealing, eliminating variation and giving consistent sealing every time.

Features

- Variable time and temperature control
- Simple, three button operation
- Real-time temperature display
- Auto stand-by function and switch-off mode
- Compatible with all SBS microplates
- Small footprint
- RS-232 serial port
- 12 month warranty

Benefits

- Enables optimization of any seal / plate type
- Ease of use
- Rapid heating element enables fast start-up, saving time
- Conserves energy by reducing temperature of heating element down to 60°C or switching off after defined time
- Flexibility for wide range of consumables, eg competitor plates
- Minimal bench space
- RS-232 port enables unit integration within robotic systems
- Instrument reliability



Technical Specifications

Parameter	Value
Dimensions (W x L x H)	181 x 275 x 310 mm
Sealing Temperature Range	60-200°C
Sealing Time Range	0-10 sec
Weight	9 kg
Power Supply	110/230V
Power Consumption	300 W (max)
Connection	RS-232 serial port

The semi-automated sheet heat sealer is compatible with a wide range of seals & plates of differing designs & heights




Reproducible
Sealing

Adapters

The Semi-Automated Sheet Heat Sealer employs a unique and cost saving adapter system

- Deep well plates and most other skirted plates can be sealed on the standard plate support adapter (59-2001)
- For all other 96 well PCR plates, an additional plate support adapter (59-2003, optional extra) must be used
- Both adapters can also be combined for easier handling of shallow plates such as 384 well PCR plates.
- For best sealing results with Individual Access seals, 8 different adapters are available to offer the perfect plate support and seal alignment



Highest
Flexibility

Sealing Aids

- For better sealing of films which have a tendency to curl, Azenta offers sealing aids
- A sealing frame for use with all other plate designs is supplied when ordering the optional 59-2003 plate support adapter

Plate Formats

The following plate formats can be used with the appropriate plate adapter:

- Standard SBS footprint PP & PS plates & deepwell blocks
- PCR plates (skirted, semi-skirted and non-skirted formats)

Ordering Information

Semi-Automated Sheet Heat Sealer

Includes: Power cord, manual, standard plate support adapter (59-2001), 12 months parts and labour warranty

59-2000

Semi-Automated Sheet Heat Sealer,
includes adapter (59-2001)

Accessories

59-2001

Semi-Automated Sheet Heat Sealer Adapter, for skirted 96 and 384 well plates, 1 adapter per case

59-2002

Semi-Automated Sheet Heat Sealer Adapter, for Roche 1536 well PCR plates, 1 adapter per case

59-2003

Semi-Automated Sheet Heat Sealer Adapter, for 96 well PCR plates, includes Semi-Automated Sheet Heat Sealer Sealing Frame (59-2009) 1 adapter and frame per case

59-2004

Semi-Automated Sheet Heat Sealer Adapter, for 384 well PCR plates, 1 adapter per case

59-2005

Semi-Automated Sheet Heat Sealer Adapter, for Individual Access plates, 1 adapter per case

59-2006

Semi-Automated Sheet Heat Sealer Adapter, for 96 and 384 well storage plates, 1 adapter per case

59-2007

Semi-Automated Sheet Heat Sealer Adapter, for 4TI-LB0109 and 96 well PCR plates, 1 adapter per case

59-2008

Semi-Automated Sheet Heat Sealer Weighted Sealing Platen, 1 platen per case

59-2009

Semi-Automated Sheet Heat Sealer Sealing Frame, for use with Semi-Automated Sheet Heat Sealer Adapter, for use with 59-2003, 1 frame per case



59-2001
Plate Support
Adapter, Standard



59-2003
Plate Support
Adapter, PCR 96



59-2005
Plate Support
Adapter, Individual
Access, Low Profile



59-2008
Weighted Sealing
Platen



59-2009
Sealing Frame



AZENTA
LIFE SCIENCES

Automated Plate Seal Remover

Automatically removes seals from a wide range of microplate types with the single touch of a button.

A robust and elegantly-simple automated system, the Automated Plate Seal Remover eliminates the need for repetitive, manual removal of plate seals and enables the adoption of the gold-standard operating model (sealed plates, no lids).

The Automated Plate Seal Remover Tape eliminates the need for mechanical removal mechanisms which are often prone to failure. Automated Plate Seal Remover is highly reliable and can be used manually or integrated into automated systems with external robotics.

Key Features

Compatible with Virtually All Plate Types and Seal Types


- Can be used with a wide variety of microplates, including full-skirted PCR plates, low-base microplates and deep-well (up to 2ml) plates
- Compatible with a variety of full-plate seals, including heat and pressure applied seals
- Uses proprietary Automated Plate Seal Remover Tape to de-seal microplates, up to 400 seal removals per roll

Preserves Sample Integrity

- Eliminates cross contamination common with manual seal removal techniques
- Supports Quality Control measures requiring samples to be sealed until their moment of use
- Holds the plate down whilst the seal is peeled away from the plate, eliminating another contamination issue
- Operating mode minimizes plate or seal damage
- When used as part of an integrated workflow, seal removal verification feature reduces plate handling errors

Easy to Use, Easy to Integrate

- Can be used as a standalone system, or integrated into automated and robotic workflows via Serial RS232 remote interface
- One-touch, push-button operation to de-seal plates makes the Automated Plate Seal Remover an ideal standalone device for busy laboratories
- Capacity to remove up to 200 plate seals per hour
- Robust, time-proven device with hundreds of units placed globally in a range of manual and automated environments



Specifications

Parameter	Value
Seal Removal Capacity	Up to 400 per Tape Roll
Seal Verification Sensor	Reflective with Sensitivity Adjustment
Communication	Serial RS232
Motion Parameters	Tape Adhere Time, De-seal Speed, Plate Output Orientation, Begin Peel Location
General Parameters	Auto Tape Advance, Plate Verification Menu Language
Weight	35kg (76lbs)
Power Requirements	115VAC, 4A, 60Hz 230VAC, 2A, 50Hz
Throughput	Up to 200 Plate Seal Removals per Hour

Ordering Information

Automated Plate Seal Remover Includes: Power cord, manual, plate support adapters, 12 months parts and warranty.	
XP-A_100V	Automated Plate Seal Remover, 100V, version for Asia; 1 unit
XP-A	Automated Plate Seal Remover, 115V, version for North America; 1 unit
XP-A_230V	Automated Plate Seal Remover, 230V, version for Europe; 1 unit
Compatible Tapes	
X-Tape_2000	Automated Plate Seal Remover Tape, 5 spools per case



AZENTA
LIFE SCIENCES

Adhesive Sealing Consumables & Sealing Accessories



Adhesive Sealing Consumables & Sealing Accessories


Adhesive seals provide a versatile, yet robust, method for protecting samples. The seals are supplied as sheets and some are also available in roll format. Most adhesive seals are supplied with convenient tabs on both ends for easy application and removal. These tabs also enable easy peeling to remove a seal without leaving adhesive residue on the sealing surface.

All our seals are produced and processed under strictly controlled environmental conditions and certified free from DNase, RNase and human genomic DNA.

To choose the most suitable seal, please refer to the comparison table on page 232 which describes the recommended applications and technical features of each seal.

To obtain the best sealing results with adhesive seals, we strongly recommend the use of high quality plates with raised rim sealing rings for optimal sealing integrity and guaranteed flatness. All Azenta PCR plates are designed with these features.

To improve seal application by ensuring even pressure is applied, we offer a seal roller and a seal applicator, for complete and secure application of all our adhesive seals. We also offer supplementary products like the compression pad, for details see page 230.



4titude PCR Seal

Clear adhesive film, strong adhesive, peelable; suitable for PCR and optical applications; available as a full sheet, perforated for tearing into 8 well strips or 12 well strips


- Our PCR Seal is a durable transparent polyester film with a strong adhesive layer
- The seal enables a high integrity and efficiently prevents sample evaporation
- Recommended for PCR, qPCR, and other optical applications (e.g. fluorescence or colorimetric measurements) as well as sample storage
- The PCR Seal is also available in two flexible formats with perforated sheets, to enable tearing into 8 well and 12 well strips, respectively
- Allows for sealing of complete 96 well plates, but also individual or multiple Breakable Horizontally or Vertically or Breakable Vertically strips, perfectly complementing these products
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features

- Non-pierceable, please refer to our PCR Foil Seal Strong for a pierceable variant
- Peelable
- Seal integrity range: -20°C to 110°C
- Non-gamma treated
- Free from DNase, RNase, and human genomic DNA

Use

- Applications: PCR, qPCR
- Removal: features adhesive-free end tabs for easy removal, will not leave a sticky residue on the plate surface following removal



Seal Integrity
Temperature Range

-20°C

110°C

*110°C with pressurized heated lid

Peelable



Options

- Sheet format: 135 x 80 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates
- Perforated sheet format: 115 x 100 mm, for tearing into 8 well strips
- Perforated sheet format: 137 x 71 mm, for tearing into 12 well strips

Ordering Information

4ti-0500	PCR Seal, clear adhesive film, strong adhesive, 100 sheets (135 x 80mm) per case
4ti-0500/8	PCR Seal 96/8, clear adhesive film, vertically perforated, strong adhesive, 100 sheets (115 x 100mm) per case
4ti-0500/12	PCR Seal 96/12, clear adhesive film, horizontally perforated, strong adhesive, 100 sheets (137 x 71mm) per case



AZENTA
LIFE SCIENCES

4titude Optically Clear Windowed qPCR Seal

Adhesive film with 96 optically clear windows, peelable, suitable for qPCR and optical applications


- These unique seals combine the strong sealing integrity of our PCR Seal with improved optical properties, thanks to the 96 adhesive-free windows
- The seal is made of a durably transparent polyester film, and a strong adhesive is applied across the seal, apart from the 96 round windows
- The Optically Clear Windowed qPCR Seal is recommended for qPCR and other optical applications, such as fluorescence or colorimetric measurements
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features

- Non-pierceable
- Peelable
- Seal integrity range: -20°C to 110°C
- Free from DNase, RNase, and human genomic DNA

Use

- Applications: qPCR, plate readers
- For 96-well microplates only
- Removal: will not leave a sticky residue on the plate surface following removal
- Suitable for cutting to fit part plates



Seal Integrity
Temperature Range

-200°C

110°C

*110°C with pressurized heated lid

Peelable



Options

- Sheet format: 133 x 76 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates with 96 round wells
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0565*

Optically Clear Windowed qPCR Seal,
adhesive film with 96 optically clear
windows, 100 sheets (133 x 76mm)
per case

**Not available for purchase in or onwards distribution to the USA*

qPCR Adhesive Seal

Optically clear adhesive film, pressure activated adhesive, peelable; suitable for qPCR and other imaging techniques including crystallization

- Optically clear seal specifically developed for optical applications, particularly qPCR
- It is non sticky when removed from the packaging; this aids handling when wearing gloves
- The adhesive is contained within small capsules, allowing light to pass through to ensure the optical clarity of the seal
- When the seal is in position, pressure can be applied to burst the capsules, releasing a strong adhesive only where the seal touches the raised well rims of the plate - the rest of the seal area above the wells remains optically clear

Key Features

- Non-pierceable
- Peelable
- Seal integrity range: -80°C to 110°C
- Free from DNase, RNase, and human genomic DNA

Use

- Applications: qPCR and other fluorescent applications, plate readers, microscopy and protein crystallization (96 or 384 well)
- Removal: will not leave a sticky residue on the plate surface following removal
- Replacement for: ABI® MicroAmp Optical Adhesive film, Absolute qPCR plate seals and Roche LightCycler® sealing foils

Options

- Sheet format: 140 x 77 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates with 96 round wells
- Roll format: 100 m x 80 mm, approx. 700 seals
- Non-gamma treated as standard; gamma treatment available upon request



Seal Integrity
Temperature Range

-80°C

110°C

*110°C with pressurized heated lid

Peelable



Specifications

- This is a pressure-activated seal: the adhesive is released when pressure is applied firmly and evenly to the seal
- Our Adhesive Seal Roller and Adhesive Seal Applicator are ideal for use with this product
- We also recommend the use of our Optical Film Compression Pad during PCR with this product

Ordering Information

4ti-0560	qPCR Adhesive Seal, optically clear film, pressure activated adhesive, 100 sheets (140 x 77mm) per case
4ti-0561	qPCR Adhesive Seal, optically clear film, pressure activated adhesive, 1 roll (100m x 80mm)
4ti-0561S	qPCR Adhesive Seal, optically clear film, pressure activated adhesive, 1 sample roll (5m x 80mm)



AZENTA
LIFE SCIENCES

PCR Foil Adhesive Seal

Pierceable adhesive aluminium foil, strong adhesive, peelable; suitable for high temperature applications

- This aluminium foil seal has a strong acrylic adhesive which produces a seal of high integrity
- It was developed for PCR and other high temperature applications due to its effectiveness in preventing sample evaporation
- The PCR Foil Seal is pierceable; when pierced, the foil tears in an irregular manner which prevents the formation of a vacuum
- Perforated end tabs for easy application and removal by peeling
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features


- Pierceable
- Peelable
- Seal integrity range: -40°C to 120°C
- Free from DNase, RNase, and human genomic DNA

Use

- Applications: PCR, incubation, storage
- Suitable for cutting to fit part plates



Seal Integrity Temperature Range	-40°C 120°C
-------------------------------------	----------------



Options

- Sheet format: 130 x 80 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0550

PCR Foil Seal, pierceable adhesive aluminium foil, strong adhesive, 100 sheets (130 x 80mm) per case



AZENTA
LIFE SCIENCES

PCR Foil Adhesive Seal Strong

Adhesive aluminium foil, strong adhesive, peelable, pierceable; suitable for high temperature incubations and low temperature storage

- Our PCR Foil Seal Strong is a pierceable aluminium foil seal with a strong acrylic adhesive, recommended for PCR and low temperature storage
- This seal features all of the advantages of our PCR Seal, but in a pierceable format
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features

- Pierceable
- Peelable
- Seal integrity range: -80°C to 110°C
- Non-gamma treated
- Free from DNase, RNase, and human genomic DNA


Use

- Applications: PCR, high temperature incubations, low temperature sample storage



Seal Integrity
Temperature Range

-80°C
110°C



Options

- Sheet format: 137 x 80 mm
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0500FL

PCR Foil Seal Strong, adhesive aluminium foil, strong adhesive, 100 sheets (137 x 80mm) per case



AZENTA
LIFE SCIENCES

DMSO Resistant Foil Adhesive Seal

Peelable adhesive foil, strong adhesive, high solvent resistance; suitable for long term storage

- This aluminium foil seal has a chemically resistant silicone adhesive layer to produce a seal with high levels of solvent resistance, including to Dimethyl Sulfoxide (DMSO)
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features

- Non-pierceable
- Peelable
- Seal integrity range: -20°C to 80°C
- Free from DNase, RNase, and human genomic DNA

Use

- Applications: long-term plate storage, samples containing <80% Dimethyl Sulfoxide (DMSO) can be stored for up to 5 years
- Removal: this seal will not leave a sticky residue on the plate surface following removal
- Suitable for cutting to fit part plates



Options

- Sheet format: 122 x 80 mm, to fit all standard SBS footprint plates, PCR and qPCR plates and part plates, microplates, assay and storage plates
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0512

DMSO Resistant Foil Seal, peelable adhesive foil, strong adhesive, high solvent resistance, 100 sheets (122 x 80mm) per case

Pierceable Film Strong Adhesive, 96 Cross-Cut Windows

Strong adhesive seal cross-cut windows, peelable, pierceable; suitable for 96 well plates (auto samplers, HPLC, sequencers)

- This strong adhesive seal was developed to facilitate easy sample removal with a manual or automated system
- The seal is optically clear, being made from a transparent polyester film, and has a strong adhesive applied across the underside of the seal, except for 96 round windows which align to the 96 wells of a plate
- The optical windows are cross-cut, allowing for easy access to the sample wells with a tip or probe with minimal pressure
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features


- Pierceable
- Peelable
- Seal integrity range: -20°C to 110°C
- Non-gamma treated
- Free from DNase, RNase, and human genomic DNA

Use

- Applications: sample application or retrieval, for use in automated systems and sample analyzers such as HPLC and DNA sequencers
- Removal: will not leave a sticky residue on the plate surface following removal
- Suitable for cutting to fit part plates
- Replacement for: ABI® septa mats on capillary sequencers and to replace silicone and EVA storage plate cap mats



Seal Integrity Temperature Range	-20°C 110°C
-------------------------------------	----------------



Options

- Sheet format: 135 x 77 mm, to fit all standard SBS footprint PCR and qPCR plates, microplates, assay and storage plates with 96 round wells
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0566*	Pierceable Film, strong adhesive seal with 96 cross-cut windows, 100 sheets (135 x 77mm) per case
-----------	---

**Not available for purchase in or onwards distribution to the US*



AZENTA
LIFE SCIENCES

Pierceable Film Strong Adhesive, 384 Cross-Cut Windows

Strong adhesive seal with cross-cut windows, red adhesive, peelable, pierceable; suitable for 384 square well plates


- This strong adhesive seal was developed to facilitate easy sample removal with a manual or automated system
- A red adhesive is applied across the underside of the seal, except for 384 windows which align to the 384 wells of a plate
- The optical windows are cross-cut, allowing for easy access to the sample wells with a tip or probe with minimal pressure
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features


- Pierceable
- Peelable
- Seal integrity range: -80°C to 120°C
- Adhesive-free windows avoid clogging of tips or needles during piercing
- Pre-slit well centres for access to samples with no force needed for piercing
- No cross contamination
- Good solvent and DMSO resistance
- Non-gamma treated
- Free from DNase, RNase, and human genomic DNA

Use

- Applications: sample application or retrieval, for use in automated systems and sample analyzers such as HPLC and DNA sequencers



Seal Integrity Temperature Range	-80°C
	80°C



Options

- Sheet format: 117 x 80 mm, to fit all standard SBS footprint plates, microplates, assay and storage plates with 384 wells with rounded square wells
- Non-gamma treated as standard; gamma treatment available upon request

Ordering Information

4ti-0566/384

Pierceable Film, strong adhesive seal with 384 cross-cut windows, 100 sheets (117 x 80mm) per case

Moisture Barrier Adhesive Seal 24, 96, 384

Gas permeable adhesive film, optically clear, with adhesive free windows, peelable, pierceable, gamma treated; suitable for cell culture

- This unique seal has optically clear adhesive free windows; these windows can be imaged through, and also allow for gas exchange
- The clear film is made of a porous material which allows for a uniform air and CO₂ exchange, whilst acting as a moisture barrier and preventing evaporation
- The seal is coated with a strong adhesive (except for the optical windows) which has a high sealing integrity
- This ensures a reliable seal and prevention of well-to-well contamination
- The gas permeable seal enables long term incubations without intervention, whilst allowing for imaging through the optically clear windows
- Can reduce the risk of sample contamination, evaporation and can improve your experimental workflow
- Use of this seal within plate readers can prevent moisture release into the sensitive equipment
- A number of plate reader manufacturers recommend the use of the Azenta Moisture Barrier Seal 24,96,384 (Gas Permeable Moisture Barrier Seal)
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features

- Pierceable
- Peelable
- Seal integrity range: -20°C to 80°C
- Gas permeability rate: 0.6 m³/m²/day
- Moisture vapor transmission rate: 1 g/m²/day
- Optically clear for imaging
- Gamma treated
- Free from RNase and DNase

Use

- Applications: plate readers, eukaryotic cell culture, bacterial cell culture, long-term incubation, live cell assays, confocal microscopy



Seal Integrity
Temperature Range

-20°C

80°C

Peelable



Pierceable




Options

- Sheet format: 140 x 80 mm (4ti-0516/24) and 137 x 80 mm (4ti-0516/96 and 4ti-0516/384)
- Available with 24/96/384 adhesive free windows for use with 24/96/384 well plates

Ordering Information

4ti-0516/24*	Moisture Barrier Seal, gas permeable adhesive film, optically clear, with 24 adhesive free windows, gamma treated, 5 x 10 sheets (137 x 77mm) per case
4ti-0516/96*	Moisture Barrier Seal, gas permeable adhesive film, optically clear, with 96 adhesive free windows, gamma treated, 5 x 10 sheets (137 x 77mm) per case
4ti-0516/384*	Moisture Barrier Seal, gas permeable adhesive film, optically clear, with 384 adhesive free windows, gamma treated, 5 x 10 sheets (137 x 77mm) per case

**Not available for purchase in or onwards distribution to the USA*



AZENTA
LIFE SCIENCES

Cell Culture Adhesive Seal

Gas permeable adhesive seal, peelable; suitable for cell culture

- Our Cell Culture Adhesive Seal is a gas permeable adhesive seal which seals assay and tissue culture plates, microplates and storage plates
- The Cell Culture Adhesive Seal prevents evaporation and contamination whilst enabling cells to breathe. It is made of a non-woven fiber with an acrylate adhesive layer for effective sealing
- The seal has a low moisture transfer rate and a porosity enabling gas exchange
- Due to its paper-based material it should not be used in wet conditions
- Suitable for cell culture and enables long term culture with significantly reduced evaporation
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features

- Peelable
- Seal integrity range: -20°C to 40°C
- Gas permeability rate: 8,900 m³/m²/day
- Moisture vapor transmission rate (very low): 4,200 g/m²/day
- Air porosity: 10 sec/100 cc/in²
- RNase/DNase free

Use

- Application: bacterial and eukaryotic cell culture



Seal Integrity
Temperature Range

-20°C
40°C



Options

- Sheet format: 135 x 80 mm
- Available gamma treated (4ti-0517/ST)

Ordering Information

4ti-0517	Cell culture adhesive seal, gas permeable, 100 sheets (135 x 80mm) per case
4ti-0517/ST	Cell culture adhesive seal, gas permeable, 10x10 sheets (135 x 80mm) per case



AZENTA
LIFE SCIENCES

Microplate Seal Low Strength Adhesive Film

Low strength adhesive film, transparent, peelable; suitable for short term storage

- This transparent polyester-based film has a low strength adhesive
- It is designed as a low-cost sealing option, and useful for temporary storage and as a cover for applications such as centrifugation
- End tabs allow for easy application and removal
- This seal is removable without residue on the plate
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features

- Non-pierceable
- Peelable
- Seal integrity range: -20°C to 80°C
- Free from DNase, RNase, and human genomic DNA

Use

- Applications: sample storage (aqueous samples only)
- Suitable for all plate types



Seal Integrity
Temperature Range

-20°C
80°C

Peelable



Options

- Sheet format: 130 x 80 mm

Ordering Information

4ti-0510

Microplate Seal, low strength adhesive film, transparent, 100 sheets (130 x 80mm) per case



AZENTA
LIFE SCIENCES

Double Sided Re-Sealable Adhesive Film

Double sided adhesive film, black, with 96 holes and a protective liner, peelable; suitable for re-sealing without the need for a heat sealer

- A double sided black adhesive microplate film to facilitate the sealing, accessing (piercing) and resealing of 8 Well PCR Tube Strips (4ti-0753) prior to PCR
- The strips or plates are filled with reagents and sealed with Pierce Heat Seal (4ti-0530) resulting in lowest evaporation and best possible long-term storage
- With the use of a pierceable seal, samples can easily be added
- The protective Double Sided Re-Sealable Adhesive Film makes sure there will be no damage to the Pierce Heat Seal during transportation
- The film is overlaid and the strips or plates processed with a cutter into individually sealed strips
- The strips can then be transported and opened at a customer site by peeling off the protective layer of the Double Sided Re-Sealable Adhesive Film and accessing the sample through the Pierce Heat Seal
- Can then be resealed with a foil or film which is applied to the exposed black adhesive without the need for a heat sealer
- For all adhesive seals, the best sealing results are achieved using the hand held Adhesive Seal Roller (4ti-0502)

Key Features

- Peelable
- Resealable
- Seal integrity range: -20°C to 110°C
- Free from RNase and DNase

Use

- Applications: ideal for kit manufactures to allow flexible access to individual wells


Options

- Sheet format: 148 x 98 mm, to fit all SBS footprint 96 well plates, designed to work with 96 Well Skirted PCR Plate for Removable 8 Well Tube Strips, Breakable Vertically PCR Plates, and Breakable Horizontally or Vertically plate ranges

Ordering Information

4ti-0519

Double sided, re-sealable, adhesive film, with 96 holes and a protective liner, 100 sheets (148 x 98mm) per case



Seal Integrity
Temperature Range

-20°C

110°C

Peelable






Typical Molecular Diagnostic Workflows – Example

At the kit manufacturer

- The plate (e.g. 96 Well Skirted PCR Plate with Removable 8 Well Strips Frame loaded with 8 Well PCR Tube Strips, e.g. 4ti-0753/C/757) is filled with reagents on a robotic system
- The filled plate is sealed with a pierceable heat seal (e.g. Pierce Seal, 4ti-0530) resulting in lowest evaporation and best possible long-term storage (A)
- The protective Double Sided Re-Sealable Adhesive Film is placed on top of the pierceable seal to make sure there will be no damage of the seal during transportation (B)
- The sealed plate is cut into strips - the pre-filled strips are ready for transport, e.g. to the point of diagnostic use (C)

At the point of diagnostic use

- The Double Sided Re-Sealable Adhesive Film protective layer is removed making the pierceable seal accessible (D)
- The sample is added to the pre-filled strip by piercing the seal (E)
- The strips can easily be resealed with foil or film strips by applying to the exposed black adhesive without the need for a heat sealer (F)
- After resealing, the strips are ready for PCR analysis



AZENTA
LIFE SCIENCES

Adhesive Seal Roller and Applicator

Our Adhesive Seal Roller & Adhesive Seal Applicator ensure even pressure is applied across the adhesive seal for a complete application to the plate

- To obtain the best sealing results when using our adhesive seals, we strongly recommend the use of the Adhesive Seal Roller (4ti-0502) or Adhesive Seal Applicator (4ti-0503)
- Both application tools ensure even pressure is applied across the adhesive seal for complete and secure application to your plate, across every well
- The handle of the Adhesive Seal Roller is made of a durable plastic, with a semi-hard padded rolling wheel
- The straight rigid sides of the small-sized Adhesive Seal Applicator allow for even pressure application
- When applying adhesive seals to 384 well plates, we recommend using the application tools in conjunction with our FrameStar 384 Holder (4ti-0391) to support the 384 well plates during seal applications, and to give a level base



4ti-0503

Ordering Information

4ti-0502	Adhesive Seal Roller; ensures consistent seal application across all wells, 1 roller
4ti-0503	Adhesive Seal Applicator; ensures consistent seal application across all wells, 1 applicator

Support Adapters to Improve Performance

Optical Film Compression Pad

A silicone foam mat laminated to a non-stick PTFE film, to be used with adhesive seals, compatible with heated lid cyclers

- When used in conjunction with an adhesive seal - for instance the qPCR Seal (4ti-0560) - and a thermal cycler heated lid, the pad enhances the adhesion between the seal and the PCR plate
- This in turn improves results by reducing sample evaporation
- The 96 holes align with the wells of the PCR plate, ensuring the mat is compatible with qPCR instrumentation which image through the top of the well



Ordering Information

4ti-0563

Optical Film Compression Pad,
5 pads per case

Pierce Plate

Metal block with 96 pins, suitable for piercing every well of a heat or adhesive sealed 96 well plate

- This useful tool is a machine-engineered metal block with 96 pins aligned central to each well of a 96 well plate
- The pierce plate's 96 pins pierce every well of a heat or adhesive-sealed 96 well PCR or microplate (pierceable seals only)
- Enables instant access to samples with a single or multichannel pipette or automated system
- The Pierce Plate can be cleaned between uses to avoid contamination by using most cleaning agents that are suitable for use on aluminium e.g. RNase removal solutions, bleach or UV treatment



Ordering Information

4ti-0398

Pierce Plate, metal block with 96 pins,
suitable for piercing every well of a heat
or adhesive sealed 96 well plate, 1 plate



AZENTA
LIFE SCIENCES

Adhesive Sealing Consumables Comparison Table

	Clear Seals			Foil Seals			Cross-cut Seals
Name	PCR Seal	Optically Clear Windowed qPCR Seal	qPCR Adhesive Seal	PCR Foil Seal	PCR Foil Seal Strong	DMSO Resistant Foil Seal	Pierceable Film Strong Adhesive, 96 Cross-Cut Windows
Specifications							
Application	PCR	qPCR, fluorescence 96-well microplates only	qPCR & other fluorescent applications Imaging techniques incl. crystallisation Plate readers, microscopy	PCR & sample storage Incubations	High temperature incubations & low temperature storage	Microplate sealing containing solvents incl. DMSO	Sample application or retrieval 96 well plates only
Special Properties	Good optical clarity	Discreet optical windows for 96-well plates	Good optical clarity	Irregular tearing when pierced prevents formation of vacuum	Strong adhesive	High solvent resistance	Cross-cut reduces tip or probe becoming clogged
Seal Integrity Min Temperature	-20°C	-20°C	-80°C	-40°C	-80°C	-20°C	-20°C
Seal Integrity Max Temperature	110°C	110°C	110°C	120°C	110°C	80°C	110°C
Gas Permeability Rate	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Moisture Vapor Transmission Rate	N/A	N/A	N/A	N/A	N/A	N/A	N/A
Gamma Irradiated							
Pierceable				✓	✓		✓
Peelable	✓	✓	✓	✓	✓	✓	✓
RNase/DNase free	✓	✓	✓	✓	✓	✓	✓
Product Codes							
Code	4ti-0500	4ti-0565	4ti-0560	4ti-0550	4ti-0500FL	4ti-0512	4ti-0566
Format	Sheets	Sheets	Sheets	Sheets	Sheets	Sheets	Sheets
Dimensions	135 mm x 80 mm	133 mm x 76 mm	140 mm x 77 mm	130 mm x 80 mm	137 mm x 80 mm	122 mm x 80 mm	135 mm x 77 mm
Code	4ti-0500/8		4ti-0561				
Format	Perforated sheets		Roll				
Dimensions	115 mm x 100 mm		100 m x 80 mm				
Code	4ti-0500/12		4ti-0561/S				
Format	Perforated sheets		Sample roll				
Dimensions	137 mm x 71 mm		5 m x 80 mm				
Code	4ti-0500/HP		4ti-0560/HP				
Format	Half plate sheets		Half plate sheets				
Dimensions	70 mm x 80 mm		70 mm x 77 mm				




Cross-cut Seals	Permeable Seals				
Pierceable Film Strong Adhesive, 384 Cross-Cut Windows	Moisture Barrier Seal 24, 96, 384	Cell Culture Adhesive Seal	Microplate Seal Low Strength Adhesive Film	Double Sided Re-Sealable Adhesive Film	Name
					Specifications
Sample application or retrieval 384 well plates only	Eukaryotic cell culture, bacterial culture and long-term live assays Suitable for use on plate readers	Bacterial or cell culture	Aqueous sample storage	Assay kit production	Application
Cross-cut reduces tip or probe becoming clogged Good solvent resistance, including DMSO	Gas permeable that allows air and CO ₂ exchange, but prevents moisture evaporation Optically clear	Very low moisture transfer rate Suitable for bacterial or cell culture Air porosity: 10 sec/100 cc/in ²	Medium strength transparent seal	Two sealing surfaces Optical windows	Special Properties
-80°C	-20°C	-20°C	-20°C	-20°C	Seal Integrity Min Temperature
120°C	80°C	80°C	80°C	110°C	Seal Integrity Max Temperature
N/A	0.6 m ³ /m ² /day	8,900 m ³ /m ² /day	N/A	N/A	Gas Permeability Rate
N/A	1 g/m ² /day	4,200 g/m ² /day	N/A	N/A	Moisture Vapor Transmission Rate
	✓	✓			Gamma Irradiated
✓	✓				Pierceable
✓	✓	✓	✓	✓	Peelable
✓	✓	✓	✓	✓	RNase/DNase free
					Product Codes
4ti-0566/384	4ti-0516/24	4ti-0517	4ti-0510	4ti-0519	Code
Sheets	Sheets	Sheets	Sheets	Sheets	Format
117 mm x 80 mm	140 mm x 80 mm	135 mm x 80 mm	130 mm x 80 mm	148 mm x 98 mm	Dimensions
	4ti-0516/96	4ti-0517/ST			Code
	Sheets	Sheets, gamma treated			Format
	137 mm x 80 mm	135 mm x 80 mm			Dimensions
	4ti-0516/384				Code
	Sheets				Format
	137 mm x 80 mm				Dimensions
					Code
					Format
					Dimensions



Plate Lids, Caps & Mats



Plate Lids, Caps & Mats



As an alternative to sealing films, Azenta offers multiple types of cap strips for sealing both plates and tubes - domed, flat, strips of 8, strips of 12, and our new optically superior caps.

A variety of rigid polystyrene lids are available for PCR plates and microplates, including lids compatible with our FrameStar®, Optically Clear Tissue Culture Plate and Ultra Optically Clear Plate ranges. Azenta also stock silicone sealing mats for use with our storage plate ranges and in a variety of formats depending on the well size, number and shape.



Strips of 8 & 12 Flat Sealing Caps

Clear polypropylene sealing caps, available as strips of 8 caps (domed or flat optical) and strips of 12 caps (flat optical)

- Compatible with our tube strips and 96 well PCR plates
- These strips are molded from virgin polypropylene in our UK-based Class 7 ISO certified clean-room production facility, and comply to the same stringent requirements as our FrameStar range

Key Features

- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- Flat caps are optically clear for fluorescence detection (e.g. qPCR)
- Easy to apply
- Large end tabs for easy removal
- Labeled for orientation

Options

- Strips of 8 flat optical caps (4ti-0751/4ti-0783)*: for use with our FrameStar PCR plates (including FrameStar Breakable Vertically PCR Plate), 8 Well PCR Tube Strip with PC Frame, 8 Well Removeable Tube Strips and general PCR plates, as well as with PCR tube strips (4ti-0781)
- Strips of 8 domed caps (4ti-0752/4ti-0782)*: for use with our FrameStar PCR plates (including FrameStar Breakable Vertically PCR Plate), 8 Well PCR Tube Strip with PC Frame, 8 Well Removeable Tube Strips and general PCR plates, as well as with PCR tube strips (4ti-0781)
- Strips of 12 flat optical caps (4ti-0788): recommended for use with our 96 Well Non-Skirted PCR Plate Breakable Horizontally or Vertically and FrameStar Breakable Horizontally and Vertically PCR Plates to allow for flexible sample usage; not compatible with low profile plates



Ordering Information

4ti-0751	Strips of 8 Flat Optical Caps, 300 strips per case
4ti-0783	Strips of 8 Flat Optical Caps, 125 strips per case
4ti-0752	Strips of 8 Domed Caps, 300 strips per case
4ti-0782	Strips of 8 Domed Caps, 125 strips per case
4ti-0788	Strips of 12 Flat Optical Caps, 200 strips per case ¹

¹ Recommended for use with our 96 Well Non-Skirted PCR Plate Breakable Horizontally or Vertically and FrameStar Breakable Horizontally and Vertically PCR Plates to allow for flexible sample usage; not compatible with low profile plates

Strips of 8 Flat Optical Caps Crystal Clear

Strips of 8 flat optical caps, crystal clear; designed for low volume applications such as low volume qPCR


- With the Strips of 8 Flat Optical Caps Crystal Clear Azenta offers the perfect supplement to the existing range of PCR Cap Strips
- Due to their improved optical properties and evaporation-safe fit, the strips are ideally suited for applications where small volumes are used, e.g. low volume qPCR

Key Features

- Made of a special polymer with improved optical properties leading to high transmission rates; ideally suited for small samples with low signal intensity
- Reduced shrinking during heating and cooling phases; very tight sealing
- Highest flexibility between the pitch of the individual caps
- Variations in the pitch of any 96 well plate are accommodated, the caps can easily be aligned to the corresponding tubes as they allow for more independent movement
- Free from DNase, RNase, human genomic DNA, and endotoxin/pyrogen

Use

- For use with our FrameStar PCR plates (including FrameStar Breakable Vertically PCR Plate), 8 Well PCR Tube Strip with PC Frame and general PCR plates



Ordering Information


4ti-0755	Strips of 8 Flat Optical Caps, crystal clear, 300 strips per case
4ti-0755/120	Strips of 8 Flat Optical Caps, crystal clear, 120 strips per case

For use with our our FrameStar PCR plates (including FrameStar Breakable Vertically PCR Plate), 8 Well PCR Tube Strip with PC Frame and general PCR plates

PCR Plate Lids, FrameStar Lids and Microplate Lids

Rigid polystyrene lids for PCR plates, optically clear tissue culture plates, and assay plates

- Designed to give a quick and easy sealing solution to protect samples from contamination and evaporation




Ordering Information

PCR Plate Lids & FrameStar Lids

4ti-0285	Ultra-Low Universal Lid, without condensation rings, clear, ultra-low profile, no cut corner, 100 lids per case
4ti-0288	PCR Plate Lid, without condensation rings, clear, low profile, cut corner H1, 50 lids per case
4ti-0287	FrameStar 96 Next Generation Sequencing Lid, with condensation rings, clear, low profile, cut corner H12, for use with 4ti-0960/RIG, 50 lids per case
4ti-0289	FrameStar 96 Lid, without condensation rings, clear, low profile, cut corner A12, for use with 4ti-0770, 50 lids per case

Microplate Lids

4ti-0280	Microplate 384 Lid, without condensation rings, clear, low profile, cut corner A1/P1, 100 lids per case
4ti-0281	Microplate 384 Lid, gamma treated, without condensation rings, clear, low profile, cut corner A1/P1, 100 lids per case
4ti-0282	Microplate 96 Lid, with condensation rings, clear, low profile, cut corner A1/H1, 80 lids per case
4ti-0283	Microplate 96 Lid, gamma treated, with condensation rings, clear, low profile, cut corner A1/H1, 80 lids per case
4ti-0284	Microplate 24 Lid, with condensation rings, clear, low profile, cut corner A1, 80 lids per case
4ti-0286	Microplate 24 Lid, gamma treated, with condensation rings, clear, low profile, cut corner A1, 80 lids per case
4ti-0290	Universal Microplate Lid, without condensation rings, clear, low profile, no cut corner, 50 lids per case



Auto-Sealing PCR Plate Lid

Low profile, with integrated compression pad, white, no cut corner, for PCR applications on integrated instruments

- Developed to support fully automated sealing in the absence of specific instrumentation
- Meets the needs of customers using integrated platforms that often lack suitable heat sealing instrumentation especially when low throughput is needed
- Specialized alternative to standard lids: while standard lids protect reagents, the Auto-Sealing PCR Lid helps minimize reagent evaporation during longer incubations
- Universal fit: due to the lack of cut corners, the lid can be applied to most plates, showing good sealing results when pushed down onto the plate

Key Features

- Rigid polycarbonate frame padded with an elastic foam
- No cut corners
- Stackable

Use

- Suitable for automation
- Universal fit to PCR plates
- Recommended when full sealing automation is required in the absence of dedicated instrumentation
- Suitable for low throughput workflows
- Alternative to heat sealing, when heat sealing materials and instrumentation are not an option
- Alternative to standard lids for longer incubations



Specifications

Parameter	Value
Lid length	128.10 ± 0.10 mm
Lid width	85.80 ± 0.10 mm
Lid height	8.20 ± 0.05 mm



Ordering Information

4ti-0291	Auto-Sealing PCR Plate Lid, with integrated compression pad, white, low profile, no cut corner, 20 lids per case
-----------------	--

Sealing Cap Mats

Silicone rubber or TPE mats, durable to high temperatures; for sealing storage plates to protect samples from evaporation

- The clear cap mats are made of silicone rubber, a material that is highly durable to high temperatures, and so can be used to seal storage plates being used for high-temperature storage to protect samples from evaporation
- Azenia silicone mats are DNase/RNase and pyrogen-free to meet the highest standard of both laboratory experiment
- All our clear silicone cap mats are pierceable



Ordering Information

4ti-0124	96 Round well Sealing Cap Mat , clear silicone, for use with 4ti-0125, 50 mats per case
4ti-0137	96 Square Well Sealing Cap Mat , clear silicone, for use with square 96 well microplates and deep well storage microplates, 50 mats per case
4ti-0138	96 Round Well Sealing Cap Mat , clear silicone, for use with round 96 well microplates and deep well storage microplates (not for use with 4ti-0120 and 4ti-0110), 50 mats per case
4ti-0139	384 Square Well Sealing Cap Mat , clear silicone, for use with square 384 well microplates and deep well storage microplates, 50 mats per case
4ti-0135	96 Round Well Sealing Cap Mat , white silicone, for use with 4ti-0120 only, 100 mats per case

Cap Mat for PCR Plates

96 individual caps in sheet format, blue TPE, pierceable; suitable for sealing all of our 96 well PCR plates

- These caps in sheet format are universally compatible with our 96 well PCR plates
- The caps can be individually applied and removed once detached from the backing liner, making the mats ideally suited for use with our flexible PCR consumables, including Individual Access and divisible plates
- The mats offer an alternative to adhesive and heat sealing, in particular as a temporary solution when samples need to be repeatedly accessed
- They are easily pierceable with pipette tips to access samples, and they are easily removable using 1- and 8-way decappers or, alternatively, using Azenta Automated Plate Seal Remover if a seal is overlaid on top of the caps

Key Features

- Caps made from TPE mounted on an easily removable backing liner
- Pierceable with pipette tips

Use

- Universal fit to 96 well PCR plates
- Ideally suitable for use with Individual Access and breakable PCR plates
- Applications: endpoint PCR, storage
- Not recommended for qPCR
- Alternative to adhesive and heat sealing, as temporary sealing solution when sample access must be carried out multiple times
- Caps perform equivalently to polypropylene cap strips during 25 cycle PCR
- The 96 caps can be applied all at once or individually
- The 96 caps can be removed all at once using Azenta Automated Plate Seal Remover (a seal must be placed on top of the caps), or individually using 1- or 8-way decappers



Ordering Information

4ti-0778

Cap Mat for PCR Plates, 96 individual caps in sheet format, blue TPE, pierceable, 50 mats per case



AZENTA
LIFE SCIENCES

Custom Capabilities



Custom Capabilities for OEM Diagnostic Kit and Medical Device Manufacturing

Our team brings decades of development expertise to any custom project. Our experience in working with partners worldwide has placed us as the company of choice for introducing practical solutions to changing workflows, using innovative design and bespoke manufacture.

Our commitment to quality is reflected in our ISO certified management system which is applied at all levels, from manufacture, to technical support, to packaging and delivery of plates and seals.

Azenta has an integrated quality management system where plate and seal products undergo a wide range of QC inspections. We constantly perform visual, physical and biological tests to ensure both the absence of contaminants, as well as the integrity of the products. Our ISO 13485:2016 certification is an endorsement of our excellent manufacturing practices.

The ultimate quality of any product is dependent not only on the design of the component to be produced, but also the accuracy, construction and precision of the tooling and manufacturing processes.

Azenta's highly skilled engineers have extensive knowledge and experience in the design and manufacture of precision components for the life science industry.

By working with our in-house engineers we can offer a completely integrated project management service to ensure there is a smooth transition from initial project idea through to finished product. We aim to help enhance end products, reduce overall part and production costs, and streamline assembly within manufacturing. Whatever your custom requirement, you will receive ongoing support and advice from a designated sales contact, our QC department and Customer Services team.

We understand the costs and complexities involved with OEM products and would be happy to discuss a range of solutions for your project. Please contact us to discuss your specific requirements in detail and complete confidence.



Custom and OEM Services Overview

- Product design and manufacturing for injection-molded parts
- Instrumentation design and manufacture: Heat sealers, press tools, cutters, liquid handling instrumentation
- Heat and adhesive sealing: Custom material design and manufacture
- Sample tracking solutions: Custom applications and specifications for linear and 2D coding
- Surface treatment options
- Gamma treatment options
- Packaging choices
- Specified QC procedures

Customizable PCR Plates, Microplates and Seals

We partner with our customers to provide flexible solutions that improve reproducibility and achieve consistent results. We can provide the complete end-to-end development of custom, premium plastic consumables including PCR plates, strips and tubes, sealing materials and instrumentation. Our completely tailored solutions include custom designs, tool making and contract manufacturing from concept to completion.

Plates, Strips, and Tubes



Customizable options:

- Non-SBS format
- Plate color
- Well color, including clear wells and white wells
- Application specific direct marking and branding options
- Coating
- Treatments
- 2D coding and barcoding for tracking of samples

Plate Seals



Customizable options:

- Non-standard size seals, in rolls or sheets
- Application specific direct marking and branding options
- Custom instrumentation to support your workflows
- Material
- Chemical Compatibility



AZENTA
LIFE SCIENCES

Sample Cooling and Heating



Sample Cooling and Heating Standardization



Prepare.

Ice-free sample cooling and freezing

- Consistent and reproducible
- Ideal for working in a hood



CoolBox™ Thermoconductive Tube Rack and Ice-Free Cooling Workstation systems



Protect.

Controlled-rate cell freezing

- No alcohol
- High post-thaw recovery and viability
- Proven for stem cells, primary cells, PBMC, cell lines and more



CoolCell® Alcohol Free Cell Freezing Containers



Preserve.

Archival storage integrity

- Hinged lid helps box and lid stay together
- LN₂ drain holes and water proof fiberboard



TruCool™ Hinged CryoBoxes



AZENTA
LIFE SCIENCES

CoolCell® Alcohol-Free Cell Freezing Containers



AZENTA
LIFE SCIENCES

CoolCell Alcohol-Free Cell Freezing Containers

Alcohol-free cell freezing containers ensure standardized controlled-rate $-1^{\circ}\text{C}/\text{minute}$ cell freezing in a -80°C freezer - without alcohol or any fluids. Proven for use with a variety of cell types including stem cells, primary cells, PBMC cell lines, insect cells, yeast and others. The Alcohol-Free Cell Freezing technology utilizes a thermo-conductive alloy core and highly-insulative outer material to control the rate of heat removal and provide reproducible cell cryopreservation. Alcohol-Free Cell Freezing units are easy to use and deliver comparable results to expensive programmable freezers.

Alcohol-free cell freezing containers are proven to work with many cell types including:

Stem Cells


- Human Embryonic Stem
- Preadipocytes
- Breast Cancer Stem
- Colon Cancer Stem
- Glioblastoma Stem
- Mouse Embryonic Stem
- Human Endothelial
- Progenitor

Primary Cells

- Neonatal Keratinocytes
- Human WBCs
- Mouse
- WBCs
- Human CD34+
- Muscle
- Human Tendon
- Fibroblasts
- Melanoma Tumor
- Human Cardiac
- Ventricular
- Human Cardiac Atrial

Cell Lines


- CHO
- LnCap
- HTB77
- A549
- HeLa



Cell Freezing Containers, in combination with a -80°C freezer, will provide the freezing rate of $-1^{\circ}\text{C}/\text{minute}$ that is ideal for cryopreservation of most cells and cell lines. Using a combination of uniform-density cross-linked polyethylene foam, a solid state core, and radial vial symmetry, freezing profiles are consistent and reproducible. It is important to fully load Cell Freezing Containers prior to freezing. Foam is non-absorbent and will impose negligible change in the freezer environment; thereby protecting nearby frozen samples. The low heat content also ensures that Cell Freezing Containers will rapidly return to room temperature when removed from the freezer.



Alcohol-free cell freezing containers	Isopropanol (IPA) Container
No alcohol <ul style="list-style-type: none"> • No fluids • No pre-cooling • Saves 12L/unit of IPA per year 	Requires isopropanol <ul style="list-style-type: none"> • Replace alcohol every 5 uses • Track number of uses • Pre-cool alcohol in refrigerator
No variability <ul style="list-style-type: none"> • All vials have uniform freeze rate • Radially symmetric design ensures vial consistency 	Inconsistent freeze rate <ul style="list-style-type: none"> • Alcohol degradation induces variability • Two circles of wells; two freeze rates
No on-going cost <ul style="list-style-type: none"> • No alcohol purchase or disposal 	Approximately \$350/year <ul style="list-style-type: none"> • Change alcohol weekly • Dispose of hazardous waste
No stuck lids <ul style="list-style-type: none"> • Ergonomic lid comes off easily when frozen • Not cold to the touch when removing from the -80°C freezer 	Difficult to handle and open <ul style="list-style-type: none"> • Screw cap difficult to remove when frozen • Frozen unit is slippery and cold to touch
Quick re-use time <ul style="list-style-type: none"> • Ready to use again after five minutes 	Wait between runs <ul style="list-style-type: none"> • Takes >1 hr for alcohol to warm-up




Cell Freezing Containers for 30 x 1mL or 2mL Cryo Tubes utilize a solid state core and controlled micro-convection technology to evenly draw in -80°C freezer air through a bottom base vent, uniformly disperse the cold air around each vial in the central chamber and then release the thermal load from the vials through a top vent choke. The inner vial module holds 30 cryogenic vials and can be removed in one step. Each vial achieves a uniform and reproducible $-1^{\circ}\text{C}/\text{minute}$ freezing profile and thermal profiles are highly reproducible. Due to the low thermal mass of the uniform-density cross-linked polyethylene foam container, freezing can be conducted without a rise in local freezer temperature, thereby protecting nearby samples.



AZENTA
LIFE SCIENCES

CoolCell Alcohol-Free Cell Freezing Containers

Alcohol-Free Cell Freezing Container Reproducibility




Performance test: A temperature probe was placed into a 2.0 mL cryogenic vial containing 1.0 mL of water and the tube was inserted into a room temperature Alcohol-Free Cell Freezing Container. The container was placed into a -80°C freezer and the temperature rate and profile was recorded over a 3 hour period. The test was repeated 5 consecutive times. **Conclusion:** The Alcohol-Free Cell Freezing Container generated identical fusion time and cooling profiles over five consecutive freeze cycles.

Alcohol-Free Cell Freezing Container Protocols




Cryopreservation and Thawing of Cells

Wayne M. Yokoyama, Maria L. Thompson, Rolf O. Ehrhardt
University of California School of Medicine, San Francisco, CA BioCision LLC, Larkspur, California
Curr. Protoc. Immunology. 2012 Nov; 99 Appendix 3G



Standardized Cryopreservation of Pluripotent Stem Cells


Rick I. Cohen, Maria L. Thompson, Brian Schryver, Rolf O. Ehrhardt
Rutgers University, Piscataway, New Jersey
BioCision LLC, San Rafael, California
Curr. Protoc. Stem Cell Biol. 28:1C.14.1-1C.14.10




Standardized Cryopreservation of Human Primary Cells

Thomas V. Ramos, Aby J. Mathew, Maria L. Thompson, Rolf O. Ehrhardt
HemaCare Corporation, Van Nuys, California, BioLife Solutions, Bothell, Washington, BioCision, Larkspur, California
Curr. Protoc. Cell Biology. 2014 Sept; 64 Appendix 3I.


Alcohol-Free Cell Freezing Container Performance vs. IPA Container




Human embryonic stem cells, RC-10 were frozen using the technique indicated, thawed after 2 weeks in LN₂, and counted immediately (Day 1) or after 3 days of growth (Day 3).



Alamar blue reduction assay for proliferation assessment showed cells frozen in an Alcohol-Free Cell Freezing Container grew more quickly, leading to more total cells.



HeLa, CHO-K, K562, NIH3T3. 12-well Alcohol-Free 30-well Cell Freezing Containers, Cell Freezing Containers or “Mr. Frosty” freezing containers were used to freeze all four cell lines. Identical transfection efficiencies and viabilities were observed after thawing.




Identical growth of cells was observed 24 hours post-thaw.




AZENTA
LIFE SCIENCES

CoolCell Alcohol-Free Cell Freezing Containers

Alcohol-Free Cell Freezing Container Performance vs. Programmable Freezer



Ag-Tregs. Effects of freezing on antigen-specific Treg (Ag-Treg) cell therapy products; Ag-Tregs (n = 6) were frozen at concentration of 1 to 10 x 10⁶ cells/mL using the Alcohol-Free Cell freezing Container or controlled-rate freezer (CRF) with a freezing rate of -1°C/min. Viability and absolute viable cell count of thawed Ag-Treg cell therapy products were evaluated by flow cytometry. - Data generated by TxCell SA



Comparison of freezing methods. Graph comparing % of cell recovery after freezing with the Alcohol-Free Cell Freezing Container (blue) versus freezing using a programmable step-down freezer (gray) in 3 different samples at high cell concentration. There was no significant difference between the two freezing methods. - Data performed by UCSF Diabetes Center

For 1.0 mL or 2.0 mL Cryogenic Vials



Cell Freezing Containers for 12 x 1mL 96-format Sample Tubes

For 12 standard 1mL storage tubes. 0.5mL to 1mL fill per vial. Optimized for freezing 1mL 96-format sample tubes. Radially symmetric for uniform freezing. Numbered wells for easy identification. Beveled lid for secure gripping and easy opening.



Cell Freezing Containers for 12 x 1mL or 2mL Cryo Tubes

For 12 standard 1.0 mL to 2.0 mL cryogenic vials, 1.0 mL fill per vial. Radially symmetric for uniform vial freezing. Numbered wells for easy sample identification. Beveled lid for secure gripping and easy opening. Exposed vial tops when lid is open for quick, organized removal of frozen samples.

Ordering Information

BCS-407P	Cell Freezing Container, for 12 x 1ml 96-format sample tubes, purple
BCS-407O	Cell Freezing Container, for 12 x 1ml 96-format sample tubes, orange

Ordering Information

BCS-405	Cell Freezing Container, for 12 x 1ml or 2ml cryo tubes, purple
BCS-405G	Cell Freezing Container, for 12 x 1ml or 2ml cryo tubes, green
BCS-405O	Cell Freezing Container, for 12 x 1ml or 2ml cryo tubes, orange
BCS-405PK	Cell Freezing Container, for 12 x 1ml or 2ml cryo tubes, pink
BCS-405MC	Cell Freezing Container, for 12 x 1ml or 2ml cryo tubes, multipack with 4 colors - purple, green, orange and pink



AZENTA
LIFE SCIENCES

For 3.5 mL to 5.0 mL Cryogenic Vials



Cell Freezing Containers for 12 x 3.5mL to 5mL Cryo Tubes

For 12 standard 3.5 mL to 5.0 mL fill cryogenic vials, 3.5 to 5.0 mL fill per vial. Radially symmetric for uniform vial freezing. Numbered wells for easy sample identification. Beveled lid for secure gripping and easy opening. Exposed vial tops when lid is open for quick, organized removal of frozen samples.



Cell Freezing Containers for 30 x 1mL or 2mL Cryo Tubes

For 30 standard 1.0 mL to 2.0 mL cryogenic vials, 1.0 mL fill per vial. Controlled micro-convection for uniform freezing of 30 vials. Removable vial tray for one-step transfer of samples into and out of freezing chamber.

Ordering Information

BCS-406	Cell Freezing Container, for 12 x 3.5mL to 5mL cryo tubes, purple
---------	---

Ordering Information

BCS-170	Cell Freezing Container, for 30 x 1ml or 2ml cryo tubes, purple
BCS-170G	Cell Freezing Container, for 30 x 1ml or 2ml cryo tubes, green
BCS-170O	Cell Freezing Container, for 30 x 1ml or 2ml cryo tubes, orange
BCS-170PK	Cell Freezing Container, for 30 x 1ml or 2ml cryo tubes, pink

For Injectable Cell Therapy Ampules

Cell Freezing Containers for 12 x 2mL Injectable Cell Therapy Ampules and Cell Freezing Containers for 6 x 10mL Injectable Cell Therapy Ampules

For 12 standard 2.0 mL injectable ampules, 1.0 mL fill per ampule (Cell Freezing Containers for 12 x 2mL Injectable Cell Therapy Ampules). For 6 standard 10.0 mL injectable ampules, 5.0 mL fill per ampule (Cell Freezing Containers for 6 x 10mL Injectable Cell Therapy Ampules). Radially symmetric for uniform freezing of injectable ampules. Easy open lid. Exposed vial tops when lid is open for quick, organized removal of frozen samples.



Ordering Information

BCS-172	Cell Freezing Container, for 12 x 2ml injectable cell therapy ampules, purple
BCS-262	Cell Freezing Container, for 6 x 10ml injectable cell therapy ampules, purple

Cell Cryopreservation Systems



Cell Freezing Containers for 12 x 2mL Injectable Cell Therapy Ampules Stem Cell Cryopreservation System

Ordering Information

BCS-172CS	Stem Cell Cryopreservation System, containing 1 x Cell Freezing Container, for 12 x 2ml injectable cell therapy ampules, purple and 1 x Thermoconductive Tube Rack, holds 12 x 2ml injectable cell therapy ampules, cylindrical wells, gray
-----------	---



Cell Freezing Containers for 6 x 10mL Injectable Cell Therapy Ampules Stem Cell Cryopreservation System

Ordering Information

BCS-262CS	Stem Cell Cryopreservation System, containing 1 x Cell Freezing Container, for 6 x 10ml injectable cell therapy ampules, purple and 1 x Thermoconductive Tube Rack, holds 12 x 10ml injectable cell therapy ampules, cylindrical wells, gray
-----------	--



Note: For optimal freezing it is important to fully load each Cell Freezing Container prior to freezing. Cell Freezing Container Filler Vials are recommended for filling any empty wells.

Cell Freezing Container Filler Vials

To ensure cell freezing rate consistency and uniform results when using Azenta containers, insert a Cell Freezing Container Filler Vial into empty wells when freezing less than a full load. Suitable for repeated use and compatible with Cell Freezing Containers for 12 x 1mL or 2mL Cryo Tubes, Cell Freezing Containers for 30 x 1mL or 2mL Cryo Tubes and Cell Freezing Containers for 12 x 3.5mL to 5mL Cryo Tubes containers. 6 per pack.

Ordering Information

BCS-3105	Cell Freezing Container Filler Vials, 6 x 2ml
BCS-3106	Cell Freezing Container Filler Vials, 6 x 5ml
BCS-3107	Cell Freezing Container Filler Vials, 6 x 1mL



Cell Freezing Container Vial Module for 30 x 1ml or 2ml Cryo Tubes

Cell Freezing Container Vial Module for 30 x 1ml or 2ml Cryo Tubes is a holder for 30 1.0 mL or 2.0 mL cryogenic vials that allows one-step insertion and removal of all 30 vials at once. Fits into a standard 5.0 x 5.0 x 2.0 inch cryostorage box. Compatible with dry ice and liquid nitrogen.

Ordering Information

BCS-210	Removable Cryo Tube Module for use with the Cell Freezing Container for 30 x 1mL or 2ml Cryo Tubes
---------	--




CoolBox™ Ice-Free Cooling Workstations



AZENTA
LIFE SCIENCES

CoolBox Ice-Free Cooling Workstations







Ice-Free Cooling Workstations are bench top cooling workstations that provide sample cooling or freezing without ice, electricity or batteries. Ice-Free Cooling Workstations are versatile and accommodate a variety of sample formats and temperatures. The modular design enables the use of Thermoconductive Tube Rack and Thermoconductive Sink sample modules to hold microfuge tubes, cryogenic vials, PCR tubes or plates, assay plates and more.



How It Works

- Freeze the Cooling Core in -20°C freezer
- Place frozen core in base
- Place Thermoconductive Tube Rack or Thermoconductive Sink module on core
- Module will equilibrate and maintain temperature via thermo-conductivity



How to Configure an Ice-Free Cooling Workstation System

	1. Identify Tube or Plate	2. Choose Thermoconductive Tube Rack or Thermoconductive Sink Module	3. Choose Ice-Free Cooling Workstation capacity and color
Example 1	up to 24 microcentrifuge tubes	 Thermoconductive Tube Rack for 24 Microcentrifuge Tubes	 Cooling Workstation, Single Capacity
Example 2	up to 48 microcentrifuge tubes	 2 x Thermoconductive Tube Racks for 24 Microcentrifuge Tubes	 Cooling Workstation, Double Capacity
Example 3	24 microcentrifuge tubes and one PCR plate, 12 PCR strips, or 0.2 mL PCR tubes	 Thermoconductive Tube Rack for 24 Microcentrifuge Tubes + Thermoconductive Tube Rack for 96-Well PCR Plates	 Cooling Workstation, Double Capacity



AZENTA
LIFE SCIENCES

Choose Your Thermoconductive Tube Rack Module



Open	4 hr	10 hr	10 hr
Closed	--	16 hr	16 hr

For Use With	Wells	Well Shape	Dimensions (L x W x H)	Well Dia.	Well Depth	Description	Item No.	Capacity		
Microfuge tube modules: Thermoconductive Tube Racks for Microcentrifuge Tubes										
1.5 mL or 2.0 mL tubes	6	Cylindrical	6.0 x 4.3 x 3.8 cm	11.1 mm	32.7 mm	Thermoconductive Tube Rack for 6 Microcentrifuge Tubes	BCS-163	1	up to 3	up to 8
1.5 mL or 2.0 mL tubes	15	Cylindrical	10.2 x 6.4 x 3.8 cm	11.1 mm	32.7 mm	Thermoconductive Tube Rack for 15 Microcentrifuge Tubes	BCS-125	1	1	up to 4
1.5 mL conical tubes	15	Conical	10.2 x 6.4 x 3.8 cm	11.1 mm	35.3 mm	Thermoconductive Tube Rack for 15 Microcentrifuge Tubes, Conical Wells	BCS-127	1	1	up to 4
1.5 mL Or 2.0 mL tubes	24	Cylindrical	12.8 x 8.5 x 3.8 cm	11.1 mm	32.7 mm	Thermoconductive Tube Rack for 24 Microcentrifuge Tubes*	BCS-535	1	1	up to 2
5.0 mL centrifuge tubes	12	Conical	12.7 x 8.6 x 5.0 cm	16.5 mm	48.7 mm	Thermoconductive Tube Rack for 12 x 5mL Microcentrifuge Tubes*	BCS-539	1	1	up to 2
1.5 mL or 2.0 mL tubes	30	Cylindrical	12.0 x 10.2 x 3.8 cm	11.1 mm	32.7 mm	Thermoconductive Tube Rack for 30 Microcentrifuge Tubes	BCS-108	1	--	up to 2
1.5 mL conical tubes	30	Conical	12.0 x 10.2 x 3.8 cm	11.1 mm	35.3 mm	Thermoconductive Tube Rack for 30 Microcentrifuge Tubes, Conical Wells	BCS-128	1	--	up to 2
500 uL conical tubes	30	Conical	12.0 x 10.2 x 3.8 cm	11.1 mm	35.3 mm	Thermoconductive Tube Rack for 30 Microcentrifuge Tubes, 500µl	BCS-137	1	--	up to 2
Cryogenic vial and FACS tube modules: Thermoconductive Tube Racks for Cryo or FACS Tubes										
cryogenic vials or FACS tubes	15	Cylindrical	10.2 x 6.4 x 3.8 cm	12.7 mm	32.7 mm	Thermoconductive Tube Rack for 15 Cryo or FACS Tubes	BCS-126	1	1	up to 4
cryogenic vials or FACS tubes	24	Cylindrical	12.8 x 8.5 x 3.8 cm	12.7 mm	32.7 mm	Thermoconductive Tube Rack for 24 Cryo or FACS Tubes*	BCS-534	1	1	up to 2
cryogenic vials or FACS tubes	30	Cylindrical	12.0 x 10.2 x 3.8 cm	12.7 mm	32.7 mm	Thermoconductive Tube Rack for 30 Cryo or FACS Tubes®	BCS-138	1	--	up to 2
cryogenic vials or FACS tubes	45	Cylindrical	17.3 x 9.7 x 3.8 cm	12.7 mm	32.7 mm	Thermoconductive Tube Rack for 45 Cryo or FACS Tubes	BCS-105	--	--	1
PCR plate, strip well or tube modules: Thermoconductive Tube Racks for PCR Plates										
One 96-well PCR plate, strip wells, 0.2mL tubes	96	Tapered	12.7 x 8.6 x 2.5 cm	-	13.2 mm	Thermoconductive Tube for 96-Well PCR Plates*	BCS-529	1	1	up to 2
6 strip wells and 12 x 1.5 or 2.0 mL microfuge tubes	48(PCR) 12(M)	Tapered(PCR) Cylindrical(M)	12.7 x 8.6 x 3.8 cm	- 11.1 mm	13.2 mm 32.7 mm	Thermoconductive Tube Rack for Microcentrifuge Tubes Plus Strip Wells*	BCS-523	1	1	up to 2
One 384-well PCR plate	384	Tapered	12.7 x 8.6 x 1.9 cm	-	7.6 mm	Thermoconductive Tube Rack for 384-Well PCR Plates*	BCS-538	1	1	up to 2
2D coded storage tube modules: Thermoconductive Tube Racks for 96-Well Coded Tubes										
0.5 mL 2D storage tubes	96	Cylindrical	13.1 x 8.9 x 3.6 cm	8.4 mm	24.6 mm	Thermoconductive Tube Rack for 96 x 0.5mL Barcoded Tubes	BCS-231	1	1	up to 2
1.4 mL 2D storage tubes	96	Cylindrical	13.2 x 8.9 x 3.6 cm	8.3 mm	32.7 mm	Thermoconductive Tube Rack for 96 x 1mL Barcoded Tubes	BCS-149	1	1	up to 2
Cell therapy injectable ampule modules: Thermoconductive Tube Racks for Injectable Cell Therapy Ampules										
2.0 mL injectable cell therapy ampules	12	Cylindrical	12.7 x 8.6 x 3.8 cm	16.0 mm	24.0 mm	Thermoconductive Tube Rack for 12 x 2mL Injectable Cell Therapy Ampules	BCS-266	1	1	up to 2
10.0 mL injectable cell therapy ampules	12	Cylindrical	12.7 x 8.6 x 3.8 cm	23.6 mm	27.9 mm	Thermoconductive Tube Rack for 12 x 10mL Injectable Cell Therapy Ampules	BCS-265	1	1	up to 2

* SBS-compatible ø "gripping" wells for one-hand vial opening/closing



CoolBox Ice-Free Cooling Workstations

For Use With	Wells	Well Shape	Dimensions (L x W x H)	Well Dia.	Well Depth	Description	Item No.	Capacity		
Tall tube modules: Thermoconductive Tube Racks for 15 mL, 50mL and 250 mL Centrifuge Tubes										
15 mL centrifuge tubes	12	Cylindrical	13.7 x 9.5 x 11.8 cm	17.5 mm	105.4 mm	Thermoconductive Tube Rack for 12 x 15mL Centrifuge Tubes, with insulative exterior†	BCS-232	1	1 ^Δ	up to 2 ^Δ
15 mL centrifuge tubes	9	Cylindrical	8.9 x 8.9 x 10.7 cm	17.1 mm	106.7 mm	Thermoconductive Tube Rack for 9 x 15mL Centrifuge Tubes	BCS-153	1	1 ^Δ	up to 2 ^Δ
50 mL centrifuge tubes	4	Cylindrical	8.9 x 8.9 x 10.7 cm	29.5 mm	101.6 mm	Thermoconductive Tube Rack for 4 x 50mL Centrifuge Tubes	BCS-154	1	1 ^Δ	up to 2 ^Δ
250 mL centrifuge tube	1	Conical	8.9 x 8.9 x 14.0 cm	60.4 mm	133.3 mm	Thermoconductive Tube Rack for 1 x 250mL Centrifuge Tube	BCS-532	1	1 ^{**}	up to 2 ^{**}
250 mL centrifuge tube	1	Cylindrical	8.9 x 8.9 x 7.2 cm	73.6 mm	66. mm	n/a	BCS-533	1	1	up to 2
Blood collection tube modules: Thermoconductive Tube Racks for Blood Tubes										
13 mm or 16 mm blood tubes	12	Cylindrical	13.7 x 9.5 x 9.6 cm	16.6 mm	83.3 mm	Thermoconductive Tube Rack for 12 x 13mm or 16mm Blood Tubes, with insulative exterior†	BCS-235	1	1	up to 2 ^Δ
13x75 mm blood tubes	9	Cylindrical	8.9 x 8.9 x 6.1 cm	13.0 mm	61.0 mm	Thermoconductive Tube Rack for 9 13x75mm Blood Tubes	BCS-157	1	1 ^Δ	up to 2 ^Δ
13x100 mm blood tubes or 5 mL cryogenic vials	9	Cylindrical	8.9 x 8.9 x 8.4 cm	13.0 mm	83.8 mm	Thermoconductive Tube Rack for 9 13x100mm Blood Tubes	BCS-155	1	1 ^Δ	up to 2 ^Δ
16x100 mm blood tubes	9	Cylindrical	8.9 x 8.9 x 8.4 cm	16.0 mm	83.8 mm	Thermoconductive Tube Rack for 9 16x100mm Blood Tubes	BCS-156	1	1 ^Δ	up to 2 ^Δ

† Thermo-conductive base and insulative exterior Δ Requires extension collar accessory for closed lid cooling ** Lid closure not possible even with the addition of extension collar

Choose Your CoolBox Ice-Free Cooling Workstation System



	Cooling Workstation Open Platform, Single Capacity	Cooling Workstation, Single Capacity & Cooling Workstation, Double Capacity
Holds Tubes	Yes	Yes
Holds Plates	Yes	Yes
0.5° - 4°C cooling with lid open	4 hours	10 hours
0.5° - 4°C cooling with lid closed	n/a	16 hours
<0°C freezing with lid open	n/a	5 hours
<0°C freezing with lid closed	n/a	8 hours



AZENTA
LIFE SCIENCES

Cooling Workstation Open Platform, Single Capacity



An open-platform ice-free cooler that accommodates most Thermoconductive Tube Racks and Thermoconductive Sink modules. Low profile and small footprint make it ideal for use in the hood, keeping samples cold (0.5° to 4.0°C) up to four hours. 1°C to 8°C temperature indicator provides visual assurance of temperature performance. To extend the cooling duration, keep an additional Cooling Workstation Cooling Core in the freezer and rotate the Cores as needed.

Ordering Information

BCS-504	Cooling Workstation System , single capacity open platform, cooling core included, purple , 1 system
BCS-513	Cooling Workstation , single capacity, open platform holder, purple , 1 holder
BCS-511	Cooling Workstation Cooling Core , 0.5°C to 4°C, blue

Cooling Workstation, Single Capacity or Double Capacity System

Keep sample tubes or plates cold for over 16 hours with the lid on, and over 10 hours with the lid off. Use optional Cooling Workstation Freezing Core to maintain frozen samples for over 8 hours. Dry ice may be used in place of the cores to create a compact snap freezing workstation.



Cooling Workstation, Single Capacity

Includes: Cooling Workstation, Single Capacity base, collar, lid and (1) Cooling Workstation Cooling Core for 0.5° to 4°C cooling.



Cooling Workstation, Double Capacity

Includes: Cooling Workstation, Double Capacity base, collar, lid and (2) Cooling Workstation Cooling Core for 0.5° to 4°C cooling.

Ordering Information

BCS-502	Cooling Workstation , single capacity, cooling core included, purple
BCS-502G	Cooling Workstation , single capacity, cooling core included, green
BCS-502O	Cooling Workstation , single capacity, cooling core included, orange
BCS-502PK	Cooling Workstation , single capacity, cooling core included, pink
BCS-502-F	Cooling Workstation , single capacity, freezing core included, purple

* Internal height of open space when core is in the base.

Ordering Information

BCS-503	Cooling Workstation , double capacity, cooling core included, purple
BCS-503G	Cooling Workstation , double capacity, cooling core included, green
BCS-503O	Cooling Workstation , double capacity, cooling core included, orange
BCS-503PK	Cooling Workstation , double capacity, cooling core included, pink
BCS-503-F	Cooling Workstation , double capacity, freezing core included, purple

* Internal height of open space when core is in the base.



Optional Accessories




Cooling Workstation Single Capacity and Cooling Workstation Double Capacity Extension Collar

For use with Cooling Workstation Single Capacity and Cooling Workstation Double Capacity systems to accommodate tall tube modules. The collar is magnetized and easily adheres to the unit base.

Ordering Information

BCS-502-C	Cooling Workstation Extension Collar, for Cooling Workstation, purple
BCS-502-CG	Cooling Workstation Extension Collar, for Cooling Workstation, green
BCS-502-CO	Cooling Workstation Extension Collar, for Cooling Workstation, orange
BCS-502-CPK	Cooling Workstation Extension Collar, for Cooling Workstation, pink
BCS-503-C	Cooling Workstation Extension Collar, for Cooling Workstation Large, purple
BCS-503-CG	Cooling Workstation Extension Collar, for Cooling Workstation Large, green
BCS-503-CO	Cooling Workstation Extension Collar, for Cooling Workstation Large, orange
BCS-503-CPK	Cooling Workstation Extension Collar, for Cooling Workstation Large, pink



Cooling Workstation Cores

Keep additional cooling or freezing cores in the freezer for flexibility and extended duration. Cooling Workstation Cooling Core features a 1 to 8°C temperature indicator. Both cooling and freezing cores feature a thermo-conductive surface for uniform temperature distribution.

Ordering Information

BCS-511	Cooling Workstation Cooling Core, 0.5°C to 4°C, blue
BCS-512	Cooling Workstation Freezing Core, below 0°C, blue

Popular Pre-assembled Configurations

Cooling Workstation Open Platform, Single Capacity PCR Cooling Systems, pre-assembled



Ordering Information

BCS-556	Cooling Workstation System, pre-assembled open-platform, for use with PCR plates, includes 1 x BCS-504 (Cooling Workstation System) and 1 x BCS-529 (Thermoconductive Tube Rack), purple
BCS-557	Cooling Workstation System, pre-assembled open-platform, for use with PCR strip wells, includes 1 x BCS-504 (Cooling Workstation System) and 1 x BCS-523 (Thermoconductive Tube Rack), purple

Cooling Workstation Single Capacity and Cooling Workstation Double Capacity Systems, pre-assembled



Ordering Information

BCS-576	Cooling Workstation System, pre-assembled for use with 24 microtubes, includes 1 x BCS-502 (Cooling Workstation) and 1 x BCS-535 (Thermoconductive Tube Rack), purple
BCS-575	Cooling Workstation System, pre-assembled for use with 24 cryo tubes, includes 1 x BCS-502 (Cooling Workstation) and 1 x BCS-534 (Thermoconductive Tube Rack), purple
BCS-570	Cooling Workstation System, pre-assembled for use with PCR plates, includes 1 x BCS-502 (Cooling Workstation) and 1 x BCS-529 (Thermoconductive Tube Rack), purple
BCS-572	Cooling Workstation System, pre-assembled for use with PCR strip wells, includes 1 x BCS-502 (Cooling Workstation) and 1 x BCS-523 (Thermoconductive Tube Rack), purple
BCS-573	Cooling Workstation System, pre-assembled for use with PCR plates, includes 1 x BCS-503 (Cooling Workstation), 1 x BCS-529 and 1 x BCS-535 (Thermoconductive Tube Rack), purple



Thermoconductive Tube Racks



AZENTA
LIFE SCIENCES

Thermoconductive Tube Racks



Thermoconductive tube modules eliminate variability which originates from tubes placed directly into ice, dry ice, alcohol baths, water baths and other temperature sources. Place the Thermoconductive Tube Rack module directly onto a temperature source between -196°C to $>100^{\circ}\text{C}$ and it will rapidly adapt to that temperature. Thermoconductive Tube Rack modules ensure $\pm 0.1^{\circ}\text{C}$ temperature uniformity across all tubes when cooling, snap freezing, heating or thawing. Suggested applications include cooling reagents such as restriction enzymes, dNTPs and antibodies, alcohol-free dry ice snap freezing of tissue, virus and bacteria samples and bench top cryogenic tube sorting in liquid nitrogen. All Thermoconductive Tube Rack modules may be autoclaved, high heat gamma irradiated or decontaminated with bleach, alcohol or other disinfectants or lab detergents. Certain Thermoconductive Tube Rack modules are SBS-compatible.

Problem: Samples in Ice

- Non-uniform ice contact results in variable sample temperature
- Disorganized samples, wet labels
- Shifting, sinking tubes; contamination risk
- Non-reproducible method




Solution: Samples in Thermoconductive Tube Rack Module

- All samples $<4^{\circ}\text{C}$ and uniform in temperature ($\pm 0.1^{\circ}\text{C}$)
- Samples organized, secure and dry
- All tubes upright and indexed
- Reproducible method




How It Works

Cooling




Thermoconductive Tube Rack on Ice: Heat from the relatively warmer Thermoconductive Tube Rack module is transferred to cooling source (wet or dry ice, cartridge, LN_2) until equilibrium is reached.

Heating



Thermoconductive Tube Rack in Water Bath: Heat is transferred from water bath toward relatively cooler Thermoconductive Tube Rack until equilibrium is reached.




Thermoconductive Tube Racks and Sinks are precision-engineered sample modules manufactured from a novel thermo-conductive alloy material. Thermoconductivity is the transfer of heat energy from a higher temperature region to a lower temperature region. Tube Rack modules evenly distribute the temperature across all wells providing very uniform and consistent temperature to all samples ($\pm 0.1^{\circ}\text{C}$).



AZENTA
LIFE SCIENCES

Thermoconductive Tube Racks

Thermoconductive Tube Rack Reproducibility




Performance test: A temperature probe was placed into a 2.0 mL cryogenic vial containing 1.0 mL of water. The tube was inserted into a Thermoconductive Tube Rack for 45 Cryo or FACS Tubes module. The module was placed onto a Thermoconductive Tray platform in a 37°C water bath and allowed to equilibrate. The Tube Rack for 45 Cryo or FACS Tubes module was then removed and placed onto dry ice and equilibrated to -78°C (0 - 15 minutes) and then returned to the water bath to re-equilibrate to 37°C (15 - 30 minutes). This experiment was repeated five consecutive times and temperature profiles were recorded.

Conclusion: The Thermoconductive Tube Rack for 45 Cryo or FACS Tubes module showed identical cooling profiles and phase transition (orange circles) over five consecutive freeze-thaw cycles.

Thermoconductive Tube Rack Versatility and Performance


On Ice

- Adapts from ambient (25°C) to <4°C in 60-90 seconds*
- Samples and labels stay dry, organized
- Hours of ice cooling without direct ice contact
- Reproducible method




On Dry Ice

- Adapts from ambient (25°C) to -78°C in approximately 5-7 minutes*
- Eliminates ethanol from snap freezing
- Samples are upright and organized as they freeze
- Equal or better freezing rate as compared to direct immersion into dry ice or alcohol slurry
- Reproducible method



In Liquid Nitrogen (LN₂)

- Adapts from ambient (25°C) to approximately -150°C in approximately 12-14 minutes*
- Vapor barrier protects from ambient air
- Samples are upright and organized as they freeze
- No direct contact between samples and LN₂
- Reproducible method



With Heat Sources

- Use with water baths, hot plates, incubators and other heat sources to keep samples warm



* Average cooling rate from room temperature



AZENTA
LIFE SCIENCES

Thermoconductive Tube Rack Modules

Thermoconductive Tube Racks for Microcentrifuge Tubes



Ordering Information

BCS-163	Thermoconductive Tube Rack, holds 6 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, gray
BCS-165	Thermoconductive Tube Rack, holds 6 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, orange
BCS-164	Thermoconductive Tube Rack, holds 6 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, green
BCS-125	Thermoconductive Tube Rack, holds 15 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, gray
BCS-125G	Thermoconductive Tube Rack, holds 15 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, green
BCS-125O	Thermoconductive Tube Rack, holds 15 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, orange
BCS-127	Thermoconductive Tube Rack, holds 15 x 1.5 conical tubes, conical wells, gray
BCS-535	Thermoconductive Tube Rack, holds 24 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, SBS compatible, gray*
BCS-539	Thermoconductive Tube Rack, holds 12 x 5ml microcentrifuge tubes, conical wells, SBS compatible, gray*
BCS-108	Thermoconductive Tube Rack, holds 30 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, gray
BCS-108G	Thermoconductive Tube Rack, holds 30 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, green
BCS-108O	Thermoconductive Tube Rack, holds 30 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, orange
BCS-128	Thermoconductive Tube Rack, holds 30 x 1.5 conical tubes, conical wells, gray
BCS-137	Thermoconductive Tube Rack, holds 30 x 500ul microcentrifuge tubes, conical wells, gray
BCS-102	Thermoconductive Tube Rack, holds 90 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, gray
BCS-116	Thermoconductive Tube Rack, holds 96 x 1.5 or 2ml microcentrifuge tubes, cylindrical wells, SBS compatible, row and column indexing, gray**

* SBS-compatible ** Thermoconductive Tube Rack for 96 Microcentrifuge Tubes has A-H and 1-12 row and column indexing

Thermoconductive Tube Racks for Cryo or FACS Tubes



Ordering Information

BCS-126	Thermoconductive Tube Rack, holds 15 cryo tubes or FACS tube modules, cylindrical wells, gray
BCS-534	Thermoconductive Tube Rack, holds 24 cryo tubes or FACS tube modules, cylindrical "gripping" wells for one-hand opening/closing vials, SBS compatible, gray*
BCS-138	Thermoconductive Tube Rack, holds 30 cryo tubes or FACS tube modules, cylindrical "gripping" wells for one-hand opening/closing vials, gray†
BCS-105	Thermoconductive Tube Rack, holds 45 cryo tubes or FACS tube modules, cylindrical wells, gray

* SBS-compatible † "gripping" wells for one-hand vial opening/closing



Thermoconductive Tube Racks

Thermoconductive Tube Racks for
PCR Plate, Strip Well or Tubes



Ordering Information

BCS-529	Thermoconductive Tube Rack, holds one 96-well PCR plate, 12 x strip wells or 96 tubes, tapered wells, SBS compatible, gray*
BCS-523	Thermoconductive Tube Rack, holds 6 strips wells and 12 x 1.5 or 2ml microcentrifuge tubes, 48 tapered wells for strips and 12 cylindrical wells, SBS compatible, gray*
BCS-538	Thermoconductive Tube Rack, holds one 384-well PCR plate, tapered wells, SBS compatible, gray*

* SBS-compatible

Thermoconductive Tube Racks for
96-Well 2D Coded Storage Tubes



Ordering Information

BCS-231	Thermoconductive Tube Rack, holds 96 x 0.5ml 2D storage tubes, cylindrical wells, gray
BCS-149	Thermoconductive Tube Rack, holds 96 x 1ml 2D storage tubes, cylindrical wells, gray

Thermoconductive Tube Racks for
Cell Therapy Injectable Ampules



Ordering Information

BCS-266	Thermoconductive Tube Rack, holds 12 x 10ml injectable cell therapy ampules, cylindrical wells, gray
BCS-265	Thermoconductive Tube Rack, holds 12 x 2ml injectable cell therapy ampules, cylindrical wells, gray

Tall Tube Modules

Thermoconductive Tube Racks for 15mL, 50mL and 250mL Centrifuge Tubes



Ordering Information

BCS-232	Thermoconductive Tube Rack, holds 12 x 15ml centrifuge tubes, cylindrical wells, with thermoconductive base and insulative exterior, purple*
BCS-153	Thermoconductive Tube Rack, holds 9 x 15ml centrifuge tubes, cylindrical wells, gray
BCS-154	Thermoconductive Tube Rack, holds 4 x 50ml centrifuge tubes, cylindrical wells, gray
BCS-532	Thermoconductive Tube Rack, holds one 250ml centrifuge tube, conical well, gray
BCS-533	Thermoconductive Tube Rack, holds one 250ml centrifuge tube, cylindrical well, gray

* Thermo-conductive base and insulative exterior

Thermoconductive Tube Rack for Blood Collection Tubes



Ordering Information

BCS-235	Thermoconductive Tube Rack, holds 12 x13mm or 16mm blood tubes, cylindrical wells, with thermoconductive base and insulative exterior, purple
BCS-157	Thermoconductive Tube Rack, holds 9 13x75mm blood tubes, cylindrical wells, gray
BCS-155	Thermoconductive Tube Rack, holds 9 13x100mm blood tubes, cylindrical wells, gray
BCS-156	Thermoconductive Tube Rack, holds 9 16x100mm blood tubes, cylindrical wells, gray

* Thermo-conductive base and insulative exterior

Thermoconductive Sinks




Thermoconductive Sinks

Thermoconductive plate and reservoir modules provide uniform temperature to all wells, regardless of position. When placed onto a temperature source such as ice, dry ice, liquid nitrogen or water baths, the Thermoconductive Sink module will rapidly adapt to that temperature - from -196°C to >+100°C. Sink modules ensure temperature sample uniformity when cooling, snap freezing, heating or thawing samples. All Thermoconductive Sink modules may be autoclaved, high heat gamma irradiated or decontaminated with bleach, alcohol or other disinfectants or lab detergents. All modules are compatible with all temperature sources.

Problem: Non-Uniform Plate Cooling with Crushed Ice

Final equilibrium well temperature for a 96-well flat bottom plate in direct contact with crushed ice.



Solution: Uniform Plate Cooling with Thermoconductive Sink for use with Flat Bottom Plates Module

Final equilibrium well temperature for a 96-well flat bottom plate in direct contact with crushed ice.




Plate and Reservoir Modules

Thermoconductive Sink, SBS-Compatible Plate Modules



Thermoconductive Sink for use with 55mL Reagent Reservoirs



Ordering Information

BCS-536	Thermoconductive Sink, for use with 6-, 12-, 24-, 48-, 96-, 384-well flat bottom plates, SBS compatible, gray
BCS-537	Thermoconductive Sink, for use with one 96-well U-bottom plate, SBS compatible, gray

Ordering Information

BCS-184	Thermoconductive Sink, for use with 55ml reagent reservoirs, gray
----------------	---

Accessories



Accessories





Ice Pans

Non-toxic, recyclable ethyl-vinyl acetate (EVA) foam containers for use with ice, dry ice, liquid nitrogen, alcohol slurries. Will not sweat, leak or skid on bench.



Ordering Information

	Ice Pan without Lid, Square 1L
BCS-211PL	Ice Pan, without Lid, square, 1L, purple
BCS-211GR	Ice Pan, without Lid, square, 1L, lime green
BCS-211OR	Ice Pan, without Lid, square, 1L, orange
BCS-211PK	Ice Pan, without Lid, square, 1L, pink
BCS-211B	Ice Pan, without Lid, square, 1L, blue
BCS-211	Ice Pan, without Lid, square, 1L, green
BCS-212	Ice Pan, without Lid, square, 1L, red
	Ice Pan without Lid, Rectangle 4L
BCS-113PL	Ice Pan, without Lid, rectangle, 4L, purple
BCS-113GR	Ice Pan, without Lid, rectangle, 4L, lime green
BCS-113OR	Ice Pan, without Lid, rectangle, 4L, orange
BCS-113PK	Ice Pan, without Lid, rectangle, 4L, pink
BCS-113B	Ice Pan, without Lid, rectangle, 4L, blue
BCS-113	Ice Pan, without Lid, rectangle, 4L, green
BCS-114	Ice Pan, without Lid, rectangle, 4L, red
	Ice Pan with Lid, Rectangle 4L
BCS-117PL	Ice Pan, with Lid, rectangle, 4L, purple
BCS-117GR	Ice Pan, with Lid, rectangle, 4L, lime green
BCS-117OR	Ice Pan, with Lid, rectangle, 4L, orange
BCS-117PK	Ice Pan, with Lid, rectangle, 4L, pink
BCS-117B	Ice Pan, with Lid, rectangle, 4L, blue

	Ice Pan without Lid, Rectangle 9L
BCS-111PL	Ice Pan, without Lid, rectangle, 9L, purple
BCS-111GR	Ice Pan, without Lid, rectangle, 9L, lime green
BCS-111OR	Ice Pan, without Lid, rectangle, 9L, orange
BCS-111PK	Ice Pan, without Lid, rectangle, 9L, pink
BCS-111B	Ice Pan, without Lid, rectangle, 9L, blue
BCS-111	Ice Pan, without Lid, rectangle, 9L, green
BCS-112	Ice Pan, without Lid, rectangle, 9L, red
	Ice Pan with Lid, Rectangle 9L
BCS-118PL	Ice Pan, with Lid, rectangle, 9L, purple
BCS-118GR	Ice Pan, with Lid, rectangle, 9L, lime green
BCS-118B	Ice Pan, with Lid, rectangle, 9L, blue
	Ice Pan with Lid, Round 2.5L
BCS-115-25PL	Ice Bucket, with Lid, round, 2.5L, purple
BCS-115-25GR	Ice Bucket, with Lid, round, 2.5L, lime green
BCS-115-25OR	Ice Bucket, with Lid, round, 2.5L, orange
BCS-115-25PK	Ice Bucket, with Lid, round, 2.5L, pink
BCS-115-25B	Ice Bucket, with Lid, round, 2.5L, blue
BCS-115-25G	Ice Bucket, with Lid, round, 2.5L, green
BCS-115-25R	Ice Bucket, with Lid, round, 2.5L, red
	Ice Pan with Lid, Round 4L
BCS-115PL	Ice Bucket, with Lid, round, 4L, purple
BCS-115GR	Ice Bucket, with Lid, round, 4L, lime green
BCS-115OR	Ice Bucket, with Lid, round, 4L, orange
BCS-115PK	Ice Bucket, with Lid, round, 4L, pink
BCS-115B	Ice Bucket, with Lid, round, 4L, blue
BCS-115	Ice Bucket, with Lid, round, 4L, green
BCS-115R	Ice Bucket, with Lid, round, 4L, red



Hinged CryoBoxes

Patented hinged lid offers convenience and archival integrity, ensuring markings and vials remain in sync. Lid stays attached to base minimizing risk of separation and lid contamination. Lid is easy to open when frozen. Available in 9x9, 10x10, and vapor phase LN₂ compatible formats. Plastic 81-place grid has adjustable slats to accommodate multiple vial types. 2-inch box holds 1.0 mL or 2.0 mL cryogenic vials and microcentrifuge tubes. 3.5-inch box holds 3.0 mL to 5.0 mL cryogenic vials.

Hinged cryoboxes can be customized to suit various requirements. Options include new colors, logos, designs, grid sizes and additional components.



Ordering Information

Hinged CryoBox 2 Inch, 81-Place	
BCS-206	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, white , 5 per case
BCS-206B	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, blue , 5 per case
BCS-206G	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, green , 5 per case
BCS-206O	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, orange , 5 per case
BCS-206P	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, purple , 5 per case
BCS-206PK	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, pink , 5 per case
BCS-206MC	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, multipack, no white , 5 per case
BCS-207	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, white , 50 per case
BCS-207B	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, blue , 50 per case
BCS-207G	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, green , 50 per case
BCS-207O	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, orange , 50 per case
BCS-207P	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, purple , 50 per case
BCS-207PK	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, pink , 50 per case

Hinged CryoBox 2 Inch, 100-Place	
BCS-209G	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 100-place, green , 5 per case
BCS-209P	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 100-place, purple , 5 per case
BCS-220G	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 100-place, green , 50 per case
BCS-220P	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 100-place, purple , 50 per case
Hinged CryoBox 3.5in, 81-Place	
BCS-215G	Cryobox, 3.5-inch cryobox, hinged, adjustable plastic grid, 81-place, green , 6 per case
BCS-215P	Cryobox, 3.5-inch cryobox, hinged, adjustable plastic grid, 81-place, purple , 6 per case
BCS-219G	Cryobox, 3.5-inch cryobox, hinged, adjustable plastic grid, 81-place, green , 30 per case
BCS-219P	Cryobox, 3.5-inch cryobox, hinged, adjustable plastic grid, 81-place, purple , 30 per case
Hinged CryoBox 2 Inch, 81-Place, With Drain Holes	
BCS-217G	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, LN ₂ drain holes, green , 5 per case
BCS-217P	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, LN ₂ drain holes, purple , 5 per case
BCS-221G	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, LN ₂ drain holes, green , 50 per case
BCS-221P	Cryobox, 2-inch cryobox, hinged, adjustable plastic grid, 81-place, LN ₂ drain holes, purple , 50 per case



Thermoconductive Tube Module Temperature Strips



Adhesive temperature display shows the temperature of a surface with 1°C resolution. Ideal for placement on Thermoconductive Tube Rack, Sink and Tray modules.

Ordering Information

BCS-143	Thermoconductive Tube Module Temperature Strips, 1°-8°C, 3pk.
---------	---

Thermoconductive Tube Module Sleeves



Ordering Information

BCS-205	Thermoconductive Tube Module Sleeves, 4pk
---------	---

Cryo Tube Grippers



Cryo Tube Grippers feature a unique design to grasp internal- or external-thread cryogenic vials. Easily sort or move vials while preventing contamination and protecting fingers from frozen vials, dry ice and liquid nitrogen. 5 per pack.

Ordering Information

BCS-213MC	Cryo Tube Grippers, multi-color, 5 per case
-----------	---

Cryo Tube Locking Racks



Cryo Tube Locking Racks feature a locking mechanism that allow one-hand opening for self-standing cryogenic vials. Accommodates both round bottom and self-standing vial formats. Racks have A - J and 1 - 5 row and column indexing for easy organization. Autoclavable. 5 per pack.

Ordering Information

BCS-222	Cryo Tube Locking Racks, multi-pack, 5 per case
---------	---

1D-coded Cryo Tubes

Leak-proof, auto-cap cryogenic tubes are ideal for cell culture and biobanking. The screw cap features a co-molded thermally-fused gasket which prevents leaking, slipping and risk of contamination. The gasket is 95kPa certified to provide a leak-proof seal. The star socket on cap top is compatible with auto-decapping equipment. Each vial is individually barcoded with a unique identifier that can be read with common barcode readers. Recommended for storage down to vapor phase liquid nitrogen but not suitable for use directly in LN₂. 500 per case.



Internal Threads

External Threads

Ordering Information

1.0ml - 5ml 1D-coded Cryo Tube, Internal Thread	
BCS-2510	1ml 1D-coded Cryo Tube, Internal Thread, self-standing, 500 tubes per case
BCS-2511	2ml 1D-coded Cryo Tube, Internal Thread, self-standing, 500 tubes per case
BCS-2512	2ml 1D-coded Cryo Tube, Internal Thread, round-bottom, 500 tubes per case
BCS-2513	4ml 1D-coded Cryo Tube, Internal Thread, round-bottom, 500 tubes per case
BCS-2514	4ml 1D-coded Cryo Tube, Internal Thread, self-standing, 500 tubes per case
BCS-2515	5ml 1D-coded Cryo Tube, Internal Thread, round-bottom, 500 tubes per case
BCS-2516	5ml 1D-coded Cryo Tube, Internal Thread, self-standing, 500 tubes per case

1.0ml - 5ml 1D-coded Cryo Tube, External Thread	
BCS-2517	1ml 1D-coded Cryo Tube, Internal Thread, self-standing, 500 tubes per case
BCS-2501	2ml 1D-coded Cryo Tube, External Thread, round-bottom, 500 tubes per case
BCS-2502	2ml 1D-coded Cryo Tube, External Thread, self-standing, 500 tubes per case
BCS-2503	3ml 1D-coded Cryo Tube, External Thread, self-standing, 500 tubes per case
BCS-2504	4ml 1D-coded Cryo Tube, External Thread, self-standing, 500 tubes per case
BCS-2505	5ml 1D-coded Cryo Tube, External Thread, self-standing, 500 tubes per case

Cryo Tubes Cap Inserts

Inserts for auto caps. 1,000 per pack.

Caps designed to color code tubes. Ideal for labeling different specimen tubes and cataloging sample inventory.

Ordering Information

BCS-2436	Cryo Tube Cap Insert, violet, 1000 inserts per case
BCS-2432	Cryo Tube Cap Insert, pink, 1000 inserts per case
BCS-2431	Cryo Tube Cap Insert, green, 1000 inserts per case
BCS-2434	Cryo Tube Cap Insert, yellow, 1000 inserts per case
BCS-2435	Cryo Tube Cap Insert, white, 1000 inserts per case
BCS-2433	Cryo Tube Cap Insert, red, 1000 inserts per case
BCS-2438	Cryo Tube Cap Insert, gray, 1000 inserts per case
BCS-2430	Cryo Tube Cap Insert, blue, 1000 inserts per case
BCS-2437	Cryo Tube Cap Insert, orange, 1000 inserts per case



Part Number Index

Part No.	Page
4ti-0110	169, 171, 240
4ti-0116	171
4ti-0117	171
4ti-0120	169, 171, 240
4ti-0124	170, 240
4ti-0125	170, 240
4ti-0126	168
4ti-0130	169
4ti-0131	172
4ti-0132	168
4ti-0133	172
4ti-0135	169, 240
4ti-0136	168
4ti-0137	168, 240
4ti-0138	169, 171, 240
4ti-0139	167, 240
4ti-0147	167
4ti-0150	174
4ti-0151	173
4ti-0152	175
4ti-0201	177
4ti-0203	177
4ti-0204	177
4ti-0205	177
4ti-0206	177
4ti-0214	180
4ti-0221	178
4ti-0223	178
4ti-0224	178
4ti-0225	178
4ti-0226	178
4ti-0234	181
4ti-0241	179
4ti-0243	179
4ti-0244	179

Part No.	Page
4ti-0245	179
4ti-0246	179
4ti-0254	185
4ti-0262	182
4ti-0263	183
4ti-0264	185
4ti-0273	183
4ti-0274	185
4ti-0280	177, 180, 185, 238
4ti-0281	177, 185, 238
4ti-0282	178, 183, 238
4ti-0283	178, 183, 238
4ti-0284	179, 182, 238
4ti-0285	238
4ti-0286	179, 182, 238
4ti-0287	110, 238
4ti-0288	238
4ti-0289	118, 238
4ti-0290	140, 181, 238
4ti-0291	239
4ti-0292	139
4ti-0370	139
4ti-0371	139
4ti-0372	140
4ti-0373	132, 140
4ti-0374	156
4ti-0380	109
4ti-0380/C	109
4ti-0381	109
4ti-0382	109
4ti-0383	109
4ti-0384	107, 162, 164
4ti-0384/RIG	107
4ti-0385	107
4ti-0386	107

Part No.	Page
4ti-0387	107
4ti-0391	230
4ti-0398	188, 189, 195, 201, 202, 203, 204, 231
4ti-0500	218, 232
4ti-0500FL	222, 232
4ti-0502	218, 219, 221, 222, 223, 224, 225, 226, 227, 228, 229, 230
4ti-0503	230
4ti-0510	228, 232
4ti-0512	223, 232
4ti-0516	177, 178, 179, 182, 183, 184, 226, 232
4ti-0517	227, 232
4ti-0519	229, 232
4ti-0520	192, 206
4ti-0520S	192, 206
4ti-0521	135, 192, 194, 206
4ti-0522	135, 192, 206
4ti-0522S	192, 206
4ti-0523	193, 206
4ti-0523S	193, 206
4ti-0524	193, 206
4ti-0524S	193, 206
4ti-0530	196, 206, 229
4ti-0530S	196, 206
4ti-0531	135, 196, 197, 206
4ti-0532	135, 196, 206
4ti-0532S	196, 206
4ti-0535	200, 206
4ti-0535S	200, 206
4ti-0536	200, 206
4ti-0537	200, 206
4ti-0537S	200, 206
4ti-0538	198, 206



Part Number Index

Part No.	Page
4ti-0538S.....	198, 206
4ti-0539.....	135, 198, 199, 206
4ti-0539S.....	198, 206
4ti-0540.....	188, 204, 206
4ti-0540S.....	188, 206
4ti-0541.....	188, 204, 206
4ti-0542.....	188, 206
4ti-0542S.....	188, 206
4ti-0545.....	201, 206
4ti-0545S.....	201, 206
4ti-0546.....	201, 206
4ti-0547.....	201, 206
4ti-0548.....	190, 206
4ti-0549.....	190, 206
4ti-0550.....	221, 232
4ti-0560.....	105, 108, 109, 113, 114, 220, 231, 232
4TI-0560.....	113, 114, 148
4ti-0561.....	220, 232
4ti-0561S.....	220
4ti-0563.....	231
4ti-0565.....	110, 115, 116, 117, 219, 232
4TI-0565.....	111, 113
4ti-0566.....	224, 225, 232
4ti-0573.....	189, 206
4ti-0573S.....	189, 206
4ti-0574.....	189, 206
4ti-0574S.....	189, 206
4ti-0575.....	189, 206
4ti-0580.....	191, 206
4ti-0580S.....	191, 206
4ti-0581.....	191, 206
4ti-0582.....	191, 206
4ti-0582S.....	191, 206
4ti-0585.....	195, 206
4ti-0585S.....	195, 206
4ti-0586.....	195, 206
4ti-0586S.....	195, 206
4ti-0587.....	195, 206

Part No.	Page
4ti-0590.....	202, 206
4ti-0590S.....	202, 206
4ti-0591.....	202, 206
4ti-0592.....	202, 206
4ti-0592S.....	202, 206
4ti-0597.....	203, 206
4ti-0598.....	203, 206
4ti-0598S.....	203, 206
4ti-0599.....	203, 206
4ti-0599S.....	203, 206
4ti-0640.....	208
4ti-0641.....	208
4ti-0642.....	208
4TI-0656.....	115
4ti-0665.....	210
4ti-0665-1.....	210
4ti-0665-2.....	210
4ti-0665-3.....	210
4ti-0665-4.....	210
4ti-0665-5.....	210
4ti-0665-6.....	210
4ti-0665-7.....	210
4ti-0665-8.....	210
4ti-0665-41.....	210
4ti-0665-56.....	210
4ti-0665-57.....	210
4ti-0665-58.....	210
4ti-0665-59.....	210
4ti-0665-60.....	210
4ti-0680.....	92
4ti-0680-1.....	92
4ti-0681.....	92
4ti-0683.....	92
4ti-0684.....	92
4ti-0685.....	92
4ti-0686.....	92
4ti-0688.....	92
4ti-0689.....	92
4ti-0710.....	105, 121, 162, 164
4ti-0711.....	121

Part No.	Page
4ti-0720.....	122, 132, 162, 164
4ti-0721.....	122
4ti-0730.....	116, 117, 162, 164
4ti-0735.....	149, 162, 164
4ti-0736.....	149
4ti-0740.....	146, 162, 164
4ti-0741.....	146
4ti-0750.....	150, 152, 158, 159, 162, 164
4ti-0750/8/B.....	152
4ti-0750/8/G.....	152
4ti-0750/8/P.....	152
4ti-0750/8/R.....	152
4ti-0750/8/Y.....	152
4ti-0750/16/B.....	152
4ti-0750/16/G.....	152
4ti-0750/16/P.....	152
4ti-0750/16/R.....	152
4ti-0750/16/Y.....	152
4ti-0750/24/B.....	152
4ti-0750/24/G.....	152
4ti-0750/24/P.....	152
4ti-0750/24/R.....	152
4ti-0750/24/Y.....	152
4ti-0750/32/B.....	152
4ti-0750/32/G.....	152
4ti-0750/32/P.....	152
4ti-0750/32/R.....	152
4ti-0750/32/Y.....	152
4ti-0750/48/B.....	152
4ti-0750/48/G.....	152
4ti-0750/48/P.....	152
4ti-0750/48/R.....	152
4ti-0750/48/Y.....	152
4ti-0751.....	141, 236, 237
4ti-0752.....	236
4ti-0753.....	135, 141, 142, 162, 164, 229
4ti-0754.....	141
4ti-0755.....	237



Part Number Index

Part No.	Page
4ti-0757	142
4ti-0760	147, 162, 164
4ti-0761	147
4ti-0770	102, 119, 162, 164
4ti-0771	102, 119
4ti-0772	102, 119
4ti-0775	138
4ti-0778	241
4ti-0780	153
4ti-0781	153, 162, 164, 236
4ti-0782	153, 236
4ti-0783	153, 236
4ti-0784	153
4ti-0785	137, 162, 164
4ti-0786	138, 162, 164
4ti-0788	236
4ti-0789	138
4ti-0790	157
4ti-0790/2D	157
4ti-0791	156
4ti-0792	154, 162, 164
4ti-0793	154, 162, 164
4ti-0794	154, 162, 164
4ti-0795	157
4ti-0796	155, 162, 164
4ti-0900	120, 162, 164
4ti-0901	120
4ti-0910	115
4TI-0910/C	115
4ti-0911	115
4ti-0912	115
4ti-0950	113, 162, 164
4TI-0950/C	113
4ti-0950W-F	142
4ti-0951	113
4TI-0951	113
4ti-0952	113
4ti-0953	113
4ti-0954	114, 162, 164
4TI-0954	114

Part No.	Page
4ti-0955	148, 162, 164
4TI-0955	148
4ti-0960	110, 111, 128, 131, 135, 162, 164, 199
4TI-0960	111
4ti-0961	111
4ti-0966	111
4ti-0970	112, 129, 131
4ti-0975	131
4ti-1000	125, 162, 164
4ti-1001	125
4ti-1200	125, 135, 162, 164
4ti-1201	125
4ti-1300	126, 162, 164
4ti-1381	145, 162, 164
4ti-1384	144, 162, 164
4ti-1385	144
4ti-1387	144
4ti-1400	126, 135, 162, 164
4ti-05231	193, 206
4ti-05381	135, 198, 199, 206
4ti-05481	190, 206
4ti-LB0109	171
4ti-LB0125	170
4ti-LB0147	167
4ti-LB0384/C	107
4ti-LB0384/RIG	107
4ti-LB0770	118, 119
4ti-LB0960	110, 111
4ti-0X730	117
4ti-0X770C/SBC	118, 119
4ti-0X960	110, 111
6.09.661	86
6.09.663	86
6.09.664	86
10-5010	88
10-5020	88
42-1001	86
42-1003	86
46-2004-115V	86

Part No.	Page
46-2004-230V	86
46-5030	85
46-5051	85
46-5061	85
46-6001	78
46-6002-1	78
46-6002-2	78
46-6002-3	78
46-6002-4	78
46-6002-5	78
46-6002-6	78
46-6002-7	78
46-6002-8	78
46-6002-9	78
46-6002-10	78
46-6002-11	78
46-6002-12	78
46-6002-13	78
46-6002-14	78
46-6002-15	78
46-6002-16	78
46-6002-17	78
46-6002-18	78
46-6501	81
46-6502	81
46-6511	81
46-6512	81
46-6513	81
46-6521	81
46-6601	81
46-6602	81
46-6604	81
46-6605	81
46-6606	81
46-8010	83
46-8011	83
46-8012	83
46-8014	83, 156
46-8018	84
46-8112	83, 84



Part Number Index

Part No.	Page
46-9001.....	79
46-9008.....	79
46-9012.....	79
48-8013-18.....	84
48-9013-01.....	79
48-9013-02.....	79
48-9013-03.....	79
48-9013-04.....	79
48-9013-05.....	79
48-9013-06.....	79
59-1000.....	135, 211
59-1001.....	211
59-1002.....	211
59-1003.....	211
59-1004.....	135, 211
59-2000.....	134, 135, 214
59-2001.....	135, 214
59-2002.....	214
59-2003.....	214
59-2004.....	214
59-2005.....	134, 135, 214
59-2006.....	214
59-2007.....	214
59-2008.....	214
59-2009.....	214
65-7514.....	20, 31
65-7515.....	20, 31
65-7516.....	20, 31
65-7517.....	20, 31
65-7572.....	50
65-7573.....	50
65-7574.....	50
65-7575.....	50
65-7576.....	50
65-7577.....	50
65-7640.....	20, 30
65-7641.....	20, 30
65-7642.....	20, 30
65-7643.....	20, 30
65-7644.....	30

Part No.	Page
65-7645.....	30
65-7646.....	30
65-7647.....	30
65-7650.....	28, 29
65-7651.....	28, 29
65-7652.....	28, 29
65-7653.....	28, 29
65-7660.....	20, 28
65-7661.....	20, 28
65-7662.....	20, 28
65-7663.....	20, 28
65-7664.....	28
65-7665.....	28
65-7666.....	28
65-7667.....	28
65-9303.....	20, 32
65-9451.....	28, 30
65-9460.....	31
65-9801.....	53
65-54000.....	88
65-54001.....	88
65-54004.....	88
65-73000.....	52
65-73001.....	52
65-73002.....	52
65-73003.....	52
65-73004.....	52
65-74000.....	52
66-0196-01.....	13, 42
66-0196-02.....	13
66-0196-03.....	13
66-0700-00.....	46
66-0700-01.....	46
66-0700-02.....	46
66-0700-10.....	46
66-0700-11.....	46
66-0700-12.....	46
66-1000-00.....	46
66-1000-01.....	46
66-1000-02.....	46

Part No.	Page
66-1000-10.....	46
66-1000-11.....	46
66-1000-12.....	46
66-1001.....	205
66-1021.....	205
66-1800.....	13, 30
66-1801.....	13, 28
66-1802.....	13
66-1803.....	13
66-9302.....	20, 32
66-9401.....	32, 51, 53
66-9402.....	51
66-9455.....	32
66-9951.....	53
66-32033.....	47
66-32033-Y6.....	47
66-32034.....	47
66-32034-L.....	47
66-32034-Y6.....	47
66-32034-Y6-L.....	47
66-32040.....	47
66-32040-Y6.....	47
66-32041.....	47
66-32041-Y6.....	47
66-32041-Y6-L.....	47
66-32042.....	47
66-32042-L.....	47
66-32043.....	47
66-32043-L.....	47
66-32043-Y6.....	47
66-32043-Y6-L.....	47
66-32062.....	47
66-32062-Y6.....	47
66-32141.....	47
66-51003.....	47
66-51004.....	25, 46, 47
66-51014.....	47
66-51016.....	42, 46, 47
66-51017.....	47
66-51020.....	27



Part Number Index

Part No.	Page
66-51021	26
66-51022	34
66-51023	35
66-51025	33
66-51026	25, 41, 46
66-61002	39, 42, 46, 47
66-61003	40, 47
66-62318	20, 39
66-62318-Y6	20, 39
66-62319	20, 39
66-62319-Y6	20, 39
66-62325	20, 38
66-62325-Y6	20, 38
66-62326	20, 38
66-62326-Y6	20, 38
66-62330	20, 40
66-62330-Y6	20, 40
66-62345	20, 40
66-62345-Y6	20, 40
66-63100-Y1	51
66-63100-Y2	51
66-63100-Y3	51
66-63100-Y4	51
66-63100-Y5	51
66-63100-Y6	51
66-63100-Y8	51
66-63100-Y10	51
66-63100-Y11	51
66-63100-Y12	51
66-63100-Y13	51
67-0200-00	45
67-0203-00	45
67-0203-01	20, 45
67-0203-02	20, 45
67-0203-10	20, 45
67-0203-11	20, 45
67-0203-51	45
67-0753-00	20, 33
67-0753-02	20, 33
67-0753-10	20, 33

Part No.	Page
67-0753-12	20, 33
67-0755-00	20, 34
67-0755-01	20, 34
67-0755-10	20, 34
67-0755-11	20, 34
67-0757-00	20, 35
67-0757-01	20, 35
67-0757-10	20, 35
67-0757-11	20, 35
67-63111-10	53
67-63111-50	53
68-0300-20	37
68-0301-00	20
68-0301-01	20
68-0301-10	20
68-0301-11	20
68-0303-00	20, 37
68-0303-01	20, 37
68-0303-10	20, 37
68-0303-11	20, 37
68-0701-00	20, 41
68-0701-02	20, 41
68-0701-10	20, 41
68-0701-11	41
68-0701-12	20
68-0703-00	20, 25
68-0703-02	20, 25, 57
68-0703-10	20, 25
68-0703-11	25
68-0703-12	20, 25
68-0704-00	25
68-0704-02	25
68-0704-10	25
68-0704-12	25
68-0801-00	20, 26
68-0801-01	20, 26
68-0801-10	20, 26
68-0801-11	20, 26
68-0802-00	26
68-0802-01	26

Part No.	Page
68-0802-10	26
68-0802-11	26
68-1001-00	20, 42
68-1001-01	20, 42
68-1001-10	20, 42
68-1001-11	20, 42
68-1003-00	20, 27
68-1003-01	20, 26, 27
68-1003-10	20, 27
68-1003-11	20, 27
68-1004-00	27
68-1004-01	27
68-1004-10	27
68-1004-11	27
68-53100-Z1N	50
68-53100-Z2N	50
68-53100-Z3N	50
68-53100-Z4N	50
68-53100-Z5N	50
68-53100-Z6N	50
68-53100-Z8N	50
68-53100-Z10N	50
68-53100-Z11N	50
68-53100-Z12N	50
68-53100-Z13N	50
68-53111-10N	50, 53
68-53111-10X	53
68-53111-50N	50, 53
68-53111-50X	53
69-0200-11	20
96-0001	93
96-0002	93
96-0003	93
96-0004	93
97-0001	93
98-0001	93
98-0002	93
98-0003	93
98-0004	93
243354-001	98



Part Number Index

Part No.	Page	Part No.	Page	Part No.	Page
252885	98	BCS-116	264	BCS-206P	271
252886	98	BCS-117B	270	BCS-206PK	271
252888-001	98	BCS-117GR	270	BCS-207	271
252888-002	98	BCS-117OR	270	BCS-207B	271
252888-003	98	BCS-117PK	270	BCS-207G	271
252888-004	98	BCS-117PL	270	BCS-207O	271
252888-005	98	BCS-118B	270	BCS-207P	271
BCS-102	264	BCS-118GR	270	BCS-207PK	271
BCS-105	256, 264	BCS-118PL	270	BCS-209G	271
BCS-108	256, 264	BCS-125	256, 264	BCS-209P	271
BCS-108G	264	BCS-125G	264	BCS-210	252
BCS-108O	264	BCS-125O	264	BCS-211	270
BCS-111	270	BCS-126	256, 264	BCS-211B	270
BCS-111B	270	BCS-127	256, 264	BCS-211GR	270
BCS-111GR	270	BCS-128	256, 264	BCS-211OR	270
BCS-111OR	270	BCS-137	256, 264	BCS-211PK	270
BCS-111PK	270	BCS-138	256, 264	BCS-211PL	270
BCS-111PL	270	BCS-143	272	BCS-212	270
BCS-112	270	BCS-149	256, 265	BCS-213MC	272
BCS-113	270	BCS-153	257, 266	BCS-215G	271
BCS-113B	270	BCS-154	257, 266	BCS-215P	271
BCS-113GR	270	BCS-155	257, 266	BCS-217G	271
BCS-113OR	270	BCS-156	257, 266	BCS-217P	271
BCS-113PK	270	BCS-157	257, 266	BCS-219G	271
BCS-113PL	270	BCS-163	256, 264	BCS-219P	271
BCS-114	270	BCS-164	264	BCS-220G	271
BCS-115	270	BCS-165	264	BCS-220P	271
BCS-115-25B	270	BCS-170	251	BCS-221G	271
BCS-115-25G	270	BCS-170G	251	BCS-221P	271
BCS-115-25GR	270	BCS-170O	251	BCS-222	272
BCS-115-25OR	270	BCS-170PK	251	BCS-231	256, 265
BCS-115-25PK	270	BCS-172	251	BCS-232	257, 266
BCS-115-25PL	270	BCS-172CS	252	BCS-235	257, 266
BCS-115-25R	270	BCS-184	268	BCS-262	251
BCS-115B	270	BCS-205	272	BCS-262CS	252
BCS-115GR	270	BCS-206	271	BCS-265	256, 265
BCS-115OR	270	BCS-206B	271	BCS-266	256, 265
BCS-115PK	270	BCS-206G	271	BCS-405	250
BCS-115PL	270	BCS-206MC	271	BCS-405G	250
BCS-115R	270	BCS-206O	271	BCS-405MC	250



Part Number Index

Part No.	Page
BCS-4050	250
BCS-405PK	250
BCS-406	251
BCS-4070	250
BCS-407P	250
BCS-502	258, 260
BCS-502-C	259
BCS-502-CG	259
BCS-502-CO	259
BCS-502-CPK	259
BCS-502-F	258
BCS-502G	258
BCS-502O	258
BCS-502PK	258
BCS-503	258, 260
BCS-503-C	259
BCS-503-CG	259
BCS-503-CO	259
BCS-503-CPK	259
BCS-503-F	258
BCS-503G	258
BCS-503O	258
BCS-503PK	258
BCS-504	258, 260
BCS-511	258, 259
BCS-512	259
BCS-513	258
BCS-523	256, 260, 265
BCS-529	256, 260, 265
BCS-532	257, 266
BCS-533	257, 266
BCS-534	256, 260, 264
BCS-535	256, 260, 264
BCS-536	268

Part No.	Page
BCS-537	268
BCS-538	256, 265
BCS-539	256, 264
BCS-556	260
BCS-557	260
BCS-570	260
BCS-572	260
BCS-573	260
BCS-575	260
BCS-576	260
BCS-2430	273
BCS-2431	273
BCS-2432	273
BCS-2433	273
BCS-2434	273
BCS-2435	273
BCS-2436	273
BCS-2437	273
BCS-2438	273
BCS-2501	273
BCS-2502	273
BCS-2503	273
BCS-2504	273
BCS-2505	273
BCS-2510	273
BCS-2511	273
BCS-2512	273
BCS-2513	273
BCS-2514	273
BCS-2515	273
BCS-2516	273
BCS-2517	273
BCS-3105	252
BCS-3106	252

Part No.	Page
BCS-3107	252
DP5-C-STD	69
DP5-C-STD-80	69
DP5-CYS-C	69, 73
DP5-CYS-E	73
DP5-CYS-M	65, 73
DP5-E-NET	67
DP5-E-NET-80	67
DP5-E-STD	67
DP5-E-STD-80	67
DP5-M-LUX-80	70
DP5-M-LUX-80-EB	70
DP5-M-NET	65
DP5-M-NET-80	65
DP5-M-RIT-80	71
DP5-M-STD	65
DP5-M-STD-80	65
DP5-NET2	71
DP5-RIT-APL	71
DP5-RIT-TAG	71
XP-A	215
XP-A_100V	215
XP-A_230V	215
X-Tape_2000	215
ZTS-DUO-80	63
ZTS-HND2	73
ZTS-MHK-48	90
ZTS-MHK-96	90
ZTS-MHK-LOW	90
ZTS-UNO-80	63





AZENTA
LIFE SCIENCES

azenta.com

© 2025 Azenta US, Inc. All rights reserved. All trademarks are property of Azenta US, Inc. unless otherwise specified. 40001-CAT 0125