

Deltalab Sterile Conical High Resistance Tubes (50ml)

Product Code: 409922

Brand: Deltalab

Made of transparent and copolymer polypropylene, Deltalab's Sterile Conical High Resistance Tubes (50ml) are autoclavable and feature a solvent resistant white panel and black graduations for both dark or clear samples.

Their blue caps are made of polyethylene with an internal liner (the liner is not autoclavable), which ensures leakproofness. Each tube is sterilised by radiation and withstands temperatures down to -80°C.



Specifications:

Capacity (ml)	50
Dimensions (mm)	29.6 x 114.6
Resistance to Centrifugation	20,000xg
Material(s)	Polypropylene (tube) and Polyethylene (cap)
Sterility	Sterile
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com