

Deltalab Autoclave Bags (50 x 75cm)

Product Code: 200100

Brand: Deltalab

Deltalab's Autoclave Bags (50 x 75cm) are printed with a biohazard graphic symbol and instructions of use in five languages. Do not close the bags hermetically when placing them in the autoclave at 121°C.

Specifications:

Dimensions (cm)	50 x 75
Material	Polypropylene
Approx. Capacity (L)	34
Thickness (microns)	30
Sterility	Non-Sterile
Quantity	100



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Autoclave Bags (60 x 75cm)

Product Code: 200312

Brand: Deltalab

Deltalab's Autoclave Bags (60 x 75cm) are printed with a biohazard graphic symbol and instructions of use in five languages. Do not close the bags hermetically when placing them in the autoclave at 121°C.

Specifications:

Dimensions (cm)	60 x 75
Material	Polypropylene
Approx. Capacity (L)	49
Thickness (microns)	37
Sterility	Non-Sterile
Quantity	200



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Autoclave Bags (40 x 75cm)

Product Code: 200311

Brand: Deltalab

Deltalab's Autoclave Bags (40 x 75cm) are printed with a biohazard graphic symbol and instructions of use in five languages. Do not close the bags hermetically when placing them in the autoclave at 121°C.

Specifications:

Dimensions (cm)	40 x 75
Material	Polypropylene
Approx. Capacity (L)	24
Thickness (microns)	37
Sterility	Non-Sterile
Quantity	300



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Autoclave Bags (30 x 60cm)

Product Code: 200310

Brand: Deltalab

Deltalab's Autoclave Bags (30 x 60cm) are printed with a biohazard graphic symbol and instructions of use in five languages. Do not close the bags hermetically when placing them in the autoclave at 121°C.

Specifications:

Dimensions (cm)	30 x 60
Material	Polypropylene
Approx. Capacity (L)	12
Thickness (microns)	37.5
Sterility	Non-Sterile
Quantity	500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com