

Deltalab Natural Colour Filter Pipette Tip Rack (100 - 1000µl)

Product Code: 327-30R

Brand: Deltalab

Made of polypropylene and autoclavable up to 121°C, Deltalab's Natural Colour Filter Pipette Tip Rack (100 - 1000µl) has a hinged transparent lid and base. It can hold up to 96 tips.

Specifications:

Description	96 tips rack filter
Sterility	Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Yellow Eppendorf Pipette Tip Rack (5 - 200µl)

Product Code: 200009R

Brand: Deltalab

Deltalab's Yellow Eppendorf Pipette Tip Rack (5 - 200µl) features a transparent lid and blue base, which can hold up to 96 yellow pipette tips.

It is made of autoclavable polypropylene up to 121°C.

Specifications:

Description	96 tips rack
Sterility	Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Yellow Gilson Pipette Tip Rack (2 - 200µl)

Product Code: 200016R

Brand: Deltalab

Deltalab's Yellow Gilson Pipette Tip Rack (2 - 200µl) features a transparent lid and blue base and it can hold up to 96 yellow pipette tips.

The rack is made of autoclavable (up to 121°C) polypropylene.

Specifications:

Description	96 tips rack
Sterility	Non-Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Natural Colour Eppendorf Pipette Tip Rack (10 - 300µl)

Product Code: 200077FR

Brand: Deltalab

Deltalab's Natural Colour Eppendorf Pipette Tip Rack (10 - 300µl) has a hinged transparent lid and base. It can hold up to 96 tips.

This rack is made of autoclavable (up to 121°C) polypropylene.

Specifications:

Description	96 tips rack filter
Material	Polypropylene
Sterility	Sterile
Quantity	10



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com