

## Deltalab Sterile Water Sampling Bottles with Sodium Thiosulfate (1000ml)

**Product Code:** 282324.BU

**Brand:** Deltalab

Deltalab's 1000ml Water Sampling Bottles are specifically designed for examining water that is used for drinking, swimming pools, along with other types of water whereby chlorine may alter the arrangement of the sample whilst being transported.

Each water sampling bottle is sterile and has a printed label specifying description, code, lot number, and expiry date.

### Individually wrapped

### Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

### Specifications:

Capacity (ml)	1000
Dimensions (mm)	83 x 65 x 235
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	61g
Material	Polyethylene
Sterility	Sterile
Quantity	68



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)

## Deltalab Sterile Water Sampling Bottles with Sodium Thiosulfate (1000ml)

**Product Code:** 282324

**Brand:** Deltalab

Deltalab's 1000ml Water Sampling Bottles are specifically designed for examining water that is used for drinking, swimming pools, along with other types of water whereby chlorine may alter the arrangement of the sample whilst being transported.

Each water sampling bottle is sterile and has a printed label specifying description, code, lot number, and expiry date.

### Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

### Specifications:

Capacity (ml)	1000
Dimensions (mm)	83 x 65 x 235
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	61g
Material	Polyethylene
Sterility	Sterile
Quantity	68



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)

## Deltalab Sterile Water Sampling Bottles with Sodium Thiosulfate (500ml) - Packs of 111

**Product Code:** 282323

**Brand:** Deltalab

Deltalab's Sterile Water Sampling Bottles with Sodium Thiosulfate (500ml) - Packs of 111 are specifically designed for examining water that is used for drinking and swimming pools, along with other types of water where chlorine may alter the arrangement of the sample while being transported.

Each water sampling bottle is labelled specifying description, code, lot, and expiry date.

### Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

### Specifications:

Capacity (ml)	500
Dimensions (mm)	83 x 65 x 135
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	40g
Material	Polyethylene
Sterility	Sterile
Quantity	111



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)

## Deltalab Sterile Water Sampling Bottles with Sodium Thiosulfate (1000ml)

**Product Code:** 282321

**Brand:** Deltalab

Deltalab's 1000ml Water Sampling Bottles are specifically designed for examining water that is used for drinking, swimming pools, along with other types of water whereby chlorine may alter the arrangement of the sample whilst being transported.

Each water sampling bottle is sterile and has a printed label specifying description, code, lot number, and expiry date.

### Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

### Specifications:

Capacity (ml)	1000
Dimensions (mm)	83 x 65 x 235
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	61g
Material	Polyethylene
Sterility	Sterile
Quantity	20



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)

## Deltalab Sterile Water Sampling Bottles with Sodium Thiosulfate (500ml) - Packs of 24

**Product Code:** 282320

**Brand:** Deltalab

Deltalab's Sterile Water Sampling Bottles with Sodium Thiosulfate (500ml) - Packs of 24 are specifically designed for examining water that is used for drinking and swimming pools, along with other types of water where chlorine may alter the arrangement of the sample while being transported.

Each water sampling bottle is labelled specifying description, code, lot, and expiry date.

### Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

### Specifications:

Capacity (ml)	500
Dimensions (mm)	83 x 65 x 135
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	40g
Material	Polyethylene
Sterility	Sterile
Quantity	24



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)

## Deltalab Sterile Water Sampling Bottles with Sodium Thiosulfate (500ml)

**Product Code:** 282323.BU

**Brand:** Deltalab

Deltalab's 500ml Water Sampling Bottles are specifically designed for examining water that is used for drinking, swimming pools, along with other types of water whereby chlorine may alter the arrangement of the sample whilst being transported.

Each water sampling bottle is sterile and has a printed label specifying description, code, lot number, and expiry date.

### Individually wrapped

### Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

### Specifications:

Capacity (ml)	500
Dimensions (mm)	83 x 65 x 135
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	40g
Material	Polyethylene
Sterility	Sterile
Quantity	111



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)