

Deltalab Sterile Water Sampling Jars - Without Sodium Thiosulfate (1000ml)

Product Code: 282351

Brand: Deltalab

Deltalab's Sterile Water Sampling Jars - Without Sodium Thiosulfate (1000ml) can be used to sample water for physico-chemical analysis as well as for other sampling requiring sterile bottles.

Features:

- Transparent. Made of PET (Polyethylene Terephthalate).
- Wide mouth and grips at two sides, easing the sample collection.
- Red cap made of polyethylene with an internal liner made of pollexan with a label that performs as tamper evident.

Specifications:

Internal Mouth Diameter (mm)	55
Capacity (ml)	1000
Jar Weight (g)	65
Sterility	Sterilised by radiation
Material(s)	Polyethylene Terephthalate (jar) and Polyethylene (cap)
Quantity	48

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Sterile Water Sampling Jars - Without Sodium Thiosulfate (500ml)

Product Code: 282350

Brand: Deltalab

Deltalab's Sterile Water Sampling Jars - Without Sodium Thiosulfate (500ml) can be used to sample water for physico-chemical analysis as well as for other sampling requiring sterile bottles.

Features:

- Transparent. Made of PET (Polyethylene Terephthalate).
- Wide mouth and grips at two sides, easing the sample collection.
- Red cap made of polyethylene with an internal liner made of pollexan with a label that performs as tamper evident.

Specifications:

Internal Mouth Diameter (mm)	55
Capacity (ml)	500
Jar Weight (g)	44
Sterility	Sterilised by radiation
Material(s)	Polyethylene Terephthalate (jar) and Polyethylene (cap)
Quantity	44



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Sterile Water Sampling Bottles Without Sodium Thiosulfate (1000ml) - Packs of 68

Product Code: 282334

Brand: Deltalab

Deltalab's Sterile Water Sampling Bottles Without Sodium Thiosulfate (1000ml) - Packs of 68 feature tamper-evident red screw caps, which assures the bottle's sterility before use. These laboratory containers can be used to sample water for physico-chemical analysis, along with other sampling requiring sterile bottles.

Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

Specifications:

Size (ml)	1000
Dimensions (mm)	83 x 65 x 235
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	61g
Material	Polyethylene
Sterility	Sterilised by radiation
Quantity	68



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Sterile Water Sampling Bottles Without Sodium Thiosulfate (500ml) - Packs of 111

Product Code: 282333

Brand: Deltalab

Deltalab's Sterile Water Sampling Bottles Without Sodium Thiosulfate (500ml) - Packs of 111 feature tamper-evident red screw caps, which assures the bottle's sterility before use. These laboratory containers can be used to sample water for physico-chemical analysis, along with other sampling requiring sterile bottles.

Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

Specifications:

Size (ml)	500
Dimensions (mm)	83 x 65 x 135
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	40g
Material	Polyethylene
Sterility	Sterilised by radiation
Quantity	111



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Sterile Water Sampling Bottles Without Sodium Thiosulfate (1000ml) - Packs of 20

Product Code: 282331

Brand: Deltalab

Deltalab's Sterile Water Sampling Bottles Without Sodium Thiosulfate (1000ml) - Packs of 20 feature tamper-evident red screw caps, which assures the bottle's sterility before use. These laboratory containers can be used to sample water for physico-chemical analysis, along with other sampling requiring sterile bottles.

Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

Specifications:

Capacity (ml)	1000
Dimensions (mm)	83 x 65 x 235
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	61g
Material	Polyethylene
Sterility	Sterile
Quantity	20



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Deltalab Sterile Water Sampling Bottles Without Sodium Thiosulfate (500ml) - Packs of 24

Product Code: 282330

Brand: Deltalab

Deltalab's Sterile Water Sampling Bottles Without Sodium Thiosulfate (500ml) - Packs of 24 feature tamper-evident red screw caps, which assures the bottle's sterility before use. These laboratory containers can be used to sample water for physico-chemical analysis, along with other sampling requiring sterile bottles.

Features:

- Minor sides are gripped for a better handling.
- Larger sides are flat for labelling.
- Tamper evident red cap with an internal joint.
- Hermetical closure.

Specifications:

Capacity (ml)	500
Dimensions (mm)	83 x 65 x 135
Mouth Diameter (mm)	28
Unitary Weight (empty bottle)	40g
Material	Polyethylene
Sterility	Sterile
Quantity	24



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com