

Medline High Temperature Autoclave Bags (610 x 810mm)

Product Code: 200312R

Brand: Medline Scientific

Medline's High Temperature Autoclave Bags (610 x 810mm) are manufactured from heavy gauge polypropylene material and are suitable for autoclaving up to a temperature of 136°C.

Each autoclave bag has blue biohazard labelling, designed for contaminated waste disposal.

They also have a write on panel to enable the user to mark each bag if required.

Specifications

Dimensions (mm)	610 x 810
Material	Polypropylene
Sterility	Non-Sterile
Quantity	200



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Medline High Temperature Autoclave Bags (400 x 700mm)

Product Code: 200311R

Brand: Medline Scientific

Medline's High Temperature Autoclave Bags (400 x 700mm) are manufactured from heavy gauge polypropylene material and are suitable for autoclaving up to a temperature of 136°C.

Each autoclave bag has blue biohazard labelling, designed for contaminated waste disposal.

They also have a write on panel to enable the user to mark each bag if required.

Specifications

Dimensions (mm)	400 x 700
Material	Polypropylene
Sterility	Non-Sterile
Quantity	300



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

Medline High Temperature Autoclave Bags (310 x 660mm)

Product Code: 200310R

Brand: Medline Scientific

Medline's High Temperature Autoclave Bags (310 x 660mm) are manufactured from heavy gauge polypropylene material and are suitable for autoclaving up to a temperature of 136°C.

Each autoclave bag has blue biohazard labelling, designed for contaminated waste disposal.

They also have a write on panel to enable the user to mark each bag if required.

Specifications

Dimensions (mm)	310 x 660
Material	Polypropylene
Sterility	Non-Sterile
Quantity	500



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com