



Non-contact FOREHEAD THERMOMETER

No-touch infrared technology
For safety and hygiene.

Forehead use
Displays "oral equivalent"
temperature.

Simple to use
Point the thermometer at
the forehead and press START
It's that easy!

High temperature alert
Audible alert when
temperatures
are > 37.5 °C (99.5 °F)

Memory
Stores the 25 most recent
readings for easy recall.

FDA-cleared
For use in adults and children.

Quick one second readings
Instant results - no waiting.
Beeps when done.

Mute setting
Disable audible alert
that could disturb a
sleeping child.

Bright LED backlight
Makes temperatures
easy to read, even in
the dark.



NON-CONTACT, ONE-SECOND READINGS FOR QUICK, SAFE SCREENING.

*FDA-cleared for safety & accuracy for kids and adults. ETI non-contact forehead thermometer measures infrared energy radiated from the skin at the centre of the forehead area. This captured energy is collected through the lens and converted to a body temperature value (displays oral equivalent temperature).

It's helpful to know each individual's normal temperature when they are well. This is the best way to accurately diagnose a 'high temperature'. To get a reference value, we suggest taking multiple readings (5 to 10) over one minute, taking the highest number to determine "normal temperature".

Premier Scientific Ltd

Unit 3 Westbank Drive | Belfast | BT3 9LA | Tel: 02890 773030 | Email: sales@premierscientific.co.uk