

## KF Flake Ice Machines

Flake ice can be used in many areas both in laboratories and medical institutions from physiotherapy to instrument calibration or transportation of donated organs. The ice remains free flowing, even during long storage periods.

We offer a comprehensive range from counter top to under bench models. Each machine is supplied with a high quality stainless steel body work. They are easy to sanitise and fully automatic in operation.

### Features

- Stainless steel bodywork with one piece insulated storage bin
- Spray system (KL)  
Belt system (KF)
- Automatic production cycle (KL)  
Electronic controlled production (KF)
- Built-in cleaning system (KL)
- Low power and water consumption
- Crystal clear, hygienic and compact ice cubes with the size of  $\varnothing 27 \times 32$  mm (kl)
- Hygienic and dry flake ice with a temperature of  $-0,5^{\circ}\text{C}$
- Production rates are based on water entering the ice machine at  $+10^{\circ}\text{C}$  and air on condenser at  $+10^{\circ}\text{C}$ .



## KF Flake Ice Machines

### KF 45

Air cooled, production 40 Kg / 24 h  
Water cooled, production 42 Kg / 24 h  
Exterior dimensions 345 x 520 x 610 + 25 mm (width x depth x height + feet)  
Bin storage 10 kg  
Connection power 0,41 kW / 230 V / 50 Hz  
Gross weight 41 Kg

### KF 85

Air cooled, production 80 Kg / 24 h  
Water cooled, production 83 Kg / 24 h  
Exterior dimensions 485 x 586 x 720 + 140 mm (width x depth x height + feet)  
Bin storage 25 kg  
Connection power 0,52 kW / 230 V / 50 Hz  
Gross weight 62 Kg

### KF 145

Air cooled, production 135 Kg / 24 h  
Water cooled, production 141 Kg / 24 h  
Exterior dimensions 680 x 588 x 890 + 140 mm (width x depth x height + feet)  
Bin storage 42 kg  
Connection power 0,59 kW / 230 V / 50 Hz  
Gross weight 72 Kg

### KF 185A

Air cooled, production 175 Kg / 24 h  
Water cooled, production 183 Kg / 24 h  
Exterior dimensions 953 x 590 x 781 + 140 mm (width x depth x height + feet)  
Bin storage 60 kg  
Gross weight 89 Kg



**KF 45**



**KF 85**



**KF 145**

For pricing, please ask your local dealer for the price list

Porkka reserves the right to make product alterations and or price adjustments at any time and cannot be held responsible for any printing errors or omissions.