



Your Source for Laboratory Equipment and Chemicals



PREMIER SCIENTIFIC CATALOGUE SPRING/SUMMER

www.premierscientific.co.uk
028 9077 3030



WELCOME

At Premier Scientific, we are proud to continue supporting you with reliable, high-quality scientific products you can trust.

This brochure showcases our current product range, highlighting the equipment and solutions we provide to help you work efficiently and achieve accurate results.

As your needs evolve, we remain committed to offering dependable products, expert knowledge, and consistent service. Inside, you'll find an overview of supplies available, carefully selected to meet the standards you expect from us.

Thank you for choosing Premier Scientific; we look forward to continuing to support your work.

Professionals trust Premier Scientific for dependable products and consistent service. Every item in this catalogue has been selected to meet the standards required in laboratories, education, and industry environments.

Our commitment to you:


- ✓ Quality-tested scientific products
- ✓ Trusted suppliers and manufacturers
- ✓ Reliable availability and delivery
- ✓ Competitive pricing
- ✓ Knowledgeable customer support
- ✓ After Sales Service and Support



INTEGRA PIPETBOY

ACU 2

- Autoclavable sterile module, easy to replace filter
- LED provides information about charging status
- Continuous preset of speed via thumb wheel
- Sensitive control by applying various finger pressure
- Lightweight and well-balanced
- Over 6.5 hours of non-stop pipetting



Code	Description	List Price	Offer Price
155019	INTEGRA PIPETBOY ACU 2	£297.00	£280.36



SARTORIUS

SARTORIUS BIOHIT MIDI PLUS CONTROLLER, 1-100 ML

TOP FEATURES OF THE MIDI PLUS CONTROLLER

- LIGHTWEIGHT, CORDLESS PIPETTING CONTROLLER
- FITS ALL COMMONLY USED GLASS OR PLASTIC PIPETTES
- FINE-TUNE SPEED
- EASY TO MAINTAIN/CLEAN

POPULAR APPLICATIONS FOR THE MIDI PLUS CONTROLLER

- CELL CULTURE
- LIQUID PREPARATION IN LARGE VOLUMES
- WORKING WITH T-FLASKS



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
710931	Biohit Midi Plus Controller, 1-100 ml	£463.00	£333.36



View Our Website



ABDOS PIPETTE FILLER

- Fast, efficient pipetting for higher productivity
- Ergonomic, lightweight design for comfortable use
- Precise speed control for accurate handling
- Built-in filter for contamination protection
- Long-lasting rechargeable battery
- Compatible with a wide range of pipettes

Code	Description	Offer Price
E11114	Rave Plus Pipette Filler 1un	£269.00
E11120	Rave Eco Pipette Filler 1un	£191.00



ABDOS FILTERED TIPS

CODE	DESCRIPTION	OFFER PRICE
P11060	Pipette Filter tips, 0.1-10 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£53.74
P11061	Pipette Filter tips, 0.1-10 uL Extended, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£58.12
P11062	Pipette Filter tips, 1-20 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£53.74
P11063	Pipette Filter tips, 1-50 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£58.12
P11064	Pipette Filter tips, 1-100 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£53.74
P11065	Pipette Filter tips, 1-200 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£58.12
P11066	Pipette Filter tips, 1-300 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£53.74
P11067	Pipette Filter tips, 100 -1000 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£58.12
P11068	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£53.74
P11069	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£58.12
P11072	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£53.75
P11073	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£ 58.12
P11074	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£ 53.75
P11075	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£58.12
P11076	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£60.66
P11077	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£65.05
P11078	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£70.46
P11079	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£74.93
P11080	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£70.46
P11081	Pipette Filter tips, 100 -1250 uL, Low Retention, Racked, Natural, Sterile, DNase/RNase free - 10 x 96 un.	£74.93

GILSON STARTER KITS

- Autoclavable sterile module, easy to replace filter
- LED provides information about charging status
- Continuous preset of speed via thumb wheel
- Sensitive control by applying various finger pressure
- Lightweight and well-balanced
- Over 6.5 hours of non-stop pipetting



Code	Description	List Price	Offer Price
F167900	PIPETMAN STARTER KIT EACH PIPETMAN STARTER KIT INCLUDES: - THREE PIPETMAN PIPETTES: P20, P200, AND P1000 EQUIPPED WITH METAL EJECTOR - TWO PIPETMAN® DIAMOND TIPS TIPACKS (MODELS: D200 AND D1000) - THREE SINGLE PIPETTE HOLDERS - THREE COMFORT HANDLES - ONE SHEET OF 12 COMFORT HANDLE LABELS	£825.00	£600.60
F167800	PIPETMAN MICRO-VOLUME KIT EACH PIPETMAN MICRO-VOLUME KIT INCLUDES: - THREE PIPETMAN PIPETTES: P2, P10, AND P100 EQUIPPED WITH METAL EJECTOR - TWO PIPETMAN® DIAMOND TIPS TIPACKS (MODELS: DL10 AND D200) - THREE SINGLE PIPETTE HOLDERS - THREE PIPETMAN COMFORT HANDLES - ONE SHEET OF 12 PIPETMAN COMFORT HANDLE LABELS	£825.00	£600.60
F167920	PIPETMAN LARGE-VOLUME KIT INCLUDES P5000 AND P10ML Each PIPETMAN Large-volume kit includes: - Two PIPETMAN pipettes: P5000 and P10mL models - Two PIPETMAN® DIAMOND Tips TIPACKS (models: D5000 and D10mL) - Two SINGLE Pipette Holders - Two PIPETMAN Comfort Handles - One sheet of 12 PIPETMAN Comfort Handle labels - Two bags of ten polyethylene filters	£552.00	£401.86



ABDOS TIPS

CODE	DESCRIPTION	OFFER PRICE
P10115	Pipette tips, 0.1-10 uL, Natural, Bulk, DNase/RNase free - 1000 un. - ABDOS	£6.96
P10137	Pipette tips, 0.1-10 uL Extended, Natural, Bulk, DNase/RNase free - 1000 un. - ABDOS	£6.96
P10129	Pipette tips, PP, 1-200 uL, Yellow Bulk, DNase/RNase free - 1000 un. - ABDOS	£6.62
P10130	Pipette tips, PP, 1-200 uL, Natural Bulk, DNase/RNase free - 1000 un. - ABDOS	£6.62
P10101	Pipette tips, PP, 1-200 uL, Yellow, Ultra Point Tip, Bulk, DNase/RNase free - 1000 un. - ABDOS	£7.26
P10105	Pipette tips, PP, 1-200 uL, Natural, Ultra Point Tip, Bulk, DNase/RNase free - 1000 un. - ABDOS	£7.26
P10138	Pipette tips, PP, 1-300 uL, Natural, Bulk DNase/RNase free - 1000 un. - ABDOS	£8.22
P10102	1000 µl (1-1000µl) Bevelled Pipette Tip Blue, Graduated, 500 in resealable bags-ABDOS	£5.13
P10174	Pipette tips, 100-1000 XL uL, Natural, Bulk, DNase/RNase free - 500 un. - ABDOS	£5.13
P10106	Pipette tips, 100-1000 uL, Natural, Bulk, DNase/RNase free - 500 un. - ABDOS	£5.70
P10175	Pipette tips, 100-1250 uL, Natural, Bulk, DNase/RNase free - 500 un. - ABDOS	£5.70



NOVASBIO FILTERED TIPS

CODE	DESCRIPTION	OFFER PRICE
Filter Tips		
FT-10-R-S	10µl, Aerosol Barrier Pipette Tips, Racked, Sterile (FT-300-R-S)	£316.33
FT-10XL-R-S	10µl Extended, Aerosol Barrier Pipette Tips, Racked, Sterile (FT-400-R-S)	£316.33
FT-20-R-S	20µl, Aerosol Barrier Pipette Tips, Racked, Sterile	£316.33
FT-100-R-S	100µl, Aerosol Barrier Pipette Tips, Racked, Sterile	£316.33
FT-200-R-S	200µl, Aerosol Barrier Pipette Tips, Racked, Sterile	£316.33
FT-1000XL-R-S	1000µl, Extended, Aerosol Barrier Pipette Tips, Racked, Sterile	£388.37
FT-1250-R-S	1250µl, Aerosol Barrier Pipette Tips, Racked, Sterile	£388.37
Filter tips, Low Retention		
FT-10-L-R-S	10µl, Aerosol Barrier Pipette Tips, Low Retention, Racked, Sterile (FT-300-L-R-S)	£350.78
FT-10XL-L-R-S	10µl Extended, Aerosol Barrier Pipette Tips, Low Retention, Racked, Sterile (FT-400-L-R-S)	£350.78
FT-20-L-R-S	20µl, Aerosol Barrier Pipette Tips, Low Retention, Racked, Sterile	£360.18
FT-100-L-R-S	100µl, Aerosol Barrier Pipette Tips, Low Retention, Racked, Sterile	£360.18
FT-200-L-R-S	200µl, Aerosol Barrier Pipette Tips, Low Retention, Racked, Sterile	£360.18
FT-300-L-R-S	300µl, Aerosol Barrier Pipette Tips, Graduations, Low Retention, Racked, Sterile	£371.14
FT-1000-L-R-S	1000µl, Aerosol Barrier Pipette Tips, Low Retention, Racked, Sterile	£471.37
FT-1000XL-L-R-S	1000µl, Extended, Aerosol Barrier Pipette Tips, Low Retention, Racked, Sterile	£458.84
FT-1250-L-R-S	1250µl, Aerosol Barrier Pipette Tips, Low Retention, Racked, Sterile	£433.78



ABDOS SPECIALS

Centrifuge tubes

P10202	Micro Centrifuge Tubes, PP, 1.5 mL, DNase/RNase free - 500	£5.56
P10203	Micro Centrifuge Tubes, PP, 2.0 mL, DNase/RNase free - 500	£8.34
P10209	Last Drop, Low Retention Micro Centrifuge Tubes, PP, 1.5 mL, DNase/RNase free - 500	£8.34
P10210	Last Drop, Low Retention Micro Centrifuge Tubes, PP, 2.0 mL, DNase/RNase free - 500	£10.01
P10402	Centrifuge Tube, Conical, PP, 15 mL, Sterile, DNase/RNase free - 20 x 25	£48.40
P10404	Centrifuge Tube, Conical, PP, 50 mL, Sterile, DNase/RNase free - 20 x 25	£50.00
P10432	Centrifuge Tube, Conical, PP, 50 mL, Self Standing, Bulk, DNase/RNase free - 500	£78.62
P10430	Centrifuge Tube, Conical, PP, 15 mL, Bulk, DNase/RNase free - 500	£48.97
P10431	Centrifuge Tube, Conical, PP, 50 mL, Bulk, DNase/RNase free - 500	£56.69
P10432	Centrifuge Tube, Conical, PP, 50 mL, Self Standing, Bulk, DNase/RNase free - 500	£65.36

Cryo Vials & boxes

P60116	Cryo Vial, External Thread, Star Foot, PP, 1.8 mL, Sterile, DNase/RNase free - 1000 un. - ABDOS	£94.93
P60109	Cryo Vial, Internal Thread, Star Foot, Silicone Seal, PP, 1.8 mL, Sterile, DNase/RNase free - 1000 un. - ABDOS	£158.27
P20602W	Cryo Box, PC, 1.0 & 1.8 mL Cryo Vials, 81 Place, White - 4 un. - ABDOS	£18.65
P20602Y	Cryo Box, PC, 1.0 & 1.8 mL Cryo Vials, 81 Place, Yellow - 4 un. - ABDOS	£18.65
P20602B	Cryo Box, PC, 1.0 & 1.8 mL Cryo Vials, 81 Place, Blue - 4 un. - ABDOS	£18.65
P20602 ASSORTED	Cryo Box, PC, 1.0 & 1.8 mL Cryo Vials, 81 Place, Assorted Colours - 4 un. - ABDOS	£18.65

PCR

P10701	Pcr tubes, pp clear flat cap 0.2ml , DNase/RNase free - 1000 un. - ABDOS	£18.95
P10726	Combo - Clear Strip tubes & Optically Clear Flat cap - 8 strip low profile tubes 0.1ml, pp, DNase/RNase free - 125un - ABDOS	£39.49
P10727	Combo - Clear Strip tubes & Optically Clear Flat cap - 8 strip high profile tubes 0.2ml, pp, DNase/RNase free - 125un - ABDOS	£35.53
P10734	Combo - Clear Strip tubes & Clear Domed caps - 8 strip high profile tubes 0.2ml, pp, DNase/RNase free - 125un - ABDOS	£35.53
P10735	Combo - Clear Strip tubes & Clear Domed cap - 8 strip low profile tubes 0.1ml, pp, DNase/RNase free - 125un - ABDOS	£39.49
P10713	Pcr plates , pp - 96 unskirted standard well, DNase/RNase free - 25 un. - ABDOS	£37.30
P10714	Pcr plates , pp - 96 unskirted low profile well, DNase/RNase free - 25 un. - ABDOS	£39.49
P10722	Pcr plates , pp - 96 semi-skirted standard well, DNase/RNase free - 25 un. - ABDOS	£45.90
P10716	Pcr plates , pp - 96 semi-skirted standard well, DNase/RNase free - 10 un. - ABDOS	£32.67



View Our
Website

ABDOS SPECIALS

Liquid Handling

P10603	Serological Pipettes, Sterile, PS, 5 mL, Ind. Wrapped, DNase/RNase free - 200 un. - ABDOS	£25.92
P10604	Serological Pipettes, Sterile, PS, 10 mL, Ind. Wrapped, DNase/RNase free - 200 un. - ABDOS	£27.69
P10605	Serological Pipettes, Sterile, PS, 25 mL, Ind. Wrapped, DNase/RNase free - 200 un. - ABDOS	£51.85
P10606	Serological Pipettes, Sterile, PS, 50 mL, Ind. Wrapped, DNase/RNase free - 100 un. - ABDOS	£52.19
P31201	Pasteur Pipettes, LDPE, 1.0 mL, Low-Affinity, DNase/RNase free - 500 un. - ABDOS	£8.62
P31205	Pasteur Pipettes, LDPE, 3.0 mL, Low-Affinity, DNase/RNase free - 500 un. - ABDOS	£8.19

Multipurpose Racks & Boxes

P11104	Frosty Mate™ Freezer Storage Rack, PC, -20°C, 32 Place (8 x 4 Array) - ABDOS	£41.28
P11105	Frosty Mate™ Freezer Storage Rack, PC, 0°C, 32 Place (8 x 4 Array) - ABDOS	£41.28
P11106	Frosty Mate™ Colour Change Freezer Racks, PP, 24 Micro Centr. Tubes (6 x 4 Array) - ABDOS	£31.99
P11107	Frosty Mate™ Colour Change Freezer Racks, PP, 96 PCR Tubes (12 x 8 Array) - ABDOS	£31.99


Cell Culture

P21011	25cm ² Cell Culture Flask, Vent Cap, TC, Sterile 10/pk, 200 un.- ABDOS	£118.20
P21013	75cm ² Cell Culture Flask, Vent Cap, TC, Sterile 5/pk, 100 un.- ABDOS	£118.20
P21015	175cm ² Cell Culture Flask, Vent Cap, TC, Sterile 5/pk, 40 un.- ABDOS	£90.48
P21017	225cm ² Cell Culture Flask, Vent Cap, TC, Sterile 5/pk, 25 un.- ABDOS	£81.40
P21020	6 Well Cell Culture Plate, Flat, TC, Sterile 1/pk, 50 un.- ABDOS	£64.03
P21021	12 Well Cell Culture Plate, Flat, TC, Sterile 1/pk, 50 un.- ABDOS	£66.55
P21022	24 Well Cell Culture Plate, Flat, TC, Sterile 1/pk, 50 un.- ABDOS	£66.55
P21023	48 Well Cell Culture Plate, Flat, TC, Sterile 1/pk, 50 un.- ABDOS	£67.93
P21024	96 Well Cell Culture Plate, Flat, TC, Sterile 1/pk, 100 un.- ABDOS	£125.42
P21025	96 Well Cell Culture Plate, U-bottom, TC, Sterile 1/pk, 100 un.- ABDOS	£133.95
P21031	60mm Cell Culture Dish, TC, Sterile 20/bag, 500 un.- ABDOS	£98.02
P21032	100mm Cell Culture Dish, TC, Sterile 20/bag, 300 un.- ABDOS	£90.44
P21033	150mm Cell Culture Dish, TC, Sterile 5/pk, 100 un.- ABDOS	£105.30



View Our
Website

HAIER REFRIGERATED MICROCENTRIFUGE



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
LX-165T2R	Benchtop Refrigerated Centrifuge (Needs Rotor)	£4,699.00	£2,884.00
ACCESSORIES			
BE12C0Z00	Standard Fixed-angle rotor (Plastic Cover) 24*1.5/2ml for refrigerated centrifuge	£489.00	£300.00
BE11D0Z00	Standard Fixed-angle rotor (Sealed) 24*1.5/2ml for refrigerated centrifuge	£539.00	£330.00
BE12E0Z00	Optional fixed-angle rotor 4*8* 0.2ml PCR tube for refrigerated centrifuge	£649.00	£398.00
BE12D0Z00	Optional fixed-angle rotor (Sealed) 12 * 5ml tubes for refrigerated Centrifuge	£699.00	£428.00
BE12DR000	Fixed-angle Rotor 12*10ml Sealed -Large Blood Collection Tube	POA	POA
BE12J2003	0.2ML Adaptors	£49.00	£30.00
BE12JN002	0.5ML Adaptors	£49.00	£30.00

CAMLAB MINI CENTRIFUGE



THE 8-PLACE CAMLAB CHOICE D1008 MINI MICRO CENTRIFUGE IS IDEAL FOR QUICK SPIN DOWNS OF MICROTUBES AND PCR TUBES. THE MICROFUGE IS SUPPLIED WITH A ROTOR FOR 8X MICRO TUBES AND A ROTOR FOR 16X PCR TUBES.

KEY FEATURES:

- MAXIMUM SPEED RPM: 7000RPM
- MAXIMUM RCF[XG]: 2680
- SUPPLIED WITH AN 8 PLACE 1.5/2.0ML MICROTUBE ROTOR AND 16 X 0.2ML PCR ROTOR
- QUIET RUNNING, NOISE LEVEL 45 DB
- ROTOR CHANGE WITHOUT TOOLS
- COMPACT DESIGN
- ROTOR CHANGE WITHOUT TOOLS
- AUTO LID OPENS AND CLOSES AUTOMATICALLY
- UNIVERSAL VOLTAGE
- 2 YEAR WARRANTY
- WEIGHT 2.1 LBS
- DIMENSIONS: (LXWXH) 150X150X125MM
- AVAILABLE IN OTHER COLOURS

Code	Description
1190944	Camlab D1008 Mini Centrifuge with Blue Lid



HAIER NON REFRIGERATED MICROCENTRIFUGE

KEY FEATURES

- THE MACHINE HAS A HIGH-SENSITIVITY, SIX AXIS VIBRATION DETECTION MECHANISM. IN CASE OF ABNORMAL VIBRATION OR ROTOR IMBALANCE CAUSED BY NON BALANCING, ASYMMETRY, TUBE BREAKAGE AND OR ROTOR DAMAGE, AN ALARM WILL SOUND AND THE MACHINE WILL STOP IMMEDIATELY
- ANTI-VIBRATION LID LOCK ENSURES LID REMAINS LOCKED DURING EXCESS VIBRATION FOR USER SAFETY
- AUTOMATIC ROTOR IDENTIFICATION LIMITS MAXIMUM SPEED FOR SAFETY
- CABINET CONSTRUCTION INCLUDES METAL LINER, FULLY COVERED STEEL PLATE PROTECTIVE SLEEVE AND 3MM STEEL SHELL
- MULTIPLE DAMPING AND NOISE REDUCTION TECHNOLOGY
- MULTIPLE ANTI-VIBRATION TECHNOLOGIES ALLOW FOR LOW NOISE AND SMOOTH OPERATION EVEN AT MAXIMUM SPEED (16,500RPM 58DB)



Code	Description	List Price	Offer Price
LX-165T2-J	DESKTOP NON-REFRIGERATED HIGH-SPEED MICRO CENTRIFUGE	£1,499.00	£919.00
BE12C4Z00	LX-165T2-J fixed-angled rotor 24×1.5/2.0ml (sealed)	£539.00	£330.00
BE12D4Z00	LX-165T2-J fixed-angled rotor 12×5ml (sealed)	£699.00	£428.00
BE12E4Z00	LX-165T2-J fixed-angled rotor-4-PCR	£649.00	£398.00
BE12C5Z00	LX-165T2-J fixed-angled rotor 24×1.5/2.0ml (plastic cover)	£489.00	£300.00

ABDOS MICROCENTRIFUGE TUBES

Code	Description	Offer Price
P10208	Last Drop, Low Retention Micro Centrifuge Tubes, PP, 0.5 mL, DNase/RNase free - 1000 un	£14.87
323P10209	Last Drop, Low Retention Micro Centrifuge Tubes, PP, 1.5 mL, DNase/RNase free - 500 un.)	£9.50
P10210	Last Drop, Low Retention Micro Centrifuge Tubes, PP, 1.5 mL, DNase/RNase free - 500 un.)	£11.52



STANDARD MICRO TUBES

Code	Description	Offer Price
200400P	DELTALAB 1.5ML GRADUATED MICROTUBES PK500	£7.15
P10201	MICRO CENTRIFUGE TUBE 0.5 ML-ABDOS PK1000	£10.77
P10202	MICRO CENTRIFUGE TUBES,PP,1.5ML,DNASE/RNASE FREE PK500	£6.95
P10203	MICRO CENTRIFUGE TUBES,PP,2.0ML,DNASE,RNASE FREE PK500	£8.34



WANT TO SEE MORE MICRO TUBES?
HEAD OVER TO OUR WEBSITE TO SEE
MORE

CAMLAB CHOICE MX-S VORTEX MIXER 0-2500RPM WITH UK PLUG



CODE

1189983

DESCRIPTION

Choice MX-S Vortex Mixer 0-2500rpm With UK plug

MX-S VORTEX MIXER IS A SIMPLE LABORATORY MIXER DESIGNED FOR WORKING WITH A RANGE OF TUBE SIZES AND SAMPLES.

FEATURES:

- WIDE SPEED RANGE: 0-2500RPM, STEPLESS SPEED REGULATION OF MOTOR
- SUPPLIED WITH 1X 1189984 STANDARD TOP FOR \varnothing30MM TUBES
- CAN BE USED FOR VARIOUS SHAKING APPLICATION WITH OPTIONAL ADAPTERS FOR EPPENDORF AND OTHER MICROCENTRIFUGE TUBES
- ADAPTED SECURELY CLOKED ONTO APPLIANCE IN ANY POSITION
- ROBUST ALUMINIUM-CAST CONSTRUCTION
- SPECIALLY DESIGNED SILICONE FEET AND STEEL BASE FOR EXCELLENT VIBRATION

SPECIFICATIONS:

- POWER: 60W
- SPEED RANGE: 0-2500RPM
- ORBITAL DIAMETER: 4MM
- WEIGHT: 3.5KG
- PERMISSIBLE AMBIENT TEMPERATURE 5-40 °C



SARTORIUS TRADE IN DEAL TACTA PIPETTES

35% OFF LIST PRICE

Code	Description	List Price	Offer Price
Single Channel Pipettes			
LH-729010	Tacta Mech. Pipette, 1-ch, 0.3-3 µl	£291.00	£189.15
LH-729070	Tacta Mech. Pipette, 1-ch, 100-1000 µl	£291.00	£189.15
LH-729050	Tacta Mech. Pipette, 1-ch, 10-100 µl	£291.00	£189.15
LH-729020	Tacta Mech. Pipette, 1-ch, 0.3-3 µl	£291.00	£189.15
LH-729090	Tacta Mech. Pipette, 1-ch, 1-10 ml	£291.00	£189.15
LH-729060	Tacta Mech. Pipette, 1-ch, 20-200 µl	£291.00	£189.15
LH-729030	Tacta Mech. Pipette, 1-ch, 2-20 µl	£291.00	£189.15
LH-729080	Tacta Mech. Pipette, 1-ch, 500-5000µl	£291.00	£189.15
Eight Channel Pipettes			
LH-729130	Tacta Mech. Pipette, 8-ch, 10-100 µl	£777.00	£505.00
LH-729120	Tacta Mech. Pipette, 8-ch, 1-10 µl	£777.00	£505.00
LH-729140	Tacta Mech. Pipette, 8-ch, 30-300 µl	£777.00	£505.00
725130	Mline Mech. Pipette, 8-ch, 10-100 µl	£847.00	£550.55
725120	Mline Mech. Pipette, 8-ch, 1-10 µl	£847.00	£550.55
725140	Mline Mech. Pipette, 8-ch, 30-300 µl	£847.00	£550.55



SARTORIUS TRADE IN DEAL PICUS PIPETTES

17% OFF LIST PRICE

Code	Description	List Price	Offer Price
Eight Channel Pipettes			
LH-747321	Picus 2, 8-ch, 0.5 - 10 µL	£1,130.00	£937.90
LH-747361	Picus 2, 8-ch, 10 - 300 µL	£1,130.00	£937.90
LH-747341	Picus 2, 8-ch, 5 - 120 µL	£1,130.00	£937.90
LH-747391	Picus 2, 8-ch, 50 - 1,200 µL	£1,130.00	£937.90
Twelve Channel Pipettes			
LH-747421	Picus 2, 12-ch, 0.5 - 10 µL	£1,330.00	£1,103.90
LH-747461	Picus 2, 12-ch, 10 - 300 µL	£1,330.00	£1,103.90
LH-747441	Picus 2, 12-ch, 5 - 120 µL	£1,330.00	£1,103.90
LH-747491	Picus 2, 12-ch, 50 - 1,200 µL	£1,330.00	£1,103.90
Single Channel Pipettes			
LH-747021	Picus 2, 1-ch, 0.5 - 10 µL	£711.00	£590.13
LH-747061	Picus 2, 1-ch, 10 - 300 µL	£711.00	£590.13
LH-747101	Picus 2, 1-ch, 100 - 5,000 µL	£711.00	£590.13
LH-747041	Picus 2, 1-ch, 5 - 120 µL	£711.00	£590.13
LH-747081	Picus 2, 1-ch, 50 - 1,000 µL	£711.00	£590.13
LH-747111	Picus 2, 1-ch, 500 - 10,000 µL	£711.00	£590.13



SARTORIUS



View Our
Website

STIRRING HOTPLATES



Cole-Parmer SHP-200-C Undergrad Analog Stirring Hot Plate

Designed with microprocessor-controlled dual thermocouples, these stirring hot plates offer you accurate temperature control and protection from overheating. Units feature a colored LED indicator which surrounds the heat control knob illuminating the actual surface temperature in real time. Powerful magnets and motor provide you with stirring speeds of up to 2000 rpm and mixing capability of up to 15 L of liquid. The square plate design allows you flexibility of vessel sizes and combinations.

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
04807-58	Undergrad Analog Stirring Hot Plate, 15 x 15 cm, Ceramic; 230 VAC	£466.19	£328.94




Cole-Parmer®

Cole-Parmer SHP-200D-C Undergrad Digital Stirring Hot Plate

Designed with microprocessor-controlled dual thermocouples, this stirring hot plate provides highly accurate temperature regulation and reliable protection against overheating. The unit features a clear digital display that shows the actual surface temperature in real time, allowing for precise monitoring and control during experiments. A powerful motor and magnetic drive system enable stirring speeds of up to 2000 rpm and efficient mixing of volumes up to 15 L. The square plate design accommodates a wide range of vessel sizes, offering flexibility for different laboratory applications.

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
04807-62	Cole-Parmer SHP-200D-C Undergrad Digital Stirring Hot Plate, 15 x 15 cm, Ceramic; 230V	£493.91	£348.51



Cole-Parmer SHP-200-S Undergrad Analog Stirring Hot Plate, 15 x 15cm, Aluminum; 230 VAC

Designed with microprocessor-controlled dual thermocouples, these stirring hot plates provide precise temperature regulation and reliable overheating protection. A coloured LED ring surrounding the heat control knob displays the surface temperature in real time, offering clear visual feedback during operation. The powerful motor and magnetic system deliver stirring speeds of up to 2000 rpm and can effectively mix volumes of up to 15 L. The square plate design accommodates a variety of vessel sizes, ensuring flexibility for different laboratory setups.

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
04807-60	Undergrad Analog Stirring Hot Plate, 15 x 15 cm, Aluminum; 230 VAC	£416.79	£294.08

Cole-Parmer SHP-200D-S Undergrad Digital Stirring Hot Plate, 15 x 15 cm, Aluminum; 230V

Designed with microprocessor-controlled dual thermocouples, these stirring hot plates provide precise temperature control along with reliable protection against overheating. The actual surface temperature is clearly displayed in real time on an easy-to-read digital screen, allowing for accurate monitoring throughout use. A powerful motor and magnetic system enable stirring speeds of up to 2000 rpm and efficient mixing of volumes up to 15 L. The square plate design accommodates a wide range of vessel sizes, offering flexibility for various laboratory applications.

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
04807-64	Cole-Parmer SHP-200D-S Undergrad Digital Stirring Hot Plate, 15 x 15 cm, Aluminum; 230V	£480.93	£339.35



View Our Website

HOTPLATES



OHAUS HOTPLATE STIRRERS

The OHAUS Guardian hotplate stirrer range—models 30680270, 30680279, 30541642 and 30541648—provides a scalable solution for laboratory heating and stirring, progressing from basic to advanced performance. The Guardian 2000 offers simple, reliable analogue control for everyday tasks, while the Guardian 3000 introduces digital control, improved accuracy and enhanced safety features such as a hot-top indicator and optional probe compatibility. At the top end, the Guardian 5000 models deliver high-performance heating and stirring with advanced control modes, superior temperature stability and additional functions including timers and enhanced safety systems, with the 30541648 also offering a larger plate for higher-volume applications.

Key Features

- Scalable range (2000–5000 series) to suit basic, routine, and advanced laboratory applications
- Powerful heating up to ~500–550 °C with rapid heat-up performance
- Magnetic stirring up to 15–22 L capacity for a wide range of sample volumes



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
30680270	Hotplate-Stirrer e-G21HS07C 230V 50Hz EU	£392.00	£289.92
30680279	Hotplate-Stirrer e-G31HS07C 230V EU	£510.00	£376.99
30541642	Hotplate-Stirrer e-G52HS07C 230V EU	£680.00	£502.65
30541648	Hotplate-Stirrer e-G52HS10C 230V EU	£925.00	£683.76



Cole-Parmer HP-200-C Undergrad Analog Hot Plate, 15 x 15 cm, Ceramic; 230 VAC

Designed with microprocessor-controlled dual thermocouples, these stirring hot plates provide precise temperature control along with reliable protection against overheating. The actual surface temperature is clearly displayed in real time on an easy-to-read digital screen, allowing for accurate monitoring throughout use. A powerful motor and magnetic system enable stirring speeds of up to 2000 rpm and efficient mixing of volumes up to 15 L. The square plate design accommodates a wide range of vessel sizes, offering flexibility for various laboratory applications.

Key Features

- Compact, stackable design takes up less benchtop and storage space
- Hot warning light remains on even when unplugged
- BioCote® antimicrobial protection is built in for a cleaner, more hygienic product
- Built with tough aluminum casing for longevity
- Designed to prevent accidental spillage on control panel
- 3 year warranty
- New, improved version of UC150

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
04807-50	Undergrad Analog Hot Plate, 15 x 15 cm, Ceramic; 230 VAC	£361.32	£254.96
04807-52	Undergrad Analog Hot Plate, 15 x 15 cm, Aluminum; 230 VAC	£334.55	£236.06

Cole-Parmer®



View Our Website

STIRRERS

Cole-Parmer®

UNDERGRAD ANALOG STIRRER, 15X15 CM, CERAMIC, 230VAC

KEY FEATURES

- BUILT WITH TOUGH ALUMINUM CASING FOR LONGEVITY
- DESIGNED TO PREVENT ACCIDENTAL SPILLAGE ON CONTROL PANEL
- IP32 RATED FOR WATER AND DUST RESISTANCE
- BIOCOTE® ANTIMICROBIAL PROTECTION IS BUILT IN FOR A CLEANER, MORE HYGIENIC PRODUCT
- 3 YEAR WARRANTY
- NEW, IMPROVED VERSION OF UC151



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
4807-54	Undergrad Analog Stirrer, 15 x 15 cm, Ceramic; 230 VAC	£303.25	£213.97

COLE-PARMER ST-200-S UNDERGRAD ANALOG STIRRER, 15 X 15 CM, STAINLESS STEEL; 230 VAC

CODE	4807-56
DESCRIPTION	Undergrad Analog Stirrer, 15 x 15 cm, Stainless Steel; 230 VAC
LIST PRICE	£289.00
OFFER PRICE	£210.96



KEY FEATURES

- BUILT WITH TOUGH ALUMINUM CASING FOR LONGEVITY
- DESIGNED TO PREVENT ACCIDENTAL SPILLAGE ON CONTROL PANEL
- IP32 RATED FOR WATER AND DUST RESISTANCE
- BIOCOTE® ANTIMICROBIAL PROTECTION IS BUILT IN FOR A CLEANER, MORE HYGIENIC PRODUCT
- 3 YEAR WARRANTY
- NEW, IMPROVED VERSION OF US151



View Our Website

SARTORIUS BCE AND BCA BALANCES ENTRIS® II

EVERY SARTORIUS BALANCE OFFERS QUALITY, VALUE, AND CONSISTENCY.

NO MATTER WHAT YOU'RE WEIGHING, THE NEW ENTRIS® II BALANCE IS ALWAYS THE RIGHT CHOICE.

OFFERING UNRIVALLED VALUE AND BACKED BY ALMOST 150 YEARS OF GERMAN ENGINEERING EXPERTISE, THE ENTRIS® II COMES IN TWO PRODUCT LINES, SO YOU CAN FIND THE BALANCE THAT MEETS YOUR SPECIFIC WEIGHING NEEDS.



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
BCA124-1S	Analytical balance 120g 0.1mg	£2,890.00	£2,184.00
BCA124i-1S	Analytical balance int.cal. 120g 0.1mg	£3,120.00	£2,352.00
BCA224-1S	Analytical balance 220g 0.1mg	£3,170.00	£2,402.00
BCA224i-1S	Analytical balance int.cal. 220g 0.1mg	£3,400.00	£2,570.40

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
BCE124-1S	Analytical balance 120g 0.1mg	£2,090.00	£1,375.92
BCE124i-1S	Analytical balance int.cal. 120g 0.1mg	£2,310.00	£1,514.24
BCE224-1S	Analytical balance 220g 0.1mg	£2,230.00	£1,470.56
BCE224i-1S	Analytical balance int.cal. 220g 0.1mg	£2,420.00	£1,587.04
BCE2201-1S	Precision balance 2200g 0.1g	£1,200.00	£786.24
BCE2201i-1S	Precision balance int.cal. 2200g 0.1g	£1,450.00	£953.68
BCE2202-1S	Precision balance 2200g 10mg	£1,460.00	£960.96
BCE2202i-1S	Precision balance int.cal. 2200g 10mg	£1,700.00	£1,106.56



View Our Website

OHAUS ADVENTURER BALANCE RANGE



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
30852308	Semi-Micro Balance AX85	4,380.00	3120.94
30852307	Semi-Micro Balance AX125D	4,515.00	3217.13
30852309	Semi-Micro Balance AX225D	5,345.00	3808.54
30122610	Analytical Balance AX124	2,955.00	2105.56
30122657	Analytical Balance AX124/E	2,620.00	1866.86
30122611	Analytical Balance AX224	3,110.00	2216.01
30122612	Analytical Balance AX224/E	2,770.00	1973.74
30122613	Analytical Balance AX324	3,145.00	2240.95
30127960	Precision Balance AX223	2,090.00	1489.21
30127961	Precision Balance AX223/E	1,735.00	1236.26
30127962	Precision Balance AX423	2,310.00	1645.97
30127963	Precision Balance AX423/E	1,950.00	1389.46
30122614	Precision Balance AX523	2,745.00	1955.93
30122615	Precision Balance AX523/E	2,395.00	1706.54
30805890	Precision Balance AX623	2,745.00	1955.93
30805891	Precision Balance AX623/E	2,395.00	1706.54
30206079	Precision Balance AX422	1,700.00	1211.32
30206080	Precision Balance AX422/E	1,320.00	940.56
30206081	Precision Balance AX822	1,875.00	1336.02
30206082	Precision Balance AX822/E	1,515.00	1079.50
30127966	Precision Balance AX1502	2,070.00	1474.96
30127967	Precision Balance AX1502/E	1,710.00	1218.45
30127968	Precision Balance AX2202	2,105.00	1499.90



OHAUS ADVENTURER BALANCE RANGE



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
30127969	Precision Balance AX2202/E	£1,755.00	£1,250.51
30122616	Precision Balance AX4202	£2,285.00	£1,628.16
30122617	Precision Balance AX4202/E	£1,920.00	£1,368.08
30122618	Precision Balance AX5202	£2,745.00	£1,955.93
30805893	Precision Balance AX6202/E	£2,395.00	£1,706.54
30805892	Precision Balance AX6202	£2,745.00	£1,955.93
30206083	Precision Balance AX2201	£1,630.00	£1,161.44
30206084	Precision Balance AX2201/E	£1,270.00	£904.93
30127970	Precision Balance AX4201	£1,870.00	£1,332.46
30127971	Precision Balance AX4201/E	£1,510.00	£1,075.94
30122619	Precision Balance AX8201	£2,050.00	£1,460.71
30122620	Precision Balance AX8201/E	£1,685.00	£1,200.63
30805894	Precision Balance AX12001	£2,150.00	£1,531.97
30805895	Precision Balance AX12001/E	£1,770.00	£1,261.20





Haier

HAIER LOW ENERGY ULT FREEZER - LCD TOUCHSCREEN

KEY FEATURES:

- 10-INCH HIGH-PERFORMANCE LCD CAPACITIVE SCREEN, SENSITIVE TOUCH OPERATION; INTUITIVE DISPLAY OF INSIDE TEMPERATURE, AMBIENT TEMPERATURE, INPUT VOLTAGE AND OTHER DATA AND TEMPERATURE CURVES.
- THE SINGLE VERSION OF THE INVENTORY MANAGEMENT ALLOWS USERS WHO DO NOT HAVE A RACKING SYSTEM WITHIN THE FREEZER TO EASILY RECORD ITEM LOCATIONS, ENTRY AND EXIT RECORDS.
- STANDARD USB INTERFACE, CAPABLE OF STORING DATA FOR 10 YEARS
- STANDARD PASSWORD LOCK, OPTIONAL FINGERPRINT MODULE, PUNCH CARD MODULE, REALIZING SAFE AND SECURE MULTI-USER MANAGEMENT.
- MULTIPLE ALARM FUNCTIONS INCLUDE HIGH TEMPERATURE ALARM, LOW TEMPERATURE ALARM, SENSOR ERROR ALARM, POWER FAILURE ALARM, LOW BATTERY ALARM, DOOR AJAR ALARM AND HIGH AMBIENT TEMPERATURE ALARM.
- ADJUSTABLE TEMPERATURE SET-POINT BETWEEN -40°C ~ -86°C AND THE ABILITY TO MONITOR TEMPERATURE CURVES IN REAL TIME, EQUIPPED WITH SOPHISTICATED TEMPERATURE SAFETY ALARM SYSTEM TO ENSURE PRODUCT SAFETY.

CODE

DW-86L579

DESCRIPTION

LT FREEZER, UPRIGHT, ENERGY EFFICIENT, TOUCH SCREEN, -86 DEGREES CELSIUS, 579L CAPACITY

LIST PRICE

£13,388.73

OFFER PRICE

£9,164.00



View Our Website

HAIER ULTRAECO ULTRA LOW IMPACT ULT FREEZER

LOW IMPACT, HIGH PERFORMANCE

- NEW GREEN INSULATION WITH OPTIMISED ULTRA-LOW IMPACT COOLING SYSTEM
- FAST RE-ENTRY FOR SEAMLESS SAMPLE ACCESS
- ULTRA-LOW ENERGY



CODE	DW-86L730BPT
DESCRIPTION	ULT FREEZER, UPRIGHT, ULTRAECO, ULTRA ENERGY EFFICIENT, TOUCH SCREEN, -86 DEGREES CELSIUS, 730L CAPACITY
LIST PRICE	£18,282.55
OFFER PRICE	£11,866.00

ULTRAECO ULTRA LOW IMPACT ULT FREEZER

THE ULTRAECO IS OUR NEXT GENERATION ULTRA LOW IMPACT ULT FREEZER. THE COMBINATION OF CUTTING EDGE REFRIGERATION TECHNOLOGY AND NEW SUSTAINABLE INSULATION PROVIDES MARKET LEADING PERFORMANCE WITH WORLD LEADING ENERGY EFFICIENCY. DELIVERING RESULTS WITHOUT COMPROMISE.


CODE	DW-86L830BPT
DESCRIPTION	ULT FREEZER, UPRIGHT, ULTRAECO, ULTRA ENERGY EFFICIENT, TOUCH SCREEN, -86 DEGREES CELSIUS, 830L CAPACITY
LIST PRICE	£19,171.64
OFFER PRICE	£12,443.00



**OVER
35% OFF**



View Our Website



LEC UNDERCOUNTER REFRIGERATOR

1158L UNDERCOUNTER REFRIGERATOR WITH SOLID DOOR

CODE LSFRT110

DESCRIPTION LEC 158L Undercounter Refrigerator with solid door, spark free, lockable

OFFER PRICE £472.00



- INTEGRATED LOCK - SUPPLIED WITH 2 KEYS FOR ADDED SECURITY AND ACCESS CONTROL TO KEEP CONTENTS SAFE
- SPARK FREE - REMOVED ITEMS SUCH AS FANS AND LIGHTS SO NORISK OF IGNITING FLAMMABLE LIQUIDS
- 2 YEAR PARTS & LABOUR WARRANT

LEC UNDERCOUNTER FREEZER

- ENHANCED SAFETY FOR FLAMMABLE MATERIALS: SPARK-FREE INTERIOR WITH NO FANS OR LIGHTS, ELIMINATING IGNITION RISKS.
- TEMPERATURE ACCURACY: AIR TEMPERATURE PROBE ENSURES RELIABLE MONITORING OF INTERNAL CONDITIONS.
- SECURE STORAGE: INTEGRATED LOCK WITH TWO KEYS PREVENTS UNAUTHORISED ACCESS.
- SPECIALIST DESIGN: BUILT SPECIFICALLY FOR MEDICAL AND LABORATORY USE, PRIORITISING SAFETY AND RELIABILITY.



CODE LSFST110

DESCRIPTION LEC 98L UNDERCOUNTER FREEZER WITH SOLID DOOR, SPARK FREE, LOCKABLE

OFFER PRICE £472.00



LIEBHERR UNDERCOUNTER SPARK FREE FRIDGE



LIEBHERR

CODE	LKUexv1610
DESCRIPTION	LIEBHERR UNDERCOUNTER SPARK FREE FRIDGE +3-+16C 141L FAN ASSISTED WITH COMFORT CONTROLLER
LIST PRICE	£1,579.00
OFFER PRICE	£1,184.00

KEY FEATURES:

- **SAFETY & COMPLIANCE:** EXPLOSION-PROOF (ATEX 95) DESIGN MAKES IT SUITABLE FOR STORING FLAMMABLE AND EXPLOSIVE MATERIALS, WITH BUILT-IN ACOUSTIC AND VISUAL TEMPERATURE ALARMS FOR ADDED PROTECTION.
- **RELIABLE MONITORING:** PRECISE DIGITAL TEMPERATURE DISPLAY AND COMPATIBILITY WITH EXTERNAL MONITORING SYSTEMS ENSURE CONSTANT OVERSIGHT WITHOUT OPENING THE UNIT.
- **PERFORMANCE:** FORCED AIR COOLING PROVIDES FAST COOLING AND MAINTAINS A STABLE, EVEN INTERNAL TEMPERATURE.
- **QUALITY ASSURANCE:** TESTED FOR TEMPERATURE STABILITY (EN 60068-3) TO GUARANTEE CONSISTENT PERFORMANCE.
- **USER CONVENIENCE:** INCLUDES AN ALARM TEST FUNCTION AND REVERSIBLE DOOR HINGES FOR FLEXIBLE INSTALLATION.

LIEBHERR UNDERCOUNTER FREEZER

CONFORMS TO THE ATEX DIRECTIVE AND MEETS EXACTING REQUIREMENTS IN TERMS OF SAFETY, TEMPERATURE STABILITY, AND EASE OF USE: THIS UNDERCOUNTER FREEZER FOR LABORATORIES IS SUITABLE FOR THE RELIABLE STORAGE OF CLOSED CONTAINERS CONTAINING COMBUSTIBLE SUBSTANCES THAT COULD FORM EXPLOSIVE ATMOSPHERES.



CODE	SFTfg 1501-001
DESCRIPTION	UNDERBENCH ATEX FREEZER WITH SHELVES
LIST PRICE	£1,479.00
OFFER PRICE	£1,110.00



View Our Website

LIEBHERR FREE STANDING FRIDGE



CODE

SRFFg4001

DESCRIPTION

SCIENTIFIC ATEX TALL FRIDGE +3/+16C 394 LITRE

LIST PRICE

£2,249.00

OFFER PRICE

£1,688.00

KEY FEATURES:

- **EXPLOSION-PROOF (ATEX 95 CERTIFIED):** SAFE FOR STORING FLAMMABLE AND EXPLOSIVE MATERIALS, MEETING STRICT EU SAFETY STANDARDS.
- **PRECISE TEMPERATURE CONTROL:** EXTERNAL DIGITAL DISPLAY SHOWS EXACT INTERNAL TEMPERATURE WITHOUT OPENING THE DOOR.
- **BE CONNECTED TO EXTERNAL MONITORING SYSTEMS** FOR CONTINUOUS OVERSIGHT.
- **PROVEN RELIABILITY:** TESTED FOR TEMPERATURE STABILITY TO RECOGNISED STANDARDS (EN 60068-3).
- **FLEXIBLE INSTALLATION:** REVERSIBLE DOOR HINGES ALLOW ADAPTATION TO DIFFERENT SPACES.
- **ECO-CONSCIOUS:** USES ENVIRONMENTALLY FRIENDLY REFRIGERANT.

LIEBHERR FREE STANDING FREEZER

THIS LABORATORY FRIDGE DOES NOT JUST GUARANTEE CONSTANT TEMPERATURES FOR STORING SENSITIVE SUBSTANCES, IT ALSO BOASTS EXTENSIVE MONITORING AND ALARM FUNCTIONS. IN ADDITION, THE ROBUST POLYSTYRENE INTERIOR IS MANUFACTURED IN ONE PIECE WITHOUT ANY JOINTS, MAKING IT QUICK AND EASY TO CLEAN

KEY FEATURES:

- **EXCELLENT TEMPERATURE STABILITY (± 5 °C):** PROTECTS SENSITIVE SUBSTANCES FROM FLUCTUATIONS, PRESERVING QUALITY AND INTEGRITY.
- **ADVANCED SAFETY (ATEX COMPLIANT):** SUITABLE FOR STORING FLAMMABLE MATERIALS IN SEALED CONTAINERS, MEETING STRICT EU SAFETY STANDARDS.
- **LOW MAINTENANCE COOLING (SMARTFROST):** MINIMISES ICE BUILD-UP, REDUCES DEFROSTING FREQUENCY, AND IMPROVES ENERGY EFFICIENCY.



Code

Description

List Price

Offer Price

SFFfg4001

SCIENTIFIC FREEZER -9C/-30C 316L ATEX

£2,249.00

£1,688.00

SFFfg4001 740

SCIENTIFIC FREEZER (WITH DRAWERS) -9/-30C 316L ATEX

£2,499.00

£1,877.00



View Our Website

LABCOLD UPRIGHT FRIDGE



- CONFIGURATION FREE STANDING
- CAPACITY (LITRES) 439
- HEIGHT X WIDTH X DEPTH
• 1875MM X 600MM X 700MM
- WEIGHT 93KG
- SHELVES 6
- TEMPERATURE RANGE +0°C TO +10°C
- DOOR TYPE SOLID
- LOCK TYPE KEY LOCK WITH TWO KEYS SUPPLIED
- AUTO DEFROST
- SPARKFREE SAFETY INTERIOR
- 6 ADJUSTABLE SHELVES
- DOOR LOCK
- TEMPERATURE ALARM AND DISPLAY
- HIGH/LOW TEMPERATURE ALARM
- DOOR OPEN ALARM
- ACCESS PORT
- INTERNAL FAN
- 2 YEAR PARTS AND LABOUR WARRANTY (UK ONLY)

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
RLPR1517	Labcold Upright Fridge	£2,125.00	£1,593.00

LABCOLD UPRIGHT FREEZER

- CONFIGURATION FREE STANDING
- CAPACITY (LITRES) 406
- HEIGHT X WIDTH X DEPTH
• 1875MM X 600MM X 700MM
- WEIGHT 97KG
- SHELVES 6
- TEMPERATURE RANGE -18°C TO -25°C
- DOOR TYPE SOLID
- LOCK TYPE KEY LOCK WITH TWO KEYS
- ACCESS PORT
- SPARKFREE SAFETY INTERIOR
- MANUAL DEFROST
- DOOR LOCK
- SELF CLOSING DOOR
- TEMPERATURE DISPLAY
- HIGH/LOW TEMPERATURE ALARM
- DOOR OPEN ALARM
- 6 SHELVES
- FAN CIRCULATED
- 2 YEAR PARTS AND LABOUR WARRANTY (UK ONLY)



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
RLVF1517	Labcold Upright Freezer	£2,490.00	£1,866.00



View Our Website

CHEMGENE MEDLAB

CHEMGENE HAS ADOPTED NEW TECHNOLOGY TO ALIGN WITH UPCOMING REGULATIONS. (THE MEDLAB RANGE IS NOW COMPLIANT WITH THE BIOCIDAL PRODUCTS REGULATION (BPR) IN EUROPE AND THE UK, AS WELL AS THE CLASSIFICATION, LABELLING AND PACKAGING (CLP) AND REACH REGULATIONS).

CHEMGENE'S UNIQUE NEW FORMULATION USES A COMBINATION OF ACTIVE SUBSTANCES, RESULTING IN LESS RISK OF ANTIMICROBIAL RESISTANCE

- DENATURES DNA/RNA
- BIODEGRADABLE INGREDIENTS
- NON-CORROSIVE*
- ALCOHOL-FREE FORMULATION



NEW CODE	DESCRIPTION	OFFER PRICE
SKU077A	CHEMGENE MEDLAB Multi-Surface Disinfectant Concentrate- Eucalyptus 5 L	£47.50
SKU075A	CHEMGENE MEDLAB Multi-Surface Disinfectant Concentrate- Un-fragranced 5 L	£47.50
SKU095A	CHEMGENE MEDLAB Multi-Surface Disinfectant – Eucalyptus 750ml x 6	£60.00

COLOUR CODED PH BUFFERS 1L



CODE	DESCRIPTION	OFFER PRICE
1040C.	Buffer solution , colour coded pH 4.0 ± 0.01 @20°C	£15.50
1070C.	Buffer solution, colour coded pH 7.0 ± 0.01 @20°C	£15.50
1100C.	Buffer solution , colour coded pH 10.0 ± 0.01 @20°C	£15.50

SYRINGE FILTERS

PRODUCT CODE	DESCRIPTION	PACK SIZE	OFFER PRICE
SFPES025045SI	25mm, PES Sterile syringe filter, 0.22um	100	£44.00
SFPES025022SI	25mm, PES Sterile syringe filter, 0.45um	100	£44.00







AZENTA
LIFE SCIENCES

PCR PLATES AND SEALS

FRAMESTAR® ROCHE STYLE

CODE	Description	LIST PRICE	OFFER PRICE
4ti-0951	FrameStar 96 Well Semi-Skirted PCR Plate, Roche style, white PP wells, clear PC frame, low profile, cut corner H12, 50 plates per case	£335.00	£146.87
4ti-0560	qPCR Adhesive Seal, optically clear film, pressure activated adhesive, 100 sheets (140 x 77mm) per case	£192.00	£144.38
COMBI PACK			
4ti-0952	FrameStar 96 Well Semi-Skirted PCR Plate, Roche style, plus qPCR Seal, <i>4TI-0951 plus 4TI-0560</i> , combi pack, 50 plates and seals per case	£407.00	£211.88



FRAMESTAR® ABI® STYLE

SUITABLE FOR USE WITH QUANTSTUDIO


CODE	Description	LIST PRICE	OFFER PRICE
4ti-0730	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, purple PC frame, with upstand, high profile, cut corner A12, 50 plates per case	£310.50	£133.72
4ti-0565	Optically Clear Windowed qPCR Seal, adhesive film with 96 optically clear windows, 100 sheets (133 x 76mm) per case	£157.50	£90.00
COMBI PACK			
4ti-0730/C/0565	FrameStar 96 Well Semi-Skirted PCR Plate, ABI style, clear PP wells, clear PC frame, with upstand, high profile, cut corner A12, 50 plates per case <u>plus 4ti-0565</u>	£380.50	£167.02

AZENTA SEMI AUTOMATED HEAT SEALER

SEMI-AUTOMATED SHEET HEAT SEALER COMPATIBLE WITH A WIDE RANGE OF SEALS AND MICROPLATES, PCR AND ASSAY PLATES, AS WELL AS STORAGE PLATES OF DIFFERING DESIGNS AND HEIGHTS

CODE	DESCRIPTION	OFFER PRICE
59-2000	Semi-Automated Sheet Heat Sealer, includes adapter (59-2001)	£3,576.27





AZENTA
LIFE SCIENCES

ADHESIVE SEALS

CODE		LIST PRICE	OFFER PRICE
PCR SEAL			
4TI-0500	PCR Seal, clear adhesive film, strong adhesive, 100 sheets (135 x 80mm) per case	£114.50	£69.38
QPCR SEAL			
4TI-0560	qPCR Adhesive Seal, optically clear film, pressure activated adhesive, 100 sheets (140 x 77mm) per case	£192.50	£144.38
PCR FOIL SEAL			
4TI-0550	PCR Foil Seal, pierceable adhesive aluminium foil, strong adhesive, 100 sheets (130 x 80mm) per case	£90.00	£55.62
DMSO RESISTANT FOIL SEAL			
4TI-0512	DMSO Resistant Foil Seal, peelable adhesive foil, strong adhesive, high solvent resistance, 100 sheets (122 x 80mm) per case	£318.50	£215.61
PIERCEABLE FILM			
4TI-0566	Pierceable Film, strong adhesive seal with 96 cross-cut windows, 100 sheets (135 x 77mm) per case	£157.50	£70.00

HEAT SEALS



HEAT SEALS

4TI-0541 Clear Heat Seal, peelable heat sealing film, 100 sheets (125 x 78mm) per case

PIERCE HEAT SEAL

4TI-0531 Pierce Heat Seal, pierceable heat sealing foil, 100 sheets (125 x 78mm) per case

FOIL HEAT SEAL

4TI-0536 Foil Heat Seal, aluminium heat sealing foil, 100 sheets (125 x 78mm) per case

GAS PERMEABLE HEAT SEAL

4TI-0597 Gas Permeable Heat Seal, heat sealing membrane, 100 sheets (125 x 78mm) per case



DMSO RESISTANT PEEL HEAT SEAL

OFFER PRICE

4TI-0587 DMSO Resistant Peel Heat Seal, solvent resistant heat sealing foil, 1 perforated roll with 100 sheets (125 x 78mm)

£55.62



DNA GEL ELECTROPHORESIS

CLEAR SIGHT MSMAXIDUO



KEY FEATURES

THESE UNITS ARE PRIMARILY DESIGNED FOR SEPARATING HIGH NUMBERS OF SAMPLES FROM PCR OR CLONING:

- AVAILABLE WITH 20 X 25CM, 20 X 20CM, 20 X 15CM OR 20 X 10CM GEL TRAYS
- RUN UP TO 550 SAMPLES
- LOW BUFFER VOLUMES



CODE	DESCRIPTION	OFFER PRICE
CLV-MSMAXIDUO	multiSUB Maxi, Maxi Horizontal Electrophoresis System (with 20x10cm and 20x20cm gel trays)	£670.08

POWERPRO 300 POWER SUPPLY, 300V, 700MA, 150W


- PERFECT FOR ROUTINE ELECTROPHORESIS APPLICATIONS
- 300V, 700MA, 150W
- 5 PAIRS OF OUTLETS FOR MULTIPLE GEL TANKS



POWERPRO 300

CODE	DESCRIPTION	OFFER PRICE
CLV-POWERPRO300	PowerPRO 300 Power Supply, 300V, 700mA, 150W	£466.15

MULTISUB MINI ONE



KEY FEATURES

- ALL IN ONE HORIZONTAL ELECTROPHORESIS SYSTEM
- IN BUILT POWER SUPPLY WITH 35V, 50V AND 100V SETTINGS
- TIMER FUNCTION FOR RUNS FORM 0 - 99 MINUTES
- 2 GEL SIZES AND REVERSIBLE COMB OPTIONS

CODE	DESCRIPTION	OFFER PRICE
CLV-MSMINIONE	multiSUB MINI ONE, All in One Horizontal Electrophoresis System	£373.85

CONSUMABLES

CODE	DESCRIPTION	OFFER PRICE
CLV-CSL-MDNA-1KBPLUS	1Kb PLUS DNA ladder, 100-10Kb	£27.69
CLV-CSL-MDNA-50BP	50bp DNA ladder, 50-1500bp	£49.23
CLV-CSL-AG500	Agarose, Low EEO 500g	£215.17
CLV-CSL-RUNSAFE	runSAFE Safe DNA staining reagent	£36.92
TAE50X1L	Cleaver Buffer Tris-Acetate EDTA (TAE) Running Buffer 50 x 1L	£83.98



PROTEIN GEL ELECTROPHORESIS

KEY FEATURES

- RUN 1 – 2 PRECAST OR HAND-CAST MINI PAGE GELS
- INTUITIVE CASTING AND GEL CLAMPING SYSTEMS
- BONDED SPACERS FOR REPRODUCIBILITY IN CASTING
- INJECTION MOULDED CONSTRUCTION PREVENTS LEAKS AND EXTENDS LIFETIME.

PROPAGE 2



CODE	DESCRIPTION	OFFER PRICE
CVS82BSYS1.5	proPAGE 2 – mini PAGE system for 2 gels	£581.21

MINI VERTICAL BLOTTER



KEY FEATURES

THE PROPAGE MINIBLOT SYSTEMS ARE COMPLETE STAND-ALONE UNITS FOR WET TRANSFER, COMPRISING TANK AND LID, ELECTROBLOTTING INSERT, AND CASSETTES AND FIBRE PADS.

- BUILT-IN COOLING UNITS ABSORB HEAT EFFICIENTLY.
- COMPLETE RUN-TIME FLEXIBILITY
- INJECTION MOLDED MINI UNIT – DURABLE, LEAK-PROOF ENVIRONMENT FOR COMPLETE SAFETY AND LONG LIFE.



CODE	DESCRIPTION	OFFER PRICE
CVS8-SBU	miniBLOT – Mini Vertical Blotter	£493.75

CONSUMABLES

CODE	DESCRIPTION	OFFER PRICE
CSL-BP1010	Blot-Absorbent Filter paper, 10 x 10cm, 1-mm thick pack of 50	£37.33
CSL-RNC45	Nitrocellulose Blotting Membrane Roll 0.2 uM pore size (30 cm x 3 metres)	£241.65
CSL-BBL	Blue Wide Range Protein Ladder (10- 245kDa) (500 µL)	£107.14
CLV-ECLONE	LumiGO ECL substrate kit: 125ml Luminol/enhancer solution; 125ml Peroxide solution	£195.75
PVDF0.45R	Roll PVDF 30cm X 3m, 0.45um	£241.65
WBLUF0100	Immobilon Forte Western HRP Substrate, 100ml Western blotting, dot/slot blotting	£148.73



CODE		LIST PRICE	OFFER PRICE
F25108T	Immobilon Classico Western HRP Substrate	£105.00	£95.00
F25109T	Immobilon Crescendo Western HRP Substrate	£125.00	£113.00
F25110T	Immobilon Forte Western HRP Substrate	£160.00	£144.00
F25111T	Immobilon Western Chemiluminescent HRP Substrate	£145.00	£131.00
F25112T	Immobilon ECL Ultra Western HRP Substrate	£239.00	£216.00
F25113T	Immobilon-E PVDF Membrane, 7x8.4cm, 0.45um	£452.00	£409.00
F25114T	Immobilon-P PVDF Membrane, 7x8.4cm, 0.45um	£464.00	£419.00
F25115T	Immobilon-PSQ PVDF Membrane, 7x8.4cm, 0.2um	£457.00	£413.00
F25116T	Immobilon-FL PVDF Membrane, 7x8.4cm, 0.45um	£99.50	£90.00
F25117T	Immobilon Blotting Filter Paper, 8.5x13.5cm	£129.00	£116.00

GS BIOTECH 170L CO2 INCUBATOR

WITH HIGH TEMPERATURE DECONTAMINATION

Engineered for reliability and precision, the 170L laboratory incubator is ideal for medium to large laboratories requiring consistent environmental control. Designed to support a wide range of applications, it combines advanced monitoring, flexible configuration, and user-friendly operation to ensure optimal performance and sample protection.

Key Features

- Stable Chamber Environment – Gentle, fan-free heating ensures uniform conditions while reducing contamination risk
- Precise Control – Programmable CO₂ (0.2–20%) and temperature (ambient +4°C to 50°C) for accurate conditions
- Infrared CO₂ Sensor – Self-referencing technology for reliable, consistent readings
- Easy to Clean Design – Seamless chamber with rounded corners and removable components
- Enhanced Visibility – Sealed inner glass door allows monitoring without disrupting conditions
- Flexible Configuration – Wide range of options to tailor the unit to your application
- Optional Oxygen Control – Accurate low-oxygen environments (0.1–19%) for specialised work
- Optional Decontamination Cycle – 120°C sterilisation programme for improved hygiene
- Smart Monitoring Options – Data logging and reporting available for QA and compliance
- User-Friendly Interface – Large LCD display with clear, easy-to-read controls
- Operational Safety – Optional humidity tray alarm helps prevent downtime
- Practical Design – Non-tip shelves and 25mm access port for added convenience
- High Capacity – 170L chamber ideal for busy laboratory environments
- Warranty Included – 3-year warranty, extendable for added peace of mind



COST	DESCRIPTION	OFFER PRICE
170/GOLD/HT	170L INCUBATOR WITH HIGH TEMPERATURE DECONTAMINATION PLUS HT DECON	£4,200.00



pH METERS

HANNA HI-2002 EDGE PH METER

SINGLE PARAMETER PH METER FEATURES USER-FRIENDLY TABLET DESIGN WITH BENCH, PORTABLE OR WALL MOUNT OPTIONS FOR ULTIMATE FLEXIBILITY. SWITCH FROM PORTABLE TO BENCH IN THE SAME PROCEDURE OR CLEAR THE SURFACE FOR EASIER WORKING AND HANG ON THE WALL.

- RANGE -2.00 TO 16.00PH;
- PH
- DIGITAL SMART ELECTRODE COMPATIBLE
- ACCURACY: ± 0.01 PH



	CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
	HI-2002-02	HI-2002 EDGE PH METER	£551.30	£463.09

HANNA HI-2600 BENCHTOP PH-CONDUCTIVITY-DO METER



A compact, high-performance multiparameter meter designed for fast, accurate measurement of pH, conductivity, dissolved oxygen, and ORP using smart digital probes. Ideal for both laboratory and portable use, with intuitive operation and reliable data management.

Key Features:

- Automatic Probe Recognition – Instantly identifies connected sensors for seamless setup
- All-in-One Measurement – pH, EC, TDS, salinity, DO, and ORP in one device
- Smart Data Logging – Store up to 1000 records with full GLP data and easy USB transfer
- Portable & Versatile – Rechargeable battery with up to 8-hour use, lab or field ready



	CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
	HI-2600-02	pH, Conductivity, Disolved Oxygen and ORP Meter	£495.00	£415.80



COLE PARMER JENWAY STANDARD DIGITAL PH METER KIT

Benchtop pH meters feature up to three decimal place resolution, support up to three pH calibration points, and offer manual and automatic buffer recognition for DIN, JIS, and NIST buffers. A large LCD displays status and results of pH and temperature readings. Easily export data to a computer or printer via an RS-232 or an analog interface.

Advanced digital pH meters help ensure compliance with Good Laboratory Practice (GLP) guidelines with calibration reminder interval, alarm output, and security code setup features. Both standard and advanced pH meters are available as meter only or as a kit that includes a pH electrode, ATC probe, and buffers.

INCLUDES

GLASS PH ELECTRODE, ATC PROBE, ELECTRODE HOLDER, PH BUFFERS (4, 7, AND 10), BNC SHORT PLUG, AND POWER CORD WITH UK PLUG.

KEY FEATURES

- SUPPORTS ONE, TWO, OR THREE POINT PH CALIBRATION
- UP TO THREE DECIMAL PLACE RESOLUTION FOR ACCURATE READINGS
- KIT INCLUDES A PH ELECTRODE, ATC PROBE, AND PH BUFFERS
- STORES UP TO 32 READINGS
- 230 VAC OPERATION
- COMES WITH UK PLUG
- 3 YEAR WARRANTY



SPECIFICATION:

Specification	Details
Meter Type	Benchtop
Display Type	Backlit Custom LCD
Connections	BNC
pH Range	-2 to 19.999
pH Resolution	0.001, 0.01, 0.1
pH Accuracy	±0.003
Buffer Recognition	DIN, NIST, JIS, Custom
Calibration	1, 2 or 3 Point
mV Range	-1999.9 to 1999.9
Temperature Range	-10 to 105°C (14 to 221°F)
Temperature Accuracy	±1°F (±0.5°C)
Interface	Analog & RS-232
Data Logging	32 Data Sets
Power	230 VAC
Dimensions (W × L × H)	21 × 25 × 5.5 cm
Power	25 × 21 × 21 (W) × 5.5 × (H)
Warranty	3 Year

CODE

DESCRIPTION

LIST PRICE

OFFER PRICE

35459-02

Jenway 3510 Standard Digital pH
Meter Kit; 230VAC/UK

£735.16

£526.97



Cole-Parmer®

JENWAY

SPECTROPHOTOMETER



TOP FEATURES

- SCANNING DIODE ARRAY TECHNOLOGY
- COLOUR TOUCHSCREEN NAVIGATION
- SMALL FOOTPRINT AND LIGHTWEIGHT
- FAST SCAN SPEED
- ENGLISH, FRENCH AND GERMAN LANGUAGE OPTIONS

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
83056-02	UV/Visible 72 Series Diode Array Scanning Spectrophotometer	£2,999.77	£2150.24

VACUUBRAND®

VACUUBRAND BVC FLUID ASPIRATION SYSTEM

- POWERFUL AND LONG-LASTING CHEMISTRY DIAPHRAGM PUMP FOR UNIVERSAL ECONOMICAL USE
- EXTREMELY COMPACT, SPACE SAVING AND LOW NOISE AND THUS THE IDEAL TOOL IN THE WORKPLACE
- EASY OPERATION DUE TO ELECTRONICALLY ADJUSTABLE SUCTION POWER VIA TOUCH PANEL, FOR STRONG OR SENSITIVE SUCTION
- AUTOCLAVABLE COLLECTION BOTTLE DESIGNED TO REDUCE AEROSOL AND FOAM FORMATION; WITH HYDROPHOBIC PROTECTION FILTER 0.2 MICROMETER, AUTOCLAVABLE QUICK-COUPPLINGS AVAILABLE AS ACCESSORY
- READY TO ACCEPT A SECOND ASPIRATION HAND CONTROLLER (VHC PRO) TO SUPPORT TWO USERS



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
5206	BVC Control Fluid Aspiration System	£1,790.00	£1,478.54
5215	BVC Professional Fluid Aspiration System	£2,390.00	£1,974.14



JB ACADEMY UNSTIRRED WATER BATHS

**lids not included; available to purchase separately*

KEY FEATURES

- AMBIENT +5°C TO 95°C
- UNIQUE SET AND FORGET® TECHNOLOGY - FAST HEAT-UP, RELIABLE TEMPERATURE CONTROL
- STABILITY ±0.5°C
- SIMPLE, INTUITIVE CONTROLS - QUICK AND EASY TO SET TEMPERATURE
- PRACTICAL FRONT PANEL LOCK - DISABLES FRONT PANEL CONTROLS PREVENTING UNINTENTIONAL TEMPERATURE CHANGES



CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
JBA5	Water bath, digital, 5L ambient +5 to 95°C, includes base tray only	£516.00	£433.44
JBA12	Water bath, digital, 12L ambient +5 to 95°C, includes base tray only	£593.00	£498.12
JBA18	Water bath, digital, 18L ambient +5 to 95°C, includes base tray only	£689.00	£578.76

LIDS FOR SAP, JBN, JBA, SBB BATHS

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
STAINLESS STEEL GABLED LIDS			
LU6	Lid gabled stainless steel for SAP5, JBN5, JBA5 and SBB AQUA 5 PLUS.	£184.00	£154.56
LU14	Lid gabled stainless steel for SAP12, JBN12, JBA12, and SBB AQUA 12 PLUS.	£200.00	£168.00
LU28	Lid gabled stainless steel for SAP18/26, JBN18/26, JBA18 and SBB AQUA 26 PLUS.	£230.00	£193.20

STAINLESS STEEL FLAT LIDS WITH RING SETS FOR SAP, JBN & JBA BATHS

LF6	Lid flat stainless steel with ring set for SAP5, JBN5, JBA5 and SBB AQUA 5 PLUS. 2 ring sets	£222.00	£186.48
LF14	Lid flat stainless steel with ring set for SAP12, JBN12, JBA12 and SBB AQUA 12 PLUS. 4 ring sets	£341.00	£286.44
LF28	Lid flat stainless steel with ring set for SAP18/26, JBN18/26, JBA18 and SBB AQUA 26 PLUS. 6 ring sets	£480.00	£403.20



GRANT STIRRED WATER BATHS



THE GRANT T100 SERIES OF GENERAL PURPOSE HEATING CIRCULATING BATHS OFFER ACCURATE TEMPERATURE CONTROL FOR STRAIGHTFORWARD APPLICATIONS IN THE RANGE OF 0°C TO 100°C.

APPLICATIONS

- CLINICAL, MICROBIOLOGY AND PATHOLOGY LABS – MEDIA TEMPERING, THAWING AND INCUBATING SAMPLES
- TEACHING LABS, HIGHER EDUCATION/UNIVERSITIES – PRACTICAL DEMONSTRATION/EXPERIMENTATION, SAMPLE PREPARATION

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
T100-ST5	5l stainless steel tank & immersion thermostat, ambient +15°C to 100°C	£1,293.00	£1,086.12
T100-ST12	12l stainless steel tank & immersion thermostat, 0°C to 100°C	£1,590.00	£1,335.60
T100-ST18	18l stainless steel tank & immersion thermostat, 0°C to 100°C	£1,623.00	£1,363.62
T100-ST26	26l stainless steel tank & immersion thermostat, 0°C to 100°C	£1,701.00	£1,428.84
T100-ST38	38l stainless steel tank & immersion thermostat, 0°C to 100°C	£2,266.00	£1903.44

LIDS

CODE	DESCRIPTION	LIST PRICE	OFFER PRICE
STL5	Flat stainless steel lid for 5L bath	£88.00	£73.92
STL12	Gabled, hinged (removable) stainless steel lid for 12L bath	£310.00	£260.40
STL26	Gabled, hinged (removable) stainless steel lid for 18L and 26L bath	£378.00	£317.92
STL38	Gabled, hinged (removable) stainless steel lid for 38L bath	£565.00	£474.60



COLE-PARMER® ESSENTIALS PIPETTE REAGENT RESERVOIR



- STERILE: YES
- MAX VOLUME (ML): 50
- NUMBER OF CHANNELS: 1
- WETTED MATERIAL: POLYSTYRENE
- COLOUR CODE: WHITE
- PACKAGE OF: INDIVIDUALLY WRAPPED 50/CASE
- COMPATIBLE WITH: 8- AND 12-CHANNEL PIPETTORS
- DESCRIPTION: REAGENT RESERVOIR, STERILE, 1-CHANNEL, 50 ML, INDIVIDUALLY WRAPPED; 50/CS

£60.27
per pack

CODE	04395-26	04395-29
DESCRIPTION	Cole-Parmer Reagent Reservoir, Sterile, 1-Channel, 50 mL, Individually Wrapped; 50/Cs	Cole-Parmer Reagent Reservoir, Sterile, 1-Channel, 100 mL, Individually Wrapped; 50/Cs

CARDBOARD FREEZER BOXES

- SNUG-FITTING LID PREVENTS ACCIDENTAL SEPARATION IF DROPPED
- USE IN MECHANICAL FREEZERS OR LIQUID NITROGEN TANKS TO -196°C
- CONSTRUCTED OF STURDY CARDBOARD WITH MOISTURE-REPELLENT COATING
- DESIGNED WITH SMALL HOLES ON THE BOTTOM TO ALLOW FOR DRAINAGE
- KEEP YOUR SAMPLES ORGANIZED WITH NUMERICAL GRID



CODE	06755-72	06755-73
DESCRIPTION	Cole-Parmer Cardboard Freezer Box, 5-1/4" x 5-1/4" x 2", with 81-Place Divider	Cole-Parmer Cardboard Freezer Box, 5-1/4" x 5-1/4" x 2", with 100-Place Divider
OFFER PRICE	£4.38	£4.38

PDI SANI CLOTH 70% ALCOHOL WIPES -PACK OF 200



ONLY
£4.99



View Our Website



VIRKON DISINFECTANT

PRODUCT CODE	DESCRIPTION	PACK SIZE	OFFER PRICE
1127956	Virkon Disinfectant Tablets packed in Tubs of 10x5g Tablets, Case of 20 Tubs	20	£254.88



PRODUCT CODE	DESCRIPTION	PACK SIZE	OFFER PRICE
1127955	Virkon Disinfectant Powder 5kg Drum	5kg	£204.91



PRODUCT CODE	DESCRIPTION	PACK SIZE	OFFER PRICE
1203133	Virkon Disinfectant Powder Sachets Pack of 50 x 50g sachets	50x50g	£304.89



Terms and Conditions:

AT PREMIER SCIENTIFIC WE ARE COMMITTED TO DELIVERING HIGH-QUALITY LABORATORY SUPPLIES TO MEET YOUR SPECIFIC NEEDS. PLEASE NOTE THAT CARRAIGE MAY APPLY TO YOUR ORDERS AND THESE SHALL BE DETAILED WITHIN YOUR INVOICE.

WE UNDERSTAND THAT ALL LABORATORIES ARE UNIQUE AND WE ARE HERE TO ASSIST YOU IN FINDING THE MOST TAILORED SOLUTIONS.

FOR A DIRECT BESPOKE QUOTE, OR TO DISCUSS ANY QUERIES PLEASE DO NOT HESITATE TO CONTACT US AT SALES@PREMIERSCIENTIFIC.CO.UK.


ALL PRICES EXCLUDE VAT

ALL ORDERS SUBJECT TO AVAILABILITY

WE VALUE YOUR BUSINESS AND LOOK FORWARD TO SERVING YOU WITH EXCELLENCE.



View Our Website



YOUR SOURCE FOR LABORATORY EQUIPMENT AND AND CHEMICALS



LEADING THE WAY WITH PREMIER BRANDS