

PRODUCT CARD

CHAIR GEMINI SPECIAL CHL BLACK



High swivel laboratory chair. Non-slip seat made of soft polyurethane in black (RAL 9005). Other colors available upon request. Seat height adjustment with a pneumatic lift, up and down. Chrome-plated metal base with glides.

Seat/backrest material:
Non-slip, easily washable polyurethane, resistant to abrasion and cleaning agents.

Seat height adjustment range:
50 cm - 80 cm.

Seat dimensions:
Width 32 cm to 45 cm - depth 44 cm.

Backrest dimensions:
Width 41 cm - height 31 cm.

Wheels/glides:
Standard glides / wheels available upon request.

Footrest:
Chrome-plated, adjustable, diameter 45 cm.

Base:
Chrome-plated metal. Diameter 64 cm.

Box:
60 x 58 x 23 cm.

Requires self-assembly.

Manufacturer:
PLASTPUR Sp. z o.o. - Made in Poland.