



CONTAINERS - CLINICAL WASTE



Looking for a reliable and safe way to dispose of clinical waste in your lab or healthcare facility? Look no further than these Clinical Waste Containers, These containers are designed to meet the highest standards of safety and compliance.

They are made from high-quality materials and are puncture-resistant, ensuring the safe disposal of hazardous waste. These Clinical Waste Containers come with a convenient handle and a snap-lock lid, making them easy to transport and dispose of.

Product	Capacity (L)	Dimensions (mm)
WD380-12	30	355 x 400 x 315
WD380-22	60	336 x 406 x 319
MDSC30MONO	30	400 x 400 x 390
MDSPK601	60	400 x 400 x 575



PREMIER SCIENTIFIC
in partnership with

medline
scientific

A CALIBRE SCIENTIFIC COMPANY



WWW.PREMIERSCIENTIFIC.CO.UK
SALES@PREMIERSCIENTIFIC.CO.UK
TEL: 028 9077 3030
Unit 3, 5 West Bank Drive, Belfast BT3 9LA