

WHEATON® PLASTIC PASTEUR PIPETTES

Ideal for the accurate volume transfer of liquid
and bacteriological samples.




DWK LIFE SCIENCES
DURAN WHEATON KIMBLE

WHEATON® SINGLE-USE POLYETHYLENE PASTEUR PIPETTES ARE IDEAL FOR GENERAL LIQUID TRANSFER APPLICATIONS AND EXTRACTING SAMPLES FOR ANALYSIS.

WHEATON® polyethylene pasteur pipettes are an everyday laboratory essential. Available in a range of sizes from 1ml to jumbo-sized pipettes with fine tip variants, there is a product to suit most user applications for use in areas such as blood banking, bacteriology, general life science and industry. The integral bulb design and single-use feature eliminate the need for plugging or autoclaving, helping to reduce the potential of cross-contamination during everyday sample handling.

PASTEUR PIPETTE FEATURES

- Manufactured from soft, flexible, inert non-toxic polyethylene
- Ideal for safe and easy transfer of liquid samples
- Clear graduations on the 1ml and 3ml versions support ease of volume control
- Suitable for use as a freeze vial when the bulb is fully heat sealed
- Jumbo pipettes are ideal for larger volume transfer or accessing liquids in taller vessels
- Easy to handle with gloved hands
- Non-toxic material can be safely incinerated for disposal
- Available as both sterile and non-sterile options in easy-to-open, paper-plastic peel pouches



ORDERING INFORMATION

Code	Description	Sterility	Inner Case Qty.	Outer Case Qty.
711115	1ml Graduated Pasteur Pipette, Bulk	NS	500	3000
726129	1ml Graduated Pasteur Pipette, Individual Wrap	IRR	1	500
711117S/10	1ml Graduated Pasteur Pipette, Pack of 10	IRR	10	1500
711117S/20	1ml Graduated Pasteur Pipette, Pack of 20	IRR	20	1500
711116	3ml Graduated Pasteur Pipette, Bulk	NS	500	3000
726128	3ml Graduated Pasteur Pipette, Individual Wrap	IRR	1	500
711118S/10	3ml Graduated Pasteur Pipette, Pack of 10	IRR	10	1500
711118S/20	3ml Graduated Pasteur Pipette, Pack of 20	IRR	20	1500
711117F	1ml Fine-Tip Pasteur Pipette, Bulk	NS	500	3000
711117F1S	1ml Fine-Tip Pasteur Pipette, Individual Wrap	IRR	1	500
711117F10S	1ml Fine-Tip Pasteur Pipette, Pack of 10	IRR	10	1500
711117F20S	1ml Fine-Tip Pasteur Pipette, Pack of 20	IRR	20	1500
711230	Jumbo Pasteur Pipette	NS	100	600
711117MF/3000	Fine-Tip Micro-Pasteur Pipette, Bulk	NS	500	3000

WH.PPIP-ENG-V5-02.07.2020



DWK Life Sciences Limited
Unit 4, Riverside 2, Campbell Road,
Stoke-on-Trent, ST4 4RJ,
UK

Phone: +44 (0)1782 444406
Fax: +44 (0)1782 940436
enquiries.uk@dwk.com
www.dwkLtd.com